




# JAY<sup>™</sup>

## JAY CUSHIONS AND BACKS PARTS MANUAL

JANUARY 2017

Sunrise Medical appreciates your business. For your ordering convenience, our Customer Service Representatives can be contacted at:

**To place orders:  
(800)333-4000**

**Quickie Technical Assistance:**

You can reach our Technical Support Team by choosing prompt 1 followed by the # sign, then prompt 2 for Technical Support. Choose option 1 for Quickie Technical Support.

**For Assistance with a Web Order,  
or to sign up for [sunrisemedical.com](http://sunrisemedical.com)  
(800) 333-4000 Option 5**

**FOR YOUR CONVENIENCE YOU MAY PLACE ORDERS THROUGH  
OUR WEBSITE OR BY FAX 24 HOURS A DAY:**

**[www.sunrisemedical.com](http://www.sunrisemedical.com)**

**Fax: (800)300-7502**

For easy reference please fill in your Customer Number and Sunrise Account Manager's name.

Customer Account Number \_\_\_\_\_  
Ship to Account Number \_\_\_\_\_  
Sunrise Account Manager \_\_\_\_\_

\*Part prices and configurations or specifications of products included in this Parts Manual are subject to change without prior notice.

|                                                        |    |
|--------------------------------------------------------|----|
| JAY CUSHIONS & BACKS .....                             | 6  |
| JAY BACKS .....                                        | 6  |
| JAY BASIC BACK .....                                   | 6  |
| JAY BASIC BACK CLAMP STYLE .....                       | 6  |
| JAY BASIC BACK STRAP STYLE .....                       | 8  |
| JAY BASIC BACK Discontinued .....                      | 9  |
| JAY GO BACK .....                                      | 11 |
| JAY CARE BACK .....                                    | 12 |
| JAY FOCUS POINT BACK .....                             | 13 |
| JAY J2 BACK .....                                      | 14 |
| JAY J2 PLUS .....                                      | 14 |
| JAY J2 PLUS BACK COVER .....                           | 15 |
| JAY J2 PLUS BACK PLATE ASSM .....                      | 16 |
| HEADREST .....                                         | 18 |
| JAY J2 PLUS BACK SIDE PLATE .....                      | 19 |
| JAY J2 PLUS BACK SWING-AWAY LATERALS .....             | 21 |
| JAY J2 PLUS BACK FIXED LATERALS .....                  | 22 |
| JAY J2 BACK .....                                      | 23 |
| JAY J2 DEEP CONTOUR BACK .....                         | 26 |
| JAY J2 & J2 DEEP CONTOUR MOUNTING HARDWARE GUIDE ..... | 28 |
| JAY J2 & J2 DEEP CONTOUR MOUNTING HARDWARE CHART ..... | 30 |
| JAY J2 & J2 DEEP CONTOUR CLAMP STYLE HARDWARE .....    | 31 |
| HEADREST .....                                         | 32 |
| UNIVERSAL HEADREST MOUNTING KIT .....                  | 33 |
| JAY J3 BACK .....                                      | 33 |
| J3 LT BACKS .....                                      | 33 |
| J3 LT PA BACKREST .....                                | 34 |
| J3 LT PDC BACKREST .....                               | 36 |
| J3 CARBON LT PA BACK .....                             | 37 |
| J3 MT BACKS .....                                      | 38 |
| J3 MT PA BACKREST .....                                | 39 |
| J3 HD MT PA BACKREST .....                             | 42 |
| J3 MT PL BACKREST .....                                | 44 |
| J3 MT PD BACK .....                                    | 47 |
| J3 MT PDC BACKREST .....                               | 49 |
| J3 CARBON MT PA BACK .....                             | 51 |
| J3 UT BACKS .....                                      | 51 |
| J3 UT PA BACKREST .....                                | 52 |
| J3 HD UT PA BACKREST .....                             | 55 |
| J3 UT PL BACKREST .....                                | 56 |
| J3 UT PDL BACKREST .....                               | 59 |
| J3 UT PD BACK .....                                    | 61 |
| J3 UT PDC BACKREST .....                               | 63 |
| J3 SH BACKS .....                                      | 64 |
| J3 SH PA BACKREST .....                                | 65 |
| J3 SH PL BACKREST .....                                | 67 |
| J3 SH PDL BACKREST .....                               | 69 |
| J3 HD SH PA BACKREST .....                             | 71 |
| J3 SH PD BACK .....                                    | 72 |
| J3 SH PDC BACKREST .....                               | 74 |
| MOUNTING HARDWARE .....                                | 75 |
| J3 BACK HARDWARE .....                                 | 76 |
| J3 4 POINT BACK HARDWARE .....                         | 78 |
| J3 BACK LITE MOUNT HARDWARE .....                      | 79 |
| J3 RH & RHF BACK HARDWARE .....                        | 80 |
| J3 TH BACK HARDWARE .....                              | 81 |
| LATERALS .....                                         | 81 |
| CRADLE-U LATERALS .....                                | 82 |
| J3 NEW SA LATERAL .....                                | 83 |
| J3 SWING-AWAY LATERAL .....                            | 85 |
| J3 FIXED LATERAL .....                                 | 86 |
| HARNESSES .....                                        | 86 |
| CENTER OPENING HARNESS .....                           | 87 |
| TRUNK STRAP PAD .....                                  | 88 |
| CLASSIC HARNESS .....                                  | 89 |
| CONTOUR HARNESS .....                                  | 90 |
| J3 HARNESS ATTACHMENT HARDWARE .....                   | 91 |
| HEADRESTS .....                                        | 91 |
| AXYS HEADREST .....                                    | 92 |
| COBRA HEADREST .....                                   | 93 |

|                                                  |     |
|--------------------------------------------------|-----|
| LINX HEADREST .....                              | 94  |
| ONYX HEADREST .....                              | 95  |
| UNIVERSAL HEADREST MOUNT .....                   | 97  |
| ADJUST-A-PLUSH HEADREST PAD .....                | 98  |
| CONTOUR CRADLE HEADREST PAD .....                | 99  |
| PLUSH HEADREST PAD .....                         | 101 |
| HEADS UP .....                                   | 103 |
| VANITY FLAP .....                                | 106 |
| J3 TOOLS .....                                   | 107 |
| JAY J3 BACK SPINE ALIGN KIT .....                | 108 |
| JAY ZIP BACK .....                               | 108 |
| JAY ZIP BACK .....                               | 109 |
| JAY ZIP BACK HARDWARE .....                      | 111 |
| JAY ZIP SWING AWAY LATERALS .....                | 112 |
| JAY ZIP FIXED LATERALS .....                     | 113 |
| JAY ZIP BACK POSITIONING COMPONENTS .....        | 114 |
| JAY ZIP BACK HARNESS ASSEMBLIES .....            | 115 |
| VANITY FLAP .....                                | 116 |
| JAY ZIP BACK HEADREST .....                      | 117 |
| JAY ZIP BACK TOOLS .....                         | 118 |
| JAY LUMBAR SUPPORT .....                         | 119 |
| >>>> DISCONTINUED JAY BACKS .....                | 119 |
| JAY ACTIVE BACK .....                            | 120 |
| JAY XTREME BACK .....                            | 123 |
| JAY MODULAR BACK .....                           | 125 |
| JAY MODULAR BACK (DEFORMITY) .....               | 126 |
| JAY MODULAR BACK (NON-DEFORMITY) .....           | 128 |
| JAY MODULAR BACK BOX .....                       | 130 |
| JAY PRECISION BACK .....                         | 130 |
| JAY PRECISION BACKREST .....                     | 131 |
| FOLDING STABILIZER BAR .....                     | 134 |
| JAY CUSHIONS .....                               | 135 |
| JAY BASIC CUSHION .....                          | 136 |
| JAY BASIC PRO CUSHION .....                      | 138 |
| JAY CARE CUSHION .....                           | 140 |
| JAY SOFT COMBI P CUSHION .....                   | 141 |
| JAY DUO CUSHION .....                            | 143 |
| JAY EASY CUSHION .....                           | 145 |
| JAY FUSION .....                                 | 146 |
| JAY FUSION EVALUATION KIT .....                  | 147 |
| JAY FUSION FLUID .....                           | 148 |
| JAY FUSION AIR .....                             | 154 |
| JAY FUSION AIR - REDUCED PROFILE .....           | 160 |
| JAY FUSION FLUID - REDUCED PROFILE .....         | 166 |
| JAY FUSION AIR POSITIONING .....                 | 172 |
| JAY FUSION FLUID POSITIONING .....               | 178 |
| JAY FUSION CORE .....                            | 184 |
| JAY GO CUSHION .....                             | 188 |
| JAY GS CUSHION .....                             | 191 |
| JAY ION .....                                    | 193 |
| JAY J2 CUSHION .....                             | 196 |
| JAY J2 CUSHION .....                             | 197 |
| JAY J2 DEEP CONTOUR CUSHION .....                | 199 |
| JAY J2 PLUS CUSHION .....                        | 204 |
| JAY J3 CUSHION .....                             | 208 |
| JAY J3 STANDARD CUSHION .....                    | 209 |
| JAY J3 POSITIONING CUSHION .....                 | 213 |
| JAY J3 CUSHION DEMO KIT .....                    | 217 |
| JAY J3 POSITIONING CUSHION DEMO KIT .....        | 218 |
| JAY J3 STAMPED PART NUMBER CROSS REFERENCE ..... | 220 |
| JAY LITE CUSHION .....                           | 225 |
| JAY UNION .....                                  | 229 |
| JAY X2 CUSHION .....                             | 234 |
| JAY ZIP .....                                    | 236 |
| JAY CUSTOM DROPSEAT OPTION PARTS .....           | 239 |
| >>>> DISCONTINUED JAY CUSHIONS .....             | 240 |
| JAY ACTIVE CUSHION .....                         | 241 |
| JAY CUSHION .....                                | 243 |
| JAY TRIAD CUSHION .....                          | 245 |
| JAY TRIAD CUSHION .....                          | 246 |

|                                             |     |
|---------------------------------------------|-----|
| JAY TRIAD P-SKIN & POSITIONING CUSHION..... | 249 |
| JAY XTREME CUSHION .....                    | 251 |
| JAY XTREME CUSHION .....                    | 252 |
| JAY XTREME SKIN PROTECTION CUSHION .....    | 254 |
| JAYJ-GEL.....                               | 256 |
| JAY J-GEL CUSHION DISCONTINUED LETTER ..... | 257 |
| JAY J-GEL CUSHION .....                     | 258 |
| SOFT COMBI CUSHION.....                     | 259 |
| JAY SOFT COMBI DISCONTINUED LETTER .....    | 260 |
| SOFT COMBI CUSHION.....                     | 261 |
| JAY COMBI .....                             | 262 |
| JAY COMBI DISCONTINUED LETTER .....         | 263 |
| JAY COMBI .....                             | 264 |
| JAY ACCESSORIES .....                       | 265 |
| JAY BOX .....                               | 266 |
| JAY J3 BACK SPINE ALIGN KIT .....           | 268 |
| JAY J3 CUSHION DEMO KIT .....               | 269 |
| JAY J3 POSITIONING CUSHION DEMO KIT.....    | 270 |
| JAY CUSTOM DROPSEAT OPTION PARTS.....       | 272 |
| JAY ADJUSTABLE SEAT .....                   | 274 |
| JAY PROTECTOR .....                         | 276 |
| JAY SPINAL FLUID PAD .....                  | 277 |
| JAY FULL FOOT BOX WITH DIVIDER.....         | 278 |
| JAY FULL FOOT BOX .....                     | 280 |
| JAY SINGLE FOOT BOX .....                   | 281 |
| JAY ARM TROUGH FIXED & ROTATIONAL .....     | 282 |
| JAY LUMBAR SUPPORT.....                     | 283 |
| back cover .....                            | 283 |



| Pos. | Item Number | Description                 | Remarks |
|------|-------------|-----------------------------|---------|
| 2    | B3300       | JBASIC BACK FOAM BASE,20x20 | 20x20   |
| 2    | B3306       | JBASIC BACK FOAM BASE,20x16 | 20x16   |
| 2    | B3308       | JBASIC BACK FOAM BASE,20x18 | 20x18   |
| 2    | B3340       | JBASIC BACK FOAM BASE,14x20 | 14x20   |
| 2    | B3346       | JBASIC BACK FOAM BASE,14x16 | 14x16   |
| 2    | B3348       | JBASIC BACK FOAM BASE,14x18 | 14x18   |
| 2    | B3360       | JBASIC BACK FOAM BASE,16x20 | 16x20   |
| 2    | B3366       | JBASIC BACK FOAM BASE,16x16 | 16x16   |
| 2    | B3368       | JBASIC BACK FOAM BASE,16x18 | 16x18   |
| 2    | B3380       | JBASIC BACK FOAM BASE,18x20 | 18x20   |
| 2    | B3386       | JBASIC BACK FOAM BASE,18x16 | 18x16   |
| 2    | B3388       | JBASIC BACK FOAM BASE,18x18 | 18x18   |
| 3    | C3300FC     | JBASIC BACK COVER,20x20     | 20x20   |
| 3    | C3306FC     | JBASIC BACK COVER,20x16     | 20x16   |
| 3    | C3308FC     | JBASIC BACK COVER,20x18     | 20x18   |
| 3    | C3340FC     | JBASIC BACK COVER,14x20     | 14x20   |
| 3    | C3346FC     | JBASIC BACK COVER,14x16     | 14x16   |
| 3    | C3348FC     | JBASIC BACK COVER,14x18     | 14x18   |

| Pos. | Item Number | Description                  | Remarks                    |
|------|-------------|------------------------------|----------------------------|
| 3    | C3360FC     | JBASIC BACK COVER,16x20      | 16x20                      |
| 3    | C3366FC     | JBASIC BACK COVER,16x16      | 16x16                      |
| 3    | C3368FC     | JBASIC BACK COVER,16x18      | 16x18                      |
| 3    | C3380FC     | JBASIC BACK COVER,18x20      | 18x20                      |
| 3    | C3386FC     | JBASIC BACK COVER,18x16      | 18x16                      |
| 3    | C3388FC     | JBASIC BACK COVER,18x18      | 18x18                      |
| 4    | JBHDWUL     | JBASIC UPPER CLAMP LEFT      | COMPLETE UPPER LEFT CLAMP  |
| 4    | JBHDWUR     | JBASIC UPPER CLAMP RIGHT     | COMPLETE LOWER RIGHT CLAMP |
| 6    | JBHDWUB     | JBASIC UPPER CLAMP BOLTS     | ALL BOLTS, LEFT OR RIGHT   |
| 4,5  | JBHDWH      | JBASIC CLAMP HARDWARE (PAIR) | COMPLETE CLAMP HDW, L & R  |
| 5-7  | 101575S     | ASSM MNT HRDW RH BASIC BK    | Right                      |
| 5-7  | 101576S     | ASSM MNT HRDW LH BASIC BK    | Left                       |

**Note:** Effective August 25 2008

**Note:** Item 1 Back Shell not available as a replacement part; please order complete back from Jay Order Form if replacement shell is needed.

**Note:** effective August 25 2008


**Note:** Item 1 Back Shell not available as a replacement part; please order complete back from Jay Order Form if replacement shell is needed.



| Pos. | Item Number | Description                 | Remarks | Pos. | Item Number | Description              | Remarks |
|------|-------------|-----------------------------|---------|------|-------------|--------------------------|---------|
| 2    | B3300       | JBASIC BACK FOAM BASE,20x20 | 20x20   | 3    | C3306FS     | JBASIC BACK COVER,20x16  | 20x16   |
| 2    | B3306       | JBASIC BACK FOAM BASE,20x16 | 20x16   | 3    | C3308FS     | JBASIC BACK COVER,20x18  | 20x18   |
| 2    | B3308       | JBASIC BACK FOAM BASE,20x18 | 20x18   | 3    | C3340FS     | JBASIC BACK COVER,14x20  | 14x20   |
| 2    | B3340       | JBASIC BACK FOAM BASE,14x20 | 14x20   | 3    | C3346FS     | JBASIC BACK COVER,14x16  | 14x16   |
| 2    | B3346       | JBASIC BACK FOAM BASE,14x16 | 14x16   | 3    | C3348FS     | JBASIC BACK COVER,14x18  | 14x18   |
| 2    | B3348       | JBASIC BACK FOAM BASE,14x18 | 14x18   | 3    | C3360FS     | JBASIC BACK COVER,16x20  | 16x20   |
| 2    | B3360       | JBASIC BACK FOAM BASE,16x20 | 16x20   | 3    | C3366FS     | JBASIC BACK COVER,16x16  | 16x16   |
| 2    | B3366       | JBASIC BACK FOAM BASE,16x16 | 16x16   | 3    | C3368FS     | JBASIC BACK COVER,16x18  | 16x18   |
| 2    | B3368       | JBASIC BACK FOAM BASE,16x18 | 16x18   | 3    | C3380FS     | JBASIC BACK COVER,18x20  | 18x20   |
| 2    | B3380       | JBASIC BACK FOAM BASE,18x20 | 18x20   | 3    | C3386FS     | JBASIC BACK COVER,18x16  | 18x16   |
| 2    | B3386       | JBASIC BACK FOAM BASE,18x16 | 18x16   | 3    | C3388FS     | JBASIC BACK COVER,18x18  | 18x18   |
| 2    | B3388       | JBASIC BACK FOAM BASE,18x18 | 18x18   | 4    | JBHDWSC     | JBASIC BACK STRAP BUCKLE |         |
| 3    | C3300FS     | JBASIC BACK COVER,20x20     | 20x20   |      |             |                          |         |

**Note:** parts only; discontinued August 25, 2008


**Note:** Item 3 Back Shell not available as a replacement part; please order complete back from Jay Order Form if replacement shell is needed.




| Pos. | Item Number | Description              | Pos. | Item Number | Description              |
|------|-------------|--------------------------|------|-------------|--------------------------|
| 1    | C3346       | COVER BASIC BACK 14 X 16 | 1    | C3300       | COVER BASIC BACK 20 X 20 |
| 1    | C3348       | COVER BASIC BACK 14 X 18 | 2    | B3346       | FOAM BASIC BACK 14 X 16  |
| 1    | C3340       | COVER BASIC BACK 14 X 20 | 2    | B3348       | FOAM BASIC BACK 14 X 18  |
| 1    | C3366       | COVER BASIC BACK 16 X 16 | 2    | B3340       | FOAM BASIC BACK 14 X 20  |
| 1    | C3368       | COVER BASIC BACK 16 X 18 | 2    | B3366       | FOAM BASIC BACK 16 X 16  |
| 1    | C3360       | COVER BASIC BACK 16 X 20 | 2    | B3368       | FOAM BASIC BACK 16 X 18  |
| 1    | C3386       | COVER BASIC BACK 18 X 16 | 2    | B3360       | FOAM BASIC BACK 16 X 20  |
| 1    | C3388       | COVER BASIC BACK 18 X 18 | 2    | B3386       | FOAM BASIC BACK 18 X 16  |
| 1    | C3380       | COVER BASIC BACK 18 X 20 | 2    | B3388       | FOAM BASIC BACK 18 X 18  |
| 1    | C3306       | COVER BASIC BACK 20 X 16 | 2    | B3380       | FOAM BASIC BACK 18 X 20  |
| 1    | C3308       | COVER BASIC BACK 20 X 18 | 2    | B3306       | FOAM BASIC BACK 20 X 16  |

| Pos. | Item Number | Description               | Pos.         | Item Number | Description                   |
|------|-------------|---------------------------|--------------|-------------|-------------------------------|
| 2    | B3308       | FOAM BASIC BACK 20 X 18   | 6, 7         | M-250C      | KIT CLAMP UNIVERSAL RETRO-FIT |
| 2    | B3300       | FOAM BASIC BACK 20 X 20   | 8, 9, 9A, 10 | M-250S      | KIT STRAP UNIVERSAL RETRO-FIT |
| 4, 5 | 100197      | BASIC BACK MOUNT HARDWARE |              |             |                               |

Note: Individual replacement parts are no longer available for 12, 16, 20 wide backs. New Jay Go back must be ordered for those sizes.




| Pos. | Item Number | Description                  |
|------|-------------|------------------------------|
| 1    | SJGOB14S    | JAY GO BCK 14W SHT STRAP ASM |
| 1    | SJGOB18S    | JAY GO BCK 18W SHT STRAP ASM |
| 1    | SJGOB22S    | JAY GO BCK 22W SHT STRAP ASM |
| 1    | SJGOB14T    | JAY GO BCK 14W TAL STRAP ASM |
| 1    | SJGOB18T    | JAY GO BCK 18W TAL STRAP ASM |
| 1    | SJGOB22T    | JAY GO BCK 22W TAL STRAP ASM |
| 2,3  | CJGOB14S    | JAY GO BACK 14W SHORT COVER  |
| 2,3  | CJGOB18S    | JAY GO BACK 18W SHORT COVER  |
| 2,3  | CJGOB22S    | JAY GO BACK 22W SHORT COVER  |
| 2,3  | CJGOB14T    | JAY GO BACK 14W TALL COVER   |
| 2,3  | CJGOB18T    | JAY GO BACK 18W TALL COVER   |
| 2,3  | CJGOB22T    | JAY GO BACK 22W TALL COVER   |
| 2,3  | CJGOB14SM   | JGO BK 14W SHORT MICRO COVER |
| 2,3  | CJGOB14TM   | JGO BK 14W TALL MICRO COVER  |
| 2,3  | CJGOB18SM   | JGO BK 18W SHORT MICRO COVER |
| 2,3  | CJGOB18TM   | JGO BK 18W TALL MICRO COVER  |
| 2,3  | CJGOB22SM   | JGO BK 22W SHORT MICRO COVER |
| 2,3  | CJGOB22TM   | JGO BK 22W TALL MICRO COVER  |




| Pos. | Item Number | Description                         | Remarks                            |
|------|-------------|-------------------------------------|------------------------------------|
| 1    | C3500B      | JAYCARE COVER BACK REGULAR AD BLK   |                                    |
| 1    | C3502B      | JAYCARE COVER BACK 20 BLK           |                                    |
| 1    | C3505B      | JAYCARE COVER BACK 16 BLK           |                                    |
| 2    | B3500       | JAYCARE BASE BACK 18                |                                    |
| 2    | B3502       | JAYCARE BASE BACK 20                |                                    |
| 2    | B3505       | JAYCARE BASE BACK 16                |                                    |
| 3    | S3500       | JAYCARE SHELL BACK 18               |                                    |
| 3    | S3502       | JAYCARE SHELL BACK 20               |                                    |
| 3    | S3505       | JAYCARE SHELL BACK 16               |                                    |
| 4A   | M-300       | STANDARD HARDWARE KIT-JAYCARE SHELL |                                    |
| 4B   | M-301       | EXTERIOR HARDWARE KIT-JAYCARE SHELL |                                    |
| 5    | 680C        | BOTTOM UNIVERSAL POST BRACKETS      |                                    |
| 6    | 681C        | TOP BRACKET ASSM - "QUARTER-TURN"   | fits 5/8" to 1.0" diameter tubing  |
| 6A   | 2681C       | TOP LONG 1/4 BRACKET ASSM           | fits 1.0" to 1.25" diameter tubing |
| 7    | 880997      | LONG STRAP KIT                      | USED ON 1 1/8" DIAMETER BACK POST  |
| 8    | M-500       | HEADREST ADAPTOR KIT                |                                    |

Note: Parts only, for complete Backrest see Jay Focus Point Back order form.



| Pos. | Item Number | Description                 | Pos. | Item Number | Description                 |
|------|-------------|-----------------------------|------|-------------|-----------------------------|
| 1    | FP-LS-1415R | FP RT LATERAL SHELL 14-15   | 7    | FP-NPLT-LNG | FP NUTPLATE LONG            |
| 1    | FP-LS-1617R | FP RT LATERAL SHELL 16-17   | 8    | 106943      | 1/4 AN WASHER FLAT          |
| 1    | FP-LS-1819R | FP RT LATERAL SHELL 18-19   | 9    | 106958      | 1/4-20 X 1 HEX C/S          |
| 1    | FP-LS-2021R | FP RT LATERAL SHELL 20-21   | 10   | FP-MBRKT    | FP MAIN BRKT                |
| 2    | FP-HPLT-ML  | FP HINGE PLT MALE           | 11   | FP-HANG-V3  | FP HANGER V3                |
| 3    | FP-NPLT-SHT | FP NUTPLATE SHORT           | 12   | M-455036-QU | 1/4-28 X 9/16 SPEC ROUND HD |
| 4    | FP-CS-1415T | FP CENTER SHELL 14-15 TALL  | 13   | FP-LS-1415L | FP LFT LATERAL SHELL 14-15  |
| 4    | FP-CS-1415  | FP CENTER SHELL 14-15       | 13   | FP-LS-1617L | FP LFT LATERAL SHELL 16-17  |
| 4    | FP-CS-1617T | FP CENTER SHELL 16-17 TALL  | 13   | FP-LS-1819L | FP LFT LATERAL SHELL 18-19  |
| 4    | FP-CS-1617  | FP CENTER SHELL 16-17       | 13   | FP-LS-2021L | FP LFT LATERAL SHELL 20-21  |
| 4    | FP-CS-1819T | FP CENTER SHELL 18-19 TALL  | 14   | FP-CF-1415T | FP BCK FOAM CENTER 1415T    |
| 4    | FP-CS-1819  | FP CENTER SHELL 18-19       | 14   | FP-CF-1415  | FP BCK FOAM CENTER 1415     |
| 4    | FP-CS-2021T | FP CENTER SHELL 20-21 TALL  | 14   | FP-CF-1617T | FP BCK FOAM CENTER 1617T    |
| 4    | FP-CS-2021  | FP CENTER SHELL 20-21       | 14   | FP-CF-1617  | FP BCK FOAM CENTER 1617     |
| 5    | FP-HPLT-FL  | FP HINGE PLT FEMALE         | 14   | FP-CF-1819T | FP BCK FOAM CENTER 1819T    |
| 6    | 106935      | NUT NYLOCK THIN 1/4-28 ROHS | 14   | FP-CF-1819  | FP BCK FOAM CENTER 1819     |


| Pos. | Item Number  | Description               | Pos. | Item Number   | Description                |
|------|--------------|---------------------------|------|---------------|----------------------------|
| 14   | FP-CF-2021T  | FP BCK FOAM CENTER 2021T  | 17   | FP-LCVR-1415R | FP LATERAL COVER RT 14-15  |
| 14   | FP-CF-2021   | FP BCK FOAM CENTER 2021   | 17   | FP-LCVR-1617R | FP LATERAL COVER RT 16-17  |
| 15   | FP-CVR-1415T | FP BCK COVER CENTER 1415T | 17   | FP-LCVR-1819R | FP LATERAL COVER RT 18-19  |
| 15   | FP-CVR-1415  | FP BCK COVER CENTER 1415  | 17   | FP-LCVR-2021R | FP LATERAL COVER RT 20-21  |
| 15   | FP-CVR-1617T | FP BCK COVER CENTER 1617T | 18   | FP-LF-1415L   | FP LATERAL FOAM LFT 14-15  |
| 15   | FP-CVR-1617  | FP BCK COVER CENTER 1617  | 18   | FP-LF-1617L   | FP LATERAL FOAM LFT 16-17  |
| 15   | FP-CVR-1819T | FP BCK COVER CENTER 1819T | 18   | FP-LF-1819L   | FP LATERAL FOAM LFT 18-19  |
| 15   | FP-CVR-1819  | FP BCK COVER CENTER 1819  | 18   | FP-LF-2021L   | FP LATERAL FOAM LFT 20-21  |
| 15   | FP-CVR-2021T | FP BCK COVER CENTER 2021T | 19   | FP-LCVR-1415L | FP LATERAL COVER LFT 14-15 |
| 15   | FP-CVR-2021  | FP BCK COVER CENTER 2021  | 19   | FP-LCVR-1617L | FP LATERAL COVER LFT 16-17 |
| 16   | FP-LF-1415R  | FP LATERAL FOAM RT 14-15  | 19   | FP-LCVR-1819L | FP LATERAL COVER LFT 18-19 |
| 16   | FP-LF-1617R  | FP LATERAL FOAM RT 16-17  | 19   | FP-LCVR-2021L | FP LATERAL COVER LFT 20-21 |
| 16   | FP-LF-1819R  | FP LATERAL FOAM RT 18-19  | 20   | M-260         | KIT CLAMP FP BACK          |
| 16   | FP-LF-2021R  | FP LATERAL FOAM RT 20-21  | 21   | M-201105      | 1/4-20 NYLOCK              |



| Pos.       | Item Number | Description                    |
|------------|-------------|--------------------------------|
| 1A         | C202023A    | J2 PLUS BACK COVER 20 - 23 AIR |
| 1C         | C202326A    | J2 PLUS BACK COVER 23 - 26 AIR |
| 2A         | B202023     | J2 PLUS BACK FOAM 20 - 23      |
| 2B         | B202426     | J2 PLUS BACK FOAM 24 - 26      |
| 3          | B200007     | ASIS FOAM BARIATRIC            |
| 4A         | B200006     | SACRUM FOAM SMALL BARIATRIC    |
| 4B         | B200005     | SACRUM FOAM LARGE BARIATRIC    |
| 1A,2A,3,4A | C202023AB   | J2 PLUS COVER ASSM 20 - 23 AIR |
| 1B,2B,3,4B | C202326AB   | J2PLUS BACK CVR ASM 23-26" AIR |


Note: Parts only, for complete Backrest see J2 Plus Back orderform.

Note: Parts only, for complete Backrest see J2 Plus Back order form.




| Pos. | Item Number | Description                        |
|------|-------------|------------------------------------|
| 1    | M-200007    | FHCS .250 - 20 X 1.00 W/ PATCH BLK |
| 3    | 200010-MJ   | NUT PLATE BARIATRIC                |
| 4    | 200110-MJ   | J2 PLUS GROWTH PLATE ASSM          |
| 5    | 202023-MJ   | J2 PLUS VANITY FLAP ASSM 20 - 23   |
| 5A   | 202326-MJ   | J2 PLUS VANITY FLAP ASSM 23 - 26   |
| 6    | B200001     | GROWTH PLATE FOAM                  |
| 7    | M-200010    | BHCS .250 - 20 X .38 W/ PATCH BLK  |
| 8    | M-200043    | BHCS .250 - 20 X .38 W/ PATCH BLK  |
| 9    | M-200039    | BHCS .250 - 20 X .38 W/ PATCH BLK  |
| 10   | 200009-MJ   | UPPER SHELL CONNECTOR              |
| 12   | 200006-MJ   | LOWER SHELL CONNECTOR              |
| 13   | M-200029    | BHCS .250 - 20 X .38 W/ PATCH BLK  |
| 14   | M-200006    | J2 PLUS HINGE FRICTION ADJUST      |
| 15   | 200014-MJ   | J2 PLUS HINGE PLATE SHIM           |


| Pos.             | Item Number | Description                          |
|------------------|-------------|--------------------------------------|
| 16               | M-200036    | WASHER .25 ANSI SERIES N PLAIN BLK   |
| 17               | M-200040    | FHCS 10 - 32 X .38 W/ PATCH BLK      |
| 18               | M-200042    | BHCS .25 - 28 X 63 W/ PATCH BLK      |
| 19               | 200023-MJ   | UNIVERSAL HEADREST UPPER MOUNTING    |
| 20               | XM2025      | MNTNG PLATE/PEDIAT.HEAD REST         |
| 21               | 200022-MJ   | UNIVERSAL HEADREST LOWER MOUNTING    |
| 22               | 226018-QU   | BHCS 10 - 32 X 3/8 W/ PATCH BLK      |
| 23               | 479034      | FHCS, 250 - 28 X 1.50                |
| 9,16-18          | 200121-MJ   | J2 PLUS HINGE ASSM                   |
| 1, 13, 14        | 200102-MJ   | J2 PLUS UPPER CONNECTOR ASM          |
| 1, 14, 15        | 200100-MJ   | J2 PLUS LOWER CONNECTOR ASSM         |
| 1-2, 3-5, 6-17   | 200130H-MJ  | BCKRST ASM 20-23 J2PLUS HEADREST     |
| 1-2A, 3-5, 6-17  | 200130-MJ   | J2 PLUS BACKREST ASSM 20-23          |
| 1-2C, 3-5, 6-17  | 200131H-MJ  | J2 PLUS BACKREST ASSM 23-26 HEADREST |
| 1-2B, 3-5A, 6-21 | 200131-MJ   | J2 PLUS BACKREST ASSM 23-26          |
| 18-21            | 200206-MJ   | J2 HEADREST HARDWARE KIT             |
| 18-21, 23        | 200207-MJ   | J2 HEADREST HARDWARE KIT, WT         |
| 18-22            | 200208-MJ   | J2 HEADREST HARDWARE KIT, OT         |




| Pos.       | Item Number | Description                    | Remarks   |
|------------|-------------|--------------------------------|-----------|
| 1          | 221115      | 1/4-20 x 3/4" SOC HD C/S BLK   |           |
| 3A         | 750696      | HEADREST W/ZIPPER MEDIUM CG    | 8x4.5 pad |
| 3B         | 750697      | HEADREST W/ZIPPER LARGE CG     | 9x4.5 pad |
| 4          | 505932      | HEADREST MOUNTING TRACK        |           |
| 5          | 476304      | 1/4-28 x .31 BHCS W/PATCH      |           |
| 6          | 221010      | 10-32 x .75 SHCS BLK W/PL      |           |
| 7          | 261204      | #10 AN WASHER .063             |           |
| 8          | 331077      | HEADREST LINK 6 INCH           |           |
| 9          | 341045      | #10-24 CONE POINT SET SCREW    |           |
| 10         | 505933      | HEADREST COLLAR                |           |
| 11         | 226072      | 1/4-20 x 1/2 FH C/S BLK        |           |
| 13         | 226095      | 1/4-20 x 3/8" SOC C/S BLK ZN   |           |
| 14         | 341046      | HEADREST CLAMPING HANDLE SCREW |           |
| 15         | 340033      | 1/4-20 WELD NUTS SLAB BASE     |           |
| 4,5,13-16  | 888101      | CLAMP HEADREST RIGID BACK      |           |
| 1,2,4,6-16 | 881003      | HEADREST MOUNTING BRACKET ASSM |           |
| 1-3A,4-16  | 881907      | HEADREST ASSEMBLY MEDIUM       | 8x4.5 pad |
| 1-3B,4-16  | 881908      | HEADREST ASSEMBLY LARGE        | 9x4.5 pad |




| Pos.         | Item Number | Description                       | Remarks |
|--------------|-------------|-----------------------------------|---------|
| 2            | M-200030    | NUT, JAM NYLOCK .4375-20          |         |
| 4            | M-200034    | FHCS, .313 -18 X 1.25 W/PATCH BLK |         |
| 5            | B200003     | SIDE PLATE FOAM J2 PLUS R/H       | Right   |
| 5            | B200004     | SIDE PLATE FOAM J2 PLUS L/H       | Left    |
| 6            | 468903      | 10-32 X .50 BHCS W/PATCH          |         |
| 7            | 2702        | #10 SAE WASHER                    |         |
| 8            | 200004-MJ   | J2 PLUS PSIS PIVOT                |         |
| 9            | 200017-MJ   | PIN RETAINER L/H J2 PLUS          | Left    |
| 9            | 200016-MJ   | PIN RETAINER R/H J2 PLUS          | Right   |
| 10           | M-200035    | FHCS, .313 -18 X 1.0 W/PATCH BLK  |         |
| 11           | M-200031    | BHCS, .313 -18 X 1.0 W/PATCH BLK  |         |
| 12           | M-200003    | J2 PLUS LOCATING PIN              |         |
| 13           | 200002-MJ   | LOWER SIDE CLAMP                  |         |
| 1-13         | 200202-MJ   | SIDE PLATE ASM RH 8 DEG J2PLUS    |         |
| 1-13         | 200203-MJ   | SIDE PLATE ASM LH 8 DEG J2PLUS    |         |
| 1,2,3A, 4-13 | 200200-MJ   | SIDE PLATE ASM RH 0 DEG J2PLUS    |         |
| 1,2,3A, 4-13 | 200201-MJ   | SIDE PLATE ASM LH 0 DEG J2PLUS    |         |




| Pos. | Item Number | Description                                 | Remarks |
|------|-------------|---------------------------------------------|---------|
| 4    | 226020      | 1/4-20 X 3/8 BH C/S W/PATCH                 |         |
| 6    | XM2141S     | 1/4-28 x .31 BHSC W/ PATCH                  |         |
| 7    | 476180      | ADJUSTABLE LATERAL SWING AWAY               |         |
| 8    | 200026-MJ   | J2 PLUS S/A LATERAL BRACKET                 |         |
| 9    | 200028-MJ   | J2 PLUS LATERAL BRACKET PLATE               |         |
| 10   | M-200036    | WASHER, .25 ANSI SERIES N PLAIN BLK         |         |
| 11   | M-200037    | BHSC .250 - 20 X .63 W/ PATCH BLK           |         |
| 1-3  | 476230      | LATERAL THORACIC SUPPORT LARGE W/O HARDWARE |         |
| 1-11 | 200124-MJ   | J2 PLUS S/A LATERAL SUPPORT ASSM            | RIGHT   |
| 1-11 | 200125-MJ   | J2 PLUS S/A LATERAL SUPPORT ASSM            | LEF     |
| 1-11 | 200127-MJ   | J2 PLUS S/A LATERAL SUPPORT ASSM            | PAIR    |



| Pos.    | Item Number | Description                                 | Remarks |
|---------|-------------|---------------------------------------------|---------|
| 1       | 226018-QU   | BHCS 10 - 32 x 3/8 W/ PATCH BLK             |         |
| 2       | 200027-MJ   | J2 PLUS FIXED LATERAL BRACKET               |         |
| 3       | M-200029    | BHCS .250 - 20 X .38 W/ PATCH BLK           |         |
| 5       | M-200037    | BHCS 10 - 32 x .63 W/ PATCH BLK             |         |
| 6       | 200028-MJ   | J2 PLUS LATERAL BRACKET PLATE               |         |
| 7       | M-200036    | WASHER, .25 ANSI SERIES N PLAIN BLK         |         |
| 4, 8, 9 | 476230      | LATERAL THORACIC SUPPORT LARGE W/O HARDWARE |         |
| 1-9     | 200114-MJ   | J2 PLUS LATERAL SUPPORT ASSM                | RIGHT   |
| 1-9     | 200115-MJ   | J2 PLUS LATERAL SUPPORT ASSM                | LEFT    |
| 1-8     | 200116-MJ   | J2 PLUS LATERAL SUPPORT ASSM                | PAIR    |



4A  
4B  
4C



5



6A  
6B  
6C

| Pos. | Item Number | Description             | Remarks |
|------|-------------|-------------------------|---------|
| 1    | C2512       | COVER J2 BACK 12 X 12.5 |         |
| 1    | C2513       | COVER J2 BACK 13 X 13.2 |         |
| 1    | C2514       | COVER J2 BACK 14 X 14.3 |         |
| 1    | C2515       | COVER J2 BACK 15 X 15.3 |         |
| 1    | C2516       | COVER J2 BACK 16 X 16.3 |         |
| 1    | C2517       | COVER J2 BACK 17 X 16.3 |         |
| 1    | C2518       | COVER J2 BACK 18 X 16.3 |         |
| 1    | C2519       | COVER J2 BACK 19 X 16.8 |         |
| 1    | C2520       | COVER J2 BACK 20 X 17.3 |         |
| 1    | C2514T      | COVER J2 BACK 14 X 21   |         |
| 1    | C2515T      | COVER J2 BACK 15 X 21   |         |
| 1    | C2516T      | COVER J2 BACK 16 X 21.8 |         |
| 1    | C2517T      | COVER J2 BACK 17 X 21.8 |         |
| 1    | C2518T      | COVER J2 BACK 18 X 21.8 |         |
| 1    | C2519T      | COVER J2 BACK 19 X 22.8 |         |
| 1    | C2520T      | COVER J2 BACK 20 X 22.8 |         |
| 2    | B2512       | BACK FOAM 12 X 12.5     |         |
| 2    | B2513       | BACK FOAM 13 X 13.2     |         |
| 2    | B2514       | BACK FOAM 14 X 14.3     |         |
| 2    | B2515       | BACK FOAM 15 X 15.3     |         |
| 2    | B2516       | BACK FOAM 16 X 16.3     |         |


| Pos.  | Item Number | Description                             | Remarks                                                 |
|-------|-------------|-----------------------------------------|---------------------------------------------------------|
| 2     | B2517       | BACK FOAM 17 X 16.3                     |                                                         |
| 2     | B2518       | BACK FOAM 18 X 16.3                     |                                                         |
| 2     | B2519       | BACK FOAM 19 X 16.8                     |                                                         |
| 2     | B2520       | BACK FOAM 20 X 17.3                     |                                                         |
| 2     | B2514T      | BACK FOAM 14 X 21                       |                                                         |
| 2     | B2515T      | BACK FOAM 15 X 21                       |                                                         |
| 2     | B2516T      | BACK FOAM 16 X 21.8                     |                                                         |
| 2     | B2517T      | BACK FOAM 17 X 21.8                     |                                                         |
| 2     | B2518T      | BACK FOAM 18 X 21.8                     |                                                         |
| 2     | B2519T      | BACK FOAM 19 X 22.8                     |                                                         |
| 2     | B2520T      | BACK FOAM 20 X 22.8                     |                                                         |
| 3     | S2512K      | J2 BACK SHELL W/ UNI MOUNT 12 X 12.5    |                                                         |
| 3     | S2513K      | J2 BACK SHELL W/ UNI MOUNT 13 X 13.2    |                                                         |
| 3     | S2514K      | J2 BACK SHELL W/ UNI MOUNT 14 X 14.3    |                                                         |
| 3     | S2515K      | J2 BACK SHELL W/ UNI MOUNT 15 X 15.3    |                                                         |
| 3     | S2516K      | J2 BACK SHELL W/ UNI MOUNT 16 X 16.3    |                                                         |
| 3     | S2517K      | J2 BACK SHELL W/ UNI MOUNT 17 X 16.3    |                                                         |
| 3     | S2518K      | J2 BACK SHELL W/ UNI MOUNT 118 X 16.3   |                                                         |
| 3     | S2519K      | J2 BACK SHELL W/ UNI MOUNT 19 X 16.8    |                                                         |
| 3     | S2520K      | J2 BACK SHELL W/ UNI MOUNT 20 X 17.3    |                                                         |
| 3     | S2514TK     | J2 BACK SHELL W/ UNI MOUNT 14 X 21      |                                                         |
| 3     | S2515TK     | J2 BACK SHELL W/ UNI MOUNT 15 X 21      |                                                         |
| 3     | S2516TK     | J2 BACK SHELL W/ UNI MOUNT 16 X 21.8    |                                                         |
| 3     | S2517TK     | J2 BACK SHELL W/ UNI MOUNT 17 X 21.8    |                                                         |
| 3     | S2518TK     | J2 BACK SHELL W/ UNI MOUNT 18 X 21.8    |                                                         |
| 3     | S2519TK     | J2 BACK SHELL W/ UNI MOUNT 19 X 22.8    |                                                         |
| 3     | S2520TK     | J2 BACK SHELL W/ UNI MOUNT 20 X 22.8    |                                                         |
| 4A    | B2560SA     | SWING AWAY BRACKET ONLY                 |                                                         |
| 4B    | B2560LSA    | SWING AWAY BRACKET ONLY -LONG           |                                                         |
| 5     | 2560J       | J2 BACK LATERAL SUPPORT 12/13 SMALL     |                                                         |
| 5     | 2560-MJ     | J2 BACK LATERAL SUPPORT 4 X 4 SMALL     |                                                         |
| 5     | 2562J       | J2 BACK LATERAL SUPPORT 12/13 MEDIUM    |                                                         |
| 5     | 2562        | J2 BACK LATERAL SUPPORT 5 X 5 MEDIUM    |                                                         |
| 5     | 2562L       | J2 BACK EXT LATERAL SUPPORT MEDIUM      |                                                         |
| 5     | 2564        | J2 BACK LARGE LATERAL SUPPORT 6 X 6     |                                                         |
| 6A    | P2560S      | SMALL J2 LATERAL-PAD ONLY               | PAIR - INCLUDES SMALL PADDLE, FOAM, COVER, AND HARDWARE |
| 6A    | B662        | BASE FOR LATERAL SUPPORT SM             | EACH - INCLUDES SMALL FOAM ONLY                         |
| 6B    | P2562M      | MEDIUM J2 LATERAL-PAD ONLY              | PAIR - INCLUDES MED PADDLE, FOAM, COVER, AND HARDWARE   |
| 6B    | B672        | BASE EXTENDED LATERAL MED               | EACH - INCLUDES MED FOAM ONLY                           |
| 6C    | P2564L      | LARGE J2 LATERAL-PAD ONLY               | PAIR - INCLUDES LARGE PADDLE, FOAM, COVER, AND HARDWARE |
| 6C    | B682        | BASE EXTENDED LATERAL LARGE             | EACH - INCLUDES LARGE FOAM ONLY                         |
| 4C,6A | 2560JSA     | J2 BACK LATERAL SUPPORT SA 12/13 SMALL  |                                                         |
| 4B,6A | 2560LSA     | SMALL J2 SA LATERAL LONG                |                                                         |
| 4A,6A | 2560SA      | J2 BACK LATERAL SUPPORT SA SMALL        |                                                         |
| 4C,6B | 2562JSA     | J2 BACK LATERAL SUPPORT SA 12/13 MEDIUM |                                                         |
| 4A,6B | 2562SA      | MEDIUM SA LATERAL SUPPORT               |                                                         |
| 4B,6B | 2562LSA     | J2 BACK LATERAL LONG SA MEDIUM          |                                                         |
| 4A,6C | 2564SA      | J2 BACK LATERAL SA LARGE                |                                                         |
| 4B,6C | 2564LSA     | LARGE J2 SA LATERAL LONG                |                                                         |
|       | M-251       | MNTING ASSEMB KIT J2 BK LAT.            | Lateral hardware kit- 4 screws, washers & nuts          |

**NOTE:** SEE JAY J2 & JAY J2 DEEP CONTOUR MOUNTING HARDWARE GUIDE FOR REPLACEMENT PART NUMBERS FOR ALL ATTACHMENT HARDWARE.

**NOTE:** CARRYING HANDLE IS NOT AVAILABLE FOR BACK SHELL WIDTHS 12" THROUGH 15" WIDE.

**NOTE:** SEE JAY J2 & JAY J2 DEEP CONTOUR MOUNTING HARDWARE GUIDE FOR REPLACEMENT PART NUMBERS FOR ALL ATTACHMENT HARDWARE.


**NOTE:** CARRYING HANDLE IS NOT AVAILABLE FOR BACK SHELL WIDTHS 12" THROUGH 15" WIDE.



| Pos. | Item Number | Description                      | Remarks |
|------|-------------|----------------------------------|---------|
| 1    | S2314       | J2 DC 14 BACK SHELL ASSM 14 X 17 | 14 X 17 |
| 1    | S2315       | J2 DC 14 BACK SHELL ASSM 15 X 18 | 15 X 18 |
| 1    | S2316       | J2 DC 14 BACK SHELL ASSM 16 X 19 | 16 X 19 |
| 1    | S2317       | J2 DC 14 BACK SHELL ASSM 17 X 19 | 17 X 19 |
| 1    | S2318       | J2 DC 14 BACK SHELL ASSM 18 X 19 | 18 X 19 |
| 1    | S2319       | J2 DC 14 BACK SHELL ASSM 19 X 20 | 19 X 20 |
| 1    | S2320       | J2 DC 14 BACK SHELL ASSM 20 X 20 | 20 X 20 |
| 2    | B2314F      | J2 DC BACK BASE 14 X 17          | 14 X 17 |
| 2    | B2315F      | J2 DC BACK BASE 15 X 18          | 15 X 18 |
| 2    | B2316F      | J2 DC BACK BASE 16 X 19          |         |
| 2    | B2317F      | J2 DC BACK BASE 17 X 19          |         |

| Pos. | Item Number | Description                          | Remarks |
|------|-------------|--------------------------------------|---------|
| 2    | B2318F      | J2 DC BACK BASE 18 X 19              |         |
| 2    | B2319F      | J2 DC BACK BASE 19 X 20              |         |
| 2    | B2320F      | J2 DC BACK BASE 20 X 20              |         |
| 3    | C2314       | J2 DC BACK COVER 14 X 17             |         |
| 3    | C2315       | J2 DC BACK COVER 15 X 18             |         |
| 3    | C2316       | J2 DC BACK COVER 16 X 19             |         |
| 3    | C2317       | J2 DC BACK COVER 17 X 19             |         |
| 3    | C2318       | J2 DC BACK COVER 18 X 19             |         |
| 3    | C2319       | J2 DC BACK COVER 19 X 20             | 19 X 20 |
| 3    | C2320       | J2 DC BACK COVER 20 X 20             | 20 X 20 |
| 4    | B2320W      | J2 DC BACK VELCRO LUMBAR SUPPORT     |         |
| 5    | 2370        | J2 DC POST LATERAL LUMBAR SUPPORTS   | PAIR    |
| 6    | 2360        | J2 LATERAL SUPPORT ADJUSTABLE WEDGES | PAIR    |

**NOTE:** IF ADDING JAY J2 DEEP CONTOUR BACK TO A CHAIR THAT HAS AN ANGLE ADJUSTABLE BACK PLEASE ORDER 880997 (ITEM 13).



| Pos. | Item Number | Description               | Remarks                              |
|------|-------------|---------------------------|--------------------------------------|
| 1    | M-254       | TOP ADJUSTERS             |                                      |
| 2    | M-253       | TOP HOOK KIT (SHORT)      | SHIPS STANDARD W/ J2 & J2 DEEP BACKS |
| 3    | M-2312      | SHORT HOOKS W/ 1/2 OFFSET |                                      |
| 4    | M-2301      | SHORT HOOKS W/ 1 OFFSET   |                                      |
| 5    | M-256       | TOP HOOK KIT (LONG)       |                                      |

| Pos.    | Item Number | Description                       | Remarks                            |
|---------|-------------|-----------------------------------|------------------------------------|
| 6       | M-2322      | LONG HOOKS W/ 1/2 OFFSET          |                                    |
| 7       | M-2311      | LONG HOOKS W/ 1 OFFSET            |                                    |
| 8       | M-241       | SHOULDER BOLT 1/4 X 1/4           |                                    |
| 9       | M-240       | SHOULDER BOLT 1/4 X 10-32         |                                    |
| 10      | 681C        | TOP BRACKET ASSM - 'QUARTER-TURN' | fits 5/8" to 1.0" diameter tubing  |
| 11      | 2681C       | TOP LONG 1/4 BRACKET ASSM         | fits 1.0" to 1.25" diameter tubing |
| 12      | 680C        | BOTTOM UNIVERSAL POST BRACKETS    |                                    |
| 13      | 880997      | LONG STRAP KIT                    | USED ON 1 1/8" DIAMETER BACK POST  |
| 14      | M-258       | BOTTOM PINS STANDARD              | PIN LENGTH = 2"                    |
| 15      | M-2313      | BOTTOM PINS 1/2 EXTENDED          | PIN LENGTH = 2.5"                  |
| 16      | M-2323      | BOTTOM PINS 1 EXTENDED            | PIN LENGTH = 3"                    |
| 17      | M-500       | HEADREST ADAPTOR KIT              |                                    |
| 18      | M-509       | HEADREST ADAPTOR KIT-TALL BACK    |                                    |
| 19      | M-228       | ALLEN WRENCH                      |                                    |
| 2,10,12 | M-250       | UNIV.RETRO-FIT KIT J2 BK          |                                    |




SHORT HOOKS OFTEN USED WITH ANGLED BACKPOST,  
LONG HOOKS MORE COMMON WITH STRAIGHT BACKPOST.

| J2    | J2 TALL | J2 DC | WIDTH    | NARROWEST SETTING | WIDEST SETTING | 1/2" OFFSET HOOK | 1" OFFSET HOOK |
|-------|---------|-------|----------|-------------------|----------------|------------------|----------------|
| 2512K | 2512TK  | N/A   | 12" WIDE | 11                | 12.5           | 13.5             | 14.5           |
| 2513K | 2513TK  | N/A   | 13" WIDE | 12                | 13.5           | 14.5             | 15.5           |
| 2514K | 2514TK  | 2314  | 14" WIDE | 13                | 14.5           | 15.5             | 16.5           |
| 2515K | 2515TK  | 2315  | 15" WIDE | 14                | 15.5           | 16.5             | 17.5           |
| 2516K | 2516TK  | 2316  | 16" WIDE | 15                | 16.5           | 17.5             | 18.5           |
| 2517K | 2517TK  | 2317  | 17" WIDE | 16                | 17.5           | 18.5             | 19.5           |
| 2518K | 2518TK  | 2318  | 18" WIDE | 17                | 18.5           | 19.5             | 20.5           |
| 2519K | 2519TK  | 2319  | 19" WIDE | 18                | 19.5           | 20.5             | N/A            |
| 2520K | 2520TK  | 2320  | 20" WIDE | 19                | 20.5           | N/A              | N/A            |


| Pos. | Description     |
|------|-----------------|
| -    | SEE CHART ABOVE |

**NOTE:** Kit includes pair of upper and lower clamps with spacers to fit 1 1/8", 1" and 7/8" Back Tubes.


**NOTE:** Please refer to "Jay J2 & Deep Contour Mounting Hardware Guide" page for shell mounting bracket hardware.



| Pos. | Item Number | Description              |
|------|-------------|--------------------------|
| 1    | XM250C      | KIT CLAMP UNIV RETRO-FIT |




| Pos.       | Item Number | Description                    | Remarks   |
|------------|-------------|--------------------------------|-----------|
| 1          | 221115      | 1/4-20 x 3/4" SOC HD C/S BLK   |           |
| 3A         | 750696      | HEADREST W/ZIPPER MEDIUM CG    | 8x4.5 pad |
| 3B         | 750697      | HEADREST W/ZIPPER LARGE CG     | 9x4.5 pad |
| 4          | 505932      | HEADREST MOUNTING TRACK        |           |
| 5          | 476304      | 1/4-28 x .31 BHCS W/PATCH      |           |
| 6          | 221010      | 10-32 x .75 SHCS BLK W/PL      |           |
| 7          | 261204      | #10 AN WASHER .063             |           |
| 8          | 331077      | HEADREST LINK 6 INCH           |           |
| 9          | 341045      | #10-24 CONE POINT SET SCREW    |           |
| 10         | 505933      | HEADREST COLLAR                |           |
| 11         | 226072      | 1/4-20 x 1/2 FH C/S BLK        |           |
| 13         | 226095      | 1/4-20 x 3/8" SOC C/S BLK ZN   |           |
| 14         | 341046      | HEADREST CLAMPING HANDLE SCREW |           |
| 15         | 340033      | 1/4-20 WELD NUTS SLAB BASE     |           |
| 4,5,13-16  | 888101      | CLAMP HEADREST RIGID BACK      |           |
| 1,2,4,6-16 | 881003      | HEADREST MOUNTING BRACKET ASSM |           |
| 1-3A,4-16  | 881907      | HEADREST ASSEMBLY MEDIUM       | 8x4.5 pad |
| 1-3B,4-16  | 881908      | HEADREST ASSEMBLY LARGE        | 9x4.5 pad |



| Pos. | Item Number | Description                    |
|------|-------------|--------------------------------|
| 1    | XM2025      | MNTNG PLATE/PEDIAT.HEAD REST   |
| 2    | 226018-QU   | 10-32 X 3/8 BH C/S W/PATCH BLK |
| 3    | 476304      | 1/4-28 x .31 BHCS W/PATCH      |
| 1-3  | 476314      | HEADREST NMT KIT W/HDWR POWER  |

**NOTE:** Parts only, for complete Backrest see J3 Back order form.

Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.




| Pos. | Item Number | Description             | Remarks | Pos. | Item Number | Description              | Remarks |
|------|-------------|-------------------------|---------|------|-------------|--------------------------|---------|
| 1    | CJ3PALTM12  | J3 CVR POS LT 12W MED   |         | 1    | CJ3PALTT18  | J3 CVR, POS LT 18W TALL  |         |
| 1    | CJ3PALTT12  | J3 CVR POS LT 12W TAL   |         | 1    | CJ3PALTT20  | J3 CVR, POS LT 20W TALL  |         |
| 1    | CJ3PALTS14  | J3 CVR, POS LT 14W STD  |         | 2    | BJ3PALTM12  | J3 FOAM POS LT 12W MED   |         |
| 1    | CJ3PALTM14  | J3 CVR, POS LT 14W MED  |         | 2    | BJ3PALTT12  | J3 FOAM POS LT 12W TAL   |         |
| 1    | CJ3PALTT14  | J3 CVR, POS LT 14W TALL |         | 2    | BJ3PALTS14  | J3 FOAM, POS LT 14W STD  |         |
| 1    | CJ3PALTS16  | J3 CVR, POS LT 16W STD  |         | 2    | BJ3PALTM14  | J3 FOAM, POS LT 14W MED  |         |
| 1    | CJ3PALTM16  | J3 CVR, POS LT 16W MED  |         | 2    | BJ3PALTT14  | J3 FOAM, POS LT 14W TALL |         |
| 1    | CJ3PALTT16  | J3 CVR, POS LT 16W TALL |         | 2    | BJ3PALTS16  | J3 FOAM, POS LT 16W STD  |         |
| 1    | CJ3PALTS18  | J3 CVR, POS LT 18W STD  |         | 2    | BJ3PALTM16  | J3 FOAM, POS LT 16W MED  |         |
| 1    | CJ3PALTM18  | J3 CVR, POS LT 18W MED  |         | 2    | BJ3PALTT16  | J3 FOAM, POS LT 16W TALL |         |

| Pos. | Item Number | Description               | Remarks | Pos. | Item Number | Description               | Remarks     |
|------|-------------|---------------------------|---------|------|-------------|---------------------------|-------------|
| 2    | BJ3PALTS18  | J3 FOAM, POS LT 18W STD   |         | 3    | SJ3PALTT18  | J3 SHELL, POS LT 18W TALL |             |
| 2    | BJ3PALTM18  | J3 FOAM, POS LT 18W MED   |         | 3    | SJ3PALTT20  | J3 SHELL, POS LT 20W TALL |             |
| 2    | BJ3PALTT18  | J3 FOAM, POS LT 18W TALL  |         | 4    | ICJ3PALTM12 | J3 BK POS LT 12W MED SH   | inner cover |
| 2    | BJ3PALTT20  | J3 FOAM, POS LT 20W TALL  |         | 4    | ICJ3PALTM14 | J3 BK POS LT 14W MED SH   | inner cover |
| 3    | SJ3PALTM12  | J3 SHELL POS LT 12W MED   |         | 4    | ICJ3PALTM16 | J3 BK POS LT 16W MED SH   | inner cover |
| 3    | SJ3PALTT12  | J3 SHELL POS LT 12W TAL   |         | 4    | ICJ3PALTM18 | J3 BK POS LT 18W MED SH   | inner cover |
| 3    | SJ3PALTS14  | J3 SHELL, POS LT 14W STD  |         | 4    | ICJ3PALTS14 | J3 BK POS LT 14W STD SH   | inner cover |
| 3    | SJ3PALTM14  | J3 SHELL, POS LT 14W MED  |         | 4    | ICJ3PALTS16 | J3 BK POS LT 16W STD SH   | inner cover |
| 3    | SJ3PALTT14  | J3 SHELL, POS LT 14W TALL |         | 4    | ICJ3PALTS18 | J3 BK POS LT 18W STD SH   | inner cover |
| 3    | SJ3PALTS16  | J3 SHELL, POS LT 16W STD  |         | 4    | ICJ3PALTT12 | J3 BK POS LT 12W TALL SH  | inner cover |
| 3    | SJ3PALTM16  | J3 SHELL, POS LT 16W MED  |         | 4    | ICJ3PALTT14 | J3 BK POS LT 14W TALL SH  | inner cover |
| 3    | SJ3PALTT16  | J3 SHELL, POS LT 16W TALL |         | 4    | ICJ3PALTT16 | J3 BK POS LT 16W TALL SH  | inner cover |
| 3    | SJ3PALTS18  | J3 SHELL, POS LT 18W STD  |         | 4    | ICJ3PALTT18 | J3 BK POS LT 18W TALL SH  | inner cover |
| 3    | SJ3PALTM18  | J3 SHELL, POS LT 18W MED  |         | 4    | ICJ3PALTT20 | J3 BK POS LT 20W TALL SH  | inner cover |


**NOTE:** Parts only, for complete Backrest see J3 Back order form.

Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.



| Pos. | Item Number | Description            | Remarks | Pos. | Item Number  | Description              | Remarks     |
|------|-------------|------------------------|---------|------|--------------|--------------------------|-------------|
| 1    | CJ3PDCLTM14 | J3 CVR PDC LT 14W MED  |         | 3    | SJ3PDCLTM14  | J3 SHELL PDC LT 14W MED  |             |
| 1    | CJ3PDCLTT14 | J3 CVR PDC LT 14W TAL  |         | 3    | SJ3PDCLTT14  | J3 SHELL PDC LT 14W TAL  |             |
| 1    | CJ3PDCLTM16 | J3 CVR PDC LT 16W MED  |         | 3    | SJ3PDCLTM16  | J3 SHELL PDC LT 16W MED  |             |
| 1    | CJ3PDCLTT16 | J3 CVR PDC LT 16W TAL  |         | 3    | SJ3PDCLTT16  | J3 SHELL PDC LT 16W TAL  |             |
| 1    | CJ3PDCLTM18 | J3 CVR PDC LT 18W MED  |         | 3    | SJ3PDCLTM18  | J3 SHELL PDC LT 18W MED  |             |
| 1    | CJ3PDCLTT18 | J3 CVR PDC LT 18W TAL  |         | 3    | SJ3PDCLTT18  | J3 SHELL PDC LT 18W TAL  |             |
| 1    | CJ3PDCLTT20 | J3 CVR PDC LT 20W TAL  |         | 3    | SJ3PDCLTT20  | J3 SHELL PDC LT 20W TAL  |             |
| 2    | BJ3PDCLTM14 | J3 FOAM PDC LT 14W MED |         | 4    | ICJ3PDCLTM14 | J3 BK PDC LT 14W MED SH  | inner cover |
| 2    | BJ3PDCLTT14 | J3 FOAM PDC LT 14W TAL |         | 4    | ICJ3PDCLTM16 | J3 BK PDC LT 16W MED SH  | inner cover |
| 2    | BJ3PDCLTM16 | J3 FOAM PDC LT 16W MED |         | 4    | ICJ3PDCLTM18 | J3 BK PDC LT 18W MED SH  | inner cover |
| 2    | BJ3PDCLTT16 | J3 FOAM PDC LT 16W TAL |         | 4    | ICJ3PDCLTT14 | J3 BK PDC LT 14W TALL SH | inner cover |
| 2    | BJ3PDCLTM18 | J3 FOAM PDC LT 18W MED |         | 4    | ICJ3PDCLTT16 | J3 BK PDC LT 16W TALL SH | inner cover |
| 2    | BJ3PDCLTT18 | J3 FOAM PDC LT 18W TAL |         | 4    | ICJ3PDCLTT18 | J3 BK PDC LT 18W TALL SH | inner cover |
| 2    | BJ3PDCLTT20 | J3 FOAM PDC LT 20W TAL |         | 4    | ICJ3PDCLTT20 | J3 BK PDC LT 20W TALL SH | inner cover |

**NOTE:** Parts only, for complete Backrest see J3 Back order form.




| Pos. | Item Number | Description                     | Pos. | Item Number | Description                     |
|------|-------------|---------------------------------|------|-------------|---------------------------------|
| 1A   | CJ3CFLTS14  | COVER J3 CF BACK PA LT SHT 14W  | 2B   | BJ3PALTT20  | J3 FOAM, POS LT 20W TALL        |
| 1A   | CJ3CFLTS16  | COVER J3 CF BACK PA LT SHT 16W  | 3A   | SJ3CFLTS14  | SHELL J3 CF BACK PA LT SHT 14W  |
| 1A   | CJ3CFLTS18  | COVER J3 CF BACK PA LT SHT 18W  | 3A   | SJ3CFLTS16  | SHELL J3 CF BACK PA LT SHT 16W  |
| 1B   | CJ3CFLTT14  | COVER J3 CF BACK PA LT TALL 14W | 3A   | SJ3CFLTS18  | SHELL J3 CF BACK PA LT SHT 18W  |
| 1B   | CJ3CFLTT16  | COVER J3 CF BACK PA LT TALL 16W | 3B   | SJ3CFLTT14  | SHELL J3 CF BACK PA LT TALL 14W |
| 1B   | CJ3CFLTT18  | COVER J3 CF BACK PA LT TALL 18W | 3B   | SJ3CFLTT16  | SHELL J3 CF BACK PA LT TALL 16W |
| 1B   | CJ3CFLTT20  | COVER J3 CF BACK PA LT TALL 20W | 3B   | SJ3CFLTT18  | SHELL J3 CF BACK PA LT TALL 18W |
| 2A   | BJ3PALTS14  | J3 FOAM, POS LT 14W STD         | 3B   | SJ3CFLTT20  | SHELL J3 CF BACK PA LT TALL 20W |
| 2A   | BJ3PALTS16  | J3 FOAM, POS LT 16W STD         | 4A   | ICJ3PALTS14 | J3 BK POS LT 14W STD SH         |
| 2A   | BJ3PALTS18  | J3 FOAM, POS LT 18W STD         | 4A   | ICJ3PALTT16 | J3 BK POS LT 16W TALL SH        |
| 2B   | BJ3PALTT14  | J3 FOAM, POS LT 14W TALL        | 4A   | ICJ3PALTS18 | J3 BK POS LT 18W STD SH         |
| 2B   | BJ3PALTT16  | J3 FOAM, POS LT 16W TALL        | 4B   | ICJ3PALTT14 | J3 BK POS LT 14W TALL SH        |
| 2B   | BJ3PALTT18  | J3 FOAM, POS LT 18W TALL        | 4B   | ICJ3PALTS16 | J3 BK POS LT 16W STD SH         |

| Pos. | Item Number | Description              | Pos. | Item Number | Description              |
|------|-------------|--------------------------|------|-------------|--------------------------|
| 4B   | ICJ3PALTT18 | J3 BK POS LT 18W TALL SH | 4B   | ICJ3PALTT20 | J3 BK POS LT 20W TALL SH |

**NOTE:** Parts only, for complete Backrest see J3 Back order form.

Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.




| Pos. | Item Number | Description             | Remarks |
|------|-------------|-------------------------|---------|
| 1A   | CJ3PAMTS12  | J3 CVR POS MT 12W STD   |         |
| 1A   | CJ3PAMTM12  | J3 CVR POS MT 12W MED   |         |
| 1A   | CJ3PAMTT12  | J3 CVR POS MT 12W TAL   |         |
| 1A   | CJ3PAMTS14  | J3 CVR, POS MT 14W STD  |         |
| 1A   | CJ3PAMTM14  | J3 CVR, POS MT 14W MED  |         |
| 1A   | CJ3PAMTT14  | J3 CVR, POS MT 14W TALL |         |
| 1A   | CJ3PAMTS16  | J3 CVR, POS MT 16W STD  |         |
| 1A   | CJ3PAMTM16  | J3 CVR, POS MT 16W MED  |         |
| 1A   | CJ3PAMTT16  | J3 CVR, POS MT 16W TALL |         |
| 1A   | CJ3PAMTS18  | J3 CVR, POS MT 18W STD  |         |

| Pos. | Item Number  | Description                   | Remarks         |
|------|--------------|-------------------------------|-----------------|
| 1A   | CJ3PAMTM18   | J3 CVR, POS MT 18W MED        |                 |
| 1A   | CJ3PAMTT18   | J3 CVR, POS MT 18W TALL       |                 |
| 1A   | CJ3PAMTS20   | J3 CVR, POS MT 20W STD        |                 |
| 1A   | CJ3PAMTM20   | J3 CVR, POS MT 20W MED        |                 |
| 1A   | CJ3PAMTT20   | J3 CVR, POS MT 20W TALL       |                 |
| 1B   | CJ3PAMTS16CU | CVR J3 BK POS MT 16W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PAMTM16CU | CVR J3 BK POS MT 16W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PAMTT16CU | CVR J3 BK POS MT 16W TALL CU  | use w/ Cradle-U |
| 1B   | CJ3PAMTS18CU | CVR J3 BK POS MT 18W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PAMTM18CU | CVR J3 BK POS MT 18W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PAMTT18CU | CVR J3 BK POS MT 18W TALL CU  | use w/ Cradle-U |
| 1B   | CJ3PAMTS20CU | CVR J3 BK POS MT 20W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PAMTM20CU | CVR J3 BK POS MT 20W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PAMTT20CU | CVR J3 BK POS MT 20W TALL CU  | use w/ Cradle-U |
| 2A   | BJ3PAMTS12   | J3 FOAM POS MT 12W STD        |                 |
| 2A   | BJ3PAMTM12   | J3 FOAM POS MT 12W MED        |                 |
| 2A   | BJ3PAMTT12   | J3 FOAM POS MT 12W TAL        |                 |
| 2A   | BJ3PAMTS14   | J3 FOAM, POS MT 14W STD       |                 |
| 2A   | BJ3PAMTM14   | J3 FOAM, POS MT 14W MED       |                 |
| 2A   | BJ3PAMTT14   | J3 FOAM, POS MT 14W TALL      |                 |
| 2A   | BJ3PAMTS16   | J3 FOAM, POS MT 16W STD       |                 |
| 2A   | BJ3PAMTM16   | J3 FOAM, POS MT 16W MED       |                 |
| 2A   | BJ3PAMTT16   | J3 FOAM, POS MT 16W TALL      |                 |
| 2A   | BJ3PAMTS18   | J3 FOAM, POS MT 18W STD       |                 |
| 2A   | BJ3PAMTM18   | J3 FOAM, POS MT 18W MED       |                 |
| 2A   | BJ3PAMTT18   | J3 FOAM, POS MT 18W TALL      |                 |
| 2A   | BJ3PAMTT20   | J3 FOAM, POS MT 20W TALL      |                 |
| 2A   | BJ3PAMTS20   | J3 FOAM, POS MT 20W STD       |                 |
| 2B   | BJ3PAMTS16CU | BASE J3 BK POS MT 16W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PAMTM16CU | BASE J3 BK POS MT 16W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PAMTT16CU | BASE J3 BK POS MT 16W TALL CU | use w/ Cradle-U |
| 2B   | BJ3PAMTS18CU | BASE J3 BK POS MT 18W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PAMTM18CU | BASE J3 BK POS MT 18W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PAMTT18CU | BASE J3 BK POS MT 18W TALL CU | use w/ Cradle-U |
| 2B   | BJ3PAMTS20CU | BASE J3 BK POS MT 20W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PAMTM20CU | BASE J3 BK POS MT 20W MED CU  | use w/ Cradle-U |
| 3    | SJ3PAMTS12   | J3 SHELL POS MT 12W STD       |                 |
| 3    | SJ3PAMTM12   | J3 SHELL POS MT 12W MED       |                 |
| 3    | SJ3PAMTT12   | J3 SHELL POS MT 12W TAL       |                 |
| 3    | SJ3PAMTS14   | J3 SHELL, POS MT 14W STD      |                 |
| 3    | SJ3PAMTM14   | J3 SHELL, POS MT 14W MED      |                 |
| 3    | SJ3PAMTT14   | J3 SHELL, POS MT 14W TALL     |                 |
| 3    | SJ3PAMTS16   | J3 SHELL, POS MT 16W STD      |                 |
| 3    | SJ3PAMTM16   | J3 SHELL, POS MT 16W MED      |                 |
| 3    | SJ3PAMTT16   | J3 SHELL, POS MT 16W TALL     |                 |
| 3    | SJ3PAMTS18   | J3 SHELL, POS MT 18W STD      |                 |
| 3    | SJ3PAMTM18   | J3 SHELL, POS MT 18W MED      |                 |
| 3    | SJ3PAMTT18   | J3 SHELL, POS MT 18W TALL     |                 |
| 3    | SJ3PAMTS20   | J3 SHELL, POS MT 20W STD      |                 |
| 3    | SJ3PAMTM20   | J3 SHELL, POS MT 20W MED      |                 |
| 3    | SJ3PAMTT20   | J3 SHELL, POS MT 20W TALL     |                 |
| 4    | ICJ3PAMTS12  | J3 BK POS MT 12W STD SH       | inner cover     |
| 4    | ICJ3PAMTM12  | J3 BK POS MT 12W MED SH       | inner cover     |
| 4    | ICJ3PAMTT12  | J3 BK POS MT 12W TALL SH      | inner cover     |
| 4    | ICJ3PAMTS14  | J3 BK POS MT 14W STD SH       | inner cover     |
| 4    | ICJ3PAMTM14  | J3 BK POS MT 14W MED SH       | inner cover     |
| 4    | ICJ3PAMTT14  | J3 BK POS MT 14W TALL SH      | inner cover     |

| Pos. | Item Number | Description              | Remarks     |
|------|-------------|--------------------------|-------------|
| 4    | ICJ3PAMTS16 | J3 BK POS MT 16W STD SH  | inner cover |
| 4    | ICJ3PAMTM16 | J3 BK POS MT 16W MED SH  | inner cover |
| 4    | ICJ3PAMTT16 | J3 BK POS MT 16W TALL SH | inner cover |
| 4    | ICJ3PAMTS18 | J3 BK POS MT 18W STD SH  | inner cover |
| 4    | ICJ3PAMTM18 | J3 BK POS MT 18W MED SH  | inner cover |
| 4    | ICJ3PAMTT18 | J3 BK POS MT 18W TALL SH | inner cover |
| 4    | ICJ3PAMTM20 | J3 BK POS MT 20W MED SH  | inner cover |
| 4    | ICJ3PAMTT20 | J3 BK POS MT 20W TALL SH | inner cover |
| 4    | ICJ3PAMTS20 | J3 BK POS MT 20W STD SH  | inner cover |

**NOTE:** Parts only, for complete Backrest see J3 Back order form.

Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.




| Pos. | Item Number | Description              | Pos. | Item Number | Description                  |
|------|-------------|--------------------------|------|-------------|------------------------------|
| 1    | CJ3PAMTS22  | J3HD CVR POS MT 22W STD  | 2    | BJ3PAMTT26  | J3HD FOAM POS MT 26W TAL     |
| 1    | CJ3PAMTT22  | J3HD CVR POS MT 22W TAL  | 3    | SJ3PAMTS22  | J3HD SHELL POS MT 22W STD    |
| 1    | CJ3PAMTS24  | J3HD CVR POS MT 24W STD  | 3    | SJ3PAMTT22  | J3HD SHELL POS MT 22W TAL    |
| 1    | CJ3PAMTT24  | J3HD CVR POS MT 24W TAL  | 3    | SJ3PAMTS24  | J3HD SHELL POS MT 24W STD    |
| 1    | CJ3PAMTS26  | J3HD CVR POS MT 26W STD  | 3    | SJ3PAMTT24  | J3HD SHELL POS MT 24W TAL    |
| 1    | CJ3PAMTT26  | J3HD CVR POS MT 26W TAL  | 3    | SJ3PAMTS26  | J3HD SHELL POS MT 26W STD    |
| 2    | BJ3PAMTS22  | J3HD FOAM POS MT 22W STD | 3    | SJ3PAMTT26  | J3HD SHELL POS MT 26W TAL    |
| 2    | BJ3PAMTT22  | J3HD FOAM POS MT 22W TAL | 4    | ICJ3PAMTS22 | J3 HD BK,POS MT 22W STD, FS  |
| 2    | BJ3PAMTS24  | J3HD FOAM POS MT 24W STD | 4    | ICJ3PAMTT22 | J3 HD BK,POS MT 22W TALL, FS |
| 2    | BJ3PAMTT24  | J3HD FOAM POS MT 24W TAL | 4    | ICJ3PAMTS24 | J3 HD BK,POS MT 24W STD, FS  |
| 2    | BJ3PAMTS26  | J3HD FOAM POS MT 26W STD | 4    | ICJ3PAMTT24 | J3 HD BK,POS MT 24W TALL, FS |

| Pos. | Item Number | Description                 | Pos. | Item Number | Description                  |
|------|-------------|-----------------------------|------|-------------|------------------------------|
| 4    | ICJ3PAMTS26 | J3 HD BK,POS MT 26W STD, FS | 4    | ICJ3PAMTT26 | J3 HD BK,POS MT 26W TALL, FS |

Parts only; for complete Backrest see J3 Back order form.

Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.




| Pos. | Item Number  | Description                  | Remarks         |
|------|--------------|------------------------------|-----------------|
| 1A   | CJ3PLMTS12   | J3 CVR SHAL MT 12W STD       |                 |
| 1A   | CJ3PLMTM12   | J3 CVR SHAL MT 12W MED       |                 |
| 1A   | CJ3PLMTT12   | J3 CVR SHAL MT 12W TAL       |                 |
| 1A   | CJ3PLMTS14   | J3 CVR, SHAL MT 14W STD      |                 |
| 1A   | CJ3PLMTM14   | J3 CVR, SHAL MT 14W MED      |                 |
| 1A   | CJ3PLMTT14   | J3 CVR, SHAL MT 14W TALL     |                 |
| 1A   | CJ3PLMTS16   | J3 CVR, SHAL MT 16W STD      |                 |
| 1A   | CJ3PLMTM16   | J3 CVR, SHAL MT 16W MED      |                 |
| 1A   | CJ3PLMTT16   | J3 CVR, SHAL MT 16W TALL     |                 |
| 1A   | CJ3PLMTS18   | J3 CVR, SHAL MT 18W STD      |                 |
| 1A   | CJ3PLMTM18   | J3 CVR, SHAL MT 18W MED      |                 |
| 1A   | CJ3PLMTT18   | J3 CVR, SHAL MT 18W TALL     |                 |
| 1A   | CJ3PLMTS20   | J3 CVR, SHAL MT 20W STD      |                 |
| 1A   | CJ3PLMTM20   | J3 CVR, SHAL MT 20W MED      |                 |
| 1A   | CJ3PLMTT20   | J3 CVR, SHAL MT 20W TALL     |                 |
| 1B   | CJ3PLMTS16CU | CVR J3 BK SHAL MT 16W STD CU | use w/ Cradle-U |
| 1B   | CJ3PLMTM16CU | CVR J3 BK SHAL MT 16W MED CU | use w/ Cradle-U |

| Pos. | Item Number  | Description                    | Remarks         |
|------|--------------|--------------------------------|-----------------|
| 1B   | CJ3PLMTT16CU | CVR J3 BK SHAL MT 16W TALL CU  | use w/ Cradle-U |
| 1B   | CJ3PLMTS18CU | CVR J3 BK SHAL MT 18W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PLMTM18CU | CVR J3 BK SHAL MT 18W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PLMTT18CU | CVR J3 BK SHAL MT 18W TALL CU  | use w/ Cradle-U |
| 1B   | CJ3PLMTS20CU | CVR J3 BK SHAL MT 20W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PLMTM20CU | CVR J3 BK SHAL MT 20W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PLMTT20CU | CVR J3 BK SHAL MT 20W TALL CU  | use w/ Cradle-U |
| 2A   | BJ3PLMTS12   | J3 FOAM SHAL MT 12W STD        |                 |
| 2A   | BJ3PLMTM12   | J3 FOAM SHAL MT 12W MED        |                 |
| 2A   | BJ3PLMTT12   | J3 FOAM SHAL MT 12W TAL        |                 |
| 2A   | BJ3PLMTS14   | J3 FOAM, SHAL MT 14W STD       |                 |
| 2A   | BJ3PLMTM14   | J3 FOAM, SHAL MT 14W MED       |                 |
| 2A   | BJ3PLMTT14   | J3 FOAM, SHAL MT 14W TALL      |                 |
| 2A   | BJ3PLMTS16   | J3 FOAM, SHAL MT 16W STD       |                 |
| 2A   | BJ3PLMTM16   | J3 FOAM, SHAL MT 16W MED       |                 |
| 2A   | BJ3PLMTT16   | J3 FOAM, SHAL MT 16W TALL      |                 |
| 2A   | BJ3PLMTS18   | J3 FOAM, SHAL MT 18W STD       |                 |
| 2A   | BJ3PLMTM18   | J3 FOAM, SHAL MT 18W MED       |                 |
| 2A   | BJ3PLMTT18   | J3 FOAM, SHAL MT 18W TALL      |                 |
| 2A   | BJ3PLMTS20   | J3 FOAM, SHAL MT 20W STD       |                 |
| 2A   | BJ3PLMTM20   | J3 FOAM, SHAL MT 20W MED       |                 |
| 2A   | BJ3PLMTT20   | J3 FOAM, SHAL MT 20W TALL      |                 |
| 2B   | BJ3PLMTS16CU | BASE J3 BK SHAL MT 16W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PLMTM16CU | BASE J3 BK SHAL MT 16W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PLMTT16CU | BASE J3 BK SHAL MT 16W TALL CU | use w/ Cradle-U |
| 2B   | BJ3PLMTS18CU | BASE J3 BK SHAL MT 18W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PLMTM18CU | BASE J3 BK SHAL MT 18W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PLMTT18CU | BASE J3 BK SHAL MT 18W TALL CU | use w/ Cradle-U |
| 2B   | BJ3PLMTS20CU | BASE J3 BK SHAL MT 20W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PLMTM20CU | BASE J3 BK SHAL MT 20W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PLMTT20CU | BASE J3 BK SHAL MT 20W TALL CU | use w/ Cradle-U |
| 3    | SJ3PLMTS12   | J3 SHELL SHAL MT 12W STD       |                 |
| 3    | SJ3PLMTM12   | J3 SHELL SHAL MT 12W MED       |                 |
| 3    | SJ3PLMTT12   | J3 SHELL SHAL MT 12W TAL       |                 |
| 3    | SJ3PLMTS14   | J3 SHELL, SHAL MT 14W STD      |                 |
| 3    | SJ3PLMTM14   | J3 SHELL, SHAL MT 14W MED      |                 |
| 3    | SJ3PLMTT14   | J3 SHELL, SHAL MT 14W TALL     |                 |
| 3    | SJ3PLMTS16   | J3 SHELL, SHAL MT 16W STD      |                 |
| 3    | SJ3PLMTM16   | J3 SHELL, SHAL MT 16W MED      |                 |
| 3    | SJ3PLMTT16   | J3 SHELL, SHAL MT 16W TALL     |                 |
| 3    | SJ3PLMTS18   | J3 SHELL, SHAL MT 18W STD      |                 |
| 3    | SJ3PLMTM18   | J3 SHELL, SHAL MT 18W MED      |                 |
| 3    | SJ3PLMTT18   | J3 SHELL, SHAL MT 18W TALL     |                 |
| 3    | SJ3PLMTS20   | J3 SHELL, SHAL MT 20W STD      |                 |
| 3    | SJ3PLMTM20   | J3 SHELL, SHAL MT 20W MED      |                 |
| 3    | SJ3PLMTT20   | J3 SHELL, SHAL MT 20W TALL     |                 |
| 4    | ICJ3PLMTS12  | J3 BK SHAL MT 12W STD SH       |                 |
| 4    | ICJ3PLMTM12  | J3 BK SHAL MT 12W MED SH       |                 |
| 4    | ICJ3PLMTT12  | J3 BK SHAL MT 12W TALL SH      |                 |
| 4    | ICJ3PLMTS14  | J3 BK SHAL MT 14W STD SH       |                 |
| 4    | ICJ3PLMTM14  | J3 BK SHAL MT 14W MED SH       |                 |
| 4    | ICJ3PLMTT14  | J3 BK SHAL MT 14W TALL SH      |                 |
| 4    | ICJ3PLMTS16  | J3 BK SHAL MT 16W STD SH       |                 |
| 4    | ICJ3PLMTM16  | J3 BK SHAL MT 16W MED SH       |                 |
| 4    | ICJ3PLMTT16  | J3 BK SHAL MT 16W TALL SH      |                 |
| 4    | ICJ3PLMTS18  | J3 BK SHAL MT 18W STD SH       |                 |
| 4    | ICJ3PLMTM18  | J3 BK SHAL MT 18W MED SH       |                 |

| Pos. | Item Number | Description               | Remarks |
|------|-------------|---------------------------|---------|
| 4    | ICJ3PLMTT18 | J3 BK SHAL MT 18W TALL SH |         |
| 4    | ICJ3PLMTS20 | J3 BK SHAL MT 20W STD SH  |         |
| 4    | ICJ3PLMTM20 | J3 BK SHAL MT 20W MED SH  |         |
| 4    | ICJ3PLMTT20 | J3 BK SHAL MT 20W TALL SH |         |

**Note:** parts only; for complete Backrest see J3 Back order form.




| Pos. | Item Number | Description               | Pos. | Item Number | Description                    |
|------|-------------|---------------------------|------|-------------|--------------------------------|
| 1A   | CJ3PDMTS14  | J3 BK CVR PD MT 14W SHT   | 2B   | BJ3PDMTT16  | J3 BK FOAM PD MT 16W TALL      |
| 1A   | CJ3PDMTS16  | J3 BK CVR PD MT 16W SHT   | 2B   | BJ3PDMTT18  | J3 BK FOAM PD MT 18W TALL      |
| 1A   | CJ3PDMTS18  | J3 BK CVR PD MT 18W SHT   | 2B   | BJ3PDMTT20  | J3 BK FOAM PD MT 20W TALL      |
| 1A   | CJ3PDMTS20  | J3 BK CVR PD MT 20W SHT   | 3A   | ICJ3PDMTS14 | J3 BK INNER CVR PD MT 14W SHT  |
| 1B   | CJ3PDMTT14  | J3 BK CVR PD MT 14W TALL  | 3A   | ICJ3PDMTS16 | J3 BK INNER CVR PD MT 16W SHT  |
| 1B   | CJ3PDMTT16  | J3 BK CVR PD MT 16W TALL  | 3A   | ICJ3PDMTS18 | J3 BK INNER CVR PD MT 18W SHT  |
| 1B   | CJ3PDMTT18  | J3 BK CVR PD MT 18W TALL  | 3A   | ICJ3PDMTS20 | J3 BK INNER CVR PD MT 20W SHT  |
| 1B   | CJ3PDMTT20  | J3 BK CVR PD MT 20W TALL  | 3B   | ICJ3PDMTT14 | J3 BK INNER CVR PD MT 14W TALL |
| 2A   | BJ3PDMTS14  | J3 BK FOAM PD MT 14W SHT  | 3B   | ICJ3PDMTT16 | J3 BK INNER CVR PD MT 16W TALL |
| 2A   | BJ3PDMTS16  | J3 BK FOAM PD MT 16W SHT  | 3B   | ICJ3PDMTT18 | J3 BK INNER CVR PD MT 18W TALL |
| 2A   | BJ3PDMTS18  | J3 BK FOAM PD MT 18W SHT  | 3B   | ICJ3PDMTT20 | J3 BK INNER CVR PD MT 20W TALL |
| 2A   | BJ3PDMTS20  | J3 BK FOAM PD MT 20W SHT  | 4A   | SJ3PDMTS14  | J3 BK SHL PD MT 14W SHT        |
| 2B   | BJ3PDMTT14  | J3 BK FOAM PD MT 14W TALL | 4A   | SJ3PDMTS16  | J3 BK SHL PD MT 16W SHT        |

| Pos. | Item Number | Description              | Pos. | Item Number | Description              |
|------|-------------|--------------------------|------|-------------|--------------------------|
| 4A   | SJ3PDMTS18  | J3 BK SHL PD MT 18W SHT  | 4B   | SJ3PDMTT16  | J3 BK SHL PD MT 16W TALL |
| 4A   | SJ3PDMTS20  | J3 BK SHL PD MT 20W SHT  | 4B   | SJ3PDMTT18  | J3 BK SHL PD MT 18W TALL |
| 4B   | SJ3PDMTT14  | J3 BK SHL PD MT 14W TALL | 4B   | SJ3PDMTT20  | J3 BK SHL PD MT 20W TALL |

**NOTE:** Parts only, for complete Backrest see J3 Back order form.


Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.



| Pos. | Item Number   | Description                  | Remarks         |
|------|---------------|------------------------------|-----------------|
| 1A   | CJ3PDCMTS14   | J3 CVR PDC MT 14W STD        |                 |
| 1A   | CJ3PDCMTT14   | J3 CVR PDC MT 14W TAL        |                 |
| 1A   | CJ3PDCMTS16   | J3 CVR PDC MT 16W STD        |                 |
| 1A   | CJ3PDCMTT16   | J3 CVR PDC MT 16W TAL        |                 |
| 1A   | CJ3PDCMTS18   | J3 CVR PDC MT 18W STD        |                 |
| 1A   | CJ3PDCMTT18   | J3 CVR PDC MT 18W TAL        |                 |
| 1A   | CJ3PDCMTS20   | J3 CVR PDC MT 20W STD        |                 |
| 1A   | CJ3PDCMTT20   | J3 CVR PDC MT 20W TAL        |                 |
| 1B   | CJ3PDCMTS16CU | CVR J3 BK PDC MT 16W STD CU  | use w/ Cradle-U |
| 1B   | CJ3PDCMTT16CU | CVR J3 BK PDC MT 16W TALL CU | use w/ Cradle-U |

| Pos. | Item Number   | Description                   | Remarks         |
|------|---------------|-------------------------------|-----------------|
| 1B   | CJ3PDCMTS18CU | CVR J3 BK PDC MT 18W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PDCMTT18CU | CVR J3 BK PDC MT 18W TALL CU  | use w/ Cradle-U |
| 1B   | CJ3PDCMTS20CU | CVR J3 BK PDC MT 20W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PDCMTT20CU | CVR J3 BK PDC MT 20W TALL CU  | use w/ Cradle-U |
| 2A   | BJ3PDCMTS14   | J3 FOAM PDC MT 14W STD        |                 |
| 2A   | BJ3PDCMTT14   | J3 FOAM PDC MT 14W TAL        |                 |
| 2A   | BJ3PDCMTS16   | J3 FOAM PDC MT 16W STD        |                 |
| 2A   | BJ3PDCMTT16   | J3 FOAM PDC MT 16W TAL        |                 |
| 2A   | BJ3PDCMTS18   | J3 FOAM PDC MT 18W STD        |                 |
| 2A   | BJ3PDCMTT18   | J3 FOAM PDC MT 18W TAL        |                 |
| 2A   | BJ3PDCMTS20   | J3 FOAM PDC MT 20W STD        |                 |
| 2A   | BJ3PDCMTT20   | J3 FOAM PDC MT 20W TAL        |                 |
| 2B   | BJ3PDCMTS16CU | BASE J3 BK PDC MT 16W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PDCMTT16CU | BASE J3 BK PDC MT 16W TALL CU | use w/ Cradle-U |
| 2B   | BJ3PDCMTS18CU | BASE J3 BK PDC MT 18W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PDCMTT18CU | BASE J3 BK PDC MT 18W TALL CU | use w/ Cradle-U |
| 2B   | BJ3PDCMTS20CU | BASE J3 BK PDC MT 20W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PDCMTT20CU | BASE J3 BK PDC MT 20W TALL CU | use w/ Cradle-U |
| 3    | SJ3PDCMTS14   | J3 SHELL PDC MT 14W STD       |                 |
| 3    | SJ3PDCMTT14   | J3 SHELL PDC MT 14W TAL       |                 |
| 3    | SJ3PDCMTS16   | J3 SHELL PDC MT 16W STD       |                 |
| 3    | SJ3PDCMTT16   | J3 SHELL PDC MT 16W TAL       |                 |
| 3    | SJ3PDCMTS18   | J3 SHELL PDC MT 18W STD       |                 |
| 3    | SJ3PDCMTT18   | J3 SHELL PDC MT 18W TAL       |                 |
| 3    | SJ3PDCMTS20   | J3 SHELL PDC MT 20W STD       |                 |
| 3    | SJ3PDCMTT20   | J3 SHELL PDC MT 20W TAL       |                 |
| 4    | ICJ3PDCMTS14  | J3 BK PDC MT 14W STD SH       |                 |
| 4    | ICJ3PDCMTT14  | J3 BK PDC MT 14W TALL SH      |                 |
| 4    | ICJ3PDCMTS16  | J3 BK PDC MT 16W STD SH       |                 |
| 4    | ICJ3PDCMTT16  | J3 BK PDC MT 16W TALL SH      |                 |
| 4    | ICJ3PDCMTS18  | J3 BK PDC MT 18W STD SH       |                 |
| 4    | ICJ3PDCMTT18  | J3 BK PDC MT 18W TALL SH      |                 |
| 4    | ICJ3PDCMTS20  | J3 BK PDC MT 20W STD SH       |                 |
| 4    | ICJ3PDCMTT20  | J3 BK PDC MT 20W TALL SH      |                 |


Parts only, for complete Backrest see J3 Back orderform.



| Pos. | Item Number | Description                     | Pos. | Item Number | Description                     |
|------|-------------|---------------------------------|------|-------------|---------------------------------|
| 1A   | CJ3CFMSTS14 | COVER J3 CF BACK PA MT SHT 14W  | 3A   | SJ3CFMSTS14 | SHELL J3 CF BACK PA MT SHT 14W  |
| 1A   | CJ3CFMSTS16 | COVER J3 CF BACK PA MT SHT 16W  | 3A   | SJ3CFMSTS16 | SHELL J3 CF BACK PA MT SHT 16W  |
| 1A   | CJ3CFMSTS18 | COVER J3 CF BACK PA MT SHT 18W  | 3A   | SJ3CFMSTS18 | SHELL J3 CF BACK PA MT SHT 18W  |
| 1A   | CJ3CFMSTS20 | COVER J3 CF BACK PA MT SHT 20W  | 3A   | SJ3CFMSTS20 | SHELL J3 CF BACK PA MT SHT 20W  |
| 1B   | CJ3CFMTT14  | COVER J3 CF BACK PA MT TALL 14W | 3B   | SJ3CFMTT14  | SHELL J3 CF BACK PA MT TALL 14W |
| 1B   | CJ3CFMTT16  | COVER J3 CF BACK PA MT TALL 16W | 3B   | SJ3CFMTT16  | SHELL J3 CF BACK PA MT TALL 16W |
| 1B   | CJ3CFMTT18  | COVER J3 CF BACK PA MT TALL 18W | 3B   | SJ3CFMTT18  | SHELL J3 CF BACK PA MT TALL 18W |
| 1B   | CJ3CFMTT20  | COVER J3 CF BACK PA MT TALL 20W | 3B   | SJ3CFMTT20  | SHELL J3 CF BACK PA MT TALL 20W |
| 2A   | BJ3PAMTS14  | J3 FOAM, POS MT 14W STD         | 4A   | ICJ3PAMTS14 | J3BK INNER CVR PA MT 14W SHT    |
| 2A   | BJ3PAMTS16  | J3 FOAM, POS MT 16W STD         | 4A   | ICJ3PAMTS16 | J3BK INNER CVR PA MT 16W SHT    |
| 2A   | BJ3PAMTS18  | J3 FOAM, POS MT 18W STD         | 4A   | ICJ3PAMTS18 | J3BK INNER CVR PA MT 18W SHT    |
| 2A   | BJ3PAMTS20  | J3 FOAM, POS MT 20W STD         | 4A   | ICJ3PAMTS20 | J3BK INNER CVR PA MT 20W SHT    |
| 2B   | BJ3PAMTT14  | J3 FOAM, POS MT 14W TALL        | 4B   | ICJ3PAMTT14 | J3BK INNER CVR PA MT 14W TALL   |
| 2B   | BJ3PAMTT16  | J3 FOAM, POS MT 16W TALL        | 4B   | ICJ3PAMTT16 | J3BK INNER CVR PA MT 16W TALL   |
| 2B   | BJ3PAMTT18  | J3 FOAM, POS MT 18W TALL        | 4B   | ICJ3PAMTT18 | J3BK INNER CVR PA MT 18W TALL   |
| 2B   | BJ3PAMTT20  | J3 FOAM, POS MT 20W TALL        | 4B   | ICJ3PAMTT20 | J3BK INNER CVR PA MT 20W TALL   |

Parts only; for complete Backrest see J3 Back order form.

Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.




| Pos. | Item Number | Description             | Remarks |
|------|-------------|-------------------------|---------|
| 1A   | CJ3PAUTS12  | J3 CVR POS UT 12W STD   |         |
| 1A   | CJ3PAUTS14  | J3 CVR, POS UT 14W STD  |         |
| 1A   | CJ3PAUTM14  | J3 CVR, POS UT 14W MED  |         |
| 1A   | CJ3PAUTT14  | J3 CVR, POS UT 14W TALL |         |
| 1A   | CJ3PAUTS16  | J3 CVR, POS UT 16W STD  |         |
| 1A   | CJ3PAUTM16  | J3 CVR, POS UT 16W MED  |         |
| 1A   | CJ3PAUTT16  | J3 CVR, POS UT 16W TALL |         |
| 1A   | CJ3PAUTS18  | J3 CVR, POS UT 18W STD  |         |
| 1A   | CJ3PAUTM18  | J3 CVR, POS UT 18W MED  |         |
| 1A   | CJ3PAUTT18  | J3 CVR, POS UT 18W TALL |         |
| 1A   | CJ3PAUTS20  | J3 CVR, POS UT 20W STD  |         |

| Pos. | Item Number  | Description                   | Remarks         |
|------|--------------|-------------------------------|-----------------|
| 1A   | CJ3PAUTM20   | J3 CVR, POS UT 20W MED        |                 |
| 1A   | CJ3PAUTT20   | J3 CVR, POS UT 20W TALL       |                 |
| 1B   | CJ3PAUTS16CU | CVR J3 BK POS UT 16W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PAUTM16CU | CVR J3 BK POS UT 16W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PAUTT16CU | CVR J3 BK POS UT 16W TALL CU  | use w/ Cradle-U |
| 1B   | CJ3PAUTS18CU | CVR J3 BK POS UT 18W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PAUTM18CU | CVR J3 BK POS UT 18W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PAUTT18CU | CVR J3 BK POS UT 18W TALL CU  | use w/ Cradle-U |
| 1B   | CJ3PAUTS20CU | CVR J3 BK POS UT 20W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PAUTM20CU | CVR J3 BK POS UT 20W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PAUTT20CU | CVR J3 BK POS UT 20W TALL CU  | use w/ Cradle-U |
| 2A   | BJ3PAUTS12   | J3 FOAM POS UT 12W STD        |                 |
| 2A   | BJ3PAUTS14   | J3 FOAM, POS UT 14W STD       |                 |
| 2A   | BJ3PAUTM14   | J3 FOAM, POS UT 14W MED       |                 |
| 2A   | BJ3PAUTT14   | J3 FOAM, POS UT 14W TALL      |                 |
| 2A   | BJ3PAUTS16   | J3 FOAM, POS UT 16W STD       |                 |
| 2A   | BJ3PAUTM16   | J3 FOAM, POS UT 16W MED       |                 |
| 2A   | BJ3PAUTT16   | J3 FOAM, POS UT 16W TALL      |                 |
| 2A   | BJ3PAUTS18   | J3 FOAM, POS UT 18W STD       |                 |
| 2A   | BJ3PAUTM18   | J3 FOAM, POS UT 18W MED       |                 |
| 2A   | BJ3PAUTT18   | J3 FOAM, POS UT 18W TALL      |                 |
| 2A   | BJ3PAUTS20   | J3 FOAM, POS UT 20W STD       |                 |
| 2A   | BJ3PAUTM20   | J3 FOAM, POS UT 20W MED       |                 |
| 2A   | BJ3PAUTT20   | J3 FOAM, POS UT 20W TALL      |                 |
| 2B   | BJ3PAUTS16CU | BASE J3 BK POS UT 16W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PAUTM16CU | BASE J3 BK POS UT 16W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PAUTT16CU | BASE J3 BK POS UT 16W TALL CU | use w/ Cradle-U |
| 2B   | BJ3PAUTS18CU | BASE J3 BK POS UT 18W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PAUTM18CU | BASE J3 BK POS UT 18W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PAUTT18CU | BASE J3 BK POS UT 18W TALL CU | use w/ Cradle-U |
| 2B   | BJ3PAUTS20CU | BASE J3 BK POS UT 20W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PAUTM20CU | BASE J3 BK POS UT 20W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PAUTT20CU | BASE J3 BK POS UT 20W TALL CU | use w/ Cradle-U |
| 3    | SJ3PAUTS12   | J3 SHELL POS UT 12W STD       |                 |
| 3    | SJ3PAUTS14   | J3 SHELL, POS UT 14W STD      |                 |
| 3    | SJ3PAUTM14   | J3 SHELL, POS UT 14W MED      |                 |
| 3    | SJ3PAUTT14   | J3 SHELL, POS UT 14W TALL     |                 |
| 3    | SJ3PAUTS16   | J3 SHELL, POS UT 16W STD      |                 |
| 3    | SJ3PAUTM16   | J3 SHELL, POS UT 16W MED      |                 |
| 3    | SJ3PAUTT16   | J3 SHELL, POS UT 16W TALL     |                 |
| 3    | SJ3PAUTS18   | J3 SHELL, POS UT 18W STD      |                 |
| 3    | SJ3PAUTM18   | J3 SHELL, POS UT 18W MED      |                 |
| 3    | SJ3PAUTT18   | J3 SHELL, POS UT 18W TALL     |                 |
| 3    | SJ3PAUTS20   | J3 SHELL, POS UT 20W STD      |                 |
| 3    | SJ3PAUTM20   | J3 SHELL, POS UT 20W MED      |                 |
| 3    | SJ3PAUTT20   | J3 SHELL, POS UT 20W TALL     |                 |
| 4    | ICJ3PAUTS12  | J3 BK POS UT 12W STD SH       |                 |
| 4    | ICJ3PAUTS14  | J3 BK POS UT 14W STD SH       |                 |
| 4    | ICJ3PAUTM14  | J3 BK POS UT 14W MED SH       |                 |
| 4    | ICJ3PAUTT14  | J3 BK POS UT 14W TALL SH      |                 |
| 4    | ICJ3PAUTS16  | J3 BK POS UT 16W STD SH       |                 |
| 4    | ICJ3PAUTM16  | J3 BK POS UT 16W MED SH       |                 |
| 4    | ICJ3PAUTT16  | J3 BK POS UT 16W TALL SH      |                 |
| 4    | ICJ3PAUTS18  | J3 BK POS UT 18W STD SH       |                 |
| 4    | ICJ3PAUTM18  | J3 BK POS UT 18W MED SH       |                 |
| 4    | ICJ3PAUTT18  | J3 BK POS UT 18W TALL SH      |                 |
| 4    | ICJ3PAUTS20  | J3 BK POS UT 20W STD SH       |                 |

| Pos. | Item Number | Description              | Remarks |
|------|-------------|--------------------------|---------|
| 4    | ICJ3PAUTM20 | J3 BK POS UT 20W MED SH  |         |
| 4    | ICJ3PAUTT20 | J3 BK POS UT 20W TALL SH |         |

**NOTE:** Parts only, for complete Backrest see J3 Back order form.


Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.



| Pos. | Item Number | Description                 |
|------|-------------|-----------------------------|
| 1    | CJ3PAUTM22  | J3HD CVR POS UT 22W MED     |
| 1    | CJ3PAUTM24  | J3HD CVR POS UT 24W MED     |
| 2    | BJ3PAUTM22  | J3HD FOAM POS UT 22W MED    |
| 2    | BJ3PAUTM24  | J3HD FOAM POS UT 24W MED    |
| 3    | SJ3PAUTM22  | J3HD SHELL POS UT 22W MED   |
| 3    | SJ3PAUTM24  | J3HD SHELL POS UT 24W MED   |
| 4    | ICJ3PAUTM22 | J3 HD BK,POS UT 22W MED, FS |
| 4    | ICJ3PAUTM24 | J3 HD BK,POS UT 24W MED, FS |

Parts only; for complete Backrest see J3 Back order form.

Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.




| Pos. | Item Number | Description              | Remarks |
|------|-------------|--------------------------|---------|
| 1A   | CJ3PLUTS12  | J3 CVR SHAL UT 12W STD   |         |
| 1A   | CJ3PLUTS14  | J3 CVR, SHAL UT 14W STD  |         |
| 1A   | CJ3PLUTM14  | J3 CVR, SHAL UT 14W MED  |         |
| 1A   | CJ3PLUTT14  | J3 CVR, SHAL UT 14W TALL |         |
| 1A   | CJ3PLUTS16  | J3 CVR, SHAL UT 16W STD  |         |
| 1A   | CJ3PLUTM16  | J3 CVR, SHAL UT 16W MED  |         |
| 1A   | CJ3PLUTT16  | J3 CVR, SHAL UT 16W TALL |         |
| 1A   | CJ3PLUTS18  | J3 CVR, SHAL UT 18W STD  |         |
| 1A   | CJ3PLUTM18  | J3 CVR, SHAL UT 18W MED  |         |
| 1A   | CJ3PLUTT18  | J3 CVR, SHAL UT 18W TALL |         |
| 1A   | CJ3PLUTS20  | J3 CVR, SHAL UT 20W STD  |         |

| Pos. | Item Number  | Description                    | Remarks         |
|------|--------------|--------------------------------|-----------------|
| 1A   | CJ3PLUTM20   | J3 CVR, SHAL UT 20W MED        |                 |
| 1A   | CJ3PLUTT20   | J3 CVR, SHAL UT 20W TALL       |                 |
| 1B   | CJ3PLUTS16CU | CVR J3 BK SHAL UT 16W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PLUTM16CU | CVR J3 BK SHAL UT 16W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PLUTT16CU | CVR J3 BK SHAL UT 16W TALL CU  | use w/ Cradle-U |
| 1B   | CJ3PLUTS18CU | CVR J3 BK SHAL UT 18W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PLUTM18CU | CVR J3 BK SHAL UT 18W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PLUTT18CU | CVR J3 BK SHAL UT 18W TALL CU  | use w/ Cradle-U |
| 1B   | CJ3PLUTS20CU | CVR J3 BK SHAL UT 20W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PLUTM20CU | CVR J3 BK SHAL UT 20W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PLUTT20CU | CVR J3 BK SHAL UT 20W TALL CU  | use w/ Cradle-U |
| 2A   | BJ3PLUTS12   | J3 FOAM SHAL UT 12W STD        |                 |
| 2A   | BJ3PLUTS14   | J3 FOAM, SHAL UT 14W STD       |                 |
| 2A   | BJ3PLUTM14   | J3 FOAM, SHAL UT 14W MED       |                 |
| 2A   | BJ3PLUTT14   | J3 FOAM, SHAL UT 14W TALL      |                 |
| 2A   | BJ3PLUTS16   | J3 FOAM, SHAL UT 16W STD       |                 |
| 2A   | BJ3PLUTM16   | J3 FOAM, SHAL UT 16W MED       |                 |
| 2A   | BJ3PLUTT16   | J3 FOAM, SHAL UT 16W TALL      |                 |
| 2A   | BJ3PLUTS18   | J3 FOAM, SHAL UT 18W STD       |                 |
| 2A   | BJ3PLUTM18   | J3 FOAM, SHAL UT 18W MED       |                 |
| 2A   | BJ3PLUTT18   | J3 FOAM, SHAL UT 18W TALL      |                 |
| 2A   | BJ3PLUTS20   | J3 FOAM, SHAL UT 20W STD       |                 |
| 2A   | BJ3PLUTM20   | J3 FOAM, SHAL UT 20W MED       |                 |
| 2A   | BJ3PLUTT20   | J3 FOAM, SHAL UT 20W TALL      |                 |
| 2B   | BJ3PLUTS16CU | BASE J3 BK SHAL UT 16W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PLUTM16CU | BASE J3 BK SHAL UT 16W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PLUTT16CU | BASE J3 BK SHAL UT 16W TALL CU | use w/ Cradle-U |
| 2B   | BJ3PLUTS18CU | BASE J3 BK SHAL UT 18W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PLUTM18CU | BASE J3 BK SHAL UT 18W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PLUTT18CU | BASE J3 BK SHAL UT 18W TALL CU | use w/ Cradle-U |
| 2B   | BJ3PLUTS20CU | BASE J3 BK SHAL UT 20W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PLUTM20CU | BASE J3 BK SHAL UT 20W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PLUTT20CU | BASE J3 BK SHAL UT 20W TALL CU | use w/ Cradle-U |
| 3    | SJ3PLUTS12   | J3 SHELL SHAL UT 12W STD       |                 |
| 3    | SJ3PLUTS14   | J3 SHELL, SHAL UT 14W STD      |                 |
| 3    | SJ3PLUTM14   | J3 SHELL, SHAL UT 14W MED      |                 |
| 3    | SJ3PLUTT14   | J3 SHELL, SHAL UT 14W TALL     |                 |
| 3    | SJ3PLUTS16   | J3 SHELL, SHAL UT 16W STD      |                 |
| 3    | SJ3PLUTM16   | J3 SHELL, SHAL UT 16W MED      |                 |
| 3    | SJ3PLUTT16   | J3 SHELL, SHAL UT 16W TALL     |                 |
| 3    | SJ3PLUTS18   | J3 SHELL, SHAL UT 18W STD      |                 |
| 3    | SJ3PLUTM18   | J3 SHELL, SHAL UT 18W MED      |                 |
| 3    | SJ3PLUTT18   | J3 SHELL, SHAL UT 18W TALL     |                 |
| 3    | SJ3PLUTS20   | J3 SHELL, SHAL UT 20W STD      |                 |
| 3    | SJ3PLUTM20   | J3 SHELL, SHAL UT 20W MED      |                 |
| 3    | SJ3PLUTT20   | J3 SHELL, SHAL UT 20W TALL     |                 |
| 4    | ICJ3PLUTS12  | J3 BK SHAL UT 12W STD SH       |                 |
| 4    | ICJ3PLUTS14  | J3 BK SHAL UT 14W STD SH       |                 |
| 4    | ICJ3PLUTM14  | J3 BK SHAL UT 14W MED SH       |                 |
| 4    | ICJ3PLUTT14  | J3 BK SHAL UT 14W TALL SH      |                 |
| 4    | ICJ3PLUTS16  | J3 BK SHAL UT 16W STD SH       |                 |
| 4    | ICJ3PLUTM16  | J3 BK SHAL UT 16W MED SH       |                 |
| 4    | ICJ3PLUTT16  | J3 BK SHAL UT 16W TALL SH      |                 |
| 4    | ICJ3PLUTS18  | J3 BK SHAL UT 18W STD SH       |                 |
| 4    | ICJ3PLUTM18  | J3 BK SHAL UT 18W MED SH       |                 |
| 4    | ICJ3PLUTT18  | J3 BK SHAL UT 18W TALL SH      |                 |
| 4    | ICJ3PLUTS20  | J3 BK SHAL UT 20W STD SH       |                 |

| Pos. | Item Number | Description               | Remarks |
|------|-------------|---------------------------|---------|
| 4    | ICJ3PLUTM20 | J3 BK SHAL UT 20W MED SH  |         |
| 4    | ICJ3PLUTT20 | J3 BK SHAL UT 20W TALL SH |         |

Parts only; for complete Backrest see J3 Back order form.


Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.



| Pos. | Item Number | Description              | Pos. | Item Number | Description               |
|------|-------------|--------------------------|------|-------------|---------------------------|
| 1    | CJ3PDLUTS14 | J3 CVR, DEEP UT 14W STD  | 1    | CJ3PDLUTT20 | J3 CVR, DEEP UT 20W TALL  |
| 1    | CJ3PDLUTM14 | J3 CVR, DEEP UT 14W MED  | 2    | BJ3PDLUTS14 | J3 FOAM, DEEP UT 14W STD  |
| 1    | CJ3PDLUTT14 | J3 CVR, DEEP UT 14W TALL | 2    | BJ3PDLUTM14 | J3 FOAM, DEEP UT 14W MED  |
| 1    | CJ3PDLUTS16 | J3 CVR, DEEP UT 16W STD  | 2    | BJ3PDLUTT14 | J3 FOAM, DEEP UT 14W TALL |
| 1    | CJ3PDLUTM16 | J3 CVR, DEEP UT 16W MED  | 2    | BJ3PDLUTS16 | J3 FOAM, DEEP UT 16W STD  |
| 1    | CJ3PDLUTT16 | J3 CVR, DEEP UT 16W TALL | 2    | BJ3PDLUTM16 | J3 FOAM, DEEP UT 16W MED  |
| 1    | CJ3PDLUTS18 | J3 CVR, DEEP UT 18W STD  | 2    | BJ3PDLUTT16 | J3 FOAM, DEEP UT 16W TALL |
| 1    | CJ3PDLUTM18 | J3 CVR, DEEP UT 18W MED  | 2    | BJ3PDLUTS18 | J3 FOAM, DEEP UT 18W STD  |
| 1    | CJ3PDLUTT18 | J3 CVR, DEEP UT 18W TALL | 2    | BJ3PDLUTM18 | J3 FOAM, DEEP UT 18W MED  |
| 1    | CJ3PDLUTS20 | J3 CVR, DEEP UT 20W STD  | 2    | BJ3PDLUTT18 | J3 FOAM, DEEP UT 18W TALL |
| 1    | CJ3PDLUTM20 | J3 CVR, DEEP UT 20W MED  | 2    | BJ3PDLUTS20 | J3 FOAM, DEEP UT 20W STD  |

| Pos. | Item Number | Description                | Pos. | Item Number  | Description                |
|------|-------------|----------------------------|------|--------------|----------------------------|
| 2    | BJ3PDLUTM20 | J3 FOAM, DEEP UT 20W MED   | 3    | SJ3PDLUTT20  | J3 SHELL, DEEP UT 20W TALL |
| 2    | BJ3PDLUTT20 | J3 FOAM, DEEP UT 20W TALL  | 4    | ICJ3PDLUTS14 | J3 BK DEEP UT 14W STD SH   |
| 3    | SJ3PDLUTS14 | J3 SHELL, DEEP UT 14W STD  | 4    | ICJ3PDLUTM14 | J3 BK DEEP UT 14W MED SH   |
| 3    | SJ3PDLUTM14 | J3 SHELL, DEEP UT 14W MED  | 4    | ICJ3PDLUTT14 | J3 BK DEEP UT 14W TALL SH  |
| 3    | SJ3PDLUTT14 | J3 SHELL, DEEP UT 14W TALL | 4    | ICJ3PDLUTS16 | J3 BK DEEP UT 16W STD SH   |
| 3    | SJ3PDLUTS16 | J3 SHELL, DEEP UT 16W STD  | 4    | ICJ3PDLUTM16 | J3 BK DEEP UT 16W MED SH   |
| 3    | SJ3PDLUTM16 | J3 SHELL, DEEP UT 16W MED  | 4    | ICJ3PDLUTT16 | J3 BK DEEP UT 16W TALL SH  |
| 3    | SJ3PDLUTT16 | J3 SHELL, DEEP UT 16W TALL | 4    | ICJ3PDLUTS18 | J3 BK DEEP UT 18W STD SH   |
| 3    | SJ3PDLUTS18 | J3 SHELL, DEEP UT 18W STD  | 4    | ICJ3PDLUTM18 | J3 BK DEEP UT 18W MED SH   |
| 3    | SJ3PDLUTM18 | J3 SHELL, DEEP UT 18W MED  | 4    | ICJ3PDLUTT18 | J3 BK DEEP UT 18W TALL SH  |
| 3    | SJ3PDLUTT18 | J3 SHELL, DEEP UT 18W TALL | 4    | ICJ3PDLUTS20 | J3 BK DEEP UT 20W STD SH   |
| 3    | SJ3PDLUTS20 | J3 SHELL, DEEP UT 20W STD  | 4    | ICJ3PDLUTM20 | J3 BK DEEP UT 20W MED SH   |
| 3    | SJ3PDLUTM20 | J3 SHELL, DEEP UT 20W MED  | 4    | ICJ3PDLUTT20 | J3 BK DEEP UT 20W TALL SH  |


**Note:** parts only; for complete Backrest see J3 Back order form.



| Pos. | Item Number | Description                   |
|------|-------------|-------------------------------|
| 1    | CJ3PDUTM14  | J3 BK CVR PD UT 14W MED       |
| 1    | CJ3PDUTM16  | J3 BK CVR PD UT 16W MED       |
| 1    | CJ3PDUTM18  | J3 BK CVR PD UT 18W MED       |
| 1    | CJ3PDUTM20  | J3 BK CVR PD UT 20W MED       |
| 2    | BJ3PDUTM14  | J3 BK FOAM PD UT 14W MED      |
| 2    | BJ3PDUTM16  | J3 BK FOAM PD UT 16W MED      |
| 2    | BJ3PDUTM18  | J3 BK FOAM PD UT 18W MED      |
| 2    | BJ3PDUTM20  | J3 BK FOAM PD UT 20W MED      |
| 3    | ICJ3PDUTM14 | J3 BK INNER CVR PD UT 14W MED |
| 3    | ICJ3PDUTM16 | J3 BK INNER CVR PD UT 16W MED |
| 3    | ICJ3PDUTM18 | J3 BK INNER CVR PD UT 18W MED |
| 3    | ICJ3PDUTM20 | J3 BK INNER CVR PD UT 20W MED |
| 4    | SJ3PDUTM14  | J3 BK SHL PD UT 14W MED       |
| 4    | SJ3PDUTM16  | J3 BK SHL PD UT 16W MED       |
| 4    | SJ3PDUTM18  | J3 BK SHL PD UT 18W MED       |
| 4    | SJ3PDUTM20  | J3 BK SHL PD UT 20W MED       |

**NOTE:** Parts only, for complete Backrest see J3 Back order form.

Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.




| Pos. | Item Number   | Description                 | Remarks         |
|------|---------------|-----------------------------|-----------------|
| 1A   | CJ3PDCUTM14   | J3 CVR PDC UT 14W MED       |                 |
| 1A   | CJ3PDCUTM16   | J3 CVR PDC UT 16W MED       |                 |
| 1A   | CJ3PDCUTM18   | J3 CVR PDC UT 18W MED       |                 |
| 1A   | CJ3PDCUTM20   | J3 CVR PDC UT 20W MED       |                 |
| 1B   | CJ3PDCUTM16CU | CVR J3 BK PDC UT 16W MED CU | use w/ Cradle-U |
| 1B   | CJ3PDCUTM18CU | CVR J3 BK PDC UT 18W MED CU | use w/ Cradle-U |
| 1B   | CJ3PDCUTM20CU | CVR J3 BK PDC UT 20W MED CU | use w/ Cradle-U |
| 2A   | BJ3PDCUTM14   | J3 FOAM PDC UT 14W MED      |                 |
| 2A   | BJ3PDCUTM16   | J3 FOAM PDC UT 16W MED      |                 |
| 2A   | BJ3PDCUTM18   | J3 FOAM PDC UT 18W MED      |                 |
| 2A   | BJ3PDCUTM20   | J3 FOAM PDC UT 20W MED      |                 |

| Pos. | Item Number   | Description                  | Remarks         |
|------|---------------|------------------------------|-----------------|
| 2B   | BJ3PDCUTM16CU | BASE J3 BK PDC UT 16W MED CU | use w/ Cradle-U |
| 2B   | BJ3PDCUTM18CU | BASE J3 BK PDC UT 18W MED CU | use w/ Cradle-U |
| 2B   | BJ3PDCUTM20CU | BASE J3 BK PDC UT 20W MED CU | use w/ Cradle-U |
| 3    | SJ3PDCUTM14   | J3 SHELL PDC UT 14W MED      |                 |
| 3    | SJ3PDCUTM16   | J3 SHELL PDC UT 16W MED      |                 |
| 3    | SJ3PDCUTM18   | J3 SHELL PDC UT 18W MED      |                 |
| 3    | SJ3PDCUTM20   | J3 SHELL PDC UT 20W MED      |                 |
| 4    | ICJ3PDCUTM14  | J3 BK PDC UT 14W MED SH      |                 |
| 4    | ICJ3PDCUTM16  | J3 BK PDC UT 16W MED SH      |                 |
| 4    | ICJ3PDCUTM18  | J3 BK PDC UT 18W MED SH      |                 |
| 4    | ICJ3PDCUTM20  | J3 BK PDC UT 20W MED SH      |                 |

Parts only; for complete Backrest see J3 Back order form.

Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.




| Pos. | Item Number  | Description                 | Remarks         |
|------|--------------|-----------------------------|-----------------|
| 1A   | CJ3PASHS14   | J3 CVR, POS SH 14W STD      |                 |
| 1A   | CJ3PASHS16   | J3 CVR, POS SH 16W STD      |                 |
| 1A   | CJ3PASHM16   | J3 CVR, POS SH 16W MED      |                 |
| 1A   | CJ3PASHS18   | J3 CVR, POS SH 18W STD      |                 |
| 1A   | CJ3PASHM18   | J3 CVR, POS SH 18W MED      |                 |
| 1A   | CJ3PASHT18   | J3 CVR, POS SH 18W TALL     |                 |
| 1A   | CJ3PASHS20   | J3 CVR, POS SH 20W STD      |                 |
| 1A   | CJ3PASHM20   | J3 CVR, POS SH 20W MED      |                 |
| 1A   | CJ3PASHT20   | J3 CVR, POS SH 20W TALL     |                 |
| 1B   | CJ3PASHS16CU | CVR J3 BK POS SH 16W STD CU | use w/ Cradle-U |

| Pos. | Item Number  | Description                   | Remarks         |
|------|--------------|-------------------------------|-----------------|
| 1B   | CJ3PASHM16CU | CVR J3 BK POS SH 16W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PASHS18CU | CVR J3 BK POS SH 18W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PASHM18CU | CVR J3 BK POS SH 18W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PASHT18CU | CVR J3 BK POS SH 18W TALL CU  | use w/ Cradle-U |
| 1B   | CJ3PASHS20CU | CVR J3 BK POS SH 20W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PASHM20CU | CVR J3 BK POS SH 20W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PASHT20CU | CVR J3 BK POS SH 20W TALL CU  | use w/ Cradle-U |
| 2A   | BJ3PASHS14   | J3 FOAM, POS SH 14W STD       |                 |
| 2A   | BJ3PASHS16   | J3 FOAM, POS SH 16W STD       |                 |
| 2A   | BJ3PASHM16   | J3 FOAM, POS SH 16W MED       |                 |
| 2A   | BJ3PASHS18   | J3 FOAM, POS SH 18W STD       |                 |
| 2A   | BJ3PASHM18   | J3 FOAM, POS SH 18W MED       |                 |
| 2A   | BJ3PASHT18   | J3 FOAM, POS SH 18W TALL      |                 |
| 2A   | BJ3PASHS20   | J3 FOAM, POS SH 20W STD       |                 |
| 2A   | BJ3PASHM20   | J3 FOAM, POS SH 20W MED       |                 |
| 2A   | BJ3PASHT20   | J3 FOAM, POS SH 20W TALL      |                 |
| 2B   | BJ3PASHS16CU | BASE J3 BK POS SH 16W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PASHM16CU | BASE J3 BK POS SH 16W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PASHS18CU | BASE J3 BK POS SH 18W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PASHM18CU | BASE J3 BK POS SH 18W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PASHT18CU | BASE J3 BK POS SH 18W TALL CU | use w/ Cradle-U |
| 2B   | BJ3PASHS20CU | BASE J3 BK POS SH 20W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PASHM20CU | BASE J3 BK POS SH 20W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PASHT20CU | BASE J3 BK POS SH 20W TALL CU | use w/ Cradle-U |
| 3    | SJ3PASHS14   | J3 SHELL, POS SH 14W STD      |                 |
| 3    | SJ3PASHS16   | J3 SHELL, POS SH 16W STD      |                 |
| 3    | SJ3PASHM16   | J3 SHELL, POS SH 16W MED      |                 |
| 3    | SJ3PASHS18   | J3 SHELL, POS SH 18W STD      |                 |
| 3    | SJ3PASHM18   | J3 SHELL, POS SH 18W MED      |                 |
| 3    | SJ3PASHT18   | J3 SHELL, POS SH 18W TALL     |                 |
| 3    | SJ3PASHS20   | J3 SHELL, POS SH 20W STD      |                 |
| 3    | SJ3PASHM20   | J3 SHELL, POS SH 20W MED      |                 |
| 3    | SJ3PASHT20   | J3 SHELL, POS SH 20W TALL     |                 |
| 4    | ICJ3PASHS14  | J3 BK POS SH 14W STD SH       |                 |
| 4    | ICJ3PASHS16  | J3 BK POS SH 16W STD SH       |                 |
| 4    | ICJ3PASHM16  | J3 BK POS SH 16W MED SH       |                 |
| 4    | ICJ3PASHS18  | J3 BK POS SH 18W STD SH       |                 |
| 4    | ICJ3PASHM18  | J3 BK POS SH 18W MED SH       |                 |
| 4    | ICJ3PASHT18  | J3 BK POS SH 18W TALL SH      |                 |
| 4    | ICJ3PASHS20  | J3 BK POS SH 20W STD SH       |                 |
| 4    | ICJ3PASHM20  | J3 BK POS SH 20W MED SH       |                 |
| 4    | ICJ3PASHT20  | J3 BK POS SH 20W TALL SH      |                 |

Parts only; for complete Backrest see J3 Back order form.

Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.




| Pos. | Item Number  | Description                  | Remarks         |
|------|--------------|------------------------------|-----------------|
| 1A   | CJ3PLSHS14   | J3 CVR, SHAL SH 14W STD      |                 |
| 1A   | CJ3PLSHS16   | J3 CVR, SHAL SH 16W STD      |                 |
| 1A   | CJ3PLSHM16   | J3 CVR, SHAL SH 16W MED      |                 |
| 1A   | CJ3PLSHS18   | J3 CVR, SHAL SH 18W STD      |                 |
| 1A   | CJ3PLSHM18   | J3 CVR, SHAL SH 18W MED      |                 |
| 1A   | CJ3PLSHT18   | J3 CVR, SHAL SH 18W TALL     |                 |
| 1A   | CJ3PLSHS20   | J3 CVR, SHAL SH 20W TALL     |                 |
| 1A   | CJ3PLSHM20   | J3 CVR, SHAL SH 20W MED      |                 |
| 1A   | CJ3PLSHT20   | J3 CVR, SHAL SH 20W TALL     |                 |
| 1B   | CJ3PLSHS16CU | CVR J3 BK SHAL SH 16W STD CU | use w/ Cradle-U |

| Pos. | Item Number  | Description                    | Remarks         |
|------|--------------|--------------------------------|-----------------|
| 1B   | CJ3PLSHM16CU | CVR J3 BK SHAL SH 16W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PLSHS18CU | CVR J3 BK SHAL SH 18W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PLSHM18CU | CVR J3 BK SHAL SH 18W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PLSHT18CU | CVR J3 BK SHAL SH 18W TALL CU  | use w/ Cradle-U |
| 1B   | CJ3PLSHS20CU | CVR J3 BK SHAL SH 20W STD CU   | use w/ Cradle-U |
| 1B   | CJ3PLSHM20CU | CVR J3 BK SHAL SH 20W MED CU   | use w/ Cradle-U |
| 1B   | CJ3PLSHT20CU | CVR J3 BK SHAL SH 20W TALL CU  | use w/ Cradle-U |
| 2A   | BJ3PLSHS14   | J3 FOAM, SHAL SH 14W STD       |                 |
| 2A   | BJ3PLSHS16   | J3 FOAM, SHAL SH 16W STD       |                 |
| 2A   | BJ3PLSHM16   | J3 FOAM, SHAL SH 16W MED       |                 |
| 2A   | BJ3PLSHS18   | J3 FOAM, SHAL SH 18W STD       |                 |
| 2A   | BJ3PLSHM18   | J3 FOAM, SHAL SH 18W MED       |                 |
| 2A   | BJ3PLSHT18   | J3 FOAM, SHAL SH 18W TALL      |                 |
| 2A   | BJ3PLSHS20   | J3 FOAM, SHAL SH 20W STD       |                 |
| 2A   | BJ3PLSHM20   | J3 FOAM, SHAL SH 20W MED       |                 |
| 2A   | BJ3PLSHT20   | J3 FOAM, SHAL SH 20W TALL      |                 |
| 2B   | BJ3PLSHS16CU | BASE J3 BK SHAL SH 16W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PLSHM16CU | BASE J3 BK SHAL SH 16W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PLSHS18CU | BASE J3 BK SHAL SH 18W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PLSHM18CU | BASE J3 BK SHAL SH 18W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PLSHT18CU | BASE J3 BK SHAL SH 18W TALL CU | use w/ Cradle-U |
| 2B   | BJ3PLSHS20CU | BASE J3 BK SHAL SH 20W STD CU  | use w/ Cradle-U |
| 2B   | BJ3PLSHM20CU | BASE J3 BK SHAL SH 20W MED CU  | use w/ Cradle-U |
| 2B   | BJ3PLSHT20CU | BASE J3 BK SHAL SH 20W TALL CU | use w/ Cradle-U |
| 3    | SJ3PLSHS14   | J3 SHELL, SHAL SH 14W STD      |                 |
| 3    | SJ3PLSHS16   | J3 SHELL, SHAL SH 16W STD      |                 |
| 3    | SJ3PLSHM16   | J3 SHELL, SHAL SH 16W MED      |                 |
| 3    | SJ3PLSHS18   | J3 SHELL, SHAL SH 18W STD      |                 |
| 3    | SJ3PLSHM18   | J3 SHELL, SHAL SH 18W MED      |                 |
| 3    | SJ3PLSHT18   | J3 SHELL, SHAL SH 18W TALL     |                 |
| 3    | SJ3PLSHS20   | J3 SHELL, SHAL SH 20W STD      |                 |
| 3    | SJ3PLSHM20   | J3 SHELL, SHAL SH 20W MED      |                 |
| 3    | SJ3PLSHT20   | J3 SHELL, SHAL SH 20W TALL     |                 |
| 4    | ICJ3PLSHS14  | J3 BK SHAL SH 14W STD SH       |                 |
| 4    | ICJ3PLSHS16  | J3 BK SHAL SH 16W STD SH       |                 |
| 4    | ICJ3PLSHM16  | J3 BK SHAL SH 16W MED SH       |                 |
| 4    | ICJ3PLSHS18  | J3 BK SHAL SH 18W STD SH       |                 |
| 4    | ICJ3PLSHM18  | J3 BK SHAL SH 18W MED SH       |                 |
| 4    | ICJ3PLSHT18  | J3 BK SHAL SH 18W TALL SH      |                 |
| 4    | ICJ3PLSHS20  | J3 BK SHAL SH 20W STD SH       |                 |
| 4    | ICJ3PLSHM20  | J3 BK SHAL SH 20W MED SH       |                 |
| 4    | ICJ3PLSHT20  | J3 BK SHAL SH 20W TALL SH      |                 |

Parts only; for complete Backrest see J3 Back order form.

Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.




| Pos. | Item Number | Description              | Pos. | Item Number | Description               |
|------|-------------|--------------------------|------|-------------|---------------------------|
| 1    | CJ3PDLSHS14 | J3 CVR, DEEP SH 14W STD  | 1    | CJ3PDLSHT20 | J3 CVR, DEEP SH 20W TALL  |
| 1    | CJ3PDLSHS16 | J3 CVR, DEEP SH 16W STD  | 2    | BJ3PDLSHS14 | J3 FOAM, DEEP SH 14W STD  |
| 1    | CJ3PDLSHM16 | J3 CVR, DEEP SH 16W MED  | 2    | BJ3PDLSHS16 | J3 FOAM, DEEP SH 16W STD  |
| 1    | CJ3PDLSHS18 | J3 CVR, DEEP SH 18W STD  | 2    | BJ3PDLSHM16 | J3 FOAM, DEEP SH 16W MED  |
| 1    | CJ3PDLSHM18 | J3 CVR, DEEP SH 18W MED  | 2    | BJ3PDLSHS18 | J3 FOAM, DEEP SH 18W STD  |
| 1    | CJ3PDLSHT18 | J3 CVR, DEEP SH 18W TALL | 2    | BJ3PDLSHM18 | J3 FOAM, DEEP SH 18W MED  |
| 1    | CJ3PDLSHS20 | J3 CVR, DEEP SH 20W STD  | 2    | BJ3PDLSHT18 | J3 FOAM, DEEP SH 18W TALL |
| 1    | CJ3PDLSHM20 | J3 CVR, DEEP SH 20W MED  | 2    | BJ3PDLSHS20 | J3 FOAM, DEEP SH 20W STD  |

| Pos. | Item Number | Description                | Pos. | Item Number  | Description                |
|------|-------------|----------------------------|------|--------------|----------------------------|
| 2    | BJ3PDLSHM20 | J3 FOAM, DEEP SH 20W MED   | 3    | SJ3PDLSHT20  | J3 SHELL, DEEP SH 20W TALL |
| 2    | BJ3PDLSHT20 | J3 FOAM, DEEP SH 20W TALL  | 4    | ICJ3PDLSHS14 | J3 BK DEEP SH 14W STD SH   |
| 3    | SJ3PDLSHS14 | J3 SHELL, DEEP SH 14W STD  | 4    | ICJ3PDLSHS16 | J3 BK DEEP SH 16W STD SH   |
| 3    | SJ3PDLSHS16 | J3 SHELL, DEEP SH 16W STD  | 4    | ICJ3PDLSHM16 | J3 BK DEEP SH 16W MED SH   |
| 3    | SJ3PDLSHM16 | J3 SHELL, DEEP SH 16W MED  | 4    | ICJ3PDLSHS18 | J3 BK DEEP SH 18W STD SH   |
| 3    | SJ3PDLSHS18 | J3 SHELL, DEEP SH 18W STD  | 4    | ICJ3PDLSHM18 | J3 BK DEEP SH 18W MED SH   |
| 3    | SJ3PDLSHM18 | J3 SHELL, DEEP SH 18W MED  | 4    | ICJ3PDLSHT18 | J3 BK DEEP SH 18W TALL SH  |
| 3    | SJ3PDLSHT18 | J3 SHELL, DEEP SH 18W TALL | 4    | ICJ3PDLSHS20 | J3 BK DEEP SH 20W STD SH   |
| 3    | SJ3PDLSHS20 | J3 SHELL, DEEP SH 20W STD  | 4    | ICJ3PDLSHM20 | J3 BK DEEP SH 20W MED SH   |
| 3    | SJ3PDLSHM20 | J3 SHELL, DEEP SH 20W MED  | 4    | ICJ3PDLSHT20 | J3 BK DEEP SH 20W TALL SH  |


**NOTE:** Parts only; for complete Backrest see J3 Back order form.

Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.



| Pos. | Item Number | Description                 |
|------|-------------|-----------------------------|
| 1    | CJ3PASHS22  | J3HD CVR POS SH 22W STD     |
| 2    | BJ3PASHS22  | J3HD FOAM POS SH 22W STD    |
| 3    | SJ3PASHS22  | J3HD SHELL POS SH 22W STD   |
| 4    | ICJ3PASHS22 | J3 HD BK,POS SH 22W STD, FS |


**Note:** parts only, for complete Backrest see J3 Back order form.



| Pos. | Item Number | Description                   |
|------|-------------|-------------------------------|
| 1    | CJ3PDSHS14  | J3 BK CVR PD SH 14W SHT       |
| 1    | CJ3PDSHS16  | J3 BK CVR PD SH 16W SHT       |
| 1    | CJ3PDSHS18  | J3 BK CVR PD SH 18W SHT       |
| 1    | CJ3PDSHS20  | J3 BK CVR PD SH 20W SHT       |
| 2    | BJ3PDSHS14  | J3 BK FOAM PD SH 14W SHT      |
| 2    | BJ3PDSHS16  | J3 BK FOAM PD SH 16W SHT      |
| 2    | BJ3PDSHS18  | J3 BK FOAM PD SH 18W SHT      |
| 2    | BJ3PDSHS20  | J3 BK FOAM PD SH 20W SHT      |
| 3    | ICJ3PDSHS14 | J3 BK INNER CVR PD SH 14W SHT |
| 3    | ICJ3PDSHS16 | J3 BK INNER CVR PD SH 16W SHT |
| 3    | ICJ3PDSHS18 | J3 BK INNER CVR PD SH 18W SHT |
| 3    | ICJ3PDSHS20 | J3 BK INNER CVR PD SH 20W SHT |
| 4    | SJ3PDSHS14  | J3 BK SHL PD SH 14W SHT       |
| 4    | SJ3PDSHS16  | J3 BK SHL PD SH 16W SHT       |
| 4    | SJ3PDSHS18  | J3 BK SHL PD SH 18W SHT       |
| 4    | SJ3PDSHS20  | J3 BK SHL PD SH 20W SHT       |

**NOTE:** Parts only; for complete Backrest see J3 Back order form.


Note: If replacing Foam in Jay 3 Backs purchased before August 6, 2013, please order new Cover.



| Pos. | Item Number   | Description                  | Remarks         |
|------|---------------|------------------------------|-----------------|
| 1A   | CJ3PDCSHS16   | J3 CVR PDC SH 16W STD        |                 |
| 1A   | CJ3PDCSHS18   | J3 CVR PDC SH 18W STD        |                 |
| 1A   | CJ3PDCSHS20   | J3 CVR PDC SH 20W STD        |                 |
| 1B   | CJ3PDCSHS16CU | CVR J3 BK PDC SH 16W STD CU  | use w/ Cradle-U |
| 1B   | CJ3PDCSHS18CU | CVR J3 BK PDC SH 18W STD CU  | use w/ Cradle-U |
| 1B   | CJ3PDCSHS20CU | CVR J3 BK PDC SH 20W STD CU  | use w/ Cradle-U |
| 2A   | BJ3PDCSHS16   | J3 FOAM PDC SH 16W STD       |                 |
| 2A   | BJ3PDCSHS18   | J3 FOAM PDC SH 18W STD       |                 |
| 2A   | BJ3PDCSHS20   | J3 FOAM PDC SH 20W STD       |                 |
| 2B   | BJ3PDCSHS16CU | BASE J3 BK PDC SH 16W STD CU | use w/ Cradle-U |
| 2B   | BJ3PDCSHS18CU | BASE J3 BK PDC SH 18W STD CU | use w/ Cradle-U |
| 2B   | BJ3PDCSHS20CU | BASE J3 BK PDC SH 20W STD CU | use w/ Cradle-U |
| 3    | SJ3PDCSHS16   | J3 SHELL PDC SH 16W STD      |                 |
| 3    | SJ3PDCSHS18   | J3 SHELL PDC SH 18W STD      |                 |
| 3    | SJ3PDCSHS20   | J3 SHELL PDC SH 20W STD      |                 |
| 4    | ICJ3PDCSHS16  | J3 BK PDC SH 16W STD SH      |                 |
| 4    | ICJ3PDCSHS18  | J3 BK PDC SH 18W STD SH      |                 |
| 4    | ICJ3PDCSHS20  | J3 BK PDC SH 20W STD SH      |                 |

**J3 FASTENER KIT GUIDE**

| Item # | Illustration | Description                | Qty |
|--------|--------------|----------------------------|-----|
| 1      |              | WASHER FLAT M6             | 4   |
| 2      |              | SWIVEL PIN THREADED        | 2   |
| 3      |              | WASHER FLAT M5 SS          | 4   |
| 4      |              | SHCS M5 x 0.8 x 25 FULL NP | 4   |
|        |              | SHCS M5 x 0.8 x 20 FULL NP | 4   |
| 5      |              | SHCS M6 x 1 x 22 BLACK     | 4   |
| 6      |              | NUT HEX NYLOCK M6 BLACK    | 8   |
| 7      |              | TWO HOLE WASHER            | 2   |
| 15     |              | HHCS M6 x 1M16 PL ROHS     | 4   |




| Pos.                  | Item Number | Description                    | Remarks                        |
|-----------------------|-------------|--------------------------------|--------------------------------|
| 5                     | J3HDWRM     | J3 BK RECEIVER MNTG HDW        | Pack of 4                      |
| 1-7,15                | J3HDWF      | J3 FASTENER KIT                |                                |
| 8B                    | J3HDWNT     | J3B NUT PLATE TALL             | Tall                           |
| 8A                    | J3HDWN      | J3B NUT PLATE                  | Short-Lower Thoracic Back only |
| 10                    | J3LOCKPIN   | TRIGGER REL PIN TS/QT/P5       |                                |
| 11                    | J3HDWLRR    | LATCH RECEIVER ASSM, RIGHT     |                                |
| 11                    | J3HDWLRL    | LATCH RECEIVER ASSM, LEFT      |                                |
| 12B,13                | J3HDWCL     | VARIABLE CLAMP, LONG           |                                |
| 7,9,15                | J3HDWPR     | J3B-MOUNTING PIN ASSEMBLY R    |                                |
| 7,9,15                | J3HDWPL     | J3B-MOUNTING PIN ASSEMBLY L    |                                |
| 2-5,10,11,12A,13      | J3SH        | J3 STANDARD HARDWARE           |                                |
| 2-5,10,11,12B,13      | J3AH        | J3 ADJUSTABLE HARDWARE         |                                |
| 1,8B,14B,6            | J3HDWTR     | J3B-MOUNTING BRACKET, TALL, R  |                                |
| 1,8B,14B,6            | J3HDWTL     | J3B-MOUNTING BRACKET, TALL, L  |                                |
| 1,8A,14A,6            | J3HDWSR     | J3B-MOUNTING BRACKET, SHORT, R | Lower Thoracic Back only       |
| 1,8A,14A,6            | J3HDWSL     | J3B-MOUNTING BRACKET, SHORT, L | Lower Thoracic Back only       |
| 1-7 ,8B-11,12A,13,14A | 117871      | J3 MNT STD TAL RETRO KIT       |                                |
| 1-7 ,8B-11,12B,13,14A | 117872      | J3 MNT ADJ TAL RETRO KIT       |                                |


| Pos.                  | Item Number | Description              | Remarks |
|-----------------------|-------------|--------------------------|---------|
| 1-7 ,8A-11,12A,13,14B | 117873      | J3 MNT STD SHT RETRO KIT |         |
| 1-7 ,8A-11,12B,13,14B | 117874      | J3 MNT ADJ SHT RETRO KIT |         |

**J3 FASTENER KIT GUIDE**

| Item # | Illustration | Description                | Qty |
|--------|--------------|----------------------------|-----|
| 1      |              | WASHER FLAT M6             | 4   |
| 2      |              | SWIVEL PIN THREADED        | 2   |
| 3      |              | WASHER FLAT M5 SS          | 4   |
| 4      |              | SHCS M5 x 0.8 x 25 FULL NP | 4   |
| 5      |              | SHCS M6 x 1 x 22 BLACK     | 4   |
| 6      |              | NUT HEX NYLOCK M6 BLACK    | 8   |
| 7      |              | TWO HOLE WASHER            | 2   |
| 15     |              | HHCS M6 x 1M16 PL ROHS     | 4   |




| Pos.          | Item Number | Description                 | Remarks   |
|---------------|-------------|-----------------------------|-----------|
| 5             | J3HDWRM     | J3 BK RECIEVER MNTG HDW     | PACK OF 4 |
| 8             | J3HDWNT     | J3B NUT PLATE TALL          |           |
| 10            | J3LOCKPIN   | TRIGGER REL PIN TS/QT/P5    |           |
| 11            | J3HDWLRR    | LATCH RECEIVER ASSM, RIGHT  | RIGHT     |
| 11            | J3HDWLRL    | LATCH RECEIVER ASSM, LEFT   | LEFT      |
| 19            | 115116      | J3 HD BK LOWER RECEIVER LFT | LEFT      |
| 19            | 115117      | J3 HD BK LOWER RECEIVER RT  | RIGHT     |
| 1-7,15-18     | 117864      | 4 POINT FASTNER KIT         |           |
| 7,9,15        | J3HDWPR     | J3B-Mounting Pin AssblyR    |           |
| 7,9,15        | J3HDWPL     | J3B-Mounting Pin AssblyL    |           |
| 2-4,12A,13,19 | 117865      | J3 LWR CLMP SHT KIT (PR)    | PAIR      |
| 2-4,12B,13,19 | 117866      | J3 LWR CLMP ADJ KIT (PR)    | PAIR      |
| 1,6,8,14      | 117867      | 4 POINT MNT BRKT TAL R      | RIGHT     |
| 1,6,8,14      | 117868      | 4 POINT MNT BRKT TAL L      | LEFT      |
| 1-12A,13-19   | 117869      | 4 PT MNT STD RETRO KIT      |           |
| 1-12B,13-19   | 117870      | 4 PT MNT ADJ RETRO KIT      |           |



| Pos.           | Item Number | Description            | Remarks               |
|----------------|-------------|------------------------|-----------------------|
| 1A             | 117856      | LT MNT 3/4 BUSHING KIT | To fit 3/4" tubing.   |
| 1B             | 117857      | LT MNT 7/8 BUSHING KIT | To fit 7/8" tubing.   |
| 2-4            | 117861      | LT MNT CLAMP ASSM      | Fits 1" tubing.       |
| 3,4            | 117860      | LT MNT CLAMP HARDWARE  |                       |
| 5A-7           | 117858      | LT MNT HARDWARE SHT    | Used w/Lower Thoracic |
| 5B-7           | 117859      | LT MNT HARDWARE LNG    | Used w/Mid Thoracic   |
| 1A,1B,2-4,5A-7 | 117862      | LT MNT SHT RETRO KIT   | Used w/Lower Thoracic |
| 1A,1B,2-4,5B-7 | 117863      | LT MNT TAL RETRO KIT   | Used w/Mid Thoracic   |

**J3 FASTENER KIT GUIDE**


| Item # | Illustration | Description                 | Qty |
|--------|--------------|-----------------------------|-----|
| 1      |              | WASHER FLAT M6              | 4   |
| 2      |              | BARREL NUT RH J3 BACK       | 2   |
| 3      |              | WASHER FLAT M5 SS           | 4   |
| 4A     |              | SHCS M5 x 0.8 x 30 PART NP  | 4   |
| 4B     |              | SHCS M5 x .8 x 20 INORGANIC | 4   |
| 5      |              | NUT HEX NYLOCK M6 BLACK     | 8   |
| 6      |              | TWO HOLE WASHER             | 2   |
| 11     |              | HHCS M6 x 1 x 12 PL ROHS    | 4   |
| 14     |              | LW SHCS M6x1x22 PL BLK IO   | 4   |
| 17     |              | HHCS M6 x 1 x 16 PL ROHS    | 4   |



| Pos.             | Item Number | Description                     | Remarks     |
|------------------|-------------|---------------------------------|-------------|
| 12               | J3LOCKPIN   | TRIGGER REL PIN TS/QT/P5        |             |
| 13,14            | 159398      | J3 RH LATCH RECEIVER ASSM RIGHT | RIGHT       |
| 13,14            | 159399      | J3 RH LATCH RECEIVER ASSM LEFT  | LEFT        |
| 1,5,10           | 159397      | J3 RH AND RHF NUT PLATE         | 1 side      |
| 1,5,9,10         | 159404      | J3 RHF MOUNTING BRACKET         | 1 side      |
| 1,5,10,16        | 159402      | J3 RH MOUNTING BRACKET          | 1 side      |
| 6,15,17          | 159400      | J3 MOUNTING PIN ASSM RIGHT      | RIGHT       |
| 6,15,17          | 159401      | J3 MOUNTING PIN ASSM LEFT       | LEFT        |
| 1-6,7,8,10,12-17 | 159395      | J3 RIGID HARDWARE (RH)          | PAIR        |
| 1-6,14,17        | 159396      | J3 RH FASTENER KIT              | for 2 sides |
| 1-11             | 159403      | RIGID HARDWARE FIXED (RHF)      | Pair        |


**J3 FASTENER KIT GUIDE**

| Item # | Illustration | Description               | Qty |
|--------|--------------|---------------------------|-----|
| 1      |              | WASHER FLAT M6            | 4   |
| 2      |              | BARREL NUT TH 13 BLACK    | 2   |
| 3      |              | WASHER FLAT M5 SS         | 4   |
| 4      |              | SHCS M5X0.8X30 PART NP    | 4   |
| 5      |              | SHCS M5X 8X20 INORGANIC   | 4   |
| 6      |              | NUT HEX NYLOCK M6 BLACK   | 8   |
| 7      |              | WASHER PLATE TILT J3 BACK | 2   |
| 15     |              | HHCS M6X1M16 PL ROHS      | 4   |




| Pos.       | Item Number | Description                          | Remarks                    |
|------------|-------------|--------------------------------------|----------------------------|
| 10         | J3LOCKPIN   | TRIGGER REL PIN TS/QT/P5             |                            |
| 1-7,15,16  | 159386      | J3 TH FASTENER KIT                   | for 2 sides                |
| 1,6,8      | 159387      | J3 TH NUT PLATE                      | for 1 side                 |
| 5,11       | 159389      | J3 TH LATCH RECEIVER ASSM RIGHT      | RIGHT                      |
| 5,11       | 159390      | J3 TH LATCH RECEIVER ASSM LEFT       | LEFT                       |
| 7,9,15     | 159391      | J3 TH MOUNTING PIN ASSM              | for 1 side                 |
| 1,6,8,14   | 159392      | J3 TH MOUNTING BRACKET               | for 1 side                 |
| 3-5,12B    | 159393      | J3 TH EXTENDED MOUNTING BRACKET      | for 2 sides                |
| 3-5,12B,13 | 159394      | J3 TH EXTENDED MOUNTING BRACKET ASSM | for 2 sides                |
| 1-16       | 159385      | J3 TILT HARDWARE (TH)                | Complete assm left & right |

Note: When retro-fitting the Jay Cradle U on an existing J3 Back please replace the back cover and back foam as the upholstery is wider with this option to ensure lateral plates are covered.




| Pos. | Item Number | Description                  | Remarks                        |
|------|-------------|------------------------------|--------------------------------|
| 1    | CRADLEU16   | CRADLE U ASSEMBLY 16 WIDE    | complete kit                   |
| 1    | CRADLEU18   | CRADLE U ASSEMBLY 18 WIDE    | complete kit                   |
| 1    | CRADLEU20   | CRADLE U ASSEMBLY 20 WIDE    | complete kit                   |
| 2    | 126254      | CRADLE U PADDLES W/HDWR PAIR | set of pads w/hardware 2 sides |
| 3    | 126255      | CRADLE U HARDWARE PAIR       | set of hardware 2 sides        |



| Pos.     | Item Number | Description               |
|----------|-------------|---------------------------|
| 1A       | CJ3LAT4x4P  | Lateral Cover 4x4         |
| 1B       | CJ3LAT4x5P  | Lateral Cover 4x5         |
| 1C       | CJ3LAT4x6P  | Lateral Cover 4x6         |
| 1D       | CJ3LAT6x4P  | Lateral Cover 6x4         |
| 1E       | CJ3LAT6x5P  | Lateral Cover 6x5         |
| 1F       | CJ3LAT6x6P  | Lateral Cover 6x6         |
| 2A       | BJ3LAT4x4P  | Lateral Base 4x4          |
| 2B       | BJ3LAT4x5P  | Lateral Base 4x5          |
| 2C       | BJ3LAT4x6P  | Lateral Base 4x6          |
| 2D       | BJ3LAT6x4P  | Lateral Base 6x4          |
| 2E       | BJ3LAT6x5P  | Lateral Base 6x5          |
| 2F       | BJ3LAT6x6P  | Lateral Base 6x6          |
| 7B       | PRT2630     | 1" OFF/ONE PC PAD MNT STD |
| 8        | 117876S     | JAY CF LAT MNT BRKT       |
| 3-7A,8,9 | 118045      | J3 LATERAL S-A BRACKET RH |
| 3-7A,8,9 | 118046      | J3 LATERAL S-A BRACKET LH |
| 3,5,6,9  | 118047      | J3 S/A-LAT HDW KIT        |

(10/2010)


Note: Effective after September 2010.




| Pos.  | Item Number | Description              | Remarks    |
|-------|-------------|--------------------------|------------|
| 1A    | CJ3LAT4x4P  | Lateral Cover 4x4        |            |
| 1B    | CJ3LAT4x5P  | Lateral Cover 4x5        |            |
| 1C    | CJ3LAT4x6P  | Lateral Cover 4x6        |            |
| 1D    | CJ3LAT6x4P  | Lateral Cover 6x4        |            |
| 1E    | CJ3LAT6x5P  | Lateral Cover 6x5        |            |
| 1F    | CJ3LAT6x6P  | Lateral Cover 6x6        |            |
| 2A    | BJ3LAT4x4P  | Lateral Base 4x4         |            |
| 2B    | BJ3LAT4x5P  | Lateral Base 4x5         |            |
| 2C    | BJ3LAT4x6P  | Lateral Base 4x6         |            |
| 2D    | BJ3LAT6x4P  | Lateral Base 6x4         |            |
| 2E    | BJ3LAT6x5P  | Lateral Base 6x5         |            |
| 2F    | BJ3LAT6x6P  | Lateral Base 6x6         |            |
| 3-6   | J3LATSABL   | J3Lateral-S/A Bracket LH | LEFT SIDE  |
| 3-6   | J3LATSABR   | J3Lateral-S/A Bracket RH | RIGHT SIDE |
| 3,5,6 | J3LATSAP    | J3 S/A-Lat. Bag-o-Bolts  |            |

Note: Available prior to September 2010.


## MEASURING GUIDE




| DESCRIPTION       | DIM A | DIM B |
|-------------------|-------|-------|
| Lateral Cover 4x4 | 4"    | 4"    |
| Lateral Cover 4x5 | 4"    | 5"    |
| Lateral Cover 4x6 | 4"    | 6"    |
| Lateral Cover 6x4 | 6"    | 4"    |
| Lateral Cover 6x5 | 6"    | 5"    |
| Lateral Cover 6x6 | 6"    | 6"    |
| Lateral Base 4x4  | 4"    | 4"    |
| Lateral Base 4x5  | 4"    | 5"    |
| Lateral Base 4x6  | 4"    | 6"    |
| Lateral Base 6x4  | 6"    | 4"    |
| Lateral Base 6x5  | 6"    | 5"    |
| Lateral Base 6x6  | 6"    | 6"    |




| Pos.  | Item Number | Description               |
|-------|-------------|---------------------------|
| 1A    | CJ3LAT4x4P  | Lateral Cover 4x4         |
| 1B    | CJ3LAT4x5P  | Lateral Cover 4x5         |
| 1C    | CJ3LAT4x6P  | Lateral Cover 4x6         |
| 1D    | CJ3LAT6x4P  | Lateral Cover 6x4         |
| 1E    | CJ3LAT6x5P  | Lateral Cover 6x5         |
| 1F    | CJ3LAT6x6P  | Lateral Cover 6x6         |
| 2A    | BJ3LAT4x4P  | Lateral Base 4x4          |
| 2B    | BJ3LAT4x5P  | Lateral Base 4x5          |
| 2C    | BJ3LAT4x6P  | Lateral Base 4x6          |
| 2D    | BJ3LAT6x4P  | Lateral Base 6x4          |
| 2E    | BJ3LAT6x5P  | Lateral Base 6x5          |
| 2F    | BJ3LAT6x6P  | Lateral Base 6x6          |
| 3-6   | J3LATFB     | J3Lateral-Fxd Bracket     |
| 3,5,6 | J3LATFP     | J3 Fixed-Lat. Bag-o-Bolts |




| Pos. | Item Number | Description                    |
|------|-------------|--------------------------------|
| 1A   | ATS9401     | COMFORFIT, CENTER OPENING, SM  |
| 1B   | ATS9402     | COMFORFIT, CENTER OPENING, MED |
| 1C   | ATS9403     | COMFORFIT, CENTER OPENING, LG  |
| 1D   | ATS9404     | COMFORFIT, CENTER OPENING, XL  |




| Pos. | Item Number | Description                    |
|------|-------------|--------------------------------|
| 1A   | ATS4000     | PADDED ANTERIOR TRUNK STRAP SM |
| 1B   | ATS4100     | PADDED ANTERIOR TRUNK STRAP MD |
| 1C   | ATS4200     | PADDED ANTERIOR TRUNK STRAP LG |
| 1D   | ATS4300     | PADDED ANTERIOR TRUNK SPRAP XL |



| Pos. | Item Number | Description                    |
|------|-------------|--------------------------------|
| 1A   | ATS8101     | COMFORFIT PULL-DOWN CLASSIC SM |
| 1B   | ATS8102     | COMFORFIT PULL-DOWN CLASSIC MD |
| 1C   | ATS8103     | COMFORFIT PULL-DOWN CLASSIC LG |
| 1D   | ATS8104     | COMFORFIT PULL-DOWN CLASSIC XL |
| 1E   | ATS8301     | COMFORFIT STRUCT CLS PD SM     |
| 1F   | ATS8302     | COMFORFIT STRUCT CLS PD MED    |
| 1G   | ATS8303     | COMFORFIT STRUCT CLS PD LG     |
| 1H   | ATS8304     | COMFORFIT STRUCT CLS PD XL     |




| Pos. | Item Number | Description                    |
|------|-------------|--------------------------------|
| 1A   | ATS9101     | COMFORFIT PULL-DOWN CONTOUR SM |
| 1B   | ATS9102     | COMFORFIT PULL-DOWN CONTOUR MD |
| 1C   | ATS9103     | COMFORFIT PULL-DOWN CONTOUR LG |
| 1D   | ATS9104     | COMFORFIT PULL-DOWN CONTOUR XL |
| 1E   | ATS9301     | COMFORFIT, STRUCT, CTR PD SM   |
| 1F   | ATS9302     | COMFORFIT, STRUCT, CTR PD MED  |
| 1G   | ATS9303     | COMFORFIT, STRUCT, CTR PD LG   |
| 1H   | ATS9304     | COMFORFIT, STRUCT, CTR PD XL   |




| Pos. | Item Number | Description           |
|------|-------------|-----------------------|
| 1    | J3GPSG      | SHOULDER STRAP GUIDES |
| 2-4  | J3HHP       | HARNESS BAG-O-BOLTS   |

**Note:** if using the Headrest Post on a Quickie Power Chair, please refer to the HEADREST CLAMPS page located in the appropriate power chair parts manual.




| Pos.   | Item Number | Description                  |
|--------|-------------|------------------------------|
| 5A     | AX5080      | BALLROD AXYS 6" SQUARE BLACK |
| 5B     | AX6055      | AXYS 12" SQ HORIZ BAR BLK    |
| 6      | MW070300    | AXYS VERT ROD BLK            |
| 3,4    | 115215      | J3 BK STD HD RST BRKT KIT    |
| 2-4    | M2500       | AXYS DETACH BRKT             |
| 1,5A,6 | AX7010      | AXYS CORE MOUNT HDWR         |

**Note:** if using the Headrest Post on a Quickie Power Chair, please refer to the HEADREST CLAMPS page located in the appropriate power chair parts manual.




| Pos. | Item Number | Description                    | Remarks                                     |
|------|-------------|--------------------------------|---------------------------------------------|
| 1A   | CBRA7010    | COBRA XTRA MOUNT HDW, ADUL     |                                             |
| 1B   | CBRA7030    | COBRA XTRA MOUNT HDW, PED      |                                             |
| 1C   | CBRA7080    | COBRA XTRA MOUNT HDW, EI       |                                             |
| 2A   | CBRA5030    | COBRA XTRA FB MOUNT HDW, PED   |                                             |
| 2B   | CBRA5010    | COBRA XTRA FB MOUNT HDW, ADUL  |                                             |
| 2C   | CBRA5080    | COBRA XTRA FB MOUNT HDW, EI    |                                             |
| 3    | M2400-WB    | COBRA XTRA DETACH BRKT         | INCLUDES QTY 2 1/4-20 X 1/2" SHCS & WASHERS |
| 4    | MW070152    | COBRA VERT STOP BLK            |                                             |
| 5    | M0030       | KNOB PLAS ADJ LVR 1/4-20 X20mm |                                             |
| 6    | M0055       | FASTENER PACK DETACH           |                                             |
| 7    | M0075       | FASTENER PACK SWINGAWAY MOUNT  |                                             |

**Note:** if using the Headrest Post on a Quickie Power Chair, please refer to the HEADREST CLAMPS page located in the appropriate power chair parts manual.




| Pos. | Item Number | Description                      | Remarks            |
|------|-------------|----------------------------------|--------------------|
| 1A   | LX7804      | LINX CORE STDHRZ/4\" VRT MOUNT   | 4\" vertical post  |
| 1B   | LX8000      | LINX CORE MNTG HDW 8\"           | 8\" vertical post  |
| 1C   | LX7000      | LINX CORE MOUNT HDW              | 10\" vertical post |
| 1D   | LX7815      | LINX CORE STDHRZ/15\" VRT MOUNT  | 15\" vertical post |
| 1E   | LX7000Q     | LINX CORE MOUNT HDW QUICK ADJ    |                    |
| 2    | MW000680    | CAP RATCHET FASTENER BLACK       |                    |
| 3    | MW070400    | LYNX SYS ASSM 2\"                |                    |
| 4    | MW006028    | COL SPLIT SHFT 1-1/8\" W/PIN BLK |                    |
| 5    | MC100-WB    | COVER 3/4\" COLLAR/VINYL BLACK   |                    |
| 6    | M0030       | KNOB PLAS ADJ LVR 1/4-20 X20mm   |                    |
| 7    | M2300       | LINX DETACH MOUNT BRKT           |                    |
| 8    | M2375       | MOUNT FIXED LINX 7/8\" SURFACE   |                    |
| 9    | M0055       | FASTENER PACK DETACH             |                    |
| 10   | M0075       | FASTENER PACK SWINGAWAY MOUNT    |                    |


**Note:** Please order Mounting Bracket separately.






| Pos.   | Item Number | Description                     |
|--------|-------------|---------------------------------|
| 1A     | OX6012      | VERTICAL ROD 12" ONYX           |
| 1B     | OX6008      | BALLROD ONYX 8" VERTICAL        |
| 1C     | OX6030      | VERT ROD, 17.5" ONYX W/OFF-SET  |
| 4      | MW004605    | BALLROD 8" FIXD 3/4 W/20DG, BLK |
| 5      | MW000800    | ENDCAP 7/16 ID VINYL BLK        |
| 12     | M2175       | MOUNT FIXED ONYX 7/16           |
| 13A    | OX6035      | MULTI-ANG VERT ROD, 12" ONYX    |
| 13B    | OX6050      | VERTICAL ROD 12"SQUARE ONYX     |
| 1A,2-6 | OX7010      | ONYX CORE MOUNT HDW, ADULT      |
| 1B,2-6 | OX7030      | ONYX CORE MOUNT HDW, PED        |
| 1C,2-6 | OX6010      | ONYX CORE OFFSET BAR MNTG HDW   |
| 7-11   | M2100-WB    | ONYX DETACH MOUNT BRKT          |



| Pos. | Item Number | Description          | Remarks                                                                                                                                                            |
|------|-------------|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 3-5  | J3UHM       | HEADREST PLATE ASSY  | Required when using any after-market headrest or old style Quickie Headrest where the Universal Headrest Mounting Plate is required to mount Headrest to a J3 Back |
| 1-4  | J3HDP       | HEADREST BAG-O-BOLTS | Use for Whitmyer Style Headrest Assemblies                                                                                                                         |




| Pos.   | Item Number | Description                    | Remarks |
|--------|-------------|--------------------------------|---------|
| 1A     | PF8044      | FOAM NARROW AJUSTAPLUSH        |         |
| 1A     | PC8044L     | COVER, NAR ADJUST-A-PL, LY     | Lycra   |
| 1A     | PC8044V     | COVER, NAR ADJUST-A-PL, DR     | Dartex  |
| 1A     | PC8044F     | COVER NARROW ADJ-PLUSH/FLEECE  | Fleece  |
| 1B     | PF8088      | FOAM STD AJUSTAPLUSH           |         |
| 1B     | PC8088L     | COVER, STD ADJUST-A-PL, LY     | Lycra   |
| 1B     | PC8088V     | COVER STD ADJ-PLUSH/RD         | Dartex  |
| 1B     | PC8088F     | COVER STD AJUSTAPLUSH/FLEECE   | Fleece  |
| 2      | MF014600    | 1/4-28X5/8" SHCS- COLLAR BOLT  |         |
| 3      | MC100-WB    | COVER 3/4"COLLAR/VINYL BLACK   |         |
| 4      | MW027100    | COLLAR HALF RAD 3/4" THD. BLCK |         |
| 1A,2-4 | P8044L      | NARROW ADJ-A-PLUSH LYCRA       | Lycra   |
| 1A,2-4 | P8044V      | NARROW ADJ-A-PLUSH DARTEX      | Dartex  |
| 1B,2-4 | P8088L      | STANDARD ADJ-A-PLUSH LYCRA     | Lycra   |
| 1B,2-4 | P8088V      | STANDARD ADJ-A-PLUSH DARTEX    | Dartex  |

**HEADREST MEASURING GUIDE**

| Description           | A    | B     | C      |
|-----------------------|------|-------|--------|
| Small Contour Cradle  | 8.5" | 5.5"  | 4"     |
| Medium Contour Cradle | 10"  | 6.25" | 4.75"  |
| Large Contour Cradle  | 13"  | 7.5"  | 6-3/8" |

| Pos. | Item Number | Description                    | Remarks    |
|------|-------------|--------------------------------|------------|
| 1    | PF200       | FOAM SMALL CONTOURED CRADLE    | Foam only  |
| 1    | PC200L      | COVER SML CONTED CRADLE/LYCRA  | Cover only |
| 1    | PC200V      | COVER SML CONTRD CRADLE/RD     | Cover only |
| 1    | PC200F      | COVER SML CONTED CRADLE/FLEECE | Cover only |
| 1    | PF300       | FOAM MED CONTOURED CRADLE      | Foam only  |
| 1    | PC300L      | COVER MED CONTED CRADLE/LYCRA  | Cover only |
| 1    | PC300V      | COVER MED CONTRD CRADLE/RD     | Cover only |
| 1    | PC300F      | COVER MED CONTED CRADLE/FLEECE | Cover only |
| 1    | PF400       | FOAM LG CONTOURED CRADLE       | Foam only  |
| 1    | PC400L      | COVER LG CONTED CRADLE/LYCRA   | Cover only |
| 1    | PC400V      | COVER LG CONTRD CRADLE/RD      | Cover only |
| 1    | PC400F      | COVER LG CONTED CRADLE/FLEECE  | Cover only |
| 2    | MF014600    | 1/4-28X5/8" SHCS- COLLAR BOLT  |            |
| 3    | MC100-WB    | COVER 3/4"COLLAR/VINYL BLACK   |            |
| 4    | MW027100    | COLLAR HALF RAD 3/4" THD. BLCK |            |
| 1-4  | P200L       | SMALL C-PLUSH LYCRA            |            |


| Pos. | Item Number | Description           | Remarks |
|------|-------------|-----------------------|---------|
| 1-4  | P200V       | SMALL C-PLUSH DARTEX  |         |
| 1-4  | P300L       | MEDIUM C-PLUSH LYCRA  |         |
| 1-4  | P300V       | MEDIUM C-PLUSH DARTEX |         |
| 1-4  | P400L       | LARGE C-PLUSH LYCRA   |         |
| 1-4  | P400V       | LARGE C-PLUSH DARTEX  |         |



| Pos. | Item Number | Description                      | Remarks |
|------|-------------|----------------------------------|---------|
| 1A   | PF605       | REFOAM KIT 6" PLUSH              |         |
| 1A   | PC605L      | COVER 6" PLUSH/LYCRA             | Lycra   |
| 1A   | PC605V      | COVER 6" PLUSH/R/VSD DARTEXX     | Dartex  |
| 1B   | PF800       | REFOAM KIT 8"PLUSH               |         |
| 1B   | PC800L      | COVER, 8" PLUSH, LYCRA           | Lycra   |
| 1B   | PC800V      | COVER, 8" PLUSH, DR              | Dartex  |
| 1B   | PC800F      | COVER 8"PLUSH/FLEECE             | Fleece  |
| 1B   | PC800P      | COVER 8"PLUSH/ALLURE NH          |         |
| 1C   | PF1000      | REFOAM KIT 10"PLUSH              |         |
| 1C   | PC1000L     | COVER, 10" PLUSH, LYCRA          | Lycra   |
| 1C   | PC1000V     | COVER, 10" PLUSH, DR             | Dartex  |
| 1C   | PC1000F     | COVER 10"PLUSH/FLEECE            | Fleece  |
| 1D   | PF1205      | REFOAM KIT 12" NARROW PLUSH      |         |
| 1D   | PC1205L     | COVER NARROW PLUSH/LYCRA         | Lycra   |
| 1D   | PC1205V     | COVER NARROW PLUSH/R/VSD DARTEXX | Dartex  |
| 1D   | PC1205F     | COVER NARROW PLUSH/FLEECE        | Fleece  |
| 2A   | PF1400      | REFOAM KIT 14"PLUSH              |         |

| Pos.   | Item Number | Description                        | Remarks |
|--------|-------------|------------------------------------|---------|
| 2A     | PC1400L     | COVER, 14" PLUSH, LYCRA            | Lycra   |
| 2A     | PC1400V     | COVER, 14" PLUSH, DR               | Dartex  |
| 2A     | PC1400F     | COVER 14"PLUSH/FLEECE              | Fleece  |
| 2B     | PF1900      | REFOAM KIT 19"PLUSH                |         |
| 2B     | PC1900L     | COVER 19"PLUSH/LYCRA               | Lycra   |
| 2B     | PC1900V     | COVER 19" PLUSH/R'VSD DARTEXX      | Dartex  |
| 2B     | PC1900F     | COVER 19"PLUSH/FLEECE              | Fleece  |
| 4      | MF014600    | 1/4-28X5/8" SHCS- COLLAR BOLT      |         |
| 5      | MC100-WB    | COVER 3/4"COLLAR/VINYL BLACK       |         |
| 6      | MW027100    | COLLAR HALF RAD 3/4" THD. BLCK     |         |
| 1A,4-6 | P605L       | PLUSH PAD, 6", LYCRA               | Lycra   |
| 1A,4-6 | P605V       | PLUSH PAD, 6", DARTEX              | Dartex  |
| 1B,4-6 | P800L       | PLUSH PAD, 8 IN LYCRA              | Lycra   |
| 1B,4-6 | P800V       | PLUSH PAD, 8 IN DARTEX             | Dartex  |
| 1C,4-6 | P1000L      | PLUSH PAD, 10 IN LYCRA             | Lycra   |
| 1C,4-6 | P1000V      | PLUSH PAD, 10 IN DARTEX            | Dartex  |
| 1D,4-6 | P1205L      | PLUSH PAD, 12 IN NARROW LYCRA      | Lycra   |
| 1D,4-6 | P1205V      | PLUSH PAD, 12 IN NARROW DARTEX     | Dartex  |
| 2A,4-6 | P1400L      | PLUSH PAD, 14 IN LYCRA             | Lycra   |
| 2A,4-6 | P1400V      | PLUSH PAD, 14 IN DARTEX            | Dartex  |
| 2B,4-6 | P1900L      | PLUSH PAD, 19 IN LYCRA             | Lycra   |
| 2B,4-6 | P1900V      | PLUSH PAD, 19 IN DARTEX            | Dartex  |
| 3A,4-6 | P1420L      | PLUSH PAD, 14 IN DUAL MOUNT LYCRA  | Lycra   |
| 3A,4-6 | P1420V      | PLUSH PAD, 14 IN DUAL MOUNT DARTEX | Dartex  |
| 3B,4-6 | P1920L      | PLUSH PAD, 19 IN DUAL MOUNT LYCRA  | Lycra   |
| 3B,4-6 | P1920V      | PLUSH PAD, 19 IN DUAL MOUNT DARTEX | Dartex  |


Note: Parts not compatible with Headrest ordered prior to 5/27/2013, please replace complete headrest if ordered prior to 5/27/2013.




| Pos. | Item Number | Description                     | Remarks                  |
|------|-------------|---------------------------------|--------------------------|
|      |             | DARTEX                          |                          |
| 1    | 125375      | COV INFANT SHT WIDE SHT DARTEX  | Short Cradle, Short Arms |
| 1    | 125376      | COV INFANT SHT WIDE LNG DARTEX  | Short Cradle, Long Arms  |
| 1    | 125379      | COV INFANT TALL THIN SHT DARTEX | Tall Cradle, Short Arms  |
| 1    | 125380      | COV INFANT TALL THIN LNG DARTEX | Tall Cradle, Long Arms   |
| 1    | 125383      | COV SMALL SHT WIDE SHT DARTEX   | Short Cradle, Short Arms |
| 1    | 125384      | COV SMALL SHT WIDE LNG DARTEX   | Short Cradle, Long Arms  |
| 1    | 125387      | COV SMALL TALL THIN SHT DARTEX  | Tall Cradle, Short Arms  |
| 1    | 125388      | COV SMALL TALL THIN LNG DARTEX  | Tall Cradle, Long Arms   |
| 1    | 125391      | COV MEDIUM SHT WIDE SHT DARTEX  | Short Cradle, Short Arms |
| 1    | 125392      | COV MEDIUM SHT WIDE LNG DARTEX  | Short Cradle, Long Arms  |
| 1    | 125395      | COV MEDIUM TALL THIN SHT DARTEX | Tall Cradle, Short Arms  |
| 1    | 125396      | COV MEDIUM TALL THIN LNG DARTEX | Tall Cradle, Long Arms   |
| 1    | 125399      | COV LARGE SHT WIDE SHT DARTEX   | Short Cradle, Short Arms |

| Pos. | Item Number | Description                    | Remarks                       |
|------|-------------|--------------------------------|-------------------------------|
| 1    | 125400      | COV LARGE SHT WIDE LNG DARTEX  | Short Cradle, Long Arms       |
| 1    | 125403      | COV LARGE TALL THIN SHT DARTEX | Tall Cradle, Short Arms       |
| 1    | 125404      | COV LARGE TALL THIN LNG DARTEX | Tall Cradle, Long Arms        |
|      |             | LYCRA                          |                               |
| 1    | 125377      | COV INFANT SHT WIDE SHT LYCRA  | Short Cradle, Short Arms      |
| 1    | 125378      | COV INFANT SHT WIDE LNG LYCRA  | Short Cradle, Long Arms       |
| 1    | 125381      | COV INFANT TALL THIN SHT LYCRA | Tall Cradle, Short Arms       |
| 1    | 125382      | COV INFANT TALL THIN LNG LYCRA | Tall Cradle, Long Arms        |
| 1    | 125385      | COV SMALL SHT WIDE SHT LYCRA   | Short Cradle, Short Arms      |
| 1    | 125386      | COV SMALL SHT WIDE LNG LYCRA   | Short Cradle, Long Arms       |
| 1    | 125389      | COV SMALL TALL THIN SHT LYCRA  | Tall Cradle, Short Arms       |
| 1    | 125390      | COV SMALL TALL THIN LNG LYCRA  | Tall Cradle, Long Arms        |
| 1    | 125393      | COV MEDIUM SHT WIDE SHT LYCRA  | Short Cradle, Short Arms      |
| 1    | 125394      | COV MEDIUM SHT WIDE LNG LYCRA  | Short Cradle, Long Arms       |
| 1    | 125397      | COV MEDIUM TALL THIN SHT LYCRA | Tall Cradle, Short Arms       |
| 1    | 125398      | COV MEDIUM TALL THIN LNG LYCRA | Tall Cradle, Long Arms        |
| 1    | 125401      | COV LARGE SHT WIDE SHT LYCRA   | Short Cradle, Short Arms      |
| 1    | 125402      | COV LARGE SHT WIDE LNG LYCRA   | Short Cradle, Long Arms       |
| 1    | 125405      | COV LARGE TALL THIN SHT LYCRA  | Tall Cradle, Short Arms       |
| 1    | 125406      | COV LARGE TALL THIN LNG LYCRA  | Tall Cradle, Long Arms        |
| 2    | 125351      | STANDARD CRADLE PAD INF SHT    |                               |
| 2    | 125352      | STANDARD CRADLE PAD INF TL     |                               |
| 2    | 125353      | STANDARD CRADLE PAD SM SHT     |                               |
| 2    | 125354      | STANDARD CRADLE PAD SM TL      |                               |
| 2    | 125355      | STANDARD CRADLE PAD MED SHT    |                               |
| 2    | 125356      | STANDARD CRADLE PAD MED TL     |                               |
| 2    | 125357      | STANDARD CRADLE PAD LRG SHT    |                               |
| 2    | 125358      | STANDARD CRADLE PAD LRG TL     |                               |
| 2    | 125359      | COMFORT CRADLE PAD INF SHT     |                               |
| 2    | 125360      | COMFORT CRADLE PAD INF TL      |                               |
| 2    | 125361      | COMFORT CRADLE PAD SM SHT      |                               |
| 2    | 125362      | COMFORT CRADLE PAD SM TL       |                               |
| 2    | 125363      | COMFORT CRADLE PAD MED SHT     |                               |
| 2    | 125364      | COMFORT CRADLE PAD MED TL      |                               |
| 2    | 125365      | COMFORT CRADLE PAD LRG SHT     |                               |
| 2    | 125366      | COMFORT CRADLE PAD LRG TL      |                               |
| 3    | 125343      | INFANT SHORT CRADLE SHELL      |                               |
| 3    | 125345      | SMALL SHORT CRADLE SHELL       |                               |
| 3    | 125347      | MEDIUM SHORT CRADLE SHELL      |                               |
| 3    | 125349      | LARGE SHORT CRADLE SHELL       |                               |
| 3    | 125344      | INFANT TALL CRADLE SHELL       |                               |
| 3    | 125346      | SMALL TALL CRADLE SHELL        |                               |
| 3    | 125348      | MEDIUM TALL CRADLE SHELL       |                               |
| 3    | 125350      | LARGE TALL CRADLE SHELL        |                               |
| 4-8  | 125415      | PUSH BUTTON ADJ MECH 1 SIDE    |                               |
| 9    | 125367      | SWING ARM WIDE EXT SHORT RT    | Used w/ Med & Large Cradle    |
| 9    | 125368      | SWING ARM WIDE EXT SHORT LT    | Used w/ Med & Large Cradle    |
| 9    | 125369      | SWING ARM WIDE EXT LONG RT     | Used w/ Med & Large Cradle    |
| 9    | 125370      | SWING ARM WIDE EXT LONG LT     | Used w/ Med & Large Cradle    |
| 9    | 125371      | SWING ARM THIN EXT SHORT RT    | Used w/ Infant & Small Cradle |
| 9    | 125372      | SWING ARM THIN EXT SHORT LT    | Used w/ Infant & Small Cradle |
| 9    | 125373      | SWING ARM THIN EXT LONG RT     | Used w/ Infant & Small Cradle |
| 9    | 125374      | SWING ARM THIN EXT LONG LT     | Used w/ Infant & Small Cradle |
| 10   | 106933      | NUT NYLOCK THIN 1/4-20 ROHS    |                               |
| 11   | MF014600    | 1/4-28x5/8" SHCS-COLLAR BOLT   |                               |
| 12   | MC100-WB    | COVER 3/4" COLLAR /VINYL BLACK |                               |
| 13   | MW027100    | COLLAR HALF RAD 3/4" THD. BLCK |                               |


| Pos. | Item Number | Description                  | Remarks |
|------|-------------|------------------------------|---------|
| 15   | 125407      | FOAM PAD CHIN PROMPT INFANT  |         |
| 15   | 125408      | FOAM PAD CHIN PROMPT SMALL   |         |
| 15   | 125409      | FOAM PAD CHIN PROMPT MEDIUM  |         |
| 15   | 125410      | FOAM PAD CHIN PROMPT LARGE   |         |
| 14   | 125411      | COVER PAD CHIN PROMPT INFANT |         |
| 14   | 125412      | COVER PAD CHIN PROMPT SMALL  |         |
| 14   | 125413      | COVER PAD CHIN PROMPT MEDIUM |         |
| 14   | 125414      | COVER PAD CHIN PROMPT LARGE  |         |
| 16   | M9275       | S.O.F.T. BPA, SINGLE SUB-OCC |         |
| 17   | P500L-WB    | OCCIP PAD, PED, LYCRA        |         |
| 17   | P600L       | OCCIP PAD, ADULT, LYCRA      |         |
| 17   | P800L       | PLUSH PAD, 8 IN LYCRA        |         |
| 17   | P1000L      | PAD 10" PLUSH/LYCRA          |         |
| 17   | P500V       | OCCIP PAD, PED, DR           |         |
| 17   | P600V       | OCCIP PAD, ADULT, DR         |         |
| 17   | P800V       | PLUSH PAD, 8", DR            |         |
| 17   | P1000V      | PLUSH PAD, 10", DR           |         |



| Pos. | Item Number | Description                   | Remarks                |
|------|-------------|-------------------------------|------------------------|
| 1    | 116527S     | VELCRO TIE,1 WRAP.75"X12",BLK |                        |
| 1,2  | 117096      | VANITY FLAP, J3 BACK 12W      | Includes 4 velcro ties |
| 1,2  | 117097      | VANITY FLAP, J3 BACK 14W      | Includes 4 velcro ties |
| 1,2  | 117098      | VANITY FLAP, J3 BACK 16W      | Includes 4 velcro ties |
| 1,2  | 117099      | VANITY FLAP, J3 BACK 18W      | Includes 4 velcro ties |
| 1,2  | 117100      | VANITY FLAP, J3 BACK 20W      | Includes 4 velcro ties |
| 1,2  | 117101      | VANITY FLAP, J3 BACK 22W      | Includes 4 velcro ties |
| 1,2  | 117102      | VANITY FLAP, J3 BACK 24W      | Includes 4 velcro ties |
| 1,2  | 117103      | VANITY FLAP, J3 BACK 26W      | Includes 4 velcro ties |
| 1,2  | VF08        | VANITY FLAP 8-9 IN            | Includes 4 velcro ties |
| 1,2  | VF10        | VANITY FLAP 10-11 IN          | Includes 4 velcro ties |
| 1,2  | VF12        | VANITY FLAP 12-13 IN          | Includes 4 velcro ties |
| 1,2  | VF14        | VANITY FLAP 14-15 IN          | Includes 4 velcro ties |




| Pos. | Item Number | Description        |
|------|-------------|--------------------|
| 1    | J3HEX3MM    | HEX KEY 3MM LONG   |
| 2    | J3HEX4MM    | HEX KEY 4MM LONG   |
| 3    | J3CW10MM    | 10 MM COMBO WRENCH |



| Pos. | Item Number | Description                 |
|------|-------------|-----------------------------|
| 1    | DJ3SPNKA    | ASSESSMENT KIT, SPINE ALIGN |


**NOTE:** Parts only, for complete Backrest see Jay Zip Back order form.




| Pos. | Item Number | Description                  | Pos. | Item Number | Description                   |
|------|-------------|------------------------------|------|-------------|-------------------------------|
| 1    | CJZIPB08M   | COVER JAY ZIP BACK 8W MED    | 2    | ICJZIPB08T  | IN COV JAY ZIP BACK 8W TALL   |
| 1    | CJZIPB08S   | COVER JAY ZIP BACK 8W SHORT  | 2    | ICJZIPB08X  | IN COV JAY ZIP BACK 8W XTALL  |
| 1    | CJZIPB08T   | COVER JAY ZIP BACK 8W TALL   | 2    | ICJZIPB10M  | IN COV JAY ZIP BACK 10W MED   |
| 1    | CJZIPB08X   | COVER JAY ZIP BACK 8W XTALL  | 2    | ICJZIPB10S  | IN COV JAY ZIP BACK 10W SHORT |
| 1    | CJZIPB10M   | COVER JAY ZIP BACK 10W MED   | 2    | ICJZIPB10T  | IN COV JAY ZIP BACK 10W TALL  |
| 1    | CJZIPB10S   | COVER JAY ZIP BACK 10W SHORT | 2    | ICJZIPB10X  | IN COV JAY ZIP BACK 10W XTALL |
| 1    | CJZIPB10T   | COVER JAY ZIP BACK 10W TALL  | 2    | ICJZIPB12M  | IN COV JAY ZIP BACK 12W MED   |
| 1    | CJZIPB10X   | COVER JAY ZIP BACK 10W XTALL | 2    | ICJZIPB12S  | IN COV JAY ZIP BACK 12W SHORT |
| 1    | CJZIPB12M   | COVER JAY ZIP BACK 12W MED   | 2    | ICJZIPB12T  | IN COV JAY ZIP BACK 12W TALL  |
| 1    | CJZIPB12S   | COVER JAY ZIP BACK 12W SHORT | 2    | ICJZIPB12X  | IN COV JAY ZIP BACK 12W XTALL |
| 1    | CJZIPB12T   | COVER JAY ZIP BACK 12W TALL  | 2    | ICJZIPB14M  | IN COV JAY ZIP BACK 14W MED   |
| 1    | CJZIPB12X   | COVER JAY ZIP BACK 12W XTALL | 2    | ICJZIPB14S  | IN COV JAY ZIP BACK 14W SHORT |
| 1    | CJZIPB14M   | COVER JAY ZIP BACK 14W MED   | 2    | ICJZIPB14T  | IN COV JAY ZIP BACK 14W TALL  |
| 1    | CJZIPB14S   | COVER JAY ZIP BACK 14W SHORT | 2    | ICJZIPB14X  | IN COV JAY ZIP BACK 14W XTALL |
| 1    | CJZIPB14T   | COVER JAY ZIP BACK 14W TALL  | 3    | BJZIPB08M   | FOAM JAY ZIP BACK 8W MED      |
| 1    | CJZIPB14X   | COVER JAY ZIP BACK 14W XTALL | 3    | BJZIPB08S   | FOAM JAY ZIP BACK 8W SHORT    |
| 2    | ICJZIPB08M  | IN COV JAY ZIP BACK 8W MED   | 3    | BJZIPB08T   | FOAM JAY ZIP BACK 8W TALL     |
| 2    | ICJZIPB08S  | IN COV JAY ZIP BACK 8W SHORT | 3    | BJZIPB08X   | FOAM JAY ZIP BACK 8W XTALL    |

| Pos. | Item Number | Description                 | Pos. | Item Number | Description                  |
|------|-------------|-----------------------------|------|-------------|------------------------------|
| 3    | BJZIPB10M   | FOAM JAY ZIP BACK 10W MED   | 4    | SJZIPB08T   | SHELL JAY ZIP BACK 8W TALL   |
| 3    | BJZIPB10S   | FOAM JAY ZIP BACK 10W SHORT | 4    | SJZIPB08X   | SHELL JAY ZIP BACK 8W XTALL  |
| 3    | BJZIPB10T   | FOAM JAY ZIP BACK 10W TALL  | 4    | SJZIPB10M   | SHELL JAY ZIP BACK 10W MED   |
| 3    | BJZIPB10X   | FOAM JAY ZIP BACK 10W XTALL | 4    | SJZIPB10S   | SHELL JAY ZIP BACK 10W SHORT |
| 3    | BJZIPB12M   | FOAM JAY ZIP BACK 12W MED   | 4    | SJZIPB10T   | SHELL JAY ZIP BACK 10W TALL  |
| 3    | BJZIPB12S   | FOAM JAY ZIP BACK 12W SHORT | 4    | SJZIPB10X   | SHELL JAY ZIP BACK 10W XTALL |
| 3    | BJZIPB12T   | FOAM JAY ZIP BACK 12W TALL  | 4    | SJZIPB12M   | SHELL JAY ZIP BACK 12W MED   |
| 3    | BJZIPB12X   | FOAM JAY ZIP BACK 12W XTALL | 4    | SJZIPB12S   | SHELL JAY ZIP BACK 12W SHORT |
| 3    | BJZIPB14M   | FOAM JAY ZIP BACK 14W MED   | 4    | SJZIPB12T   | SHELL JAY ZIP BACK 12W TALL  |
| 3    | BJZIPB14S   | FOAM JAY ZIP BACK 14W SHORT | 4    | SJZIPB12X   | SHELL JAY ZIP BACK 12W XTALL |
| 3    | BJZIPB14T   | FOAM JAY ZIP BACK 14W TALL  | 4    | SJZIPB14M   | SHELL JAY ZIP BACK 14W MED   |
| 3    | BJZIPB14X   | FOAM JAY ZIP BACK 14W XTALL | 4    | SJZIPB14S   | SHELL JAY ZIP BACK 14W SHORT |
| 4    | SJZIPB08M   | SHELL JAY ZIP BACK 8W MED   | 4    | SJZIPB14T   | SHELL JAY ZIP BACK 14W TALL  |
| 4    | SJZIPB08S   | SHELL JAY ZIP BACK 8W SHORT | 4    | SJZIPB14X   | SHELL JAY ZIP BACK 14W XTALL |


**NOTE:** Parts only, for complete Backrest see JAY ZIP Back order form.




| Pos.              | Item Number | Description                | Remarks |
|-------------------|-------------|----------------------------|---------|
| 5                 | 210084      | NUT HEX NYLOCK M6 BLACK    |         |
| 7                 | 120982      | COVER HRDWR KIT SM PR      |         |
| 7                 | 120983      | COVER HRDWR KIT LG PR      |         |
| 8                 | 250104      | WASHER FLAT M6             |         |
| 16                | 120985      | ZIP BACK LANYARD 8-9W      |         |
| 16                | 120986      | ZIP BACK LANYARD 10-11W    |         |
| 16                | 120987      | ZIP BACK LANYARD 12-13W    |         |
| 16                | 120988      | ZIP BACK LANYARD 14-15W    |         |
| 8-13              | 120979      | TUBE CLAMP ASSM 1IN PR     | Pair    |
| 14-15             | 120984      | CLAMP INSERTS KIT 2 PR     |         |
| 1-3A,5,8          | 120981      | SHELL MOUNT BRKT KIT PR    |         |
| 1-3B,5,8          | 121111      | 8W SHELL MOUNT BRKT KIT PR |         |
| 4-6,8             | 120980      | MOUNT RECEIVER KIT PR      |         |
| 1,2,5,6,8,9,11-13 | 120978      | JAY ZIP BACK FASTNER KIT   |         |




| Pos. | Item Number | Description                    | Remarks    |
|------|-------------|--------------------------------|------------|
| 1A   | 120992      | THOR SUPP PAD FLAT 3X4         | single pad |
| 1B   | 120993      | THOR SUPP PAD FLAT 4X5         | single pad |
| 1C   | 120994      | THOR SUPP PAD CURVED 3X4       | single pad |
| 1D   | 120995      | THOR SUPP PAD CURVED 4X5       | single pad |
| 4    | 250104      | WASHER FLAT M6                 |            |
| 5    | 210084      | NUT HEX NYLOCK M6 BLACK        |            |
| 6    | PRT7500     | LK, HINGE ASSEMBLY             |            |
| 7    | 226018-QU   | 10-32 X 3/8 BH C/S W/PATCH BLK |            |
| 8    | XM2141S     | 10-32X.50 BHCS W/PATCH         |            |
| 9    | 113282S     | WASHER FLAT AN #10 ROHS        |            |
| 10   | PRT7625     | LK, PAD MOUNT                  |            |
| 11   | PRT7627     | LK, PAD MOUNT PLATE, OFFSET    |            |
| 12   | PRT7626     | LK, 1" OFF ADAPT FOR E-Z       |            |
| 13   | PRT7653     | OFFSET PAD MOUNT 20 DEG, LK    |            |
| 14   | PRT7651     | 1" OFFSET MOUNT 20 DEG, LK     |            |
| 15   | COV2650     | LAT BRKT FABRIC COVER,LK,SHORT |            |
| 2-5  | 120989      | ZIP BACK SA LAT BRKT KIT PR    |            |




| Pos. | Item Number | Description                    | Remarks    |
|------|-------------|--------------------------------|------------|
| 1A   | 120992      | THOR SUPP PAD FLAT 3X4         | single pad |
| 1B   | 120993      | THOR SUPP PAD FLAT 4X5         | single pad |
| 1C   | 120994      | THOR SUPP PAD CURVED 3X4       | single pad |
| 1D   | 120995      | THOR SUPP PAD CURVED 4X5       | single pad |
| 4    | 250104      | WASHER FLAT M6                 |            |
| 5    | 210084      | NUT HEX NYLOCK M6 BLACK        |            |
| 6    | XM2141S     | 10-32X.50 BHCS W/PATCH         |            |
| 7    | 113282S     | WASHER FLAT AN #10 ROHS        |            |
| 8    | COV2650     | LAT BRKT FABRIC COVER,LK,SHORT |            |
| 2-5  | 120990      | ZIP BACK FX LAT BRKT KIT PR    |            |




| Pos. | Item Number | Description                   | Remarks            |
|------|-------------|-------------------------------|--------------------|
| 1A   | 103714-001S | LUMBAR SUPP SM BLK 3.5W W/VEL |                    |
| 1B   | 103714-005S | LUMBAR SUPP SM BLK 1.8W W/VEL |                    |
| 2    | J3SPNCR     | CURVED SHIM                   |                    |
| 3    | J3SPNFL     | FLAT SHIM                     |                    |
| 4    | J3SPNLP     | POS/LATERAL PELVIC            |                    |
| 5    | J3SPNLS     | SHIM, LATERAL                 |                    |
| 6    | J3SPNBL     | BEAN BAG, LARGE               | Approx 5.25" X 10" |
| 7    | J3SPNBS     | BEAN BAG, SMALL               | Approx 3" X 6"     |




| Pos. | Item Number | Description                         | Pos.  | Item Number | Description                           |
|------|-------------|-------------------------------------|-------|-------------|---------------------------------------|
| 1    | ATS9407     | Center Opening w/ Buckle (EI)       | 3     | ATS9100     | Contoured - Dynamic Pull Down (XS)    |
| 1    | ATS9400     | Center Opening w/ Buckle (XS)       | 3     | ATS9101     | Contoured - Dynamic Pull Down (S)     |
| 1    | ATS9401     | Center Opening w/ Buckle (S)        | 3     | ATS9102     | Contoured - Dynamic Pull Down (M)     |
| 1    | ATS9402     | Center Opening w/ Buckle (M)        | 3     | ATS9300     | Contoured - Structured Pull Down (XS) |
| 2    | ATS8107     | Classic - Dynamic Pull Down (EI)    | 3     | ATS9301     | Contoured - Structured Pull Down (S)  |
| 2    | ATS8100     | Classic - Dynamic Pull Down (XS)    | 3     | ATS9302     | Contoured - Structured Pull Down (M)  |
| 2    | ATS8101     | Classic - Dynamic Pull Down (S)     | 4     | ATS4050     | Padded Anterior Trunk Straps (XS)     |
| 2    | ATS8102     | Classic - Dynamic Pull Down (M)     | 4     | ATS4000     | Padded Anterior Trunk Straps (S)      |
| 2    | ATS8300     | Classic - Structured Pull Down (XS) | 4     | ATS4100     | Padded Anterior Trunk Straps (M)      |
| 2    | ATS8301     | Classic - Structured Pull Down (S)  | 5     | J3GPSG      | GOAL POST SHOULDER STRAP GUIDE        |
| 2    | ATS8302     | Classic - Structured Pull Down (M)  | 6,7,8 | J3HHP       | Harness Hardware Pack                 |
| 3    | ATS9107     | Contoured - Dynamic Pull Down (EI)  |       |             |                                       |



| Pos. | Item Number | Description                   | Remarks                |
|------|-------------|-------------------------------|------------------------|
| 1    | 116527S     | VELCRO TIE,1 WRAP.75"X12",BLK |                        |
| 1,2  | 117096      | VANITY FLAP, J3 BACK 12W      | Includes 4 velcro ties |
| 1,2  | 117097      | VANITY FLAP, J3 BACK 14W      | Includes 4 velcro ties |
| 1,2  | 117098      | VANITY FLAP, J3 BACK 16W      | Includes 4 velcro ties |
| 1,2  | 117099      | VANITY FLAP, J3 BACK 18W      | Includes 4 velcro ties |
| 1,2  | 117100      | VANITY FLAP, J3 BACK 20W      | Includes 4 velcro ties |
| 1,2  | 117101      | VANITY FLAP, J3 BACK 22W      | Includes 4 velcro ties |
| 1,2  | 117102      | VANITY FLAP, J3 BACK 24W      | Includes 4 velcro ties |
| 1,2  | 117103      | VANITY FLAP, J3 BACK 26W      | Includes 4 velcro ties |
| 1,2  | VF08        | VANITY FLAP 8-9 IN            | Includes 4 velcro ties |
| 1,2  | VF10        | VANITY FLAP 10-11 IN          | Includes 4 velcro ties |
| 1,2  | VF12        | VANITY FLAP 12-13 IN          | Includes 4 velcro ties |
| 1,2  | VF14        | VANITY FLAP 14-15 IN          | Includes 4 velcro ties |




| Pos.   | Item Number | Description                    |
|--------|-------------|--------------------------------|
| 1      | OX7030      | ONYX CORE MOUNT HDW, PED       |
| 5      | MW070230S   | MT .5 SQR TUBING RECIEVER BLK  |
| 6      | M0030       | KNOB PLAS ADJ LVR 1/4-20X20MM  |
| 7      | 102971-001S | NUT PLATE M6 STUD SHORT X 51   |
| 8      | 210084      | NUT HEX NYLOCK M6 BLACK        |
| 9      | 250104      | WASHER FLAT M6                 |
| 10     | 103185-006  | BHCS M6X1X16 NP BLK            |
| 1-10   | 120452S     | ZIP HEADREST ASSM              |
| 2-4    | 120991      | HDRST STOP BLK KIT             |
| 12     | MF014600    | 1/4-28X5/8" SHCS- COLLAR BOLT  |
| 13     | MC100-WB    | COVER 3/4"COLLAR/VINYL BLACK   |
| 14     | MW027100    | COLLAR HALF RAD 3/4" THD. BLCK |
| 11A-14 | P600L       | OCCIP PAD, ADULT, LYCRA        |
| 11A-14 | P800L       | PLUSH PAD, 8 IN LYCRA          |
| 11B-14 | INFT200LS   | PAD INF CONTRD CRADLE LYCRA    |
| 11B-14 | P200L       | CONT CRADLE PAD, SM, LYCRA     |



| Pos. | Item Number | Description       |
|------|-------------|-------------------|
| 1    | 104549S     | COMBO WRENCH 10MM |



## MEASURING GUIDE



| Part Number | Dim A (Width) | Dim B (Height) | Dim C (Thickness) | Dim D (Right Side Strap*) | Dim E (Left Side Strap*) |
|-------------|---------------|----------------|-------------------|---------------------------|--------------------------|
| 2250        | 15"           | 8.75"          | 2"                | 7"                        | 7"                       |
| 2251-MJ     | 11"           | 7"             | 2"                | 6.5"                      | 15"                      |

\* Strap side is located as if you are resting against the back with the lumbar support on it.

| Pos. | Item Number | Description                   | Remarks |
|------|-------------|-------------------------------|---------|
| 1A   | 2250        | LUMBAR SUPPORT SOFT COMBI     | Adult   |
| 1B   | 2251-MJ     | LUMBAR SUPPORT SFT COMBI JUNR | Junior  |



improving people's lives



Dear Valued Customer,

Effective December 12<sup>th</sup> 2013 the Jay Active Back has been made obsolete and parts will no longer be offered by Sunrise Medical. This product has been discontinued and due to the diminished demand for its parts Sunrise Medical can no longer offer replacement parts.

Sincerely,

Sunrise Medical

2842 Business Park Avenue • Fresno, CA 93727-1328 • USA

800.333.4000 [www.SunriseMedical.com](http://www.SunriseMedical.com)

| Pos. | Description      |
|------|------------------|
| 1    | Obsolete Product |



improving people's lives



Dear Valued Customer,

Effective June 20th, 2014 the Jay Xtreme Back has been made obsolete and parts will no longer be offered by Sunrise Medical. This product has been discontinued and due to the diminished demand for its parts, Sunrise Medical can no longer offer replacement parts.


Sincerely,

Sunrise Medical

2842 Business Park Avenue • Fresno, CA 93727-1328 • USA


800.333.4000 [www.SunriseMedical.com](http://www.SunriseMedical.com)

| Description |
|-------------|
| Obsolete    |




| Pos. | Item Number | Description                   | Remarks |
|------|-------------|-------------------------------|---------|
| 1    | S600        | SHELL JAY MODULAR BACK 18"    |         |
| 1    | S602        | SHELL JAY MODULAR BACK 20"    |         |
| 1    | S605        | SHELL JAY MODULAR BACK 16"    |         |
| 1    | S610        | SHELL JAY MODULAR BACK 14"    |         |
| 2    | K600        | SHIM 1/2"SEAT DEPTH FOR 18"   |         |
| 2    | K602        | SHIM 1/2"SEAT DEPTH FOR 20"   |         |
| 2    | K605        | SHIM 1/2"SEAT DEPTH FOR 16"   |         |
| 2    | K610        | SHIM 1/2"SEAT DEPTH FOR 14"   |         |
| 4    | FL600       | FLUID PAD FULL WIDTH MOD BACK | 18"     |
| 4    | FL602       | FLUID PAD FULL WIDTH MOD BACK | 20"     |
| 4    | FL605       | FLUID PAD FULL WIDTH MOD BACK | 16"     |
| 4    | FL610       | FLUID PAD FULL WIDTH MOD BACK | 14"     |
| 5A   | C600        | COVER JAY MODULAR BACK 18"    |         |

| Pos. | Item Number | Description                       | Remarks                            |
|------|-------------|-----------------------------------|------------------------------------|
| 5A   | C602        | COVER JAY MODULAR BACK 20"        |                                    |
| 5A   | C603        | COVER JAY MODULAR BACK 17"        |                                    |
| 5A   | C605        | COVER JAY MODULAR BACK 16"        |                                    |
| 5A   | C610        | COVER JAY MODULAR BACK 14"        |                                    |
| 5B   | CB600       | INSIDE PROTECTIVE COVER           | 18"                                |
| 5B   | CB602       | INSIDE PROTECTIVE COVER           | 20"                                |
| 5B   | CB605       | INSIDE PROTECTIVE COVER           | 16"                                |
| 5B   | CB610       | INSIDE PROTECTIVE COVER           | 14"                                |
| 6    | B684-MJ     | SIDE GUARDS JAY MODULAR BACK      |                                    |
| 7A   | 681C        | TOP BRACKET ASSM - "QUARTER-TURN" | fits 5/8" to 1.0" diameter tubing  |
| 7B   | 2681C       | TOP LONG 1/4 BRACKET ASSM         | fits 1.0" to 1.25" diameter tubing |
| 8    | 680C        | BOTTOM UNIVERSAL POST BRACKETS    |                                    |
| 9    | 880997      | LONG STRAP KIT                    | USED ON 1 1/8" DIAMETER BACK POST  |
| 2,3  | BL600       | BASE&SHIM CUSTOM CONTOURED 18"    |                                    |
| 2,3  | BL602       | BASE&SHIM CUSTOM CONTOURED 20"    |                                    |
| 2,3  | BL605       | BASE&SHIM CUSTOM CONTOURED 16"    |                                    |
| 2,3  | BL610       | BASE&SHIM CUSTOM CONTOURED 14"    |                                    |
| 2-5A | DI600       | CUSTOM CONTOURED INSERT 18"       |                                    |
| 2-5A | DI602       | CUSTOM CONTOURED INSERT 20"       |                                    |
| 2-5A | DI605       | CUSTOM CONTOURED INSERT 16"       |                                    |
| 2-5A | DI610       | CUSTOM CONTOURED INSERT 14"       |                                    |




| Pos. | Item Number | Description                    | Remarks |
|------|-------------|--------------------------------|---------|
| 1    | S600        | SHELL JAY MODULAR BACK 18"     |         |
| 1    | S602        | SHELL JAY MODULAR BACK 20"     |         |
| 1    | S605        | SHELL JAY MODULAR BACK 16"     |         |
| 1    | S610        | SHELL JAY MODULAR BACK 14"     |         |
| 2    | B600        | BASE MODULAR BACK STD 18"      |         |
| 2    | B602        | BASE MODULAR BACK STD 20"      |         |
| 2    | B605        | BASE MODULAR BACK STD 16"      |         |
| 2    | B610        | BASE MODULAR BACK STANDARD 14" |         |
| 3    | F605        | SPINAL FLUID PAD MODULAR BACK  |         |
| 3    | F510        | SPINAL FLUID PAD STANDARD      |         |
| 4A   | C500AB      | COVER JAY BACK 18"             |         |
| 4A   | C502AB      | COVER JAY BACK 20"             |         |
| 4A   | C505AB      | COVER JAY BACK 16"             |         |
| 4A   | C510AB      | COVER JAY BACK 14"             |         |
| 4B   | CB600       | INSIDE PROTECTIVE COVER        | 18"     |
| 4B   | CB602       | INSIDE PROTECTIVE COVER        | 20"     |
| 4B   | CB605       | INSIDE PROTECTIVE COVER        | 16"     |
| 4B   | CB610       | INSIDE PROTECTIVE COVER        | 14"     |
| 5    | B684-MJ     | SIDE GUARDS JAY MODULAR BACK   |         |

| Pos. | Item Number | Description                       | Remarks                            |
|------|-------------|-----------------------------------|------------------------------------|
| 6    | 880997      | LONG STRAP KIT                    | USED ON 1 1/8" DIAMETER BACK POST  |
| 7A   | 681C        | TOP BRACKET ASSM - "QUARTER-TURN" | fits 5/8" to 1.0" diameter tubing  |
| 7B   | 2681C       | TOP LONG 1/4 BRACKET ASSM         | fits 1.0" to 1.25" diameter tubing |
| 8    | 680C        | BOTTOM UNIVERSAL POST BRACKETS    |                                    |
| 2-4  | NI600       | STANDARD INSERT FOR S600          | 18"                                |
| 2-4  | NI602       | STANDARD INSERT FOR S602          | 20"                                |
| 2-4  | NI605       | STANDARD INSERT FOR S605          | 16"                                |
| 2-4  | NI610       | STANDARD INSERT FOR S610          | 14"                                |



| Pos. | Item Number | Description                    | Remarks |
|------|-------------|--------------------------------|---------|
| 1    | 682         | EXTENDED LATERAL SUPPORT       |         |
| 2    | B682W       | EXTENDED LATERAL WEDGE         |         |
| 3    | F682        | FLUID PAD EXTENDED LATERAL     |         |
| 4    | FD500       | FLUID PAD DEFORMITY FILL-IN    |         |
| 5    | B689        | CORNER TRANSITION WEDGE (2 EA) |         |
| 6    | B687        | TRANSITION WEDGES 1/2 WID (4)  |         |
| 7    | B688        | TRANSITION WEDGES FULL WID.(4) |         |
| 8    | 686         | LARGE FILL-IN PADS (SET OF 4)  | 5"X10"  |
| 9    | 685         | SMALL FILL-IN PADS (SET OF 4)  | 3"X6"   |
| 10   | 655         | JAY BACK ACCESSORY BOX         |         |

**NOTE:** IF MOUNTING LATERALS ONTO A SOLID BACK PLEASE REFER TO THE "LATERAL PAD ASSEMBLY" PARTS PAGE FOR MOUNTING HARDWARE.




| Pos. | Item Number | Description                                    | Remarks                 |
|------|-------------|------------------------------------------------|-------------------------|
| 1A   | 683358      | JAY BACK UPH BACK CANE COVER MEDIUM            | LEFT 10"-14" & 12"-16"  |
| 1A   | 683360      | JAY BACK UPH BACK CANE COVER MEDIUM            | RIGHT 10"-14" & 12"-16" |
| 1B   | 683357      | JAY BACK UPH BACK CANE COVER TALL              | LEFT 14"-18" & 16"-20"  |
| 1B   | 683359      | JAY BACK UPH BACK CANE COVER TALL              | RIGHT 14"-18" & 16"-20" |
| 1C   | 683362      | JAY BACK UPH BACK CANE COVER NO LATERAL MEDIUM | LEFT 10"-14" & 12"-16"  |
| 1C   | 683364      | JAY BACK UPH BACK CANE COVER NO LATERAL MEDIUM | RIGHT 10"-14" & 12"-16" |
| 1D   | 683361      | JAY BACK UPH BACK CANE COVER NO LATERAL TALL   | LEFT 14"-18" & 16"-20"  |
| 1D   | 683363      | JAY BACK UPH BACK CANE COVER NO LATERAL TALL   | RIGHT 14"-18" & 16"-20" |
| 2A   | 683305      | JAY BACK UPH 14W MEDIUM                        | 10"-14" & 12"-16"       |
| 2A   | 683307      | JAY BACK UPH 15W MEDIUM                        | 10"-14" & 12"-16"       |

| Pos.            | Item Number | Description                                       | Remarks           |
|-----------------|-------------|---------------------------------------------------|-------------------|
| 2A              | 683309      | JAY BACK UPH 16W MEDIUM                           | 10"-14" & 12"-16" |
| 2A              | 683311      | JAY BACK UPH 17W MEDIUM                           | 10"-14" & 12"-16" |
| 2A              | 683313      | JAY BACK UPH 18W MEDIUM                           | 10"-14" & 12"-16" |
| 2A              | 683315      | JAY BACK UPH 19W MEDIUM                           | 10"-14" & 12"-16" |
| 2A              | 683317      | JAY BACK UPH 20W MEDIUM                           | 10"-14" & 12"-16" |
| 2B              | 683306      | JAY BACK UPH 14W TALL                             | 14"-18" & 16"-20" |
| 2B              | 683308      | JAY BACK UPH 15W TALL                             | 14"-18" & 16"-20" |
| 2B              | 683310      | JAY BACK UPH 16W TALL                             | 14"-18" & 16"-20" |
| 2B              | 683312      | JAY BACK UPH 17W TALL                             | 14"-18" & 16"-20" |
| 2B              | 683314      | JAY BACK UPH 18W TALL                             | 14"-18" & 16"-20" |
| 2B              | 683316      | JAY BACK UPH 19W TALL                             | 14"-18" & 16"-20" |
| 2B              | 683318      | JAY BACK UPH 20W TALL                             | 14"-18" & 16"-20" |
| 3               | 481800      | UNI LATERAL BRKT JAY PRECISION                    |                   |
| 4               | 226353      | 10-32 X 1 1/2 SOC HD C/S W/PL                     |                   |
| 5               | 481810      | UNI LATERAL BRKT SCR ROLL PIN                     |                   |
| 6               | 481811      | LATERAL BRKT PREC THD-ROLL PIN                    |                   |
| 7               | 261204      | #10 AN WASHER .063                                |                   |
| 8               | 226043      | 10-32 x 1/2 BH C/S BL                             |                   |
| 9               | 541144      | LATERAL BRKT 1.5 W WLDMNT                         | RIGHT             |
| 9               | 541145      | LATERAL BRKT 1.5 W WLDMNT                         | LEFT              |
| 10              | 881738      | S/A ASSM JAY PRECISION                            | RIGHT             |
| 10              | 881739      | S/A ASSM JAY PRECISION                            | LEFT              |
| 11              | 881109      | 1 NUT PLATE W/THREADED INSERT                     |                   |
| 12              | 500914      | SYMMETRIC LATERAL PANEL                           | Size 5"X5"        |
| 13              | 751314      | 2 x 2 VELCRO LOOP                                 |                   |
| 14              | B672        | BASE EXTENDED LATERAL                             |                   |
| 17              | 501140      | ASYMMETRIC LATERAL PANEL                          | Size 10"X3.25"    |
| 18              | 697100      | ASYMMETRIC LATERAL PANEL EDGE TRIM                |                   |
| 19              | 683319      | JAY BACK UPH INT LATERAL COVER ASYMMETRIC         | LEFT              |
| 19              | 683320      | JAY BACK UPH INT LATERAL COVER ASYMMETRIC         | RIGHT             |
| 3-6             | 880959      | UNI LATERAL CLAMP ASSM                            |                   |
| 20,21           | 923300      | JAY BACK EXTENSION LATERAL COVER ASYMMETRIC 14-16 |                   |
| 20,21           | 923301      | JAY BACK EXTENSION LATERAL COVER ASYMMETRIC 17-20 |                   |
| 15,16           | 923302      | JAY BACK EXTENSION LATERAL COVER SYMMETRIC 14-16  |                   |
| 15,16           | 923303      | JAY BACK EXTENSION LATERAL COVER SYMMETRIC 17-20  |                   |
| 1C,2            | 682520      | JAY MEDIUM BACK UPH 14 W/O LATERAL                |                   |
| 1C,2            | 682521      | JAY MEDIUM BACK UPH 15 W/O LATERAL                |                   |
| 1C,2            | 682522      | JAY MEDIUM BACK UPH 16 W/O LATERAL                |                   |
| 1C,2            | 682523      | JAY MEDIUM BACK UPH 17 W/O LATERAL                |                   |
| 1C,2            | 682524      | JAY MEDIUM BACK UPH 18 W/O LATERAL                |                   |
| 1C,2            | 682525      | JAY MEDIUM BACK UPH 19 W/O LATERAL                |                   |
| 1C,2            | 682526      | JAY MEDIUM BACK UPH 20 W/O LATERAL                |                   |
| 1D,2B           | 682527      | JAY TALL BACK UPH 14 W/O LATERAL                  |                   |
| 1D,2B           | 682528      | JAY TALL BACK UPH 15 W/O LATERAL                  |                   |
| 1D,2B           | 682529      | JAY TALL BACK UPH 16 W/O LATERAL                  |                   |
| 1D,2B           | 682530      | JAY TALL BACK UPH 17 W/O LATERAL                  |                   |
| 1D,2B           | 682531      | JAY TALL BACK UPH 18 W/O LATERAL                  |                   |
| 1D,2B           | 682532      | JAY TALL BACK UPH 19 W/O LATERAL                  |                   |
| 1D,2B           | 682533      | JAY TALL BACK UPH 20 W/O LATERAL                  |                   |
| 1A,2A,3-8,9-16  | 682534      | JAY MEDIUM BACK UPH 14 W/ FIX SYMMETRIC LATERAL   |                   |
| 1A,2A,3-8,9-16  | 682535      | JAY MEDIUM BACK UPH 15 W/ FIX SYMMETRIC LATERAL   |                   |
| 1A,2A,3-8,9-16  | 682536      | JAY MEDIUM BACK UPH 16 W/ FIX SYMMETRIC LATERAL   |                   |
| 1A,2A,3-8,9-16  | 682537      | JAY MEDIUM BACK UPH 17 W/ FIX SYMMETRIC LATERAL   |                   |
| 1A,2A,3-8,9-16  | 682538      | JAY MEDIUM BACK UPH 18 W/ FIX SYMMETRIC LATERAL   |                   |
| 1A,2A,3-8,9-16  | 682539      | JAY MEDIUM BACK UPH 19 W/ FIX SYMMETRIC LATERAL   |                   |
| 1A,2A,3-8,9-16  | 682540      | JAY MEDIUM BAC KUPH 20 W/ FIX SYMMETRIC LATERAL   |                   |
| 1A,2A,3-8,10-16 | 682541      | JAY MEDIUM BACK UPH 14 W/ S/A SYMMETRIC LATERAL   |                   |

| Pos.                  | Item Number | Description                                     | Remarks |
|-----------------------|-------------|-------------------------------------------------|---------|
| 1A,2A,3-8,10-16       | 682542      | JAY MEDIUM BACK UPH 15 W/ S/A SYMMETRIC LATERAL |         |
| 1A,2A,3-8,10-16       | 682543      | JAY MEDIUM BACK UPH 16 W/ S/A SYMMETRIC LATERAL |         |
| 1A,2A,3-8,10-16       | 682544      | JAY MEDIUM BACK UPH 17 W/ S/A SYMMETRIC LATERAL |         |
| 1A,2A,3-8,10-16       | 682545      | JAY MEDIUM BACK UPH 18 W/ S/A SYMMETRIC LATERAL |         |
| 1A,2A,3-8,10-16       | 682546      | JAY MEDIUM BACK UPH 19 W/ S/A SYMMETRIC LATERAL |         |
| 1A,2A,3-8,10-16       | 682547      | JAY MEDIUM BACK UPH 20 W/ S/A SYMMETRIC LATERAL |         |
| 1B,2B,3-9,11-16       | 682562      | JAY TALL BACK UPH 14 W/ FIX SYMMETRIC LATERAL   |         |
| 1B,2B,3-9,11-16       | 682563      | JAY TALL BACK UPH 15 W/ FIX SYMMETRIC LATERAL   |         |
| 1B,2B,3-9,11-16       | 682564      | JAY TALL BACK UPH 16 W/ FIX SYMMETRIC LATERAL   |         |
| 1B,2B,3-9,11-16       | 682565      | JAY TALL BACK UPH 17 W/ FIX SYMMETRIC LATERAL   |         |
| 1B,2B,3-9,11-16       | 682566      | JAY TALL BACK UPH 18 W/ FIX SYMMETRIC LATERAL   |         |
| 1B,2B,3-9,11-16       | 682567      | JAY TALL BACK UPH 19 W/ FIX SYMMETRIC LATERAL   |         |
| 1B,2B,3-9,11-16       | 682568      | JAY TALL BACK UPH 20 W/ FIX SYMMETRIC LATERAL   |         |
| 1B,2B,3-8,10-16       | 682569      | JAY TALL BACK UPH 14 W/ S/A SYMMETRIC LATERAL   |         |
| 1B,2B,3-8,10-16       | 682570      | JAY TALL BACK UPH 15 W/ S/A SYMMETRIC LATERAL   |         |
| 1B,2B,3-8,10-16       | 682571      | JAY TALL BACK UPH 16 W/ S/A SYMMETRIC LATERAL   |         |
| 1B,2B,3-8,10-16       | 682572      | JAY TALL BACK UPH 17 W/ S/A SYMMETRIC LATERAL   |         |
| 1B,2B,3-8,10-16       | 682573      | JAY TALL BACK UPH 18 W/ S/A SYMMETRIC LATERAL   |         |
| 1B,2B,3-8,10-16       | 682574      | JAY TALL BACK UPH 19 W/ S/A SYMMETRIC LATERAL   |         |
| 1B,2B,3-8,10-16       | 682575      | JAY TALL BACK UPH 20 W/ S/A SYMMETRIC LATERAL   |         |
| 1A,2A,3-9,11,17-21    | 682576      | JAY TALL BACK UPH 14 W/ FIX ASYMMETRIC LATERAL  |         |
| 1A,2A,3-9,11,17-21    | 682577      | JAY TALL BACK UPH 15 W/ FIX ASYMMETRIC LATERAL  |         |
| 1A,2A,3-9,11,17-21    | 682578      | JAY TALL BACK UPH 16 W/ FIX ASYMMETRIC LATERAL  |         |
| 1A,2A,3-9,11,17-21    | 682579      | JAY TALL BACK UPH 17 W/ FIX ASYMMETRIC LATERAL  |         |
| 1A,2A,3-9,11,17-21    | 682580      | JAY TALL BACK UPH 18 W/ FIX ASYMMETRIC LATERAL  |         |
| 1A,2A,3-9,11,17-21    | 682581      | JAY TALL BACK UPH 19 W/ FIX ASYMMETRIC LATERAL  |         |
| 1A,2A,3-9,11,17-21    | 682582      | JAY TALL BACK UPH 20 W/ FIX ASYMMETRIC LATERAL  |         |
| 1A,2A,3-8,10,11,17-21 | 682583      | JAY TALL BACK UPH 14 W/ S/A ASYMMETRIC LATERAL  |         |
| 1A,2A,3-8,10,11,17-21 | 682584      | JAY TALL BACK UPH 15 W/ S/A ASYMMETRIC LATERAL  |         |
| 1A,2A,3-8,10,11,17-21 | 682585      | JAY TALL BACK UPH 16 W/ S/A ASYMMETRIC LATERAL  |         |
| 1A,2A,3-8,10,11,17-21 | 682586      | JAY TALL BACK UPH 17 W/ S/A ASYMMETRIC LATERAL  |         |
| 1A,2A,3-8,10,11,17-21 | 682587      | JAY TALL BACK UPH 18 W/ S/A ASYMMETRIC LATERAL  |         |
| 1A,2A,3-8,10,11,17-21 | 682588      | JAY TALL BACK UPH 19 W/ S/A ASYMMETRIC LATERAL  |         |
| 1A,2A,3-8,10,11,17-21 | 682589      | JAY TALL BACK UPH 20 W/ S/A ASYMMETRIC LATERAL  |         |

| STUD WEDGE SIZE FOR:      | SMALL         | MEDIUM        | LARGE         |
|---------------------------|---------------|---------------|---------------|
| <u>MOST OUTSIDE DIM</u>   | 3/4"          | 7/8"          | 1"            |
| <u>AVERAGE INSIDE DIM</u> | 0.575 - 0.638 | 0.700 - 0.785 | 0.825 - 0.910 |



| Pos. | Item Number | Description                       | Remarks       |
|------|-------------|-----------------------------------|---------------|
| 1A   | 750696      | HEADREST W/ ZIPPER MEDIUM CG      | 8" X 4" Pad   |
| 1B   | 750697      | HEADREST W/ ZIPPER LARGE CG       | 9" X 4.5" Pad |
| 2    | 221115      | 1/4 - 20 X 3/4 SOC HD C/S BLK     |               |
| 3    | 341045      | #10 - 24 CONE POINT SET SCREW     |               |
| 4    | 505933      | HEADREST COLLAR                   |               |
| 5    | 341046      | HEADREST CLAMPING HANDLE SCREW    |               |
| 6    | 881003      | CLAMP AND TUBE                    |               |
| 7    | 261204      | #10 AN WASHER .063                |               |
| 8    | 221010      | 10 - 32 X .75 SHCS BLK W/ PL      |               |
| 9    | 226072      | 1/4 - 20 X 1/2 FH C/S BLK         |               |
| 10A  | 520855      | NUT WEDGE SMALL                   |               |
| 10B  | 520816      | NUT WEDGE MEDIUM                  |               |
| 10C  | 520856      | NUT WEDGE LARGE                   |               |
| 11A  | 520857      | STUD WEDGE SMALL                  |               |
| 11B  | 520815      | STUD WEDGE MEDIUM                 |               |
| 11C  | 520858      | STUD WEDGE LARGE                  |               |
| 12   | 226064      | 1/4 - 28 X 3/4 BH C/S BLK W/PATCH |               |
| 13   | 106943      | WASHER FLAT AN 1/4 ROHS           |               |
| 14   | 561004      | 1 SADDLE 1/4 HOLE                 |               |

| Pos.                 | Item Number | Description                                  | Remarks |
|----------------------|-------------|----------------------------------------------|---------|
| 15                   | 349595      | BEARING IGUS 3/8 X 15/32 X 1/2               |         |
| 16A                  | 542158      | HEADREST PVT TB WELDMOUNT 14 - 17            |         |
| 16B                  | 542159      | HEADREST PVT TB WELDMOUNT 17 - 20            |         |
| 17                   | 483342      | STUD WEDGE CAP                               |         |
| 18                   | 226351      | 1/4 - 20 X 4 - 1/2 SOC HD BLK ZINC           |         |
| 19                   | 500078      | INSERT D2 EXT TUBE                           |         |
| 21                   | 106943      | WASHER FLAT AN 1/4 ROHS                      |         |
| 22                   | 226095      | 1/4 - 20 X 4 - 3/8 SOC HD BLK ZINC           |         |
| 23                   | 750552      | ZGS FOLDING STRAP                            |         |
| 24                   | 226095      | 1/4 - 20 X 4 - 3/8 SOC HD BLK ZINC           |         |
| 25                   | 226076      | 1/4 - 28 X 1/2 FH C/S W/ PATCH               |         |
| 26                   | XM2025      | MNTNG PLATE/PEDIAT.HEAD REST                 |         |
| 27                   | 342046      | MEASURING TOOL FOLD HEADREST                 |         |
| 2-9                  | 883350      | HEADREST MOUNTING BRACKET ASSM EURO          |         |
| 20-24                | 888920      | CENTER MOUNT FOLD HEADREST ASSM              |         |
| 21-23                | 750104      | FOLD STABALIZER BAR STRAP ASSM               |         |
| 10, 11, 17, 18       | 881897      | STUD WEDGE ASSM SMALL                        |         |
| 10, 11, 17, 18       | 881898      | STUD WEDGE ASSM MEDIUM                       |         |
| 10, 11, 17, 18       | 881899      | STUD WEDGE ASSM LARGE                        |         |
| 10-15, 16A, 17-20    | 820750      | STABALIZER BAR ASSM 14 - 17 SMALL            |         |
| 10-15, 16A, 17-20    | 820751      | STABALIZER BAR ASSM 14 - 17 MEDIUM           |         |
| 10-15, 16A, 17-20    | 820752      | STABALIZER BAR ASSM 14 - 17 LARGE            |         |
| 10-15, 16B, 17-20    | 820753      | STABALIZER BAR ASSM 17 - 20 SMALL            |         |
| 10-15, 16B, 17-20    | 820754      | STABALIZER BAR ASSM 17 - 20 MEDIUM           |         |
| 10-15, 16B, 17-20    | 820755      | STABALIZER BAR ASSM 17 - 20 LARGE            |         |
| 1-15, 16A, 17-20     | 820760      | STABALIZER W/ MEDIUM HEADREST 14 - 17 SMALL  |         |
| 1-15, 16A, 17-20     | 820761      | STABALIZER W/ MEDIUM HEADREST 14 - 17 MEDIUM |         |
| 1-15, 16A, 17-20     | 820762      | STABALIZER W/ MEDIUM HEADREST 14 - 17 LARGE  |         |
| 1-15, 16B, 17-20     | 820763      | STABALIZER W/ MEDIUM HEADREST 17 - 20 SMALL  |         |
| 1-15, 16B, 17-20     | 820764      | STABALIZER W/ MEDIUM HEADREST 17 - 20 MEDIUM |         |
| 1-15, 16B, 17-20     | 820765      | STABALIZER W/ MEDIUM HEADREST 17 - 20 LARGE  |         |
| 1A, 2-15, 16A, 17-20 | 820766      | STABALIZER W/ LARGE HEADREST 14 - 17 SMALL   |         |
| 1A, 2-15, 16A, 17-20 | 820767      | STABALIZER W/ LARGE HEADREST 14 - 17 MEDIUM  |         |
| 1A, 2-15, 16A, 17-20 | 820768      | STABALIZER W/ LARGE HEADREST 14 - 17 LARGE   |         |
| 1A, 2-15, 16B, 17-20 | 820769      | STABALIZER W/ LARGE HEADREST 17 - 20 SMALL   |         |
| 1A, 2-15, 16B, 17-20 | 820770      | STABALIZER W/ LARGE HEADREST 17 - 20 MEDIUM  |         |
| 1A, 2-15, 16B, 17-20 | 820771      | STABALIZER W/ LARGE HEADREST 17 - 20 LARGE   |         |



| Pos. | Item Number | Description                   | Pos. | Item Number | Description                |
|------|-------------|-------------------------------|------|-------------|----------------------------|
| 1    | C331        | COVER BASIC INC 14X14         | 1    | C349        | COVER BASIC INC 24X24      |
| 1    | C308        | COVER BASIC INCONTINENT 14X16 | 2    | B331-MJ     | BASE BASIC CUSHION (14X14) |
| 1    | C332        | COVER BASIC INC 14X18         | 2    | B308-MJ     | FOAM BASE BASIC 14X16      |
| 1    | C356        | COVER BASIC INCONTINENT 15X16 | 2    | B332-MJ     | BASE BASIC CUSHION 14X18   |
| 1    | C333        | COVER BASIC INC 15X18         | 2    | B356        | FOAM BASE BASIC 15X16      |
| 1    | C305        | COVER BASIC INCONTINENT 16X16 | 2    | B333-MJ     | BASE BASIC CUSHION (15X18) |
| 1    | C334-MJ     | COVER BASIC INC 16X17         | 2    | B305        | FOAM BASE BASIC 16X16      |
| 1    | C306-MJ     | COVER BASIC INCONTINENT 16X18 | 2    | B334-MJ     | BASE BASIC CUSHION (16X17) |
| 1    | C376        | COVER BASIC INCONTINENT 17X16 | 2    | B306        | FOAM BASE BASIC 16X18      |
| 1    | C335        | COVER BASIC INC 17X17         | 2    | B376        | FOAM BASE BASIC 17X16      |
| 1    | C336        | COVER BASIC INC 17X18         | 2    | B335-MJ     | BASE BASIC CUSHION (17X17) |
| 1    | C300-MJ     | COVER BASIC INCONTINENT 18X16 | 2    | B336        | BASE BASIC CUSHION (17X18) |
| 1    | C337        | COVER BASIC INC 18X17         | 2    | B300        | FOAM BASE BASIC 18X16      |
| 1    | C301        | COVER BASIC INCONTINENT 18X18 | 2    | B337        | BASE BASIC CUSHION (18X17) |
| 1    | C338        | COVER BASIC INC 18X20         | 2    | B301        | FOAM BASE BASIC 18X18      |
| 1    | C339        | COVER BASIC INC 19X16         | 2    | B338        | BASE BASIC CUSHION (18X20) |
| 1    | C341        | COVER BASIC INC 19X18         | 2    | B339        | BASE BASIC CUSHION (19X16) |
| 1    | C343        | COVER BASIC INC 19X20         | 2    | B341-MJ     | BASE BASIC CUSHION (19X18) |
| 1    | C344        | COVER BASIC INC 20X16         | 2    | B343-MJ     | BASE BASIC CUSHION (19X20) |
| 1    | C309        | COVER BASIC 20X18 INC         | 2    | B344-MJ     | BASE BASIC CUSHION (20X16) |
| 1    | C345        | COVER BASIC INC 20X20         | 2    | B309        | FOAM BASE BASIC 20X18      |
| 1    | C328        | COVER BASIC 22X18 INC         | 2    | B345        | BASE BASIC CUSHION 20X20   |
| 1    | C320        | COVER BASIC 22X20 INC         | 2    | B328        | FOAM BASE 22X18 BASIC      |
| 1    | C346        | COVER BASIC INC 22X22         | 2    | B320        | FOAM BASE 22X20 BASIC      |
| 1    | C348        | COVER BASIC 24X18 INC         | 2    | B346        | BASE BASIC CUSHION 22X22   |
| 1    | C340        | COVER BASIC 24X20 INC         | 2    | B348        | FOAM BASE 24X18 BASIC      |
| 1    | C347        | COVER BASIC INC 24X22         | 2    | B340        | FOAM BASE 24X20 BASIC      |

| Pos. | Item Number | Description              | Pos. | Item Number | Description                |
|------|-------------|--------------------------|------|-------------|----------------------------|
| 2    | B347        | BASE BASIC CUSHION 24X22 | 2    | B349-MJ     | BASE BASIC CUSHION (24X24) |




| Pos. | Item Number   | Description                 | Pos. | Item Number    | Description                 |
|------|---------------|-----------------------------|------|----------------|-----------------------------|
| 1    | CJBASICP1414  | OUTER CVR BASIC PRO 14X14   | 1    | CJBASICP1918   | OUTER CVR BASIC PRO 19X18   |
| 1    | CJBASICP1416  | OUTER CVR BASIC PRO 14X16   | 1    | CJBASICP1920   | OUTER CVR BASIC PRO 19X20   |
| 1    | CJBASICP1418  | OUTER CVR BASIC PRO 14X18   | 1    | CJBASICP2016   | OUTER CVR BASIC PRO 20X16   |
| 1    | CJBASICP1518  | OUTER CVR BASIC PRO 15X18   | 1    | CJBASICP2018   | OUTER CVR BASIC PRO 20X18   |
| 1    | CJBASICP15516 | OUTER CVR BASIC PRO 15.5X16 | 1    | CJBASICP2020   | OUTER CVR BASIC PRO 20X20   |
| 1    | CJBASICP15518 | OUTER CVR BASIC PRO 15.5X18 | 1    | CJBASICP2218   | OUTER CVR BASIC PRO 22X18   |
| 1    | CJBASICP1617  | OUTER CVR BASIC PRO 16X17   | 1    | CJBASICP2220   | OUTER CVR BASIC PRO 22X20   |
| 1    | CJBASICP1618  | OUTER CVR BASIC PRO 16X18   | 1    | CJBASICP2222   | OUTER CVR BASIC PRO 22X22   |
| 1    | CJBASICP1620  | OUTER CVR BASIC PRO 16X20   | 1    | CJBASICP2418   | OUTER CVR BASIC PRO 24X18   |
| 1    | CJBASICP1716  | OUTER CVR BASIC PRO 17X16   | 1    | CJBASICP2420   | OUTER CVR BASIC PRO 24X20   |
| 1    | CJBASICP1717  | OUTER CVR BASIC PRO 17X17   | 1    | CJBASICP2422   | OUTER CVR BASIC PRO 24X22   |
| 1    | CJBASICP1718  | OUTER CVR BASIC PRO 17X18   | 1    | CJBASICP2424   | OUTER CVR BASIC PRO 24X24   |
| 1    | CJBASICP1816  | OUTER CVR BASIC PRO 18X16   | 2    | ICJBASICP1414  | INNER CVR BASIC PRO 14X14   |
| 1    | CJBASICP1817  | OUTER CVR BASIC PRO 18X17   | 2    | ICJBASICP1416  | INNER CVR BASIC PRO 14X16   |
| 1    | CJBASICP1818  | OUTER CVR BASIC PRO 18X18   | 2    | ICJBASICP1418  | INNER CVR BASIC PRO 14X18   |
| 1    | CJBASICP1820  | OUTER CVR BASIC PRO 18X20   | 2    | ICJBASICP1518  | INNER CVR BASIC PRO 15X18   |
| 1    | CJBASICP1916  | OUTER CVR BASIC PRO 19X16   | 2    | ICJBASICP15516 | INNER CVR BASIC PRO 15.5X16 |
| 1    | CJBASICP1917  | OUTER CVR BASIC PRO 19X17   | 2    | ICJBASICP15518 | INNER CVR BASIC PRO 15.5X18 |

| Pos. | Item Number   | Description               | Pos. | Item Number   | Description            |
|------|---------------|---------------------------|------|---------------|------------------------|
| 2    | ICJBASICP1617 | INNER CVR BASIC PRO 16X17 | 3    | BJBASICP1518  | BASE BASIC PRO 15X18   |
| 2    | ICJBASICP1618 | INNER CVR BASIC PRO 16X18 | 3    | BJBASICP15516 | BASE BASIC PRO 15.5X16 |
| 2    | ICJBASICP1620 | INNER CVR BASIC PRO 16X20 | 3    | BJBASICP15518 | BASE BASIC PRO 15.5X18 |
| 2    | ICJBASICP1716 | INNER CVR BASIC PRO 17X16 | 3    | BJBASICP1617  | BASE BASIC PRO 16X17   |
| 2    | ICJBASICP1717 | INNER CVR BASIC PRO 17X17 | 3    | BJBASICP1618  | BASE BASIC PRO 16X18   |
| 2    | ICJBASICP1718 | INNER CVR BASIC PRO 17X18 | 3    | BJBASICP1620  | BASE BASIC PRO 16X20   |
| 2    | ICJBASICP1816 | INNER CVR BASIC PRO 18X16 | 3    | BJBASICP1716  | BASE BASIC PRO 17X16   |
| 2    | ICJBASICP1817 | INNER CVR BASIC PRO 18X17 | 3    | BJBASICP1717  | BASE BASIC PRO 17X17   |
| 2    | ICJBASICP1818 | INNER CVR BASIC PRO 18X18 | 3    | BJBASICP1718  | BASE BASIC PRO 17X18   |
| 2    | ICJBASICP1820 | INNER CVR BASIC PRO 18X20 | 3    | BJBASICP1816  | BASE BASIC PRO 18X16   |
| 2    | ICJBASICP1916 | INNER CVR BASIC PRO 19X16 | 3    | BJBASICP1817  | BASE BASIC PRO 18X17   |
| 2    | ICJBASICP1917 | INNER CVR BASIC PRO 19X17 | 3    | BJBASICP1818  | BASE BASIC PRO 18X18   |
| 2    | ICJBASICP1918 | INNER CVR BASIC PRO 19X18 | 3    | BJBASICP1820  | BASE BASIC PRO 18X20   |
| 2    | ICJBASICP1920 | INNER CVR BASIC PRO 19X20 | 3    | BJBASICP1916  | BASE BASIC PRO 19X16   |
| 2    | ICJBASICP2016 | INNER CVR BASIC PRO 20X16 | 3    | BJBASICP1917  | BASE BASIC PRO 19X17   |
| 2    | ICJBASICP2018 | INNER CVR BASIC PRO 20X18 | 3    | BJBASICP1918  | BASE BASIC PRO 19X18   |
| 2    | ICJBASICP2020 | INNER CVR BASIC PRO 20X20 | 3    | BJBASICP1920  | BASE BASIC PRO 19X20   |
| 2    | ICJBASICP2218 | INNER CVR BASIC PRO 22X18 | 3    | BJBASICP2016  | BASE BASIC PRO 20X16   |
| 2    | ICJBASICP2220 | INNER CVR BASIC PRO 22X20 | 3    | BJBASICP2018  | BASE BASIC PRO 20X18   |
| 2    | ICJBASICP2222 | INNER CVR BASIC PRO 22X22 | 3    | BJBASICP2020  | BASE BASIC PRO 20X20   |
| 2    | ICJBASICP2418 | INNER CVR BASIC PRO 24X18 | 3    | BJBASICP2218  | BASE BASIC PRO 22X18   |
| 2    | ICJBASICP2420 | INNER CVR BASIC PRO 24X20 | 3    | BJBASICP2220  | BASE BASIC PRO 22X20   |
| 2    | ICJBASICP2422 | INNER CVR BASIC PRO 24X22 | 3    | BJBASICP2222  | BASE BASIC PRO 22X22   |
| 2    | ICJBASICP2424 | INNER CVR BASIC PRO 24X24 | 3    | BJBASICP2418  | BASE BASIC PRO 24X18   |
| 3    | BJBASICP1414  | BASE BASIC PRO 14X14      | 3    | BJBASICP2420  | BASE BASIC PRO 24X20   |
| 3    | BJBASICP1416  | BASE BASIC PRO 14X16      | 3    | BJBASICP2422  | BASE BASIC PRO 24X22   |
| 3    | BJBASICP1418  | BASE BASIC PRO 14X18      | 3    | BJBASICP2424  | BASE BASIC PRO 24X24   |

| PART # | CHAIR<br>SIZE | ACTUAL<br>WIDTH | ACTUAL<br>LENGTH | FEMORAL<br>SUPPORT LENGTH | PELVIC CONTOUR<br>DEPTH | PELVIC CONTOUR<br>WIDTH | OVERALL<br>HEIGHT | UNLOADED<br>WELL DEPTH |
|--------|---------------|-----------------|------------------|---------------------------|-------------------------|-------------------------|-------------------|------------------------|
|        | A             | B               | C                | D                         | E                       | F                       | G                 |                        |
| 3100B  | 18x17.5       | 18.5            | 17.5             | 6                         | 8                       | 17                      | 3.75              | 3.1                    |
| 3101B  | 18x19.5       | 19.25           | 19.75            | 7.5                       | 9                       | 17                      | 3.75              | 2.9                    |
| 3105B  | 16x17.5       | 17              | 17.75            | 6.25                      | 8                       | 15                      | 4                 | 2.9                    |
| 3106B  | 16x19.5       | 17.5            | 19.75            | 7                         | 9                       | 15.5                    | 4                 | 2.9                    |
| 3109B  | 20x19.5       | 20.75           | 17.5             | 6.5                       | 7.75                    | 19.5                    | 3.75              | 3.1                    |

| Pos. | Item Number | Description       |
|------|-------------|-------------------|
| 1    | C3100B      | COVER JCARE 18X16 |
| 1    | C3101B      | COVER JCARE 18X18 |
| 1    | C3105B      | COVER JCARE 16X16 |
| 1    | C3106B      | COVER JCARE 16X18 |
| 1    | C3109B      | COVER JCARE 20X18 |




| Pos. | Item Number | Description                        | Remarks                |
|------|-------------|------------------------------------|------------------------|
| 1    | C2200P      | SOFT COMBI POSITIONING COVER 18x16 |                        |
| 1    | C2201P      | SOFT COMBI POSITIONING COVER 18x18 |                        |
| 1    | C2202P      | SOFT COMBI POSITIONING COVER 20x18 |                        |
| 1    | C2203P      | SOFT COMBI POSITIONING COVER 18x20 |                        |
| 1    | C2205P      | SOFT COMBI POSITIONING COVER 16x16 |                        |
| 1    | C2206P      | SOFT COMBI POSITIONING COVER 16x18 | actual size 15.5" x18" |
| 1    | C2207P      | SOFT COMBI POSITIONING COVER 16x20 |                        |
| 1    | C2210P      | SOFT COMBI POSITIONING COVER 14x14 |                        |

| Pos. | Item Number | Description                | Remarks |
|------|-------------|----------------------------|---------|
| 1    | C2228P      | COVER 22X18 SOFT COMBI POS |         |
| 1    | C2220P      | COVER 22X20 SOFT COMBI POS |         |
| 1    | C2248P      | COVER 24X18 SOFT COMBI POS |         |
| 1    | C2240P      | COVER 24X20 SOFT COMBI POS |         |
| 3    | S200        | COMBI SOLID SEAT 18x16     |         |
| 3    | S201        | COMBI SOLID SEAT 18x18     |         |
| 3    | S202        | COMBI SOLID SEAT 20x18     |         |
| 3    | S203        | COMBI SOLID SEAT 18x20     |         |
| 3    | S205        | COMBI SOLID SEAT 15.5x16   |         |
| 3    | S206        | COMBI SOLID SEAT 15.5x18   |         |
| 3    | S207        | COMBI SOLID SEAT 16x20     |         |
| 3    | S210        | COMBI SOLID SEAT 14x14     |         |
| 3    | S228        | COMBI SOLID SEAT 22X18     |         |
| 3    | S220        | COMBI SOLID SEAT 22X20     |         |
| 3    | S248        | COMBI SOLID SEAT 24X18     |         |
| 3    | S240        | COMBI SOLID SEAT 24X20     |         |

NOTE: \* FOR WHEELCHAIRS THAT ARE SLIGHTLY NARROWER THAN LISTED DIMENSION, SUCH AS WITH QUICKIE WHEELCHAIRS. THESE PANS ALSO PROVIDE ADDITIONAL CLEARANCE FOR POP-UP FOOTRESTS.


NOTE: \*\* WHEN ORDERING A DUO FOR A QUICKIE 2 WHEELCHAIR, A SEAT RAIL INSERT STRIP (ITEM 7) MUST BE ORDERED (QUANTITY 2) TO PREVENT THE CLIP FROM CATCHING IN THE UPHOLSTERY SLOT.



| Pos. | Item Number | Description           | Remarks |
|------|-------------|-----------------------|---------|
| 1    | C700        | JAY DUO COVER 18 x 16 |         |
| 1    | C701        | JAY DUO COVER 18 x 18 |         |
| 1    | C705        | JAY DUO COVER 16 x 16 |         |
| 1    | C706        | JAY DUO COVER 16 x 18 |         |
| 1    | C709        | JAY DUO COVER 20 x 18 |         |
| 1    | C726        | JAY DUO COVER 20 x 16 |         |
| 1    | C776        | JAY DUO COVER 17 x 16 |         |
| 1    | C778        | JAY DUO COVER 17 x 18 |         |
| 2    | F700        | JAY DUO PAD 18 x 16   |         |
| 2    | F701        | JAY DUO PAD 18 x 18   |         |
| 2    | F705        | JAY DUO PAD 16 x 16   |         |


| Pos.        | Item Number | Description                              | Remarks              |
|-------------|-------------|------------------------------------------|----------------------|
| 2           | F706        | JAY DUO PAD 16 x 18                      |                      |
| 2           | F709        | JAY DUO PAD 20 x 18                      |                      |
| 2           | F726        | JAY DUO PAD 20 x 16                      |                      |
| 2           | F776        | JAY DUO PAD 17 x 16                      |                      |
| 2           | F778        | JAY DUO PAD 17 x 18                      |                      |
| 3A,5A, 5B,6 | S700        | JAY DUO SEAT PAN 18 x 16                 |                      |
| 3A,5A, 5B,6 | S701        | JAY DUO SEAT PAN 18 x 18                 |                      |
| 3A,5A, 5B,6 | S705        | JAY DUO SEAT PAN 16 x 16                 |                      |
| 3A,5A, 5B,6 | S706        | JAY DUO SEAT PAN 16 x 18                 |                      |
| 3A,5A, 5B,6 | S709        | JAY DUO SEAT PAN 20 x 18                 |                      |
| 3A,5A, 5B,6 | S726        | JAY DUO SEAT PAN 20 x 16                 |                      |
| 3A,5A, 5B,6 | S776        | JAY DUO SEAT PAN 17 x 16                 |                      |
| 3A,5A, 5B,6 | S778        | JAY DUO SEAT PAN 17 x 18                 |                      |
| 3B,5A, 5B,6 | S700Q       | JAY DUO SEAT PAN 18 x 16 (FOR 18 x 16)   | * SEE NOTE           |
| 3B,5A, 5B,6 | S701Q       | JAY DUO SEAT PAN 17.5 x 18 (FOR 18 x 18) | * SEE NOTE           |
| 3B,5A, 5B,6 | S705Q       | JAY DUO SEAT PAN 15.5 x 16 (FOR 16 x 16) | * SEE NOTE           |
| 3B,5A, 5B,6 | S706Q       | JAY DUO SEAT PAN 15.5 x 18 (FOR 16 x 18) | * SEE NOTE           |
| 3B,5A, 5B,6 | S709Q       | JAY DUO SEAT PAN 19.5 x 18 (FOR 20 x 18) | * SEE NOTE           |
| 3B,5A, 5B,6 | S726Q       | JAY DUO SEAT PAN 20 x 16 (FOR 20 x 16)   | * SEE NOTE           |
| 3B,5A, 5B,6 | S776Q       | JAY DUO SEAT PAN 17 x 16 (FOR 17 x 16)   | * SEE NOTE           |
| 3B,5A, 5B,6 | S778Q       | JAY DUO SEAT PAN 17 x 18 (FOR 17 x 18)   | * SEE NOTE           |
| 4           | M-726       | 5/8" RUBBER BUMPER DUO KIT               | SHIPS AS A SET OF 4. |
| 5A,6        | M-723       | 1" NYLON RETAINER CLIP DUO               | SHIPS AS A SET OF 4. |
| 5B,6        | M-724       | 7/8" NYLON RETAINER CLIP DUO             |                      |
| 7           | 580026      | SEAT RAIL INSERT STRIP 16"               | ** SEE NOTE          |
| 7           | 580027      | SEAT RAIL INSERT STRIP 17"               | ** SEE NOTE          |
| 7           | 580028      | SEAT RAIL INSERT STRIP 18"               | ** SEE NOTE          |
| 7           | 580029      | SEAT RAIL INSERT STRIP 19"               | ** SEE NOTE          |
| 7           | 580030      | SEAT RAIL INSERT STRIP 20"               | ** SEE NOTE          |

NOTE: IN GENERAL, THE (F) FLAT BASE IS USED FOR HARD SEATING SURFACES, SUCH AS A SEAT PAN, AND THE CURVED BASE IS USED FOR SOFT SURFACES, SUCH AS A SEAT SLING.



| Pos. | Item Number | Description                 | Pos. | Item Number | Description                   |
|------|-------------|-----------------------------|------|-------------|-------------------------------|
| >>>> |             | JAY EASY COVER FLAT         | 1    | CJE2020F    | COVER JAY EASY FLAT 20 x 20   |
| 1    | CJE1416F    | COVER JAY EASY FLAT 14 x 16 | 1    | CJE2218F    | COVER JAY EASY FLAT 22 x 18   |
| 1    | CJE1418F    | COVER JAY EASY FLAT 14 x 18 | 1    | CJE2220F    | COVER JAY EASY FLAT 22 x 20   |
| 1    | CJE1516F    | COVER JAY EASY FLAT 15 x 16 | 1    | CJE2222F    | COVER JAY EASY FLAT 22 x 22   |
| 1    | CJE1517F    | COVER JAY EASY FLAT 15 x 17 | 1    | CJE2418F    | COVER JAY EASY FLAT 24 x 18   |
| 1    | CJE1616F    | COVER JAY EASY FLAT 16 x 16 | 1    | CJE2420F    | COVER JAY EASY FLAT 24 x 20   |
| 1    | CJE1617F    | COVER JAY EASY FLAT 16 x 17 | 1    | CJE2422F    | COVER JAY EASY FLAT 24 x 22   |
| 1    | CJE1618F    | COVER JAY EASY FLAT 16 x 18 | 1    | CJE2424F    | COVER JAY EASY FLAT 24 x 24   |
| 1    | CJE1816F    | COVER JAY EASY FLAT 18 x 16 | >>>> |             | JAY EASY COVER CURVED         |
| 1    | CJE1817F    | COVER JAY EASY FLAT 18 x 17 | 1A   | CJE1416C    | COVER JAY EASY CURVED 14 x 16 |
| 1    | CJE1818F    | COVER JAY EASY FLAT 18 x 18 | 1A   | CJE1418C    | COVER JAY EASY CURVED 14 x 18 |
| 1    | CJE1820F    | COVER JAY EASY FLAT 18 x 20 | 1A   | CJE1516C    | COVER JAY EASY CURVED 15 x 16 |
| 1    | CJE2016F    | COVER JAY EASY FLAT 20 x 16 | 1A   | CJE1517C    | COVER JAY EASY CURVED 15 x 17 |
| 1    | CJE2018F    | COVER JAY EASY FLAT 20 x 18 | 1A   | CJE1616C    | COVER JAY EASY CURVED 16 x 16 |

| Pos. | Item Number | Description                   | Pos. | Item Number | Description                  |
|------|-------------|-------------------------------|------|-------------|------------------------------|
| 1A   | CJE1617C    | COVER JAY EASY CURVED 16 x 17 | 3    | BJE2422F    | BASE JAY EASY FLAT 24 x 22   |
| 1A   | CJE1618C    | COVER JAY EASY CURVED 16 x 18 | 3    | BJE2424F    | BASE JAY EASY FLAT 24 x 24   |
| 1A   | CJE1816C    | COVER JAY EASY CURVED 18 x 16 | >>>> |             | JAY EASY BASE CURVED         |
| 1A   | CJE1817C    | COVER JAY EASY CURVED 18 x 17 | 3A   | BJE1416C    | BASE JAY EASY CURVED 14 x 16 |
| 1A   | CJE1818C    | COVER JAY EASY CURVED 18 x 18 | 3A   | BJE1418C    | BASE JAY EASY CURVED 14 x 18 |
| 1A   | CJE1820C    | COVER JAY EASY CURVED 18 x 20 | 3A   | BJE1516C    | BASE JAY EASY CURVED 15 x 16 |
| 1A   | CJE2016C    | COVER JAY EASY CURVED 20 x 16 | 3A   | BJE1517C    | BASE JAY EASY CURVED 15 x 17 |
| 1A   | CJE2018C    | COVER JAY EASY CURVED 20 x 18 | 3A   | BJE1616C    | BASE JAY EASY CURVED 16 x 16 |
| 1A   | CJE2020C    | COVER JAY EASY CURVED 20 x 20 | 3A   | BJE1617C    | BASE JAY EASY CURVED 16 x 17 |
| 1A   | CJE2218C    | COVER JAY EASY CURVED 22 x 18 | 3A   | BJE1618C    | BASE JAY EASY CURVED 16 x 18 |
| 1A   | CJE2220C    | COVER JAY EASY CURVED 22 x 20 | 3A   | BJE1816C    | BASE JAY EASY CURVED 18 x 16 |
| 1A   | CJE2222C    | COVER JAY EASY CURVED 22 x 22 | 3A   | BJE1817C    | BASE JAY EASY CURVED 18 x 17 |
| 1A   | CJE2418C    | COVER JAY EASY CURVED 24 x 18 | 3A   | BJE1818C    | BASE JAY EASY CURVED 18 x 18 |
| 1A   | CJE2420C    | COVER JAY EASY CURVED 24 x 20 | 3A   | BJE1820C    | BASE JAY EASY CURVED 18 x 20 |
| 1A   | CJE2422C    | COVER JAY EASY CURVED 24 x 22 | 3A   | BJE2016C    | BASE JAY EASY CURVED 20 x 16 |
| 1A   | CJE2424C    | COVER JAY EASY CURVED 24 x 24 | 3A   | BJE2018C    | BASE JAY EASY CURVED 20 x 18 |
| >>>> |             | JAY EASY FLUID PAD            | 3A   | BJE2020C    | BASE JAY EASY CURVED 20 x 20 |
| 2    | FJE14-15W   | FLUID PAD JAY EASY 14-15      | 3A   | BJE2218C    | BASE JAY EASY CURVED 22 x 18 |
| 2    | FJE16-18W   | FLUID PAD JAY EASY 16-18      | 3A   | BJE2220C    | BASE JAY EASY CURVED 22 x 20 |
| 2    | FJE20W      | FLUID PAD JAY EASY 20         | 3A   | BJE2222C    | BASE JAY EASY CURVED 22 x 22 |
| >>>> |             | JAY EASY BASE FLAT            | 3A   | BJE2418C    | BASE JAY EASY CURVED 24 x 18 |
| 3    | BJE1416F    | BASE JAY EASY FLAT 14 x 16    | 3A   | BJE2420C    | BASE JAY EASY CURVED 24 x 20 |
| 3    | BJE1418F    | BASE JAY EASY FLAT 14 x 18    | 3A   | BJE2422C    | BASE JAY EASY CURVED 24 x 22 |
| 3    | BJE1516F    | BASE JAY EASY FLAT 15 x 16    | 3A   | BJE2424C    | BASE JAY EASY CURVED 24 x 24 |
| 3    | BJE1517F    | BASE JAY EASY FLAT 15 x 17    | >>>> |             | SOLID SEAT - JAY EASY        |
| 3    | BJE1616F    | BASE JAY EASY FLAT 16 x 16    | 4    | SJE1416     | SOLID SEAT JAY EASY 14 x 16  |
| 3    | BJE1617F    | BASE JAY EASY FLAT 16 x 17    | 4    | SJE1418     | SOLID SEAT JAY EASY 14 x 18  |
| 3    | BJE1618F    | BASE JAY EASY FLAT 16 x 18    | 4    | SJE1516     | SOLID SEAT JAY EASY 15 x 16  |
| 3    | BJE1816F    | BASE JAY EASY FLAT 18 x 16    | 4    | SJE1517     | SOLID SEAT JAY EASY 15 x 17  |
| 3    | BJE1817F    | BASE JAY EASY FLAT 18 x 17    | 4    | SJE1616     | SOLID SEAT JAY EASY 16 x 16  |
| 3    | BJE1818F    | BASE JAY EASY FLAT 18 x 18    | 4    | SJE1617     | SOLID SEAT JAY EASY 16 x 17  |
| 3    | BJE1820F    | BASE JAY EASY FLAT 18 x 20    | 4    | SJE1618     | SOLID SEAT JAY EASY 16 x 18  |
| 3    | BJE2016F    | BASE JAY EASY FLAT 20 x 16    | 4    | SJE1816     | SOLID SEAT JAY EASY 18 x 16  |
| 3    | BJE2018F    | BASE JAY EASY FLAT 20 x 18    | 4    | SJE1817     | SOLID SEAT JAY EASY 18 x 17  |
| 3    | BJE2020F    | BASE JAY EASY FLAT 20 x 20    | 4    | SJE1818     | SOLID SEAT JAY EASY 18 x 18  |
| 3    | BJE2218F    | BASE JAY EASY FLAT 22 x 18    | 4    | SJE1820     | SOLID SEAT JAY EASY 18 x 20  |
| 3    | BJE2220F    | BASE JAY EASY FLAT 22 x 20    | 4    | SJE2016     | SOLID SEAT JAY EASY 20 x 16  |
| 3    | BJE2222F    | BASE JAY EASY FLAT 22 x 22    | 4    | SJE2018     | SOLID SEAT JAY EASY 20 x 18  |
| 3    | BJE2418F    | BASE JAY EASY FLAT 24 x 18    | 4    | SJE2020     | SOLID SEAT JAY EASY 20 x 20  |
| 3    | BJE2420F    | BASE JAY EASY FLAT 24 x 20    |      |             |                              |



| Pos. | Item Number | Description         | Remarks               |
|------|-------------|---------------------|-----------------------|
| 1-3  | JFUSIONEK   | JAY FUSION EVAL KIT | Order part# JFUSIONEK |

Note: Cushion widths 14" & 15" use A size Fluid Bladders, 16" - 20" use B size bladders and 21"-24" use C size bladders.



| Pos. | Item Number  | Description                 | Remarks |
|------|--------------|-----------------------------|---------|
| 1    | CJFUSION1414 | OUTER CVR JFUSION1414 FLUID |         |
| 1    | CJFUSION1415 | OUTER CVR JFUSION1415 FLUID |         |
| 1    | CJFUSION1416 | OUTER CVR JFUSION1416 FLUID |         |
| 1    | CJFUSION1417 | OUTER CVR JFUSION1417 FLUID |         |
| 1    | CJFUSION1418 | OUTER CVR JFUSION1418 FLUID |         |
| 1    | CJFUSION1514 | OUTER CVR JFUSION1514 FLUID |         |
| 1    | CJFUSION1515 | OUTER CVR JFUSION1515 FLUID |         |
| 1    | CJFUSION1516 | OUTER CVR JFUSION1516 FLUID |         |
| 1    | CJFUSION1517 | OUTER CVR JFUSION1517 FLUID |         |
| 1    | CJFUSION1518 | OUTER CVR JFUSION1518 FLUID |         |
| 1    | CJFUSION1519 | OUTER CVR JFUSION1519 FLUID |         |
| 1    | CJFUSION1520 | OUTER CVR JFUSION1520 FLUID |         |
| 1    | CJFUSION1521 | OUTER CVR JFUSION1521 FLUID |         |

| Pos. | Item Number  | Description                 | Remarks |
|------|--------------|-----------------------------|---------|
| 1    | CJFUSION1614 | OUTER CVR JFUSION1614 FLUID |         |
| 1    | CJFUSION1615 | OUTER CVR JFUSION1615 FLUID |         |
| 1    | CJFUSION1616 | OUTER CVR JFUSION1616 FLUID |         |
| 1    | CJFUSION1617 | OUTER CVR JFUSION1617 FLUID |         |
| 1    | CJFUSION1618 | OUTER CVR JFUSION1618 FLUID |         |
| 1    | CJFUSION1619 | OUTER CVR JFUSION1619 FLUID |         |
| 1    | CJFUSION1620 | OUTER CVR JFUSION1620 FLUID |         |
| 1    | CJFUSION1621 | OUTER CVR JFUSION1621 FLUID |         |
| 1    | CJFUSION1622 | OUTER CVR JFUSION1622 FLUID |         |
| 1    | CJFUSION1714 | OUTER CVR JFUSION1714 FLUID |         |
| 1    | CJFUSION1715 | OUTER CVR JFUSION1715 FLUID |         |
| 1    | CJFUSION1716 | OUTER CVR JFUSION1716 FLUID |         |
| 1    | CJFUSION1717 | OUTER CVR JFUSION1717 FLUID |         |
| 1    | CJFUSION1718 | OUTER CVR JFUSION1718 FLUID |         |
| 1    | CJFUSION1719 | OUTER CVR JFUSION1719 FLUID |         |
| 1    | CJFUSION1720 | OUTER CVR JFUSION1720 FLUID |         |
| 1    | CJFUSION1721 | OUTER CVR JFUSION1721 FLUID |         |
| 1    | CJFUSION1722 | OUTER CVR JFUSION1722 FLUID |         |
| 1    | CJFUSION1723 | OUTER CVR JFUSION1723 FLUID |         |
| 1    | CJFUSION1814 | OUTER CVR JFUSION1814 FLUID |         |
| 1    | CJFUSION1815 | OUTER CVR JFUSION1815 FLUID |         |
| 1    | CJFUSION1816 | OUTER CVR JFUSION1816 FLUID |         |
| 1    | CJFUSION1817 | OUTER CVR JFUSION1817 FLUID |         |
| 1    | CJFUSION1818 | OUTER CVR JFUSION1818 FLUID |         |
| 1    | CJFUSION1819 | OUTER CVR JFUSION1819 FLUID |         |
| 1    | CJFUSION1820 | OUTER CVR JFUSION1820 FLUID |         |
| 1    | CJFUSION1821 | OUTER CVR JFUSION1821 FLUID |         |
| 1    | CJFUSION1822 | OUTER CVR JFUSION1822 FLUID |         |
| 1    | CJFUSION1823 | OUTER CVR JFUSION1823 FLUID |         |
| 1    | CJFUSION1824 | OUTER CVR JFUSION1824 FLUID |         |
| 1    | CJFUSION1915 | OUTER CVR JFUSION1915 FLUID |         |
| 1    | CJFUSION1916 | OUTER CVR JFUSION1916 FLUID |         |
| 1    | CJFUSION1917 | OUTER CVR JFUSION1917 FLUID |         |
| 1    | CJFUSION1918 | OUTER CVR JFUSION1918 FLUID |         |
| 1    | CJFUSION1919 | OUTER CVR JFUSION1919 FLUID |         |
| 1    | CJFUSION1920 | OUTER CVR JFUSION1920 FLUID |         |
| 1    | CJFUSION1921 | OUTER CVR JFUSION1921 FLUID |         |
| 1    | CJFUSION1922 | OUTER CVR JFUSION1922 FLUID |         |
| 1    | CJFUSION1923 | OUTER CVR JFUSION1923 FLUID |         |
| 1    | CJFUSION1924 | OUTER CVR JFUSION1924 FLUID |         |
| 1    | CJFUSION2015 | OUTER CVR JFUSION2015 FLUID |         |
| 1    | CJFUSION2016 | OUTER CVR JFUSION2016 FLUID |         |
| 1    | CJFUSION2017 | OUTER CVR JFUSION2017 FLUID |         |
| 1    | CJFUSION2018 | OUTER CVR JFUSION2018 FLUID |         |
| 1    | CJFUSION2019 | OUTER CVR JFUSION2019 FLUID |         |
| 1    | CJFUSION2020 | OUTER CVR JFUSION2020 FLUID |         |
| 1    | CJFUSION2021 | OUTER CVR JFUSION2021 FLUID |         |
| 1    | CJFUSION2022 | OUTER CVR JFUSION2022 FLUID |         |
| 1    | CJFUSION2023 | OUTER CVR JFUSION2023 FLUID |         |
| 1    | CJFUSION2024 | OUTER CVR JFUSION2024 FLUID |         |
| 1    | CJFUSION2117 | OUTER CVR JFUSION2117 FLUID |         |
| 1    | CJFUSION2118 | OUTER CVR JFUSION2118 FLUID |         |
| 1    | CJFUSION2119 | OUTER CVR JFUSION2119 FLUID |         |
| 1    | CJFUSION2120 | OUTER CVR JFUSION2120 FLUID |         |
| 1    | CJFUSION2121 | OUTER CVR JFUSION2121 FLUID |         |
| 1    | CJFUSION2122 | OUTER CVR JFUSION2122 FLUID |         |
| 1    | CJFUSION2123 | OUTER CVR JFUSION2123 FLUID |         |

| Pos. | Item Number   | Description                 | Remarks              |
|------|---------------|-----------------------------|----------------------|
| 1    | CJFUSION2124  | OUTER CVR JFUSION2124 FLUID |                      |
| 1    | CJFUSION2217  | OUTER CVR JFUSION2217 FLUID |                      |
| 1    | CJFUSION2218  | OUTER CVR JFUSION2218 FLUID |                      |
| 1    | CJFUSION2219  | OUTER CVR JFUSION2219 FLUID |                      |
| 1    | CJFUSION2220  | OUTER CVR JFUSION2220 FLUID |                      |
| 1    | CJFUSION2221  | OUTER CVR JFUSION2221 FLUID |                      |
| 1    | CJFUSION2222  | OUTER CVR JFUSION2222 FLUID |                      |
| 1    | CJFUSION2223  | OUTER CVR JFUSION2223 FLUID |                      |
| 1    | CJFUSION2224  | OUTER CVR JFUSION2224 FLUID |                      |
| 1    | CJFUSION2319  | OUTER CVR JFUSION2319 FLUID |                      |
| 1    | CJFUSION2320  | OUTER CVR JFUSION2320 FLUID |                      |
| 1    | CJFUSION2321  | OUTER CVR JFUSION2321 FLUID |                      |
| 1    | CJFUSION2322  | OUTER CVR JFUSION2322 FLUID |                      |
| 1    | CJFUSION2323  | OUTER CVR JFUSION2323 FLUID |                      |
| 1    | CJFUSION2324  | OUTER CVR JFUSION2324 FLUID |                      |
| 1    | CJFUSION2419  | OUTER CVR JFUSION2419 FLUID |                      |
| 1    | CJFUSION2420  | OUTER CVR JFUSION2420 FLUID |                      |
| 1    | CJFUSION2421  | OUTER CVR JFUSION2421 FLUID |                      |
| 1    | CJFUSION2422  | OUTER CVR JFUSION2422 FLUID |                      |
| 1    | CJFUSION2423  | OUTER CVR JFUSION2423 FLUID |                      |
| 1    | CJFUSION2424  | OUTER CVR JFUSION2424 FLUID |                      |
| 2    | FAJFUSION     | FLUID INSERT SIZE A JFUSION | Used in 14"&15" wide |
| 2    | FBJFUSION     | FLUID INSERT SIZE B JFUSION | Used in 16"-20" wide |
| 2    | FCJFUSION     | FLUID INSERT SIZE C JFUSION | Used in 21"-24" wide |
| 3    | ICJFUSION1414 | INNER CVR JFUSION1414 FLUID |                      |
| 3    | ICJFUSION1415 | INNER CVR JFUSION1415 FLUID |                      |
| 3    | ICJFUSION1416 | INNER CVR JFUSION1416 FLUID |                      |
| 3    | ICJFUSION1417 | INNER CVR JFUSION1417 FLUID |                      |
| 3    | ICJFUSION1418 | INNER CVR JFUSION1418 FLUID |                      |
| 3    | ICJFUSION1514 | INNER CVR JFUSION1514 FLUID |                      |
| 3    | ICJFUSION1515 | INNER CVR JFUSION1515 FLUID |                      |
| 3    | ICJFUSION1516 | INNER CVR JFUSION1516 FLUID |                      |
| 3    | ICJFUSION1517 | INNER CVR JFUSION1517 FLUID |                      |
| 3    | ICJFUSION1518 | INNER CVR JFUSION1518 FLUID |                      |
| 3    | ICJFUSION1519 | INNER CVR JFUSION1519 FLUID |                      |
| 3    | ICJFUSION1520 | INNER CVR JFUSION1520 FLUID |                      |
| 3    | ICJFUSION1521 | INNER CVR JFUSION1521 FLUID |                      |
| 3    | ICJFUSION1614 | INNER CVR JFUSION1614 FLUID |                      |
| 3    | ICJFUSION1615 | INNER CVR JFUSION1615 FLUID |                      |
| 3    | ICJFUSION1616 | INNER CVR JFUSION1616 FLUID |                      |
| 3    | ICJFUSION1617 | INNER CVR JFUSION1617 FLUID |                      |
| 3    | ICJFUSION1618 | INNER CVR JFUSION1618 FLUID |                      |
| 3    | ICJFUSION1619 | INNER CVR JFUSION1619 FLUID |                      |
| 3    | ICJFUSION1620 | INNER CVR JFUSION1620 FLUID |                      |
| 3    | ICJFUSION1621 | INNER CVR JFUSION1621 FLUID |                      |
| 3    | ICJFUSION1622 | INNER CVR JFUSION1622 FLUID |                      |
| 3    | ICJFUSION1714 | INNER CVR JFUSION1714 FLUID |                      |
| 3    | ICJFUSION1715 | INNER CVR JFUSION1715 FLUID |                      |
| 3    | ICJFUSION1716 | INNER CVR JFUSION1716 FLUID |                      |
| 3    | ICJFUSION1717 | INNER CVR JFUSION1717 FLUID |                      |
| 3    | ICJFUSION1718 | INNER CVR JFUSION1718 FLUID |                      |
| 3    | ICJFUSION1719 | INNER CVR JFUSION1719 FLUID |                      |
| 3    | ICJFUSION1720 | INNER CVR JFUSION1720 FLUID |                      |
| 3    | ICJFUSION1721 | INNER CVR JFUSION1721 FLUID |                      |
| 3    | ICJFUSION1722 | INNER CVR JFUSION1722 FLUID |                      |
| 3    | ICJFUSION1723 | INNER CVR JFUSION1723 FLUID |                      |
| 3    | ICJFUSION1814 | INNER CVR JFUSION1814 FLUID |                      |

| Pos. | Item Number   | Description                 | Remarks |
|------|---------------|-----------------------------|---------|
| 3    | ICJFUSION1815 | INNER CVR JFUSION1815 FLUID |         |
| 3    | ICJFUSION1816 | INNER CVR JFUSION1816 FLUID |         |
| 3    | ICJFUSION1817 | INNER CVR JFUSION1817 FLUID |         |
| 3    | ICJFUSION1818 | INNER CVR JFUSION1818 FLUID |         |
| 3    | ICJFUSION1819 | INNER CVR JFUSION1819 FLUID |         |
| 3    | ICJFUSION1820 | INNER CVR JFUSION1820 FLUID |         |
| 3    | ICJFUSION1821 | INNER CVR JFUSION1821 FLUID |         |
| 3    | ICJFUSION1822 | INNER CVR JFUSION1822 FLUID |         |
| 3    | ICJFUSION1823 | INNER CVR JFUSION1823 FLUID |         |
| 3    | ICJFUSION1824 | INNER CVR JFUSION1824 FLUID |         |
| 3    | ICJFUSION1915 | INNER CVR JFUSION1915 FLUID |         |
| 3    | ICJFUSION1916 | INNER CVR JFUSION1916 FLUID |         |
| 3    | ICJFUSION1917 | INNER CVR JFUSION1917 FLUID |         |
| 3    | ICJFUSION1918 | INNER CVR JFUSION1918 FLUID |         |
| 3    | ICJFUSION1919 | INNER CVR JFUSION1919 FLUID |         |
| 3    | ICJFUSION1920 | INNER CVR JFUSION1920 FLUID |         |
| 3    | ICJFUSION1921 | INNER CVR JFUSION1921 FLUID |         |
| 3    | ICJFUSION1922 | INNER CVR JFUSION1922 FLUID |         |
| 3    | ICJFUSION1923 | INNER CVR JFUSION1923 FLUID |         |
| 3    | ICJFUSION1924 | INNER CVR JFUSION1924 FLUID |         |
| 3    | ICJFUSION2015 | INNER CVR JFUSION2015 FLUID |         |
| 3    | ICJFUSION2016 | INNER CVR JFUSION2016 FLUID |         |
| 3    | ICJFUSION2017 | INNER CVR JFUSION2017 FLUID |         |
| 3    | ICJFUSION2018 | INNER CVR JFUSION2018 FLUID |         |
| 3    | ICJFUSION2019 | INNER CVR JFUSION2019 FLUID |         |
| 3    | ICJFUSION2020 | INNER CVR JFUSION2020 FLUID |         |
| 3    | ICJFUSION2021 | INNER CVR JFUSION2021 FLUID |         |
| 3    | ICJFUSION2022 | INNER CVR JFUSION2022 FLUID |         |
| 3    | ICJFUSION2023 | INNER CVR JFUSION2023 FLUID |         |
| 3    | ICJFUSION2024 | INNER CVR JFUSION2024 FLUID |         |
| 3    | ICJFUSION2117 | INNER CVR JFUSION2117 FLUID |         |
| 3    | ICJFUSION2118 | INNER CVR JFUSION2118 FLUID |         |
| 3    | ICJFUSION2119 | INNER CVR JFUSION2119 FLUID |         |
| 3    | ICJFUSION2120 | INNER CVR JFUSION2120 FLUID |         |
| 3    | ICJFUSION2121 | INNER CVR JFUSION2121 FLUID |         |
| 3    | ICJFUSION2122 | INNER CVR JFUSION2122 FLUID |         |
| 3    | ICJFUSION2123 | INNER CVR JFUSION2123 FLUID |         |
| 3    | ICJFUSION2124 | INNER CVR JFUSION2124 FLUID |         |
| 3    | ICJFUSION2217 | INNER CVR JFUSION2217 FLUID |         |
| 3    | ICJFUSION2218 | INNER CVR JFUSION2218 FLUID |         |
| 3    | ICJFUSION2219 | INNER CVR JFUSION2219 FLUID |         |
| 3    | ICJFUSION2220 | INNER CVR JFUSION2220 FLUID |         |
| 3    | ICJFUSION2221 | INNER CVR JFUSION2221 FLUID |         |
| 3    | ICJFUSION2222 | INNER CVR JFUSION2222 FLUID |         |
| 3    | ICJFUSION2223 | INNER CVR JFUSION2223 FLUID |         |
| 3    | ICJFUSION2224 | INNER CVR JFUSION2224 FLUID |         |
| 3    | ICJFUSION2319 | INNER CVR JFUSION2319 FLUID |         |
| 3    | ICJFUSION2320 | INNER CVR JFUSION2320 FLUID |         |
| 3    | ICJFUSION2321 | INNER CVR JFUSION2321 FLUID |         |
| 3    | ICJFUSION2322 | INNER CVR JFUSION2322 FLUID |         |
| 3    | ICJFUSION2323 | INNER CVR JFUSION2323 FLUID |         |
| 3    | ICJFUSION2324 | INNER CVR JFUSION2324 FLUID |         |
| 3    | ICJFUSION2419 | INNER CVR JFUSION2419 FLUID |         |
| 3    | ICJFUSION2420 | INNER CVR JFUSION2420 FLUID |         |
| 3    | ICJFUSION2421 | INNER CVR JFUSION2421 FLUID |         |
| 3    | ICJFUSION2422 | INNER CVR JFUSION2422 FLUID |         |
| 3    | ICJFUSION2423 | INNER CVR JFUSION2423 FLUID |         |

| Pos. | Item Number   | Description                 | Remarks |
|------|---------------|-----------------------------|---------|
| 3    | ICJFUSION2424 | INNER CVR JFUSION2424 FLUID |         |
| 4    | BJFUSION1414  | BASE JFUSION1414            |         |
| 4    | BJFUSION1415  | BASE JFUSION1415            |         |
| 4    | BJFUSION1416  | BASE JFUSION1416            |         |
| 4    | BJFUSION1417  | BASE JFUSION1417            |         |
| 4    | BJFUSION1418  | BASE JFUSION1418            |         |
| 4    | BJFUSION1514  | BASE JFUSION1514            |         |
| 4    | BJFUSION1515  | BASE JFUSION1515            |         |
| 4    | BJFUSION1516  | BASE JFUSION1516            |         |
| 4    | BJFUSION1517  | BASE JFUSION1517            |         |
| 4    | BJFUSION1518  | BASE JFUSION1518            |         |
| 4    | BJFUSION1519  | BASE JFUSION1519            |         |
| 4    | BJFUSION1520  | BASE JFUSION1520            |         |
| 4    | BJFUSION1521  | BASE JFUSION1521            |         |
| 4    | BJFUSION1614  | BASE JFUSION1614            |         |
| 4    | BJFUSION1615  | BASE JFUSION1615            |         |
| 4    | BJFUSION1616  | BASE JFUSION1616            |         |
| 4    | BJFUSION1617  | BASE JFUSION1617            |         |
| 4    | BJFUSION1618  | BASE JFUSION1618            |         |
| 4    | BJFUSION1619  | BASE JFUSION1619            |         |
| 4    | BJFUSION1620  | BASE JFUSION1620            |         |
| 4    | BJFUSION1621  | BASE JFUSION1621            |         |
| 4    | BJFUSION1622  | BASE JFUSION1622            |         |
| 4    | BJFUSION1714  | BASE JFUSION1714            |         |
| 4    | BJFUSION1715  | BASE JFUSION1715            |         |
| 4    | BJFUSION1716  | BASE JFUSION1716            |         |
| 4    | BJFUSION1717  | BASE JFUSION1717            |         |
| 4    | BJFUSION1718  | BASE JFUSION1718            |         |
| 4    | BJFUSION1719  | BASE JFUSION1719            |         |
| 4    | BJFUSION1720  | BASE JFUSION1720            |         |
| 4    | BJFUSION1721  | BASE JFUSION1721            |         |
| 4    | BJFUSION1722  | BASE JFUSION1722            |         |
| 4    | BJFUSION1723  | BASE JFUSION1723            |         |
| 4    | BJFUSION1814  | BASE JFUSION1814            |         |
| 4    | BJFUSION1815  | BASE JFUSION1815            |         |
| 4    | BJFUSION1816  | BASE JFUSION1816            |         |
| 4    | BJFUSION1817  | BASE JFUSION1817            |         |
| 4    | BJFUSION1818  | BASE JFUSION1818            |         |
| 4    | BJFUSION1819  | BASE JFUSION1819            |         |
| 4    | BJFUSION1820  | BASE JFUSION1820            |         |
| 4    | BJFUSION1821  | BASE JFUSION1821            |         |
| 4    | BJFUSION1822  | BASE JFUSION1822            |         |
| 4    | BJFUSION1823  | BASE JFUSION1823            |         |
| 4    | BJFUSION1824  | BASE JFUSION1824            |         |
| 4    | BJFUSION1915  | BASE JFUSION1915            |         |
| 4    | BJFUSION1916  | BASE JFUSION1916            |         |
| 4    | BJFUSION1917  | BASE JFUSION1917            |         |
| 4    | BJFUSION1918  | BASE JFUSION1918            |         |
| 4    | BJFUSION1919  | BASE JFUSION1919            |         |
| 4    | BJFUSION1920  | BASE JFUSION1920            |         |
| 4    | BJFUSION1921  | BASE JFUSION1921            |         |
| 4    | BJFUSION1922  | BASE JFUSION1922            |         |
| 4    | BJFUSION1923  | BASE JFUSION1923            |         |
| 4    | BJFUSION1924  | BASE JFUSION1924            |         |
| 4    | BJFUSION2015  | BASE JFUSION2015            |         |
| 4    | BJFUSION2016  | BASE JFUSION2016            |         |
| 4    | BJFUSION2017  | BASE JFUSION2017            |         |

| Pos. | Item Number  | Description      | Remarks |
|------|--------------|------------------|---------|
| 4    | BJFUSION2018 | BASE JFUSION2018 |         |
| 4    | BJFUSION2019 | BASE JFUSION2019 |         |
| 4    | BJFUSION2020 | BASE JFUSION2020 |         |
| 4    | BJFUSION2021 | BASE JFUSION2021 |         |
| 4    | BJFUSION2022 | BASE JFUSION2022 |         |
| 4    | BJFUSION2023 | BASE JFUSION2023 |         |
| 4    | BJFUSION2024 | BASE JFUSION2024 |         |
| 4    | BJFUSION2117 | BASE JFUSION2117 |         |
| 4    | BJFUSION2118 | BASE JFUSION2118 |         |
| 4    | BJFUSION2119 | BASE JFUSION2119 |         |
| 4    | BJFUSION2120 | BASE JFUSION2120 |         |
| 4    | BJFUSION2121 | BASE JFUSION2121 |         |
| 4    | BJFUSION2122 | BASE JFUSION2122 |         |
| 4    | BJFUSION2123 | BASE JFUSION2123 |         |
| 4    | BJFUSION2124 | BASE JFUSION2124 |         |
| 4    | BJFUSION2217 | BASE JFUSION2217 |         |
| 4    | BJFUSION2218 | BASE JFUSION2218 |         |
| 4    | BJFUSION2219 | BASE JFUSION2219 |         |
| 4    | BJFUSION2220 | BASE JFUSION2220 |         |
| 4    | BJFUSION2221 | BASE JFUSION2221 |         |
| 4    | BJFUSION2222 | BASE JFUSION2222 |         |
| 4    | BJFUSION2223 | BASE JFUSION2223 |         |
| 4    | BJFUSION2224 | BASE JFUSION2224 |         |
| 4    | BJFUSION2319 | BASE JFUSION2319 |         |
| 4    | BJFUSION2320 | BASE JFUSION2320 |         |
| 4    | BJFUSION2321 | BASE JFUSION2321 |         |
| 4    | BJFUSION2322 | BASE JFUSION2322 |         |
| 4    | BJFUSION2323 | BASE JFUSION2323 |         |
| 4    | BJFUSION2324 | BASE JFUSION2324 |         |
| 4    | BJFUSION2419 | BASE JFUSION2419 |         |
| 4    | BJFUSION2420 | BASE JFUSION2420 |         |
| 4    | BJFUSION2421 | BASE JFUSION2421 |         |
| 4    | BJFUSION2422 | BASE JFUSION2422 |         |
| 4    | BJFUSION2423 | BASE JFUSION2423 |         |
| 4    | BJFUSION2424 | BASE JFUSION2424 |         |

Note: Cushion widths 14" & 15" use A size Air Bladders, 16" - 20" use B size bladders and 21"-24" use C size bladders



| Pos. | Item Number   | Description               | Remarks |
|------|---------------|---------------------------|---------|
| 1    | CJFUSION1414A | OUTER CVR JFUSION1414 AIR |         |
| 1    | CJFUSION1415A | OUTER CVR JFUSION1415 AIR |         |
| 1    | CJFUSION1416A | OUTER CVR JFUSION1416 AIR |         |
| 1    | CJFUSION1417A | OUTER CVR JFUSION1417 AIR |         |
| 1    | CJFUSION1418A | OUTER CVR JFUSION1418 AIR |         |
| 1    | CJFUSION1514A | OUTER CVR JFUSION1514 AIR |         |
| 1    | CJFUSION1515A | OUTER CVR JFUSION1515 AIR |         |
| 1    | CJFUSION1516A | OUTER CVR JFUSION1516 AIR |         |
| 1    | CJFUSION1517A | OUTER CVR JFUSION1517 AIR |         |
| 1    | CJFUSION1518A | OUTER CVR JFUSION1518 AIR |         |
| 1    | CJFUSION1519A | OUTER CVR JFUSION1519 AIR |         |
| 1    | CJFUSION1520A | OUTER CVR JFUSION1520 AIR |         |
| 1    | CJFUSION1521A | OUTER CVR JFUSION1521 AIR |         |
| 1    | CJFUSION1614A | OUTER CVR JFUSION1614 AIR |         |
| 1    | CJFUSION1615A | OUTER CVR JFUSION1615 AIR |         |

| Pos. | Item Number   | Description               | Remarks |
|------|---------------|---------------------------|---------|
| 1    | CJFUSION1616A | OUTER CVR JFUSION1616 AIR |         |
| 1    | CJFUSION1617A | OUTER CVR JFUSION1617 AIR |         |
| 1    | CJFUSION1618A | OUTER CVR JFUSION1618 AIR |         |
| 1    | CJFUSION1619A | OUTER CVR JFUSION1619 AIR |         |
| 1    | CJFUSION1620A | OUTER CVR JFUSION1620 AIR |         |
| 1    | CJFUSION1621A | OUTER CVR JFUSION1621 AIR |         |
| 1    | CJFUSION1622A | OUTER CVR JFUSION1622 AIR |         |
| 1    | CJFUSION1714A | OUTER CVR JFUSION1714 AIR |         |
| 1    | CJFUSION1715A | OUTER CVR JFUSION1715 AIR |         |
| 1    | CJFUSION1716A | OUTER CVR JFUSION1716 AIR |         |
| 1    | CJFUSION1717A | OUTER CVR JFUSION1717 AIR |         |
| 1    | CJFUSION1718A | OUTER CVR JFUSION1718 AIR |         |
| 1    | CJFUSION1719A | OUTER CVR JFUSION1719 AIR |         |
| 1    | CJFUSION1720A | OUTER CVR JFUSION1720 AIR |         |
| 1    | CJFUSION1721A | OUTER CVR JFUSION1721 AIR |         |
| 1    | CJFUSION1722A | OUTER CVR JFUSION1722 AIR |         |
| 1    | CJFUSION1723A | OUTER CVR JFUSION1723 AIR |         |
| 1    | CJFUSION1814A | OUTER CVR JFUSION1814 AIR |         |
| 1    | CJFUSION1815A | OUTER CVR JFUSION1815 AIR |         |
| 1    | CJFUSION1816A | OUTER CVR JFUSION1816 AIR |         |
| 1    | CJFUSION1817A | OUTER CVR JFUSION1817 AIR |         |
| 1    | CJFUSION1818A | OUTER CVR JFUSION1818 AIR |         |
| 1    | CJFUSION1819A | OUTER CVR JFUSION1819 AIR |         |
| 1    | CJFUSION1820A | OUTER CVR JFUSION1820 AIR |         |
| 1    | CJFUSION1821A | OUTER CVR JFUSION1821 AIR |         |
| 1    | CJFUSION1822A | OUTER CVR JFUSION1822 AIR |         |
| 1    | CJFUSION1823A | OUTER CVR JFUSION1823 AIR |         |
| 1    | CJFUSION1824A | OUTER CVR JFUSION1824 AIR |         |
| 1    | CJFUSION1915A | OUTER CVR JFUSION1915 AIR |         |
| 1    | CJFUSION1916A | OUTER CVR JFUSION1916 AIR |         |
| 1    | CJFUSION1917A | OUTER CVR JFUSION1917 AIR |         |
| 1    | CJFUSION1918A | OUTER CVR JFUSION1918 AIR |         |
| 1    | CJFUSION1919A | OUTER CVR JFUSION1919 AIR |         |
| 1    | CJFUSION1920A | OUTER CVR JFUSION1920 AIR |         |
| 1    | CJFUSION1921A | OUTER CVR JFUSION1921 AIR |         |
| 1    | CJFUSION1922A | OUTER CVR JFUSION1922 AIR |         |
| 1    | CJFUSION1923A | OUTER CVR JFUSION1923 AIR |         |
| 1    | CJFUSION1924A | OUTER CVR JFUSION1924 AIR |         |
| 1    | CJFUSION2015A | OUTER CVR JFUSION2015 AIR |         |
| 1    | CJFUSION2016A | OUTER CVR JFUSION2016 AIR |         |
| 1    | CJFUSION2017A | OUTER CVR JFUSION2017 AIR |         |
| 1    | CJFUSION2018A | OUTER CVR JFUSION2018 AIR |         |
| 1    | CJFUSION2019A | OUTER CVR JFUSION2019 AIR |         |
| 1    | CJFUSION2020A | OUTER CVR JFUSION2020 AIR |         |
| 1    | CJFUSION2021A | OUTER CVR JFUSION2021 AIR |         |
| 1    | CJFUSION2022A | OUTER CVR JFUSION2022 AIR |         |
| 1    | CJFUSION2023A | OUTER CVR JFUSION2023 AIR |         |
| 1    | CJFUSION2024A | OUTER CVR JFUSION2024 AIR |         |
| 1    | CJFUSION2117A | OUTER CVR JFUSION2117 AIR |         |
| 1    | CJFUSION2118A | OUTER CVR JFUSION2118 AIR |         |
| 1    | CJFUSION2119A | OUTER CVR JFUSION2119 AIR |         |
| 1    | CJFUSION2120A | OUTER CVR JFUSION2120 AIR |         |
| 1    | CJFUSION2121A | OUTER CVR JFUSION2121 AIR |         |
| 1    | CJFUSION2122A | OUTER CVR JFUSION2122 AIR |         |
| 1    | CJFUSION2123A | OUTER CVR JFUSION2123 AIR |         |
| 1    | CJFUSION2124A | OUTER CVR JFUSION2124 AIR |         |
| 1    | CJFUSION2217A | OUTER CVR JFUSION2217 AIR |         |

| Pos. | Item Number    | Description                    | Remarks              |
|------|----------------|--------------------------------|----------------------|
| 1    | CJFUSION2218A  | OUTER CVR JFUSION2218 AIR      |                      |
| 1    | CJFUSION2219A  | OUTER CVR JFUSION2219 AIR      |                      |
| 1    | CJFUSION2220A  | OUTER CVR JFUSION2220 AIR      |                      |
| 1    | CJFUSION2221A  | OUTER CVR JFUSION2221 AIR      |                      |
| 1    | CJFUSION2222A  | OUTER CVR JFUSION2222 AIR      |                      |
| 1    | CJFUSION2223A  | OUTER CVR JFUSION2223 AIR      |                      |
| 1    | CJFUSION2224A  | OUTER CVR JFUSION2224 AIR      |                      |
| 1    | CJFUSION2319A  | OUTER CVR JFUSION2319 AIR      |                      |
| 1    | CJFUSION2320A  | OUTER CVR JFUSION2320 AIR      |                      |
| 1    | CJFUSION2321A  | OUTER CVR JFUSION2321 AIR      |                      |
| 1    | CJFUSION2322A  | OUTER CVR JFUSION2322 AIR      |                      |
| 1    | CJFUSION2323A  | OUTER CVR JFUSION2323 AIR      |                      |
| 1    | CJFUSION2324A  | OUTER CVR JFUSION2324 AIR      |                      |
| 1    | CJFUSION2419A  | OUTER CVR JFUSION2419 AIR      |                      |
| 1    | CJFUSION2420A  | OUTER CVR JFUSION2420 AIR      |                      |
| 1    | CJFUSION2421A  | OUTER CVR JFUSION2421 AIR      |                      |
| 1    | CJFUSION2422A  | OUTER CVR JFUSION2422 AIR      |                      |
| 1    | CJFUSION2423A  | OUTER CVR JFUSION2423 AIR      |                      |
| 1    | CJFUSION2424A  | OUTER CVR JFUSION2424 AIR      |                      |
| 2    | ASAJFUSION     | AIR INSERT SINGLE VALVE SIZE A | Used in 14"&15" wide |
| 2    | ASBJFUSION     | AIR INSERT SINGLE VALVE SIZE B | Used in 16"-20" wide |
| 2    | ASCJFUSION     | AIR INSERT SINGLE VALVE SIZE C | Used in 21"-24" wide |
| 2    | ADAJFUSION     | AIR INSERT DUAL VALVE SIZE A   | Used in 14"&15" wide |
| 2    | ADBJFUSION     | AIR INSERT DUAL VALVE SIZE B   | Used in 16"-20" wide |
| 2    | ADCJFUSION     | AIR INSERT DUAL VALVE SIZE C   | Used in 21"-24" wide |
| 3    | ICJFUSION1414A | INNER CVR JFUSION1414 AIR      |                      |
| 3    | ICJFUSION1415A | INNER CVR JFUSION1415 AIR      |                      |
| 3    | ICJFUSION1416A | INNER CVR JFUSION1416 AIR      |                      |
| 3    | ICJFUSION1417A | INNER CVR JFUSION1417 AIR      |                      |
| 3    | ICJFUSION1418A | INNER CVR JFUSION1418 AIR      |                      |
| 3    | ICJFUSION1514A | INNER CVR JFUSION1514 AIR      |                      |
| 3    | ICJFUSION1515A | INNER CVR JFUSION1515 AIR      |                      |
| 3    | ICJFUSION1516A | INNER CVR JFUSION1516 AIR      |                      |
| 3    | ICJFUSION1517A | INNER CVR JFUSION1517 AIR      |                      |
| 3    | ICJFUSION1518A | INNER CVR JFUSION1518 AIR      |                      |
| 3    | ICJFUSION1519A | INNER CVR JFUSION1519 AIR      |                      |
| 3    | ICJFUSION1520A | INNER CVR JFUSION1520 AIR      |                      |
| 3    | ICJFUSION1521A | INNER CVR JFUSION1521 AIR      |                      |
| 3    | ICJFUSION1614A | INNER CVR JFUSION1614 AIR      |                      |
| 3    | ICJFUSION1615A | INNER CVR JFUSION1615 AIR      |                      |
| 3    | ICJFUSION1616A | INNER CVR JFUSION1616 AIR      |                      |
| 3    | ICJFUSION1617A | INNER CVR JFUSION1617 AIR      |                      |
| 3    | ICJFUSION1618A | INNER CVR JFUSION1618 AIR      |                      |
| 3    | ICJFUSION1619A | INNER CVR JFUSION1619 AIR      |                      |
| 3    | ICJFUSION1620A | INNER CVR JFUSION1620 AIR      |                      |
| 3    | ICJFUSION1621A | INNER CVR JFUSION1621 AIR      |                      |
| 3    | ICJFUSION1622A | INNER CVR JFUSION1622 AIR      |                      |
| 3    | ICJFUSION1714A | INNER CVR JFUSION1714 AIR      |                      |
| 3    | ICJFUSION1715A | INNER CVR JFUSION1715 AIR      |                      |
| 3    | ICJFUSION1716A | INNER CVR JFUSION1716 AIR      |                      |
| 3    | ICJFUSION1717A | INNER CVR JFUSION1717 AIR      |                      |
| 3    | ICJFUSION1718A | INNER CVR JFUSION1718 AIR      |                      |
| 3    | ICJFUSION1719A | INNER CVR JFUSION1719 AIR      |                      |
| 3    | ICJFUSION1720A | INNER CVR JFUSION1720 AIR      |                      |
| 3    | ICJFUSION1721A | INNER CVR JFUSION1721 AIR      |                      |
| 3    | ICJFUSION1722A | INNER CVR JFUSION1722 AIR      |                      |
| 3    | ICJFUSION1723A | INNER CVR JFUSION1723 AIR      |                      |

| Pos. | Item Number    | Description               | Remarks |
|------|----------------|---------------------------|---------|
| 3    | ICJFUSION1814A | INNER CVR JFUSION1814 AIR |         |
| 3    | ICJFUSION1815A | INNER CVR JFUSION1815 AIR |         |
| 3    | ICJFUSION1816A | INNER CVR JFUSION1816 AIR |         |
| 3    | ICJFUSION1817A | INNER CVR JFUSION1817 AIR |         |
| 3    | ICJFUSION1818A | INNER CVR JFUSION1818 AIR |         |
| 3    | ICJFUSION1819A | INNER CVR JFUSION1819 AIR |         |
| 3    | ICJFUSION1820A | INNER CVR JFUSION1820 AIR |         |
| 3    | ICJFUSION1821A | INNER CVR JFUSION1821 AIR |         |
| 3    | ICJFUSION1822A | INNER CVR JFUSION1822 AIR |         |
| 3    | ICJFUSION1823A | INNER CVR JFUSION1823 AIR |         |
| 3    | ICJFUSION1824A | INNER CVR JFUSION1824 AIR |         |
| 3    | ICJFUSION1915A | INNER CVR JFUSION1915 AIR |         |
| 3    | ICJFUSION1916A | INNER CVR JFUSION1916 AIR |         |
| 3    | ICJFUSION1917A | INNER CVR JFUSION1917 AIR |         |
| 3    | ICJFUSION1918A | INNER CVR JFUSION1918 AIR |         |
| 3    | ICJFUSION1919A | INNER CVR JFUSION1919 AIR |         |
| 3    | ICJFUSION1920A | INNER CVR JFUSION1920 AIR |         |
| 3    | ICJFUSION1921A | INNER CVR JFUSION1921 AIR |         |
| 3    | ICJFUSION1922A | INNER CVR JFUSION1922 AIR |         |
| 3    | ICJFUSION1923A | INNER CVR JFUSION1923 AIR |         |
| 3    | ICJFUSION1924A | INNER CVR JFUSION1924 AIR |         |
| 3    | ICJFUSION2015A | INNER CVR JFUSION2015 AIR |         |
| 3    | ICJFUSION2016A | INNER CVR JFUSION2016 AIR |         |
| 3    | ICJFUSION2017A | INNER CVR JFUSION2017 AIR |         |
| 3    | ICJFUSION2018A | INNER CVR JFUSION2018 AIR |         |
| 3    | ICJFUSION2019A | INNER CVR JFUSION2019 AIR |         |
| 3    | ICJFUSION2020A | INNER CVR JFUSION2020 AIR |         |
| 3    | ICJFUSION2021A | INNER CVR JFUSION2021 AIR |         |
| 3    | ICJFUSION2022A | INNER CVR JFUSION2022 AIR |         |
| 3    | ICJFUSION2023A | INNER CVR JFUSION2023 AIR |         |
| 3    | ICJFUSION2024A | INNER CVR JFUSION2024 AIR |         |
| 3    | ICJFUSION2117A | INNER CVR JFUSION2117 AIR |         |
| 3    | ICJFUSION2118A | INNER CVR JFUSION2118 AIR |         |
| 3    | ICJFUSION2119A | INNER CVR JFUSION2119 AIR |         |
| 3    | ICJFUSION2120A | INNER CVR JFUSION2120 AIR |         |
| 3    | ICJFUSION2121A | INNER CVR JFUSION2121 AIR |         |
| 3    | ICJFUSION2122A | INNER CVR JFUSION2122 AIR |         |
| 3    | ICJFUSION2123A | INNER CVR JFUSION2123 AIR |         |
| 3    | ICJFUSION2124A | INNER CVR JFUSION2124 AIR |         |
| 3    | ICJFUSION2217A | INNER CVR JFUSION2217 AIR |         |
| 3    | ICJFUSION2218A | INNER CVR JFUSION2218 AIR |         |
| 3    | ICJFUSION2219A | INNER CVR JFUSION2219 AIR |         |
| 3    | ICJFUSION2220A | INNER CVR JFUSION2220 AIR |         |
| 3    | ICJFUSION2221A | INNER CVR JFUSION2221 AIR |         |
| 3    | ICJFUSION2222A | INNER CVR JFUSION2222 AIR |         |
| 3    | ICJFUSION2223A | INNER CVR JFUSION2223 AIR |         |
| 3    | ICJFUSION2224A | INNER CVR JFUSION2224 AIR |         |
| 3    | ICJFUSION2319A | INNER CVR JFUSION2319 AIR |         |
| 3    | ICJFUSION2320A | INNER CVR JFUSION2320 AIR |         |
| 3    | ICJFUSION2321A | INNER CVR JFUSION2321 AIR |         |
| 3    | ICJFUSION2322A | INNER CVR JFUSION2322 AIR |         |
| 3    | ICJFUSION2323A | INNER CVR JFUSION2323 AIR |         |
| 3    | ICJFUSION2324A | INNER CVR JFUSION2324 AIR |         |
| 3    | ICJFUSION2419A | INNER CVR JFUSION2419 AIR |         |
| 3    | ICJFUSION2420A | INNER CVR JFUSION2420 AIR |         |
| 3    | ICJFUSION2421A | INNER CVR JFUSION2421 AIR |         |
| 3    | ICJFUSION2422A | INNER CVR JFUSION2422 AIR |         |

| Pos. | Item Number    | Description               | Remarks |
|------|----------------|---------------------------|---------|
| 3    | ICJFUSION2423A | INNER CVR JFUSION2423 AIR |         |
| 3    | ICJFUSION2424A | INNER CVR JFUSION2424 AIR |         |
| 4    | BJFUSION1414   | BASE JFUSION1414          |         |
| 4    | BJFUSION1415   | BASE JFUSION1415          |         |
| 4    | BJFUSION1416   | BASE JFUSION1416          |         |
| 4    | BJFUSION1417   | BASE JFUSION1417          |         |
| 4    | BJFUSION1418   | BASE JFUSION1418          |         |
| 4    | BJFUSION1514   | BASE JFUSION1514          |         |
| 4    | BJFUSION1515   | BASE JFUSION1515          |         |
| 4    | BJFUSION1516   | BASE JFUSION1516          |         |
| 4    | BJFUSION1517   | BASE JFUSION1517          |         |
| 4    | BJFUSION1518   | BASE JFUSION1518          |         |
| 4    | BJFUSION1519   | BASE JFUSION1519          |         |
| 4    | BJFUSION1520   | BASE JFUSION1520          |         |
| 4    | BJFUSION1521   | BASE JFUSION1521          |         |
| 4    | BJFUSION1614   | BASE JFUSION1614          |         |
| 4    | BJFUSION1615   | BASE JFUSION1615          |         |
| 4    | BJFUSION1616   | BASE JFUSION1616          |         |
| 4    | BJFUSION1617   | BASE JFUSION1617          |         |
| 4    | BJFUSION1618   | BASE JFUSION1618          |         |
| 4    | BJFUSION1619   | BASE JFUSION1619          |         |
| 4    | BJFUSION1620   | BASE JFUSION1620          |         |
| 4    | BJFUSION1621   | BASE JFUSION1621          |         |
| 4    | BJFUSION1622   | BASE JFUSION1622          |         |
| 4    | BJFUSION1714   | BASE JFUSION1714          |         |
| 4    | BJFUSION1715   | BASE JFUSION1715          |         |
| 4    | BJFUSION1716   | BASE JFUSION1716          |         |
| 4    | BJFUSION1717   | BASE JFUSION1717          |         |
| 4    | BJFUSION1718   | BASE JFUSION1718          |         |
| 4    | BJFUSION1719   | BASE JFUSION1719          |         |
| 4    | BJFUSION1720   | BASE JFUSION1720          |         |
| 4    | BJFUSION1721   | BASE JFUSION1721          |         |
| 4    | BJFUSION1722   | BASE JFUSION1722          |         |
| 4    | BJFUSION1723   | BASE JFUSION1723          |         |
| 4    | BJFUSION1814   | BASE JFUSION1814          |         |
| 4    | BJFUSION1815   | BASE JFUSION1815          |         |
| 4    | BJFUSION1816   | BASE JFUSION1816          |         |
| 4    | BJFUSION1817   | BASE JFUSION1817          |         |
| 4    | BJFUSION1818   | BASE JFUSION1818          |         |
| 4    | BJFUSION1819   | BASE JFUSION1819          |         |
| 4    | BJFUSION1820   | BASE JFUSION1820          |         |
| 4    | BJFUSION1821   | BASE JFUSION1821          |         |
| 4    | BJFUSION1822   | BASE JFUSION1822          |         |
| 4    | BJFUSION1823   | BASE JFUSION1823          |         |
| 4    | BJFUSION1824   | BASE JFUSION1824          |         |
| 4    | BJFUSION1915   | BASE JFUSION1915          |         |
| 4    | BJFUSION1916   | BASE JFUSION1916          |         |
| 4    | BJFUSION1917   | BASE JFUSION1917          |         |
| 4    | BJFUSION1918   | BASE JFUSION1918          |         |
| 4    | BJFUSION1919   | BASE JFUSION1919          |         |
| 4    | BJFUSION1920   | BASE JFUSION1920          |         |
| 4    | BJFUSION1921   | BASE JFUSION1921          |         |
| 4    | BJFUSION1922   | BASE JFUSION1922          |         |
| 4    | BJFUSION1923   | BASE JFUSION1923          |         |
| 4    | BJFUSION1924   | BASE JFUSION1924          |         |
| 4    | BJFUSION2015   | BASE JFUSION2015          |         |
| 4    | BJFUSION2016   | BASE JFUSION2016          |         |

| Pos. | Item Number  | Description       | Remarks |
|------|--------------|-------------------|---------|
| 4    | BJFUSION2017 | BASE JFUSION2017  |         |
| 4    | BJFUSION2018 | BASE JFUSION2018  |         |
| 4    | BJFUSION2019 | BASE JFUSION2019  |         |
| 4    | BJFUSION2020 | BASE JFUSION2020  |         |
| 4    | BJFUSION2021 | BASE JFUSION2021  |         |
| 4    | BJFUSION2022 | BASE JFUSION2022  |         |
| 4    | BJFUSION2023 | BASE JFUSION2023  |         |
| 4    | BJFUSION2024 | BASE JFUSION2024  |         |
| 4    | BJFUSION2117 | BASE JFUSION2117  |         |
| 4    | BJFUSION2118 | BASE JFUSION2118  |         |
| 4    | BJFUSION2119 | BASE JFUSION2119  |         |
| 4    | BJFUSION2120 | BASE JFUSION2120  |         |
| 4    | BJFUSION2121 | BASE JFUSION2121  |         |
| 4    | BJFUSION2122 | BASE JFUSION2122  |         |
| 4    | BJFUSION2123 | BASE JFUSION2123  |         |
| 4    | BJFUSION2124 | BASE JFUSION2124  |         |
| 4    | BJFUSION2217 | BASE JFUSION2217  |         |
| 4    | BJFUSION2218 | BASE JFUSION2218  |         |
| 4    | BJFUSION2219 | BASE JFUSION2219  |         |
| 4    | BJFUSION2220 | BASE JFUSION2220  |         |
| 4    | BJFUSION2221 | BASE JFUSION2221  |         |
| 4    | BJFUSION2222 | BASE JFUSION2222  |         |
| 4    | BJFUSION2223 | BASE JFUSION2223  |         |
| 4    | BJFUSION2224 | BASE JFUSION2224  |         |
| 4    | BJFUSION2319 | BASE JFUSION2319  |         |
| 4    | BJFUSION2320 | BASE JFUSION2320  |         |
| 4    | BJFUSION2321 | BASE JFUSION2321  |         |
| 4    | BJFUSION2322 | BASE JFUSION2322  |         |
| 4    | BJFUSION2323 | BASE JFUSION2323  |         |
| 4    | BJFUSION2324 | BASE JFUSION2324  |         |
| 4    | BJFUSION2419 | BASE JFUSION2419  |         |
| 4    | BJFUSION2420 | BASE JFUSION2420  |         |
| 4    | BJFUSION2421 | BASE JFUSION2421  |         |
| 4    | BJFUSION2422 | BASE JFUSION2422  |         |
| 4    | BJFUSION2423 | BASE JFUSION2423  |         |
| 4    | BJFUSION2424 | BASE JFUSION2424  |         |
| 5    | AIRPUMP      | JAY AIR PUMP      |         |
| 6    | PATCHKIT     | JAY AIR PATCH KIT |         |

Note: Cushion widths 14" & 15" use A size Air Bladders, 16" - 20" use B size bladders and 21"-24" use C size bladders



| Pos. | Item Number     | Description                 | Remarks |
|------|-----------------|-----------------------------|---------|
| 1    | CJFUSION1414RPA | COV OUT FUSION RP AIR 14X14 |         |
| 1    | CJFUSION1415RPA | COV OUT FUSION RP AIR 14X15 |         |
| 1    | CJFUSION1416RPA | COV OUT FUSION RP AIR 14X16 |         |
| 1    | CJFUSION1417RPA | COV OUT FUSION RP AIR 14X17 |         |
| 1    | CJFUSION1418RPA | COV OUT FUSION RP AIR 14X18 |         |
| 1    | CJFUSION1514RPA | COV OUT FUSION RP AIR 15X14 |         |
| 1    | CJFUSION1515RPA | COV OUT FUSION RP AIR 15X15 |         |
| 1    | CJFUSION1516RPA | COV OUT FUSION RP AIR 15X16 |         |
| 1    | CJFUSION1517RPA | COV OUT FUSION RP AIR 15X17 |         |
| 1    | CJFUSION1518RPA | COV OUT FUSION RP AIR 15X18 |         |
| 1    | CJFUSION1519RPA | COV OUT FUSION RP AIR 15X19 |         |
| 1    | CJFUSION1520RPA | COV OUT FUSION RP AIR 15X19 |         |
| 1    | CJFUSION1521RPA | COV OUT FUSION RP AIR 15X21 |         |
| 1    | CJFUSION1614RPA | COV OUT FUSION RP AIR 16X14 |         |
| 1    | CJFUSION1615RPA | COV OUT FUSION RP AIR 16X15 |         |

| Pos. | Item Number     | Description                 | Remarks |
|------|-----------------|-----------------------------|---------|
| 1    | CJFUSION1616RPA | COV OUT FUSION RP AIR 16X16 |         |
| 1    | CJFUSION1617RPA | COV OUT FUSION RP AIR 16X17 |         |
| 1    | CJFUSION1618RPA | COV OUT FUSION RP AIR 16X18 |         |
| 1    | CJFUSION1619RPA | COV OUT FUSION RP AIR 16X19 |         |
| 1    | CJFUSION1620RPA | COV OUT FUSION RP AIR 16X20 |         |
| 1    | CJFUSION1621RPA | COV OUT FUSION RP AIR 16X21 |         |
| 1    | CJFUSION1622RPA | COV OUT FUSION RP AIR 16X22 |         |
| 1    | CJFUSION1714RPA | COV OUT FUSION RP AIR 17X14 |         |
| 1    | CJFUSION1715RPA | COV OUT FUSION RP AIR 17X15 |         |
| 1    | CJFUSION1716RPA | COV OUT FUSION RP AIR 17X16 |         |
| 1    | CJFUSION1717RPA | COV OUT FUSION RP AIR 17X17 |         |
| 1    | CJFUSION1718RPA | COV OUT FUSION RP AIR 17X18 |         |
| 1    | CJFUSION1719RPA | COV OUT FUSION RP AIR 17X19 |         |
| 1    | CJFUSION1720RPA | COV OUT FUSION RP AIR 17X20 |         |
| 1    | CJFUSION1721RPA | COV OUT FUSION RP AIR 17X21 |         |
| 1    | CJFUSION1722RPA | COV OUT FUSION RP AIR 17X22 |         |
| 1    | CJFUSION1723RPA | COV OUT FUSION RP AIR 17X23 |         |
| 1    | CJFUSION1814RPA | COV OUT FUSION RP AIR 18X14 |         |
| 1    | CJFUSION1815RPA | COV OUT FUSION RP AIR 18X15 |         |
| 1    | CJFUSION1816RPA | COV OUT FUSION RP AIR 18X16 |         |
| 1    | CJFUSION1817RPA | COV OUT FUSION RP AIR 18X17 |         |
| 1    | CJFUSION1818RPA | COV OUT FUSION RP AIR 18X18 |         |
| 1    | CJFUSION1819RPA | COV OUT FUSION RP AIR 18X19 |         |
| 1    | CJFUSION1820RPA | COV OUT FUSION RP AIR 18X20 |         |
| 1    | CJFUSION1821RPA | COV OUT FUSION RP AIR 18X21 |         |
| 1    | CJFUSION1822RPA | COV OUT FUSION RP AIR 18X22 |         |
| 1    | CJFUSION1823RPA | COV OUT FUSION RP AIR 18X23 |         |
| 1    | CJFUSION1824RPA | COV OUT FUSION RP AIR 18X24 |         |
| 1    | CJFUSION1915RPA | COV OUT FUSION RP AIR 19X15 |         |
| 1    | CJFUSION1916RPA | COV OUT FUSION RP AIR 19X16 |         |
| 1    | CJFUSION1917RPA | COV OUT FUSION RP AIR 19X17 |         |
| 1    | CJFUSION1918RPA | COV OUT FUSION RP AIR 19X18 |         |
| 1    | CJFUSION1919RPA | COV OUT FUSION RP AIR 19X19 |         |
| 1    | CJFUSION1920RPA | COV OUT FUSION RP AIR 19X20 |         |
| 1    | CJFUSION1921RPA | COV OUT FUSION RP AIR 19X21 |         |
| 1    | CJFUSION1922RPA | COV OUT FUSION RP AIR 19X22 |         |
| 1    | CJFUSION1923RPA | COV OUT FUSION RP AIR 19X23 |         |
| 1    | CJFUSION1924RPA | COV OUT FUSION RP AIR 19X24 |         |
| 1    | CJFUSION2015RPA | COV OUT FUSION RP AIR 20X15 |         |
| 1    | CJFUSION2016RPA | COV OUT FUSION RP AIR 20X16 |         |
| 1    | CJFUSION2017RPA | COV OUT FUSION RP AIR 20X17 |         |
| 1    | CJFUSION2018RPA | COV OUT FUSION RP AIR 20X18 |         |
| 1    | CJFUSION2019RPA | COV OUT FUSION RP AIR 20X19 |         |
| 1    | CJFUSION2020RPA | COV OUT FUSION RP AIR 20X20 |         |
| 1    | CJFUSION2021RPA | COV OUT FUSION RP AIR 20X21 |         |
| 1    | CJFUSION2022RPA | COV OUT FUSION RP AIR 20X22 |         |
| 1    | CJFUSION2023RPA | COV OUT FUSION RP AIR 20X23 |         |
| 1    | CJFUSION2024RPA | COV OUT FUSION RP AIR 20X24 |         |
| 1    | CJFUSION2117RPA | COV OUT FUSION RP AIR 21X17 |         |
| 1    | CJFUSION2118RPA | COV OUT FUSION RP AIR 21X18 |         |
| 1    | CJFUSION2119RPA | COV OUT FUSION RP AIR 21X19 |         |
| 1    | CJFUSION2120RPA | COV OUT FUSION RP AIR 21X20 |         |
| 1    | CJFUSION2121RPA | COV OUT FUSION RP AIR 21X21 |         |
| 1    | CJFUSION2122RPA | COV OUT FUSION RP AIR 21X22 |         |
| 1    | CJFUSION2123RPA | COV OUT FUSION RP AIR 21X23 |         |
| 1    | CJFUSION2124RPA | COV OUT FUSION RP AIR 21X24 |         |
| 1    | CJFUSION2217RPA | COV OUT FUSION RP AIR 22X17 |         |

| Pos. | Item Number      | Description                    | Remarks              |
|------|------------------|--------------------------------|----------------------|
| 1    | CJFUSION2218RPA  | COV OUT FUSION RP AIR 22X18    |                      |
| 1    | CJFUSION2219RPA  | COV OUT FUSION RP AIR 22X19    |                      |
| 1    | CJFUSION2220RPA  | COV OUT FUSION RP AIR 22X20    |                      |
| 1    | CJFUSION2221RPA  | COV OUT FUSION RP AIR 22X21    |                      |
| 1    | CJFUSION2222RPA  | COV OUT FUSION RP AIR 22X22    |                      |
| 1    | CJFUSION2223RPA  | COV OUT FUSION RP AIR 22X23    |                      |
| 1    | CJFUSION2224RPA  | COV OUT FUSION RP AIR 22X24    |                      |
| 1    | CJFUSION2319RPA  | COV OUT FUSION RP AIR 23X19    |                      |
| 1    | CJFUSION2320RPA  | COV OUT FUSION RP AIR 23X20    |                      |
| 1    | CJFUSION2321RPA  | COV OUT FUSION RP AIR 23X21    |                      |
| 1    | CJFUSION2322RPA  | COV OUT FUSION RP AIR 23X22    |                      |
| 1    | CJFUSION2323RPA  | COV OUT FUSION RP AIR 23X23    |                      |
| 1    | CJFUSION2324RPA  | COV OUT FUSION RP AIR 23X24    |                      |
| 1    | CJFUSION2419RPA  | COV OUT FUSION RP AIR 24X19    |                      |
| 1    | CJFUSION2420RPA  | COV OUT FUSION RP AIR 24X20    |                      |
| 1    | CJFUSION2421RPA  | COV OUT FUSION RP AIR 24X21    |                      |
| 1    | CJFUSION2422RPA  | COV OUT FUSION RP AIR 24X22    |                      |
| 1    | CJFUSION2423RPA  | COV OUT FUSION RP AIR 24X23    |                      |
| 1    | CJFUSION2424RPA  | COV OUT FUSION RP AIR 24X24    |                      |
| 2    | ASAJFUSION       | AIR INSERT SINGLE VALVE SIZE A | Used in 14"&15" wide |
| 2    | ASBJFUSION       | AIR INSERT SINGLE VALVE SIZE B | Used in 16"-20" wide |
| 2    | ASCJFUSION       | AIR INSERT SINGLE VALVE SIZE C | Used in 21"-24" wide |
| 2    | ADAJFUSION       | AIR INSERT DUAL VALVE SIZE A   | Used in 14"&15" wide |
| 2    | ADBJFUSION       | AIR INSERT DUAL VALVE SIZE B   | Used in 16"-20" wide |
| 2    | ADCJFUSION       | AIR INSERT DUAL VALVE SIZE C   | Used in 21"-24" wide |
| 3    | ICJFUSION1414RPA | INN COV FUSION RP AIR 14X14    |                      |
| 3    | ICJFUSION1415RPA | INN COV FUSION RP AIR 14X15    |                      |
| 3    | ICJFUSION1416RPA | INN COV FUSION RP AIR 14X16    |                      |
| 3    | ICJFUSION1417RPA | INN COV FUSION RP AIR 14X17    |                      |
| 3    | ICJFUSION1418RPA | INN COV FUSION RP AIR 14X18    |                      |
| 3    | ICJFUSION1514RPA | INN COV FUSION RP AIR 15X14    |                      |
| 3    | ICJFUSION1515RPA | INN COV FUSION RP AIR 15X15    |                      |
| 3    | ICJFUSION1516RPA | INN COV FUSION RP AIR 15X16    |                      |
| 3    | ICJFUSION1517RPA | INN COV FUSION RP AIR 15X17    |                      |
| 3    | ICJFUSION1518RPA | INN COV FUSION RP AIR 15X18    |                      |
| 3    | ICJFUSION1519RPA | INN COV FUSION RP AIR 15X19    |                      |
| 3    | ICJFUSION1520RPA | INN COV FUSION RP AIR 15X19    |                      |
| 3    | ICJFUSION1521RPA | INN COV FUSION RP AIR 15X21    |                      |
| 3    | ICJFUSION1614RPA | INN COV FUSION RP AIR 16X14    |                      |
| 3    | ICJFUSION1615RPA | INN COV FUSION RP AIR 16X15    |                      |
| 3    | ICJFUSION1616RPA | INN COV FUSION RP AIR 16X16    |                      |
| 3    | ICJFUSION1617RPA | INN COV FUSION RP AIR 16X17    |                      |
| 3    | ICJFUSION1618RPA | INN COV FUSION RP AIR 16X18    |                      |
| 3    | ICJFUSION1619RPA | INN COV FUSION RP AIR 16X19    |                      |
| 3    | ICJFUSION1620RPA | INN COV FUSION RP AIR 16X20    |                      |
| 3    | ICJFUSION1621RPA | INN COV FUSION RP AIR 16X21    |                      |
| 3    | ICJFUSION1622RPA | INN COV FUSION RP AIR 16X22    |                      |
| 3    | ICJFUSION1714RPA | INN COV FUSION RP AIR 17X14    |                      |
| 3    | ICJFUSION1715RPA | INN COV FUSION RP AIR 17X15    |                      |
| 3    | ICJFUSION1716RPA | INN COV FUSION RP AIR 17X16    |                      |
| 3    | ICJFUSION1717RPA | INN COV FUSION RP AIR 17X17    |                      |
| 3    | ICJFUSION1718RPA | INN COV FUSION RP AIR 17X18    |                      |
| 3    | ICJFUSION1719RPA | INN COV FUSION RP AIR 17X19    |                      |
| 3    | ICJFUSION1720RPA | INN COV FUSION RP AIR 17X20    |                      |
| 3    | ICJFUSION1721RPA | INN COV FUSION RP AIR 17X21    |                      |
| 3    | ICJFUSION1722RPA | INN COV FUSION RP AIR 17X22    |                      |
| 3    | ICJFUSION1723RPA | INN COV FUSION RP AIR 17X23    |                      |

| Pos. | Item Number      | Description                 | Remarks |
|------|------------------|-----------------------------|---------|
| 3    | ICJFUSION1814RPA | INN COV FUSION RP AIR 18X14 |         |
| 3    | ICJFUSION1815RPA | INN COV FUSION RP AIR 18X15 |         |
| 3    | ICJFUSION1816RPA | INN COV FUSION RP AIR 18X16 |         |
| 3    | ICJFUSION1817RPA | INN COV FUSION RP AIR 18X17 |         |
| 3    | ICJFUSION1818RPA | INN COV FUSION RP AIR 18X18 |         |
| 3    | ICJFUSION1819RPA | INN COV FUSION RP AIR 18X19 |         |
| 3    | ICJFUSION1820RPA | INN COV FUSION RP AIR 18X20 |         |
| 3    | ICJFUSION1821RPA | INN COV FUSION RP AIR 18X21 |         |
| 3    | ICJFUSION1822RPA | INN COV FUSION RP AIR 18X22 |         |
| 3    | ICJFUSION1823RPA | INN COV FUSION RP AIR 18X23 |         |
| 3    | ICJFUSION1824RPA | INN COV FUSION RP AIR 18X24 |         |
| 3    | ICJFUSION1915RPA | INN COV FUSION RP AIR 19X15 |         |
| 3    | ICJFUSION1916RPA | INN COV FUSION RP AIR 19X16 |         |
| 3    | ICJFUSION1917RPA | INN COV FUSION RP AIR 19X17 |         |
| 3    | ICJFUSION1918RPA | INN COV FUSION RP AIR 19X18 |         |
| 3    | ICJFUSION1919RPA | INN COV FUSION RP AIR 19X19 |         |
| 3    | ICJFUSION1920RPA | INN COV FUSION RP AIR 19X20 |         |
| 3    | ICJFUSION1921RPA | INN COV FUSION RP AIR 19X21 |         |
| 3    | ICJFUSION1922RPA | INN COV FUSION RP AIR 19X22 |         |
| 3    | ICJFUSION1923RPA | INN COV FUSION RP AIR 19X23 |         |
| 3    | ICJFUSION1924RPA | INN COV FUSION RP AIR 19X24 |         |
| 3    | ICJFUSION2015RPA | INN COV FUSION RP AIR 20X15 |         |
| 3    | ICJFUSION2016RPA | INN COV FUSION RP AIR 20X16 |         |
| 3    | ICJFUSION2017RPA | INN COV FUSION RP AIR 20X17 |         |
| 3    | ICJFUSION2018RPA | INN COV FUSION RP AIR 20X18 |         |
| 3    | ICJFUSION2019RPA | INN COV FUSION RP AIR 20X19 |         |
| 3    | ICJFUSION2020RPA | INN COV FUSION RP AIR 20X20 |         |
| 3    | ICJFUSION2021RPA | INN COV FUSION RP AIR 20X21 |         |
| 3    | ICJFUSION2022RPA | INN COV FUSION RP AIR 20X22 |         |
| 3    | ICJFUSION2023RPA | INN COV FUSION RP AIR 20X23 |         |
| 3    | ICJFUSION2024RPA | INN COV FUSION RP AIR 20X24 |         |
| 3    | ICJFUSION2117RPA | INN COV FUSION RP AIR 21X17 |         |
| 3    | ICJFUSION2118RPA | INN COV FUSION RP AIR 21X18 |         |
| 3    | ICJFUSION2119RPA | INN COV FUSION RP AIR 21X19 |         |
| 3    | ICJFUSION2120RPA | INN COV FUSION RP AIR 21X20 |         |
| 3    | ICJFUSION2121RPA | INN COV FUSION RP AIR 21X21 |         |
| 3    | ICJFUSION2122RPA | INN COV FUSION RP AIR 21X22 |         |
| 3    | ICJFUSION2123RPA | INN COV FUSION RP AIR 21X23 |         |
| 3    | ICJFUSION2124RPA | INN COV FUSION RP AIR 21X24 |         |
| 3    | ICJFUSION2217RPA | INN COV FUSION RP AIR 22X17 |         |
| 3    | ICJFUSION2218RPA | INN COV FUSION RP AIR 22X18 |         |
| 3    | ICJFUSION2219RPA | INN COV FUSION RP AIR 22X19 |         |
| 3    | ICJFUSION2220RPA | INN COV FUSION RP AIR 22X20 |         |
| 3    | ICJFUSION2221RPA | INN COV FUSION RP AIR 22X21 |         |
| 3    | ICJFUSION2222RPA | INN COV FUSION RP AIR 22X22 |         |
| 3    | ICJFUSION2223RPA | INN COV FUSION RP AIR 22X23 |         |
| 3    | ICJFUSION2224RPA | INN COV FUSION RP AIR 22X24 |         |
| 3    | ICJFUSION2319RPA | INN COV FUSION RP AIR 23X19 |         |
| 3    | ICJFUSION2320RPA | INN COV FUSION RP AIR 23X20 |         |
| 3    | ICJFUSION2321RPA | INN COV FUSION RP AIR 23X21 |         |
| 3    | ICJFUSION2322RPA | INN COV FUSION RP AIR 23X22 |         |
| 3    | ICJFUSION2323RPA | INN COV FUSION RP AIR 23X23 |         |
| 3    | ICJFUSION2324RPA | INN COV FUSION RP AIR 23X24 |         |
| 3    | ICJFUSION2419RPA | INN COV FUSION RP AIR 24X19 |         |
| 3    | ICJFUSION2420RPA | INN COV FUSION RP AIR 24X20 |         |
| 3    | ICJFUSION2421RPA | INN COV FUSION RP AIR 24X21 |         |
| 3    | ICJFUSION2422RPA | INN COV FUSION RP AIR 24X22 |         |

| Pos. | Item Number      | Description                 | Remarks |
|------|------------------|-----------------------------|---------|
| 3    | ICJFUSION2423RPA | INN COV FUSION RP AIR 24X23 |         |
| 3    | ICJFUSION2424RPA | INN COV FUSION RP AIR 24X24 |         |
| 4    | BJFUSION1414RP   | BASE JAY FUSION RPROF 14X14 |         |
| 4    | BJFUSION1415RP   | BASE JAY FUSION RPROF 14X15 |         |
| 4    | BJFUSION1416RP   | BASE JAY FUSION RPROF 14X16 |         |
| 4    | BJFUSION1417RP   | BASE JAY FUSION RPROF 14X17 |         |
| 4    | BJFUSION1418RP   | BASE JAY FUSION RPROF 14X18 |         |
| 4    | BJFUSION1514RP   | BASE JAY FUSION RPROF 15X14 |         |
| 4    | BJFUSION1515RP   | BASE JAY FUSION RPROF 15X15 |         |
| 4    | BJFUSION1516RP   | BASE JAY FUSION RPROF 15X16 |         |
| 4    | BJFUSION1517RP   | BASE JAY FUSION RPROF 15X17 |         |
| 4    | BJFUSION1518RP   | BASE JAY FUSION RPROF 15X18 |         |
| 4    | BJFUSION1519RP   | BASE JAY FUSION RPROF 15X19 |         |
| 4    | BJFUSION1520RP   | BASE JAY FUSION RPROF 15X19 |         |
| 4    | BJFUSION1521RP   | BASE JAY FUSION RPROF 15X21 |         |
| 4    | BJFUSION1614RP   | BASE JAY FUSION RPROF 16X14 |         |
| 4    | BJFUSION1615RP   | BASE JAY FUSION RPROF 16X15 |         |
| 4    | BJFUSION1616RP   | BASE JAY FUSION RPROF 16X16 |         |
| 4    | BJFUSION1617RP   | BASE JAY FUSION RPROF 16X17 |         |
| 4    | BJFUSION1618RP   | BASE JAY FUSION RPROF 16X18 |         |
| 4    | BJFUSION1619RP   | BASE JAY FUSION RPROF 16X19 |         |
| 4    | BJFUSION1620RP   | BASE JAY FUSION RPROF 16X20 |         |
| 4    | BJFUSION1621RP   | BASE JAY FUSION RPROF 16X21 |         |
| 4    | BJFUSION1622RP   | BASE JAY FUSION RPROF 16X22 |         |
| 4    | BJFUSION1714RP   | BASE JAY FUSION RPROF 17X14 |         |
| 4    | BJFUSION1715RP   | BASE JAY FUSION RPROF 17X15 |         |
| 4    | BJFUSION1716RP   | BASE JAY FUSION RPROF 17X16 |         |
| 4    | BJFUSION1717RP   | BASE JAY FUSION RPROF 17X17 |         |
| 4    | BJFUSION1718RP   | BASE JAY FUSION RPROF 17X18 |         |
| 4    | BJFUSION1719RP   | BASE JAY FUSION RPROF 17X19 |         |
| 4    | BJFUSION1720RP   | BASE JAY FUSION RPROF 17X20 |         |
| 4    | BJFUSION1721RP   | BASE JAY FUSION RPROF 17X21 |         |
| 4    | BJFUSION1722RP   | BASE JAY FUSION RPROF 17X22 |         |
| 4    | BJFUSION1723RP   | BASE JAY FUSION RPROF 17X23 |         |
| 4    | BJFUSION1814RP   | BASE JAY FUSION RPROF 18X14 |         |
| 4    | BJFUSION1815RP   | BASE JAY FUSION RPROF 18X15 |         |
| 4    | BJFUSION1816RP   | BASE JAY FUSION RPROF 18X16 |         |
| 4    | BJFUSION1817RP   | BASE JAY FUSION RPROF 18X17 |         |
| 4    | BJFUSION1818RP   | BASE JAY FUSION RPROF 18X18 |         |
| 4    | BJFUSION1819RP   | BASE JAY FUSION RPROF 18X19 |         |
| 4    | BJFUSION1820RP   | BASE JAY FUSION RPROF 18X20 |         |
| 4    | BJFUSION1821RP   | BASE JAY FUSION RPROF 18X21 |         |
| 4    | BJFUSION1822RP   | BASE JAY FUSION RPROF 18X22 |         |
| 4    | BJFUSION1823RP   | BASE JAY FUSION RPROF 18X23 |         |
| 4    | BJFUSION1824RP   | BASE JAY FUSION RPROF 18X24 |         |
| 4    | BJFUSION1915RP   | BASE JAY FUSION RPROF 19X15 |         |
| 4    | BJFUSION1916RP   | BASE JAY FUSION RPROF 19X16 |         |
| 4    | BJFUSION1917RP   | BASE JAY FUSION RPROF 19X17 |         |
| 4    | BJFUSION1918RP   | BASE JAY FUSION RPROF 19X18 |         |
| 4    | BJFUSION1919RP   | BASE JAY FUSION RPROF 19X19 |         |
| 4    | BJFUSION1920RP   | BASE JAY FUSION RPROF 19X20 |         |
| 4    | BJFUSION1921RP   | BASE JAY FUSION RPROF 19X21 |         |
| 4    | BJFUSION1922RP   | BASE JAY FUSION RPROF 19X22 |         |
| 4    | BJFUSION1923RP   | BASE JAY FUSION RPROF 19X23 |         |
| 4    | BJFUSION1924RP   | BASE JAY FUSION RPROF 19X24 |         |
| 4    | BJFUSION2015RP   | BASE JAY FUSION RPROF 20X15 |         |
| 4    | BJFUSION2016RP   | BASE JAY FUSION RPROF 20X16 |         |

| Pos. | Item Number    | Description                 | Remarks |
|------|----------------|-----------------------------|---------|
| 4    | BJFUSION2017RP | BASE JAY FUSION RPROF 20X17 |         |
| 4    | BJFUSION2018RP | BASE JAY FUSION RPROF 20X18 |         |
| 4    | BJFUSION2019RP | BASE JAY FUSION RPROF 20X19 |         |
| 4    | BJFUSION2020RP | BASE JAY FUSION RPROF 20X20 |         |
| 4    | BJFUSION2021RP | BASE JAY FUSION RPROF 20X21 |         |
| 4    | BJFUSION2022RP | BASE JAY FUSION RPROF 20X22 |         |
| 4    | BJFUSION2023RP | BASE JAY FUSION RPROF 20X23 |         |
| 4    | BJFUSION2024RP | BASE JAY FUSION RPROF 20X24 |         |
| 4    | BJFUSION2117RP | BASE JAY FUSION RPROF 21X17 |         |
| 4    | BJFUSION2118RP | BASE JAY FUSION RPROF 21X18 |         |
| 4    | BJFUSION2119RP | BASE JAY FUSION RPROF 21X19 |         |
| 4    | BJFUSION2120RP | BASE JAY FUSION RPROF 21X20 |         |
| 4    | BJFUSION2121RP | BASE JAY FUSION RPROF 21X21 |         |
| 4    | BJFUSION2122RP | BASE JAY FUSION RPROF 21X22 |         |
| 4    | BJFUSION2123RP | BASE JAY FUSION RPROF 21X23 |         |
| 4    | BJFUSION2124RP | BASE JAY FUSION RPROF 21X24 |         |
| 4    | BJFUSION2217RP | BASE JAY FUSION RPROF 22X17 |         |
| 4    | BJFUSION2218RP | BASE JAY FUSION RPROF 22X18 |         |
| 4    | BJFUSION2219RP | BASE JAY FUSION RPROF 22X19 |         |
| 4    | BJFUSION2220RP | BASE JAY FUSION RPROF 22X20 |         |
| 4    | BJFUSION2221RP | BASE JAY FUSION RPROF 22X21 |         |
| 4    | BJFUSION2222RP | BASE JAY FUSION RPROF 22X22 |         |
| 4    | BJFUSION2223RP | BASE JAY FUSION RPROF 22X23 |         |
| 4    | BJFUSION2224RP | BASE JAY FUSION RPROF 22X24 |         |
| 4    | BJFUSION2319RP | BASE JAY FUSION RPROF 23X19 |         |
| 4    | BJFUSION2320RP | BASE JAY FUSION RPROF 23X20 |         |
| 4    | BJFUSION2321RP | BASE JAY FUSION RPROF 23X21 |         |
| 4    | BJFUSION2322RP | BASE JAY FUSION RPROF 23X22 |         |
| 4    | BJFUSION2323RP | BASE JAY FUSION RPROF 23X23 |         |
| 4    | BJFUSION2324RP | BASE JAY FUSION RPROF 23X24 |         |
| 4    | BJFUSION2419RP | BASE JAY FUSION RPROF 24X19 |         |
| 4    | BJFUSION2420RP | BASE JAY FUSION RPROF 24X20 |         |
| 4    | BJFUSION2421RP | BASE JAY FUSION RPROF 24X21 |         |
| 4    | BJFUSION2422RP | BASE JAY FUSION RPROF 24X22 |         |
| 4    | BJFUSION2423RP | BASE JAY FUSION RPROF 24X23 |         |
| 4    | BJFUSION2424   | BASE JFUSION2424            |         |
| 5    | AIRPUMP        | JAY AIR PUMP                |         |
| 6    | PATCHKIT       | JAY AIR PATCH KIT           |         |

Note: Cushion widths 14" & 15" use A size Fluid Bladders, 16" - 20" use B size bladders and 21"-24" use C size bladders.



| Pos. | Item Number    | Description                 | Remarks |
|------|----------------|-----------------------------|---------|
| 1    | CJFUSION1414RP | COV OUT FUSION RP FLD 14X14 |         |
| 1    | CJFUSION1415RP | COV OUT FUSION RP FLD 14X15 |         |
| 1    | CJFUSION1416RP | COV OUT FUSION RP FLD 14X16 |         |
| 1    | CJFUSION1417RP | COV OUT FUSION RP FLD 14X17 |         |
| 1    | CJFUSION1418RP | COV OUT FUSION RP FLD 14X18 |         |
| 1    | CJFUSION1514RP | COV OUT FUSION RP FLD 15X14 |         |
| 1    | CJFUSION1515RP | COV OUT FUSION RP FLD 15X15 |         |
| 1    | CJFUSION1516RP | COV OUT FUSION RP FLD 15X16 |         |
| 1    | CJFUSION1517RP | COV OUT FUSION RP FLD 15X17 |         |
| 1    | CJFUSION1518RP | COV OUT FUSION RP FLD 15X18 |         |
| 1    | CJFUSION1519RP | COV OUT FUSION RP FLD 15X19 |         |
| 1    | CJFUSION1520RP | COV OUT FUSION RP FLD 15X19 |         |
| 1    | CJFUSION1521RP | COV OUT FUSION RP FLD 15X21 |         |

| Pos. | Item Number    | Description                 | Remarks |
|------|----------------|-----------------------------|---------|
| 1    | CJFUSION1614RP | COV OUT FUSION RP FLD 16X14 |         |
| 1    | CJFUSION1615RP | COV OUT FUSION RP FLD 16X15 |         |
| 1    | CJFUSION1616RP | COV OUT FUSION RP FLD 16X16 |         |
| 1    | CJFUSION1617RP | COV OUT FUSION RP FLD 16X17 |         |
| 1    | CJFUSION1618RP | COV OUT FUSION RP FLD 16X18 |         |
| 1    | CJFUSION1619RP | COV OUT FUSION RP FLD 16X19 |         |
| 1    | CJFUSION1620RP | COV OUT FUSION RP FLD 16X20 |         |
| 1    | CJFUSION1621RP | COV OUT FUSION RP FLD 16X21 |         |
| 1    | CJFUSION1622RP | COV OUT FUSION RP FLD 16X22 |         |
| 1    | CJFUSION1714RP | COV OUT FUSION RP FLD 17X14 |         |
| 1    | CJFUSION1715RP | COV OUT FUSION RP FLD 17X15 |         |
| 1    | CJFUSION1716RP | COV OUT FUSION RP FLD 17X16 |         |
| 1    | CJFUSION1717RP | COV OUT FUSION RP FLD 17X17 |         |
| 1    | CJFUSION1718RP | COV OUT FUSION RP FLD 17X18 |         |
| 1    | CJFUSION1719RP | COV OUT FUSION RP FLD 17X19 |         |
| 1    | CJFUSION1720RP | COV OUT FUSION RP FLD 17X20 |         |
| 1    | CJFUSION1721RP | COV OUT FUSION RP FLD 17X21 |         |
| 1    | CJFUSION1722RP | COV OUT FUSION RP FLD 17X22 |         |
| 1    | CJFUSION1723RP | COV OUT FUSION RP FLD 17X23 |         |
| 1    | CJFUSION1814RP | COV OUT FUSION RP FLD 18X14 |         |
| 1    | CJFUSION1815RP | COV OUT FUSION RP FLD 18X15 |         |
| 1    | CJFUSION1816RP | COV OUT FUSION RP FLD 18X16 |         |
| 1    | CJFUSION1817RP | COV OUT FUSION RP FLD 18X17 |         |
| 1    | CJFUSION1818RP | COV OUT FUSION RP FLD 18X18 |         |
| 1    | CJFUSION1819RP | COV OUT FUSION RP FLD 18X19 |         |
| 1    | CJFUSION1820RP | COV OUT FUSION RP FLD 18X20 |         |
| 1    | CJFUSION1821RP | COV OUT FUSION RP FLD 18X21 |         |
| 1    | CJFUSION1822RP | COV OUT FUSION RP FLD 18X22 |         |
| 1    | CJFUSION1823RP | COV OUT FUSION RP FLD 18X23 |         |
| 1    | CJFUSION1824RP | COV OUT FUSION RP FLD 18X24 |         |
| 1    | CJFUSION1915RP | COV OUT FUSION RP FLD 19X15 |         |
| 1    | CJFUSION1916RP | COV OUT FUSION RP FLD 19X16 |         |
| 1    | CJFUSION1917RP | COV OUT FUSION RP FLD 19X17 |         |
| 1    | CJFUSION1918RP | COV OUT FUSION RP FLD 19X18 |         |
| 1    | CJFUSION1919RP | COV OUT FUSION RP FLD 19X19 |         |
| 1    | CJFUSION1920RP | COV OUT FUSION RP FLD 19X20 |         |
| 1    | CJFUSION1921RP | COV OUT FUSION RP FLD 19X21 |         |
| 1    | CJFUSION1922RP | COV OUT FUSION RP FLD 19X22 |         |
| 1    | CJFUSION1923RP | COV OUT FUSION RP FLD 19X23 |         |
| 1    | CJFUSION1924RP | COV OUT FUSION RP FLD 19X24 |         |
| 1    | CJFUSION2015RP | COV OUT FUSION RP FLD 20X15 |         |
| 1    | CJFUSION2016RP | COV OUT FUSION RP FLD 20X16 |         |
| 1    | CJFUSION2017RP | COV OUT FUSION RP FLD 20X17 |         |
| 1    | CJFUSION2018RP | COV OUT FUSION RP FLD 20X18 |         |
| 1    | CJFUSION2019RP | COV OUT FUSION RP FLD 20X19 |         |
| 1    | CJFUSION2020RP | COV OUT FUSION RP FLD 20X20 |         |
| 1    | CJFUSION2021RP | COV OUT FUSION RP FLD 20X21 |         |
| 1    | CJFUSION2022RP | COV OUT FUSION RP FLD 20X22 |         |
| 1    | CJFUSION2023RP | COV OUT FUSION RP FLD 20X23 |         |
| 1    | CJFUSION2024RP | COV OUT FUSION RP FLD 20X24 |         |
| 1    | CJFUSION2117RP | COV OUT FUSION RP FLD 21X17 |         |
| 1    | CJFUSION2118RP | COV OUT FUSION RP FLD 21X18 |         |
| 1    | CJFUSION2119RP | COV OUT FUSION RP FLD 21X19 |         |
| 1    | CJFUSION2120RP | COV OUT FUSION RP FLD 21X20 |         |
| 1    | CJFUSION2121RP | COV OUT FUSION RP FLD 21X21 |         |
| 1    | CJFUSION2122RP | COV OUT FUSION RP FLD 21X22 |         |
| 1    | CJFUSION2123RP | COV OUT FUSION RP FLD 21X23 |         |

| Pos. | Item Number     | Description                 | Remarks              |
|------|-----------------|-----------------------------|----------------------|
| 1    | CJFUSION2124RP  | COV OUT FUSION RP FLD 21X24 |                      |
| 1    | CJFUSION2217RP  | COV OUT FUSION RP FLD 22X17 |                      |
| 1    | CJFUSION2218RP  | COV OUT FUSION RP FLD 22X18 |                      |
| 1    | CJFUSION2219RP  | COV OUT FUSION RP FLD 22X19 |                      |
| 1    | CJFUSION2220RP  | COV OUT FUSION RP FLD 22X20 |                      |
| 1    | CJFUSION2221RP  | COV OUT FUSION RP FLD 22X21 |                      |
| 1    | CJFUSION2222RP  | COV OUT FUSION RP FLD 22X22 |                      |
| 1    | CJFUSION2223RP  | COV OUT FUSION RP FLD 22X23 |                      |
| 1    | CJFUSION2224RP  | COV OUT FUSION RP FLD 22X24 |                      |
| 1    | CJFUSION2319RP  | COV OUT FUSION RP FLD 23X19 |                      |
| 1    | CJFUSION2320RP  | COV OUT FUSION RP FLD 23X20 |                      |
| 1    | CJFUSION2321RP  | COV OUT FUSION RP FLD 23X21 |                      |
| 1    | CJFUSION2322RP  | COV OUT FUSION RP FLD 23X22 |                      |
| 1    | CJFUSION2323RP  | COV OUT FUSION RP FLD 23X23 |                      |
| 1    | CJFUSION2324RP  | COV OUT FUSION RP FLD 23X24 |                      |
| 1    | CJFUSION2419RP  | COV OUT FUSION RP FLD 24X19 |                      |
| 1    | CJFUSION2420RP  | COV OUT FUSION RP FLD 24X20 |                      |
| 1    | CJFUSION2421RP  | COV OUT FUSION RP FLD 24X21 |                      |
| 1    | CJFUSION2422RP  | COV OUT FUSION RP FLD 24X22 |                      |
| 1    | CJFUSION2423RP  | COV OUT FUSION RP FLD 24X23 |                      |
| 1    | CJFUSION2424RP  | COV OUT FUSION RP FLD 24X24 |                      |
| 2    | FAJFUSION       | FLUID INSERT SIZE A JFUSION | Used in 14"&15" wide |
| 2    | FBJFUSION       | FLUID INSERT SIZE B JFUSION | Used in 16"-20" wide |
| 2    | FCJFUSION       | FLUID INSERT SIZE C JFUSION | Used in 21"-24" wide |
| 3    | ICJFUSION1414RP | INN COV FUSION RP FLD 14X14 |                      |
| 3    | ICJFUSION1415RP | INN COV FUSION RP FLD 14X15 |                      |
| 3    | ICJFUSION1416RP | INN COV FUSION RP FLD 14X16 |                      |
| 3    | ICJFUSION1417RP | INN COV FUSION RP FLD 14X17 |                      |
| 3    | ICJFUSION1418RP | INN COV FUSION RP FLD 14X18 |                      |
| 3    | ICJFUSION1514RP | INN COV FUSION RP FLD 15X14 |                      |
| 3    | ICJFUSION1515RP | INN COV FUSION RP FLD 15X15 |                      |
| 3    | ICJFUSION1516RP | INN COV FUSION RP FLD 15X16 |                      |
| 3    | ICJFUSION1517RP | INN COV FUSION RP FLD 15X17 |                      |
| 3    | ICJFUSION1518RP | INN COV FUSION RP FLD 15X18 |                      |
| 3    | ICJFUSION1519RP | INN COV FUSION RP FLD 15X19 |                      |
| 3    | ICJFUSION1520RP | INN COV FUSION RP FLD 15X19 |                      |
| 3    | ICJFUSION1521RP | INN COV FUSION RP FLD 15X21 |                      |
| 3    | ICJFUSION1614RP | INN COV FUSION RP FLD 16X14 |                      |
| 3    | ICJFUSION1615RP | INN COV FUSION RP FLD 16X15 |                      |
| 3    | ICJFUSION1616RP | INN COV FUSION RP FLD 16X16 |                      |
| 3    | ICJFUSION1617RP | INN COV FUSION RP FLD 16X17 |                      |
| 3    | ICJFUSION1618RP | INN COV FUSION RP FLD 16X18 |                      |
| 3    | ICJFUSION1619RP | INN COV FUSION RP FLD 16X19 |                      |
| 3    | ICJFUSION1620RP | INN COV FUSION RP FLD 16X20 |                      |
| 3    | ICJFUSION1621RP | INN COV FUSION RP FLD 16X21 |                      |
| 3    | ICJFUSION1622RP | INN COV FUSION RP FLD 16X22 |                      |
| 3    | ICJFUSION1714RP | INN COV FUSION RP FLD 17X14 |                      |
| 3    | ICJFUSION1715RP | INN COV FUSION RP FLD 17X15 |                      |
| 3    | ICJFUSION1716RP | INN COV FUSION RP FLD 17X16 |                      |
| 3    | ICJFUSION1717RP | INN COV FUSION RP FLD 17X17 |                      |
| 3    | ICJFUSION1718RP | INN COV FUSION RP FLD 17X18 |                      |
| 3    | ICJFUSION1719RP | INN COV FUSION RP FLD 17X19 |                      |
| 3    | ICJFUSION1720RP | INN COV FUSION RP FLD 17X20 |                      |
| 3    | ICJFUSION1721RP | INN COV FUSION RP FLD 17X21 |                      |
| 3    | ICJFUSION1722RP | INN COV FUSION RP FLD 17X22 |                      |
| 3    | ICJFUSION1723RP | INN COV FUSION RP FLD 17X23 |                      |
| 3    | ICJFUSION1814RP | INN COV FUSION RP FLD 18X14 |                      |

| Pos. | Item Number     | Description                 | Remarks |
|------|-----------------|-----------------------------|---------|
| 3    | ICJFUSION1815RP | INN COV FUSION RP FLD 18X15 |         |
| 3    | ICJFUSION1816RP | INN COV FUSION RP FLD 18X16 |         |
| 3    | ICJFUSION1817RP | INN COV FUSION RP FLD 18X17 |         |
| 3    | ICJFUSION1818RP | INN COV FUSION RP FLD 18X18 |         |
| 3    | ICJFUSION1819RP | INN COV FUSION RP FLD 18X19 |         |
| 3    | ICJFUSION1820RP | INN COV FUSION RP FLD 18X20 |         |
| 3    | ICJFUSION1821RP | INN COV FUSION RP FLD 18X21 |         |
| 3    | ICJFUSION1822RP | INN COV FUSION RP FLD 18X22 |         |
| 3    | ICJFUSION1823RP | INN COV FUSION RP FLD 18X23 |         |
| 3    | ICJFUSION1824RP | INN COV FUSION RP FLD 18X24 |         |
| 3    | ICJFUSION1915RP | INN COV FUSION RP FLD 19X15 |         |
| 3    | ICJFUSION1916RP | INN COV FUSION RP FLD 19X16 |         |
| 3    | ICJFUSION1917RP | INN COV FUSION RP FLD 19X17 |         |
| 3    | ICJFUSION1918RP | INN COV FUSION RP FLD 19X18 |         |
| 3    | ICJFUSION1919RP | INN COV FUSION RP FLD 19X19 |         |
| 3    | ICJFUSION1920RP | INN COV FUSION RP FLD 19X20 |         |
| 3    | ICJFUSION1921RP | INN COV FUSION RP FLD 19X21 |         |
| 3    | ICJFUSION1922RP | INN COV FUSION RP FLD 19X22 |         |
| 3    | ICJFUSION1923RP | INN COV FUSION RP FLD 19X23 |         |
| 3    | ICJFUSION1924RP | INN COV FUSION RP FLD 19X24 |         |
| 3    | ICJFUSION2015RP | INN COV FUSION RP FLD 20X15 |         |
| 3    | ICJFUSION2016RP | INN COV FUSION RP FLD 20X16 |         |
| 3    | ICJFUSION2017RP | INN COV FUSION RP FLD 20X17 |         |
| 3    | ICJFUSION2018RP | INN COV FUSION RP FLD 20X18 |         |
| 3    | ICJFUSION2019RP | INN COV FUSION RP FLD 20X19 |         |
| 3    | ICJFUSION2020RP | INN COV FUSION RP FLD 20X20 |         |
| 3    | ICJFUSION2021RP | INN COV FUSION RP FLD 20X21 |         |
| 3    | ICJFUSION2022RP | INN COV FUSION RP FLD 20X22 |         |
| 3    | ICJFUSION2023RP | INN COV FUSION RP FLD 20X23 |         |
| 3    | ICJFUSION2024RP | INN COV FUSION RP FLD 20X24 |         |
| 3    | ICJFUSION2117RP | INN COV FUSION RP FLD 21X17 |         |
| 3    | ICJFUSION2118RP | INN COV FUSION RP FLD 21X18 |         |
| 3    | ICJFUSION2119RP | INN COV FUSION RP FLD 21X19 |         |
| 3    | ICJFUSION2120RP | INN COV FUSION RP FLD 21X20 |         |
| 3    | ICJFUSION2121RP | INN COV FUSION RP FLD 21X21 |         |
| 3    | ICJFUSION2122RP | INN COV FUSION RP FLD 21X22 |         |
| 3    | ICJFUSION2123RP | INN COV FUSION RP FLD 21X23 |         |
| 3    | ICJFUSION2124RP | INN COV FUSION RP FLD 21X24 |         |
| 3    | ICJFUSION2217RP | INN COV FUSION RP FLD 22X17 |         |
| 3    | ICJFUSION2218RP | INN COV FUSION RP FLD 22X18 |         |
| 3    | ICJFUSION2219RP | INN COV FUSION RP FLD 22X19 |         |
| 3    | ICJFUSION2220RP | INN COV FUSION RP FLD 22X20 |         |
| 3    | ICJFUSION2221RP | INN COV FUSION RP FLD 22X21 |         |
| 3    | ICJFUSION2222RP | INN COV FUSION RP FLD 22X22 |         |
| 3    | ICJFUSION2223RP | INN COV FUSION RP FLD 22X23 |         |
| 3    | ICJFUSION2224RP | INN COV FUSION RP FLD 22X24 |         |
| 3    | ICJFUSION2319RP | INN COV FUSION RP FLD 23X19 |         |
| 3    | ICJFUSION2320RP | INN COV FUSION RP FLD 23X20 |         |
| 3    | ICJFUSION2321RP | INN COV FUSION RP FLD 23X21 |         |
| 3    | ICJFUSION2322RP | INN COV FUSION RP FLD 23X22 |         |
| 3    | ICJFUSION2323RP | INN COV FUSION RP FLD 23X23 |         |
| 3    | ICJFUSION2324RP | INN COV FUSION RP FLD 23X24 |         |
| 3    | ICJFUSION2419RP | INN COV FUSION RP FLD 24X19 |         |
| 3    | ICJFUSION2420RP | INN COV FUSION RP FLD 24X20 |         |
| 3    | ICJFUSION2421RP | INN COV FUSION RP FLD 24X21 |         |
| 3    | ICJFUSION2422RP | INN COV FUSION RP FLD 24X22 |         |
| 3    | ICJFUSION2423RP | INN COV FUSION RP FLD 24X23 |         |

| Pos. | Item Number     | Description                 | Remarks |
|------|-----------------|-----------------------------|---------|
| 3    | ICJFUSION2424RP | INN COV FUSION RP FLD 24X24 |         |
| 4    | BJFUSION1414RP  | BASE JAY FUSION RPROF 14X14 |         |
| 4    | BJFUSION1415RP  | BASE JAY FUSION RPROF 14X15 |         |
| 4    | BJFUSION1416RP  | BASE JAY FUSION RPROF 14X16 |         |
| 4    | BJFUSION1417RP  | BASE JAY FUSION RPROF 14X17 |         |
| 4    | BJFUSION1418RP  | BASE JAY FUSION RPROF 14X18 |         |
| 4    | BJFUSION1514RP  | BASE JAY FUSION RPROF 15X14 |         |
| 4    | BJFUSION1515RP  | BASE JAY FUSION RPROF 15X15 |         |
| 4    | BJFUSION1516RP  | BASE JAY FUSION RPROF 15X16 |         |
| 4    | BJFUSION1517RP  | BASE JAY FUSION RPROF 15X17 |         |
| 4    | BJFUSION1518RP  | BASE JAY FUSION RPROF 15X18 |         |
| 4    | BJFUSION1519RP  | BASE JAY FUSION RPROF 15X19 |         |
| 4    | BJFUSION1520RP  | BASE JAY FUSION RPROF 15X19 |         |
| 4    | BJFUSION1521RP  | BASE JAY FUSION RPROF 15X21 |         |
| 4    | BJFUSION1614RP  | BASE JAY FUSION RPROF 16X14 |         |
| 4    | BJFUSION1615RP  | BASE JAY FUSION RPROF 16X15 |         |
| 4    | BJFUSION1616RP  | BASE JAY FUSION RPROF 16X16 |         |
| 4    | BJFUSION1617RP  | BASE JAY FUSION RPROF 16X17 |         |
| 4    | BJFUSION1618RP  | BASE JAY FUSION RPROF 16X18 |         |
| 4    | BJFUSION1619RP  | BASE JAY FUSION RPROF 16X19 |         |
| 4    | BJFUSION1620RP  | BASE JAY FUSION RPROF 16X20 |         |
| 4    | BJFUSION1621RP  | BASE JAY FUSION RPROF 16X21 |         |
| 4    | BJFUSION1622RP  | BASE JAY FUSION RPROF 16X22 |         |
| 4    | BJFUSION1714RP  | BASE JAY FUSION RPROF 17X14 |         |
| 4    | BJFUSION1715RP  | BASE JAY FUSION RPROF 17X15 |         |
| 4    | BJFUSION1716RP  | BASE JAY FUSION RPROF 17X16 |         |
| 4    | BJFUSION1717RP  | BASE JAY FUSION RPROF 17X17 |         |
| 4    | BJFUSION1718RP  | BASE JAY FUSION RPROF 17X18 |         |
| 4    | BJFUSION1719RP  | BASE JAY FUSION RPROF 17X19 |         |
| 4    | BJFUSION1720RP  | BASE JAY FUSION RPROF 17X20 |         |
| 4    | BJFUSION1721RP  | BASE JAY FUSION RPROF 17X21 |         |
| 4    | BJFUSION1722RP  | BASE JAY FUSION RPROF 17X22 |         |
| 4    | BJFUSION1723RP  | BASE JAY FUSION RPROF 17X23 |         |
| 4    | BJFUSION1814RP  | BASE JAY FUSION RPROF 18X14 |         |
| 4    | BJFUSION1815RP  | BASE JAY FUSION RPROF 18X15 |         |
| 4    | BJFUSION1816RP  | BASE JAY FUSION RPROF 18X16 |         |
| 4    | BJFUSION1817RP  | BASE JAY FUSION RPROF 18X17 |         |
| 4    | BJFUSION1818RP  | BASE JAY FUSION RPROF 18X18 |         |
| 4    | BJFUSION1819RP  | BASE JAY FUSION RPROF 18X19 |         |
| 4    | BJFUSION1820RP  | BASE JAY FUSION RPROF 18X20 |         |
| 4    | BJFUSION1821RP  | BASE JAY FUSION RPROF 18X21 |         |
| 4    | BJFUSION1822RP  | BASE JAY FUSION RPROF 18X22 |         |
| 4    | BJFUSION1823RP  | BASE JAY FUSION RPROF 18X23 |         |
| 4    | BJFUSION1824RP  | BASE JAY FUSION RPROF 18X24 |         |
| 4    | BJFUSION1915RP  | BASE JAY FUSION RPROF 19X15 |         |
| 4    | BJFUSION1916RP  | BASE JAY FUSION RPROF 19X16 |         |
| 4    | BJFUSION1917RP  | BASE JAY FUSION RPROF 19X17 |         |
| 4    | BJFUSION1918RP  | BASE JAY FUSION RPROF 19X18 |         |
| 4    | BJFUSION1919RP  | BASE JAY FUSION RPROF 19X19 |         |
| 4    | BJFUSION1920RP  | BASE JAY FUSION RPROF 19X20 |         |
| 4    | BJFUSION1921RP  | BASE JAY FUSION RPROF 19X21 |         |
| 4    | BJFUSION1922RP  | BASE JAY FUSION RPROF 19X22 |         |
| 4    | BJFUSION1923RP  | BASE JAY FUSION RPROF 19X23 |         |
| 4    | BJFUSION1924RP  | BASE JAY FUSION RPROF 19X24 |         |
| 4    | BJFUSION2015RP  | BASE JAY FUSION RPROF 20X15 |         |
| 4    | BJFUSION2016RP  | BASE JAY FUSION RPROF 20X16 |         |
| 4    | BJFUSION2017RP  | BASE JAY FUSION RPROF 20X17 |         |

| Pos. | Item Number    | Description                 | Remarks |
|------|----------------|-----------------------------|---------|
| 4    | BJFUSION2018RP | BASE JAY FUSION RPROF 20X18 |         |
| 4    | BJFUSION2019RP | BASE JAY FUSION RPROF 20X19 |         |
| 4    | BJFUSION2020RP | BASE JAY FUSION RPROF 20X20 |         |
| 4    | BJFUSION2021RP | BASE JAY FUSION RPROF 20X21 |         |
| 4    | BJFUSION2022RP | BASE JAY FUSION RPROF 20X22 |         |
| 4    | BJFUSION2023RP | BASE JAY FUSION RPROF 20X23 |         |
| 4    | BJFUSION2024RP | BASE JAY FUSION RPROF 20X24 |         |
| 4    | BJFUSION2117RP | BASE JAY FUSION RPROF 21X17 |         |
| 4    | BJFUSION2118RP | BASE JAY FUSION RPROF 21X18 |         |
| 4    | BJFUSION2119RP | BASE JAY FUSION RPROF 21X19 |         |
| 4    | BJFUSION2120RP | BASE JAY FUSION RPROF 21X20 |         |
| 4    | BJFUSION2121RP | BASE JAY FUSION RPROF 21X21 |         |
| 4    | BJFUSION2122RP | BASE JAY FUSION RPROF 21X22 |         |
| 4    | BJFUSION2123RP | BASE JAY FUSION RPROF 21X23 |         |
| 4    | BJFUSION2124RP | BASE JAY FUSION RPROF 21X24 |         |
| 4    | BJFUSION2217RP | BASE JAY FUSION RPROF 22X17 |         |
| 4    | BJFUSION2218RP | BASE JAY FUSION RPROF 22X18 |         |
| 4    | BJFUSION2219RP | BASE JAY FUSION RPROF 22X19 |         |
| 4    | BJFUSION2220RP | BASE JAY FUSION RPROF 22X20 |         |
| 4    | BJFUSION2221RP | BASE JAY FUSION RPROF 22X21 |         |
| 4    | BJFUSION2222RP | BASE JAY FUSION RPROF 22X22 |         |
| 4    | BJFUSION2223RP | BASE JAY FUSION RPROF 22X23 |         |
| 4    | BJFUSION2224RP | BASE JAY FUSION RPROF 22X24 |         |
| 4    | BJFUSION2319RP | BASE JAY FUSION RPROF 23X19 |         |
| 4    | BJFUSION2320RP | BASE JAY FUSION RPROF 23X20 |         |
| 4    | BJFUSION2321RP | BASE JAY FUSION RPROF 23X21 |         |
| 4    | BJFUSION2322RP | BASE JAY FUSION RPROF 23X22 |         |
| 4    | BJFUSION2323RP | BASE JAY FUSION RPROF 23X23 |         |
| 4    | BJFUSION2324RP | BASE JAY FUSION RPROF 23X24 |         |
| 4    | BJFUSION2419RP | BASE JAY FUSION RPROF 24X19 |         |
| 4    | BJFUSION2420RP | BASE JAY FUSION RPROF 24X20 |         |
| 4    | BJFUSION2421RP | BASE JAY FUSION RPROF 24X21 |         |
| 4    | BJFUSION2422RP | BASE JAY FUSION RPROF 24X22 |         |
| 4    | BJFUSION2423RP | BASE JAY FUSION RPROF 24X23 |         |
| 4    | BJFUSION2424   | BASE JFUSION2424            |         |

Note: Cushion widths 14" & 15" use A size Air Bladders, 16" - 20" use B size bladders and 21"-24" use C size bladders



| Pos. | Item Number    | Description                  | Remarks |
|------|----------------|------------------------------|---------|
| 1    | CJFUSION1414PA | COV OUT FUSION POS AIR 14X14 |         |
| 1    | CJFUSION1415PA | COV OUT FUSION POS AIR 14X15 |         |
| 1    | CJFUSION1416PA | COV OUT FUSION POS AIR 14X16 |         |
| 1    | CJFUSION1417PA | COV OUT FUSION POS AIR 14X17 |         |
| 1    | CJFUSION1418PA | COV OUT FUSION POS AIR 14X18 |         |
| 1    | CJFUSION1514PA | COV OUT FUSION POS AIR 15X14 |         |
| 1    | CJFUSION1515PA | COV OUT FUSION POS AIR 15X15 |         |
| 1    | CJFUSION1516PA | COV OUT FUSION POS AIR 15X16 |         |
| 1    | CJFUSION1517PA | COV OUT FUSION POS AIR 15X17 |         |
| 1    | CJFUSION1518PA | COV OUT FUSION POS AIR 15X18 |         |
| 1    | CJFUSION1519PA | COV OUT FUSION POS AIR 15X19 |         |
| 1    | CJFUSION1520PA | COV OUT FUSION POS AIR 15X19 |         |
| 1    | CJFUSION1521PA | COV OUT FUSION POS AIR 15X21 |         |
| 1    | CJFUSION1614PA | COV OUT FUSION POS AIR 16X14 |         |
| 1    | CJFUSION1615PA | COV OUT FUSION POS AIR 16X15 |         |

| Pos. | Item Number    | Description                  | Remarks |
|------|----------------|------------------------------|---------|
| 1    | CJFUSION1616PA | COV OUT FUSION POS AIR 16X16 |         |
| 1    | CJFUSION1617PA | COV OUT FUSION POS AIR 16X17 |         |
| 1    | CJFUSION1618PA | COV OUT FUSION POS AIR 16X18 |         |
| 1    | CJFUSION1619PA | COV OUT FUSION POS AIR 16X19 |         |
| 1    | CJFUSION1620PA | COV OUT FUSION POS AIR 16X20 |         |
| 1    | CJFUSION1621PA | COV OUT FUSION POS AIR 16X21 |         |
| 1    | CJFUSION1622PA | COV OUT FUSION POS AIR 16X22 |         |
| 1    | CJFUSION1714PA | COV OUT FUSION POS AIR 17X14 |         |
| 1    | CJFUSION1715PA | COV OUT FUSION POS AIR 17X15 |         |
| 1    | CJFUSION1716PA | COV OUT FUSION POS AIR 17X16 |         |
| 1    | CJFUSION1717PA | COV OUT FUSION POS AIR 17X17 |         |
| 1    | CJFUSION1718PA | COV OUT FUSION POS AIR 17X18 |         |
| 1    | CJFUSION1719PA | COV OUT FUSION POS AIR 17X19 |         |
| 1    | CJFUSION1720PA | COV OUT FUSION POS AIR 17X20 |         |
| 1    | CJFUSION1721PA | COV OUT FUSION POS AIR 17X21 |         |
| 1    | CJFUSION1722PA | COV OUT FUSION POS AIR 17X22 |         |
| 1    | CJFUSION1723PA | COV OUT FUSION POS AIR 17X23 |         |
| 1    | CJFUSION1814PA | COV OUT FUSION POS AIR 18X14 |         |
| 1    | CJFUSION1815PA | COV OUT FUSION POS AIR 18X15 |         |
| 1    | CJFUSION1816PA | COV OUT FUSION POS AIR 18X16 |         |
| 1    | CJFUSION1817PA | COV OUT FUSION POS AIR 18X17 |         |
| 1    | CJFUSION1818PA | COV OUT FUSION POS AIR 18X18 |         |
| 1    | CJFUSION1819PA | COV OUT FUSION POS AIR 18X19 |         |
| 1    | CJFUSION1820PA | COV OUT FUSION POS AIR 18X20 |         |
| 1    | CJFUSION1821PA | COV OUT FUSION POS AIR 18X21 |         |
| 1    | CJFUSION1822PA | COV OUT FUSION POS AIR 18X22 |         |
| 1    | CJFUSION1823PA | COV OUT FUSION POS AIR 18X23 |         |
| 1    | CJFUSION1824PA | COV OUT FUSION POS AIR 18X24 |         |
| 1    | CJFUSION1915PA | COV OUT FUSION POS AIR 19X15 |         |
| 1    | CJFUSION1916PA | COV OUT FUSION POS AIR 19X16 |         |
| 1    | CJFUSION1917PA | COV OUT FUSION POS AIR 19X17 |         |
| 1    | CJFUSION1918PA | COV OUT FUSION POS AIR 19X18 |         |
| 1    | CJFUSION1919PA | COV OUT FUSION POS AIR 19X19 |         |
| 1    | CJFUSION1920PA | COV OUT FUSION POS AIR 19X20 |         |
| 1    | CJFUSION1921PA | COV OUT FUSION POS AIR 19X21 |         |
| 1    | CJFUSION1922PA | COV OUT FUSION POS AIR 19X22 |         |
| 1    | CJFUSION1923PA | COV OUT FUSION POS AIR 19X23 |         |
| 1    | CJFUSION1924PA | COV OUT FUSION POS AIR 19X24 |         |
| 1    | CJFUSION2015PA | COV OUT FUSION POS AIR 20X15 |         |
| 1    | CJFUSION2016PA | COV OUT FUSION POS AIR 20X16 |         |
| 1    | CJFUSION2017PA | COV OUT FUSION POS AIR 20X17 |         |
| 1    | CJFUSION2018PA | COV OUT FUSION POS AIR 20X18 |         |
| 1    | CJFUSION2019PA | COV OUT FUSION POS AIR 20X19 |         |
| 1    | CJFUSION2020PA | COV OUT FUSION POS AIR 20X20 |         |
| 1    | CJFUSION2021PA | COV OUT FUSION POS AIR 20X21 |         |
| 1    | CJFUSION2022PA | COV OUT FUSION POS AIR 20X22 |         |
| 1    | CJFUSION2023PA | COV OUT FUSION POS AIR 20X23 |         |
| 1    | CJFUSION2024PA | COV OUT FUSION POS AIR 20X24 |         |
| 1    | CJFUSION2117PA | COV OUT FUSION POS AIR 21X17 |         |
| 1    | CJFUSION2118PA | COV OUT FUSION POS AIR 21X18 |         |
| 1    | CJFUSION2119PA | COV OUT FUSION POS AIR 21X19 |         |
| 1    | CJFUSION2120PA | COV OUT FUSION POS AIR 21X20 |         |
| 1    | CJFUSION2121PA | COV OUT FUSION POS AIR 21X21 |         |
| 1    | CJFUSION2122PA | COV OUT FUSION POS AIR 21X22 |         |
| 1    | CJFUSION2123PA | COV OUT FUSION POS AIR 21X23 |         |
| 1    | CJFUSION2124PA | COV OUT FUSION POS AIR 21X24 |         |
| 1    | CJFUSION2217PA | COV OUT FUSION POS AIR 22X17 |         |

| Pos. | Item Number     | Description                    | Remarks              |
|------|-----------------|--------------------------------|----------------------|
| 1    | CJFUSION2218PA  | COV OUT FUSION POS AIR 22X18   |                      |
| 1    | CJFUSION2219PA  | COV OUT FUSION POS AIR 22X19   |                      |
| 1    | CJFUSION2220PA  | COV OUT FUSION POS AIR 22X20   |                      |
| 1    | CJFUSION2221PA  | COV OUT FUSION POS AIR 22X21   |                      |
| 1    | CJFUSION2222PA  | COV OUT FUSION POS AIR 22X22   |                      |
| 1    | CJFUSION2223PA  | COV OUT FUSION POS AIR 22X23   |                      |
| 1    | CJFUSION2224PA  | COV OUT FUSION POS AIR 22X24   |                      |
| 1    | CJFUSION2319PA  | COV OUT FUSION POS AIR 23X19   |                      |
| 1    | CJFUSION2320PA  | COV OUT FUSION POS AIR 23X20   |                      |
| 1    | CJFUSION2321PA  | COV OUT FUSION POS AIR 23X21   |                      |
| 1    | CJFUSION2322PA  | COV OUT FUSION POS AIR 23X22   |                      |
| 1    | CJFUSION2323PA  | COV OUT FUSION POS AIR 23X23   |                      |
| 1    | CJFUSION2324PA  | COV OUT FUSION POS AIR 23X24   |                      |
| 1    | CJFUSION2419PA  | COV OUT FUSION POS AIR 24X19   |                      |
| 1    | CJFUSION2420PA  | COV OUT FUSION POS AIR 24X20   |                      |
| 1    | CJFUSION2421PA  | COV OUT FUSION POS AIR 24X21   |                      |
| 1    | CJFUSION2422PA  | COV OUT FUSION POS AIR 24X22   |                      |
| 1    | CJFUSION2423PA  | COV OUT FUSION POS AIR 24X23   |                      |
| 1    | CJFUSION2424PA  | COV OUT FUSION POS AIR 24X24   |                      |
| 2    | ASAJFUSION      | AIR INSERT SINGLE VALVE SIZE A | Used in 14"&15" wide |
| 2    | ASBJFUSION      | AIR INSERT SINGLE VALVE SIZE B | Used in 16"-20" wide |
| 2    | ASCJFUSION      | AIR INSERT SINGLE VALVE SIZE C | Used in 21"-24" wide |
| 2    | ADAJFUSION      | AIR INSERT DUAL VALVE SIZE A   | Used in 14"&15" wide |
| 2    | ADBJFUSION      | AIR INSERT DUAL VALVE SIZE B   | Used in 16"-20" wide |
| 2    | ADCJFUSION      | AIR INSERT DUAL VALVE SIZE C   | Used in 21"-24" wide |
| 3    | ICJFUSION1414PA | INN COV FUSION POS AIR 14X14   |                      |
| 3    | ICJFUSION1415PA | INN COV FUSION POS AIR 14X15   |                      |
| 3    | ICJFUSION1416PA | INN COV FUSION POS AIR 14X16   |                      |
| 3    | ICJFUSION1417PA | INN COV FUSION POS AIR 14X17   |                      |
| 3    | ICJFUSION1418PA | INN COV FUSION POS AIR 14X18   |                      |
| 3    | ICJFUSION1514PA | INN COV FUSION POS AIR 15X14   |                      |
| 3    | ICJFUSION1515PA | INN COV FUSION POS AIR 15X15   |                      |
| 3    | ICJFUSION1516PA | INN COV FUSION POS AIR 15X16   |                      |
| 3    | ICJFUSION1517PA | INN COV FUSION POS AIR 15X17   |                      |
| 3    | ICJFUSION1518PA | INN COV FUSION POS AIR 15X18   |                      |
| 3    | ICJFUSION1519PA | INN COV FUSION POS AIR 15X19   |                      |
| 3    | ICJFUSION1520PA | INN COV FUSION POS AIR 15X19   |                      |
| 3    | ICJFUSION1521PA | INN COV FUSION POS AIR 15X21   |                      |
| 3    | ICJFUSION1614PA | INN COV FUSION POS AIR 16X14   |                      |
| 3    | ICJFUSION1615PA | INN COV FUSION POS AIR 16X15   |                      |
| 3    | ICJFUSION1616PA | INN COV FUSION POS AIR 16X16   |                      |
| 3    | ICJFUSION1617PA | INN COV FUSION POS AIR 16X17   |                      |
| 3    | ICJFUSION1618PA | INN COV FUSION POS AIR 16X18   |                      |
| 3    | ICJFUSION1619PA | INN COV FUSION POS AIR 16X19   |                      |
| 3    | ICJFUSION1620PA | INN COV FUSION POS AIR 16X20   |                      |
| 3    | ICJFUSION1621PA | INN COV FUSION POS AIR 16X21   |                      |
| 3    | ICJFUSION1622PA | INN COV FUSION POS AIR 16X22   |                      |
| 3    | ICJFUSION1714PA | INN COV FUSION POS AIR 17X14   |                      |
| 3    | ICJFUSION1715PA | INN COV FUSION POS AIR 17X15   |                      |
| 3    | ICJFUSION1716PA | INN COV FUSION POS AIR 17X16   |                      |
| 3    | ICJFUSION1717PA | INN COV FUSION POS AIR 17X17   |                      |
| 3    | ICJFUSION1718PA | INN COV FUSION POS AIR 17X18   |                      |
| 3    | ICJFUSION1719PA | INN COV FUSION POS AIR 17X19   |                      |
| 3    | ICJFUSION1720PA | INN COV FUSION POS AIR 17X20   |                      |
| 3    | ICJFUSION1721PA | INN COV FUSION POS AIR 17X21   |                      |
| 3    | ICJFUSION1722PA | INN COV FUSION POS AIR 17X22   |                      |
| 3    | ICJFUSION1723PA | INN COV FUSION POS AIR 17X23   |                      |

| Pos. | Item Number     | Description                  | Remarks |
|------|-----------------|------------------------------|---------|
| 3    | ICJFUSION1814PA | INN COV FUSION POS AIR 18X14 |         |
| 3    | ICJFUSION1815PA | INN COV FUSION POS AIR 18X15 |         |
| 3    | ICJFUSION1816PA | INN COV FUSION POS AIR 18X16 |         |
| 3    | ICJFUSION1817PA | INN COV FUSION POS AIR 18X17 |         |
| 3    | ICJFUSION1818PA | INN COV FUSION POS AIR 18X18 |         |
| 3    | ICJFUSION1819PA | INN COV FUSION POS AIR 18X19 |         |
| 3    | ICJFUSION1820PA | INN COV FUSION POS AIR 18X20 |         |
| 3    | ICJFUSION1821PA | INN COV FUSION POS AIR 18X21 |         |
| 3    | ICJFUSION1822PA | INN COV FUSION POS AIR 18X22 |         |
| 3    | ICJFUSION1823PA | INN COV FUSION POS AIR 18X23 |         |
| 3    | ICJFUSION1824PA | INN COV FUSION POS AIR 18X24 |         |
| 3    | ICJFUSION1915PA | INN COV FUSION POS AIR 19X15 |         |
| 3    | ICJFUSION1916PA | INN COV FUSION POS AIR 19X16 |         |
| 3    | ICJFUSION1917PA | INN COV FUSION POS AIR 19X17 |         |
| 3    | ICJFUSION1918PA | INN COV FUSION POS AIR 19X18 |         |
| 3    | ICJFUSION1919PA | INN COV FUSION POS AIR 19X19 |         |
| 3    | ICJFUSION1920PA | INN COV FUSION POS AIR 19X20 |         |
| 3    | ICJFUSION1921PA | INN COV FUSION POS AIR 19X21 |         |
| 3    | ICJFUSION1922PA | INN COV FUSION POS AIR 19X22 |         |
| 3    | ICJFUSION1923PA | INN COV FUSION POS AIR 19X23 |         |
| 3    | ICJFUSION1924PA | INN COV FUSION POS AIR 19X24 |         |
| 3    | ICJFUSION2015PA | INN COV FUSION POS AIR 20X15 |         |
| 3    | ICJFUSION2016PA | INN COV FUSION POS AIR 20X16 |         |
| 3    | ICJFUSION2017PA | INN COV FUSION POS AIR 20X17 |         |
| 3    | ICJFUSION2018PA | INN COV FUSION POS AIR 20X18 |         |
| 3    | ICJFUSION2019PA | INN COV FUSION POS AIR 20X19 |         |
| 3    | ICJFUSION2020PA | INN COV FUSION POS AIR 20X20 |         |
| 3    | ICJFUSION2021PA | INN COV FUSION POS AIR 20X21 |         |
| 3    | ICJFUSION2022PA | INN COV FUSION POS AIR 20X22 |         |
| 3    | ICJFUSION2023PA | INN COV FUSION POS AIR 20X23 |         |
| 3    | ICJFUSION2024PA | INN COV FUSION POS AIR 20X24 |         |
| 3    | ICJFUSION2117PA | INN COV FUSION POS AIR 21X17 |         |
| 3    | ICJFUSION2118PA | INN COV FUSION POS AIR 21X18 |         |
| 3    | ICJFUSION2119PA | INN COV FUSION POS AIR 21X19 |         |
| 3    | ICJFUSION2120PA | INN COV FUSION POS AIR 21X20 |         |
| 3    | ICJFUSION2121PA | INN COV FUSION POS AIR 21X21 |         |
| 3    | ICJFUSION2122PA | INN COV FUSION POS AIR 21X22 |         |
| 3    | ICJFUSION2123PA | INN COV FUSION POS AIR 21X23 |         |
| 3    | ICJFUSION2124PA | INN COV FUSION POS AIR 21X24 |         |
| 3    | ICJFUSION2217PA | INN COV FUSION POS AIR 22X17 |         |
| 3    | ICJFUSION2218PA | INN COV FUSION POS AIR 22X18 |         |
| 3    | ICJFUSION2219PA | INN COV FUSION POS AIR 22X19 |         |
| 3    | ICJFUSION2220PA | INN COV FUSION POS AIR 22X20 |         |
| 3    | ICJFUSION2221PA | INN COV FUSION POS AIR 22X21 |         |
| 3    | ICJFUSION2222PA | INN COV FUSION POS AIR 22X22 |         |
| 3    | ICJFUSION2223PA | INN COV FUSION POS AIR 22X23 |         |
| 3    | ICJFUSION2224PA | INN COV FUSION POS AIR 22X24 |         |
| 3    | ICJFUSION2319PA | INN COV FUSION POS AIR 23X19 |         |
| 3    | ICJFUSION2320PA | INN COV FUSION POS AIR 23X20 |         |
| 3    | ICJFUSION2321PA | INN COV FUSION POS AIR 23X21 |         |
| 3    | ICJFUSION2322PA | INN COV FUSION POS AIR 23X22 |         |
| 3    | ICJFUSION2323PA | INN COV FUSION POS AIR 23X23 |         |
| 3    | ICJFUSION2324PA | INN COV FUSION POS AIR 23X24 |         |
| 3    | ICJFUSION2419PA | INN COV FUSION POS AIR 24X19 |         |
| 3    | ICJFUSION2420PA | INN COV FUSION POS AIR 24X20 |         |
| 3    | ICJFUSION2421PA | INN COV FUSION POS AIR 24X21 |         |
| 3    | ICJFUSION2422PA | INN COV FUSION POS AIR 24X22 |         |

| Pos. | Item Number     | Description                  | Remarks |
|------|-----------------|------------------------------|---------|
| 3    | ICJFUSION2423PA | INN COV FUSION POS AIR 24X23 |         |
| 3    | ICJFUSION2424PA | INN COV FUSION POS AIR 24X24 |         |
| 4    | BJFUSION1414    | BASE JFUSION1414             |         |
| 4    | BJFUSION1415    | BASE JFUSION1415             |         |
| 4    | BJFUSION1416    | BASE JFUSION1416             |         |
| 4    | BJFUSION1417    | BASE JFUSION1417             |         |
| 4    | BJFUSION1418    | BASE JFUSION1418             |         |
| 4    | BJFUSION1514    | BASE JFUSION1514             |         |
| 4    | BJFUSION1515    | BASE JFUSION1515             |         |
| 4    | BJFUSION1516    | BASE JFUSION1516             |         |
| 4    | BJFUSION1517    | BASE JFUSION1517             |         |
| 4    | BJFUSION1518    | BASE JFUSION1518             |         |
| 4    | BJFUSION1519    | BASE JFUSION1519             |         |
| 4    | BJFUSION1520    | BASE JFUSION1520             |         |
| 4    | BJFUSION1521    | BASE JFUSION1521             |         |
| 4    | BJFUSION1614    | BASE JFUSION1614             |         |
| 4    | BJFUSION1615    | BASE JFUSION1615             |         |
| 4    | BJFUSION1616    | BASE JFUSION1616             |         |
| 4    | BJFUSION1617    | BASE JFUSION1617             |         |
| 4    | BJFUSION1618    | BASE JFUSION1618             |         |
| 4    | BJFUSION1619    | BASE JFUSION1619             |         |
| 4    | BJFUSION1620    | BASE JFUSION1620             |         |
| 4    | BJFUSION1621    | BASE JFUSION1621             |         |
| 4    | BJFUSION1622    | BASE JFUSION1622             |         |
| 4    | BJFUSION1714    | BASE JFUSION1714             |         |
| 4    | BJFUSION1715    | BASE JFUSION1715             |         |
| 4    | BJFUSION1716    | BASE JFUSION1716             |         |
| 4    | BJFUSION1717    | BASE JFUSION1717             |         |
| 4    | BJFUSION1718    | BASE JFUSION1718             |         |
| 4    | BJFUSION1719    | BASE JFUSION1719             |         |
| 4    | BJFUSION1720    | BASE JFUSION1720             |         |
| 4    | BJFUSION1721    | BASE JFUSION1721             |         |
| 4    | BJFUSION1722    | BASE JFUSION1722             |         |
| 4    | BJFUSION1723    | BASE JFUSION1723             |         |
| 4    | BJFUSION1814    | BASE JFUSION1814             |         |
| 4    | BJFUSION1815    | BASE JFUSION1815             |         |
| 4    | BJFUSION1816    | BASE JFUSION1816             |         |
| 4    | BJFUSION1817    | BASE JFUSION1817             |         |
| 4    | BJFUSION1818    | BASE JFUSION1818             |         |
| 4    | BJFUSION1819    | BASE JFUSION1819             |         |
| 4    | BJFUSION1820    | BASE JFUSION1820             |         |
| 4    | BJFUSION1821    | BASE JFUSION1821             |         |
| 4    | BJFUSION1822    | BASE JFUSION1822             |         |
| 4    | BJFUSION1823    | BASE JFUSION1823             |         |
| 4    | BJFUSION1824    | BASE JFUSION1824             |         |
| 4    | BJFUSION1915    | BASE JFUSION1915             |         |
| 4    | BJFUSION1916    | BASE JFUSION1916             |         |
| 4    | BJFUSION1917    | BASE JFUSION1917             |         |
| 4    | BJFUSION1918    | BASE JFUSION1918             |         |
| 4    | BJFUSION1919    | BASE JFUSION1919             |         |
| 4    | BJFUSION1920    | BASE JFUSION1920             |         |
| 4    | BJFUSION1921    | BASE JFUSION1921             |         |
| 4    | BJFUSION1922    | BASE JFUSION1922             |         |
| 4    | BJFUSION1923    | BASE JFUSION1923             |         |
| 4    | BJFUSION1924    | BASE JFUSION1924             |         |
| 4    | BJFUSION2015    | BASE JFUSION2015             |         |
| 4    | BJFUSION2016    | BASE JFUSION2016             |         |

| Pos. | Item Number  | Description       | Remarks |
|------|--------------|-------------------|---------|
| 4    | BJFUSION2017 | BASE JFUSION2017  |         |
| 4    | BJFUSION2018 | BASE JFUSION2018  |         |
| 4    | BJFUSION2019 | BASE JFUSION2019  |         |
| 4    | BJFUSION2020 | BASE JFUSION2020  |         |
| 4    | BJFUSION2021 | BASE JFUSION2021  |         |
| 4    | BJFUSION2022 | BASE JFUSION2022  |         |
| 4    | BJFUSION2023 | BASE JFUSION2023  |         |
| 4    | BJFUSION2024 | BASE JFUSION2024  |         |
| 4    | BJFUSION2117 | BASE JFUSION2117  |         |
| 4    | BJFUSION2118 | BASE JFUSION2118  |         |
| 4    | BJFUSION2119 | BASE JFUSION2119  |         |
| 4    | BJFUSION2120 | BASE JFUSION2120  |         |
| 4    | BJFUSION2121 | BASE JFUSION2121  |         |
| 4    | BJFUSION2122 | BASE JFUSION2122  |         |
| 4    | BJFUSION2123 | BASE JFUSION2123  |         |
| 4    | BJFUSION2124 | BASE JFUSION2124  |         |
| 4    | BJFUSION2217 | BASE JFUSION2217  |         |
| 4    | BJFUSION2218 | BASE JFUSION2218  |         |
| 4    | BJFUSION2219 | BASE JFUSION2219  |         |
| 4    | BJFUSION2220 | BASE JFUSION2220  |         |
| 4    | BJFUSION2221 | BASE JFUSION2221  |         |
| 4    | BJFUSION2222 | BASE JFUSION2222  |         |
| 4    | BJFUSION2223 | BASE JFUSION2223  |         |
| 4    | BJFUSION2224 | BASE JFUSION2224  |         |
| 4    | BJFUSION2319 | BASE JFUSION2319  |         |
| 4    | BJFUSION2320 | BASE JFUSION2320  |         |
| 4    | BJFUSION2321 | BASE JFUSION2321  |         |
| 4    | BJFUSION2322 | BASE JFUSION2322  |         |
| 4    | BJFUSION2323 | BASE JFUSION2323  |         |
| 4    | BJFUSION2324 | BASE JFUSION2324  |         |
| 4    | BJFUSION2419 | BASE JFUSION2419  |         |
| 4    | BJFUSION2420 | BASE JFUSION2420  |         |
| 4    | BJFUSION2421 | BASE JFUSION2421  |         |
| 4    | BJFUSION2422 | BASE JFUSION2422  |         |
| 4    | BJFUSION2423 | BASE JFUSION2423  |         |
| 4    | BJFUSION2424 | BASE JFUSION2424  |         |
| 5    | AIRPUMP      | JAY AIR PUMP      |         |
| 6    | PATCHKIT     | JAY AIR PATCH KIT |         |



| Pos. | Item Number   | Description                  | Remarks |
|------|---------------|------------------------------|---------|
| 1    | CJFUSION1414P | COV OUT FUSION POS FLD 14X14 |         |
| 1    | CJFUSION1415P | COV OUT FUSION POS FLD 14X15 |         |
| 1    | CJFUSION1416P | COV OUT FUSION POS FLD 14X16 |         |
| 1    | CJFUSION1417P | COV OUT FUSION POS FLD 14X17 |         |
| 1    | CJFUSION1418P | COV OUT FUSION POS FLD 14X18 |         |
| 1    | CJFUSION1514P | COV OUT FUSION POS FLD 15X14 |         |
| 1    | CJFUSION1515P | COV OUT FUSION POS FLD 15X15 |         |
| 1    | CJFUSION1516P | COV OUT FUSION POS FLD 15X16 |         |
| 1    | CJFUSION1517P | COV OUT FUSION POS FLD 15X17 |         |
| 1    | CJFUSION1518P | COV OUT FUSION POS FLD 15X18 |         |
| 1    | CJFUSION1519P | COV OUT FUSION POS FLD 15X19 |         |
| 1    | CJFUSION1520P | COV OUT FUSION POS FLD 15X19 |         |
| 1    | CJFUSION1521P | COV OUT FUSION POS FLD 15X21 |         |
| 1    | CJFUSION1614P | COV OUT FUSION POS FLD 16X14 |         |
| 1    | CJFUSION1615P | COV OUT FUSION POS FLD 16X15 |         |
| 1    | CJFUSION1616P | COV OUT FUSION POS FLD 16X16 |         |
| 1    | CJFUSION1617P | COV OUT FUSION POS FLD 16X17 |         |
| 1    | CJFUSION1618P | COV OUT FUSION POS FLD 16X18 |         |
| 1    | CJFUSION1619P | COV OUT FUSION POS FLD 16X19 |         |
| 1    | CJFUSION1620P | COV OUT FUSION POS FLD 16X20 |         |
| 1    | CJFUSION1621P | COV OUT FUSION POS FLD 16X21 |         |
| 1    | CJFUSION1622P | COV OUT FUSION POS FLD 16X22 |         |
| 1    | CJFUSION1714P | COV OUT FUSION POS FLD 17X14 |         |
| 1    | CJFUSION1715P | COV OUT FUSION POS FLD 17X15 |         |
| 1    | CJFUSION1716P | COV OUT FUSION POS FLD 17X16 |         |

| Pos. | Item Number   | Description                  | Remarks |
|------|---------------|------------------------------|---------|
| 1    | CJFUSION1717P | COV OUT FUSION POS FLD 17X17 |         |
| 1    | CJFUSION1718P | COV OUT FUSION POS FLD 17X18 |         |
| 1    | CJFUSION1719P | COV OUT FUSION POS FLD 17X19 |         |
| 1    | CJFUSION1720P | COV OUT FUSION POS FLD 17X20 |         |
| 1    | CJFUSION1721P | COV OUT FUSION POS FLD 17X21 |         |
| 1    | CJFUSION1722P | COV OUT FUSION POS FLD 17X22 |         |
| 1    | CJFUSION1723P | COV OUT FUSION POS FLD 17X23 |         |
| 1    | CJFUSION1814P | COV OUT FUSION POS FLD 18X14 |         |
| 1    | CJFUSION1815P | COV OUT FUSION POS FLD 18X15 |         |
| 1    | CJFUSION1816P | COV OUT FUSION POS FLD 18X16 |         |
| 1    | CJFUSION1817P | COV OUT FUSION POS FLD 18X17 |         |
| 1    | CJFUSION1818P | COV OUT FUSION POS FLD 18X18 |         |
| 1    | CJFUSION1819P | COV OUT FUSION POS FLD 18X19 |         |
| 1    | CJFUSION1820P | COV OUT FUSION POS FLD 18X20 |         |
| 1    | CJFUSION1821P | COV OUT FUSION POS FLD 18X21 |         |
| 1    | CJFUSION1822P | COV OUT FUSION POS FLD 18X22 |         |
| 1    | CJFUSION1823P | COV OUT FUSION POS FLD 18X23 |         |
| 1    | CJFUSION1824P | COV OUT FUSION POS FLD 18X24 |         |
| 1    | CJFUSION1915P | COV OUT FUSION POS FLD 19X15 |         |
| 1    | CJFUSION1916P | COV OUT FUSION POS FLD 19X16 |         |
| 1    | CJFUSION1917P | COV OUT FUSION POS FLD 19X17 |         |
| 1    | CJFUSION1918P | COV OUT FUSION POS FLD 19X18 |         |
| 1    | CJFUSION1919P | COV OUT FUSION POS FLD 19X19 |         |
| 1    | CJFUSION1920P | COV OUT FUSION POS FLD 19X20 |         |
| 1    | CJFUSION1921P | COV OUT FUSION POS FLD 19X21 |         |
| 1    | CJFUSION1922P | COV OUT FUSION POS FLD 19X22 |         |
| 1    | CJFUSION1923P | COV OUT FUSION POS FLD 19X23 |         |
| 1    | CJFUSION1924P | COV OUT FUSION POS FLD 19X24 |         |
| 1    | CJFUSION2015P | COV OUT FUSION POS FLD 20X15 |         |
| 1    | CJFUSION2016P | COV OUT FUSION POS FLD 20X16 |         |
| 1    | CJFUSION2017P | COV OUT FUSION POS FLD 20X17 |         |
| 1    | CJFUSION2018P | COV OUT FUSION POS FLD 20X18 |         |
| 1    | CJFUSION2019P | COV OUT FUSION POS FLD 20X19 |         |
| 1    | CJFUSION2020P | COV OUT FUSION POS FLD 20X20 |         |
| 1    | CJFUSION2021P | COV OUT FUSION POS FLD 20X21 |         |
| 1    | CJFUSION2022P | COV OUT FUSION POS FLD 20X22 |         |
| 1    | CJFUSION2023P | COV OUT FUSION POS FLD 20X23 |         |
| 1    | CJFUSION2024P | COV OUT FUSION POS FLD 20X24 |         |
| 1    | CJFUSION2117P | COV OUT FUSION POS FLD 21X17 |         |
| 1    | CJFUSION2118P | COV OUT FUSION POS FLD 21X18 |         |
| 1    | CJFUSION2119P | COV OUT FUSION POS FLD 21X19 |         |
| 1    | CJFUSION2120P | COV OUT FUSION POS FLD 21X20 |         |
| 1    | CJFUSION2121P | COV OUT FUSION POS FLD 21X21 |         |
| 1    | CJFUSION2122P | COV OUT FUSION POS FLD 21X22 |         |
| 1    | CJFUSION2123P | COV OUT FUSION POS FLD 21X23 |         |
| 1    | CJFUSION2124P | COV OUT FUSION POS FLD 21X24 |         |
| 1    | CJFUSION2217P | COV OUT FUSION POS FLD 22X17 |         |
| 1    | CJFUSION2218P | COV OUT FUSION POS FLD 22X18 |         |
| 1    | CJFUSION2219P | COV OUT FUSION POS FLD 22X19 |         |
| 1    | CJFUSION2220P | COV OUT FUSION POS FLD 22X20 |         |
| 1    | CJFUSION2221P | COV OUT FUSION POS FLD 22X21 |         |
| 1    | CJFUSION2222P | COV OUT FUSION POS FLD 22X22 |         |
| 1    | CJFUSION2223P | COV OUT FUSION POS FLD 22X23 |         |
| 1    | CJFUSION2224P | COV OUT FUSION POS FLD 22X24 |         |
| 1    | CJFUSION2319P | COV OUT FUSION POS FLD 23X19 |         |
| 1    | CJFUSION2320P | COV OUT FUSION POS FLD 23X20 |         |
| 1    | CJFUSION2321P | COV OUT FUSION POS FLD 23X21 |         |

| Pos. | Item Number    | Description                  | Remarks              |
|------|----------------|------------------------------|----------------------|
| 1    | CJFUSION2322P  | COV OUT FUSION POS FLD 23X22 |                      |
| 1    | CJFUSION2323P  | COV OUT FUSION POS FLD 23X23 |                      |
| 1    | CJFUSION2324P  | COV OUT FUSION POS FLD 23X24 |                      |
| 1    | CJFUSION2419P  | COV OUT FUSION POS FLD 24X19 |                      |
| 1    | CJFUSION2420P  | COV OUT FUSION POS FLD 24X20 |                      |
| 1    | CJFUSION2421P  | COV OUT FUSION POS FLD 24X21 |                      |
| 1    | CJFUSION2422P  | COV OUT FUSION POS FLD 24X22 |                      |
| 1    | CJFUSION2423P  | COV OUT FUSION POS FLD 24X23 |                      |
| 1    | CJFUSION2424P  | COV OUT FUSION POS FLD 24X24 |                      |
| 2    | FAJFUSION      | FLUID INSERT SIZE A JFUSION  | Used in 14"&15" wide |
| 2    | FBJFUSION      | FLUID INSERT SIZE B JFUSION  | Used in 16"-20" wide |
| 2    | FCJFUSION      | FLUID INSERT SIZE C JFUSION  | Used in 21"-24" wide |
| 4    | ICJFUSION1414P | INN COV FUSION POS FLD 14X14 |                      |
| 4    | ICJFUSION1415P | INN COV FUSION POS FLD 14X15 |                      |
| 4    | ICJFUSION1416P | INN COV FUSION POS FLD 14X16 |                      |
| 4    | ICJFUSION1417P | INN COV FUSION POS FLD 14X17 |                      |
| 4    | ICJFUSION1418P | INN COV FUSION POS FLD 14X18 |                      |
| 4    | ICJFUSION1514P | INN COV FUSION POS FLD 15X14 |                      |
| 4    | ICJFUSION1515P | INN COV FUSION POS FLD 15X15 |                      |
| 4    | ICJFUSION1516P | INN COV FUSION POS FLD 15X16 |                      |
| 4    | ICJFUSION1517P | INN COV FUSION POS FLD 15X17 |                      |
| 4    | ICJFUSION1518P | INN COV FUSION POS FLD 15X18 |                      |
| 4    | ICJFUSION1519P | INN COV FUSION POS FLD 15X19 |                      |
| 4    | ICJFUSION1520P | INN COV FUSION POS FLD 15X19 |                      |
| 4    | ICJFUSION1521P | INN COV FUSION POS FLD 15X21 |                      |
| 4    | ICJFUSION1614P | INN COV FUSION POS FLD 16X14 |                      |
| 4    | ICJFUSION1615P | INN COV FUSION POS FLD 16X15 |                      |
| 4    | ICJFUSION1616P | INN COV FUSION POS FLD 16X16 |                      |
| 4    | ICJFUSION1617P | INN COV FUSION POS FLD 16X17 |                      |
| 4    | ICJFUSION1618P | INN COV FUSION POS FLD 16X18 |                      |
| 4    | ICJFUSION1619P | INN COV FUSION POS FLD 16X19 |                      |
| 4    | ICJFUSION1620P | INN COV FUSION POS FLD 16X20 |                      |
| 4    | ICJFUSION1621P | INN COV FUSION POS FLD 16X21 |                      |
| 4    | ICJFUSION1622P | INN COV FUSION POS FLD 16X22 |                      |
| 4    | ICJFUSION1714P | INN COV FUSION POS FLD 17X14 |                      |
| 4    | ICJFUSION1715P | INN COV FUSION POS FLD 17X15 |                      |
| 4    | ICJFUSION1716P | INN COV FUSION POS FLD 17X16 |                      |
| 4    | ICJFUSION1717P | INN COV FUSION POS FLD 17X17 |                      |
| 4    | ICJFUSION1718P | INN COV FUSION POS FLD 17X18 |                      |
| 4    | ICJFUSION1719P | INN COV FUSION POS FLD 17X19 |                      |
| 4    | ICJFUSION1720P | INN COV FUSION POS FLD 17X20 |                      |
| 4    | ICJFUSION1721P | INN COV FUSION POS FLD 17X21 |                      |
| 4    | ICJFUSION1722P | INN COV FUSION POS FLD 17X22 |                      |
| 4    | ICJFUSION1723P | INN COV FUSION POS FLD 17X23 |                      |
| 4    | ICJFUSION1814P | INN COV FUSION POS FLD 18X14 |                      |
| 4    | ICJFUSION1815P | INN COV FUSION POS FLD 18X15 |                      |
| 4    | ICJFUSION1816P | INN COV FUSION POS FLD 18X16 |                      |
| 4    | ICJFUSION1817P | INN COV FUSION POS FLD 18X17 |                      |
| 4    | ICJFUSION1818P | INN COV FUSION POS FLD 18X18 |                      |
| 4    | ICJFUSION1819P | INN COV FUSION POS FLD 18X19 |                      |
| 4    | ICJFUSION1820P | INN COV FUSION POS FLD 18X20 |                      |
| 4    | ICJFUSION1821P | INN COV FUSION POS FLD 18X21 |                      |
| 4    | ICJFUSION1822P | INN COV FUSION POS FLD 18X22 |                      |
| 4    | ICJFUSION1823P | INN COV FUSION POS FLD 18X23 |                      |
| 4    | ICJFUSION1824P | INN COV FUSION POS FLD 18X24 |                      |
| 4    | ICJFUSION1915P | INN COV FUSION POS FLD 19X15 |                      |
| 4    | ICJFUSION1916P | INN COV FUSION POS FLD 19X16 |                      |

| Pos. | Item Number    | Description                  | Remarks |
|------|----------------|------------------------------|---------|
| 4    | ICJFUSION1917P | INN COV FUSION POS FLD 19X17 |         |
| 4    | ICJFUSION1918P | INN COV FUSION POS FLD 19X18 |         |
| 4    | ICJFUSION1919P | INN COV FUSION POS FLD 19X19 |         |
| 4    | ICJFUSION1920P | INN COV FUSION POS FLD 19X20 |         |
| 4    | ICJFUSION1921P | INN COV FUSION POS FLD 19X21 |         |
| 4    | ICJFUSION1922P | INN COV FUSION POS FLD 19X22 |         |
| 4    | ICJFUSION1923P | INN COV FUSION POS FLD 19X23 |         |
| 4    | ICJFUSION1924P | INN COV FUSION POS FLD 19X24 |         |
| 4    | ICJFUSION2015P | INN COV FUSION POS FLD 20X15 |         |
| 4    | ICJFUSION2016P | INN COV FUSION POS FLD 20X16 |         |
| 4    | ICJFUSION2017P | INN COV FUSION POS FLD 20X17 |         |
| 4    | ICJFUSION2018P | INN COV FUSION POS FLD 20X18 |         |
| 4    | ICJFUSION2019P | INN COV FUSION POS FLD 20X19 |         |
| 4    | ICJFUSION2020P | INN COV FUSION POS FLD 20X20 |         |
| 4    | ICJFUSION2021P | INN COV FUSION POS FLD 20X21 |         |
| 4    | ICJFUSION2022P | INN COV FUSION POS FLD 20X22 |         |
| 4    | ICJFUSION2023P | INN COV FUSION POS FLD 20X23 |         |
| 4    | ICJFUSION2024P | INN COV FUSION POS FLD 20X24 |         |
| 4    | ICJFUSION2117P | INN COV FUSION POS FLD 21X17 |         |
| 4    | ICJFUSION2118P | INN COV FUSION POS FLD 21X18 |         |
| 4    | ICJFUSION2119P | INN COV FUSION POS FLD 21X19 |         |
| 4    | ICJFUSION2120P | INN COV FUSION POS FLD 21X20 |         |
| 4    | ICJFUSION2121P | INN COV FUSION POS FLD 21X21 |         |
| 4    | ICJFUSION2122P | INN COV FUSION POS FLD 21X22 |         |
| 4    | ICJFUSION2123P | INN COV FUSION POS FLD 21X23 |         |
| 4    | ICJFUSION2124P | INN COV FUSION POS FLD 21X24 |         |
| 4    | ICJFUSION2217P | INN COV FUSION POS FLD 22X17 |         |
| 4    | ICJFUSION2218P | INN COV FUSION POS FLD 22X18 |         |
| 4    | ICJFUSION2219P | INN COV FUSION POS FLD 22X19 |         |
| 4    | ICJFUSION2220P | INN COV FUSION POS FLD 22X20 |         |
| 4    | ICJFUSION2221P | INN COV FUSION POS FLD 22X21 |         |
| 4    | ICJFUSION2222P | INN COV FUSION POS FLD 22X22 |         |
| 4    | ICJFUSION2223P | INN COV FUSION POS FLD 22X23 |         |
| 4    | ICJFUSION2224P | INN COV FUSION POS FLD 22X24 |         |
| 4    | ICJFUSION2319P | INN COV FUSION POS FLD 23X19 |         |
| 4    | ICJFUSION2320P | INN COV FUSION POS FLD 23X20 |         |
| 4    | ICJFUSION2321P | INN COV FUSION POS FLD 23X21 |         |
| 4    | ICJFUSION2322P | INN COV FUSION POS FLD 23X22 |         |
| 4    | ICJFUSION2323P | INN COV FUSION POS FLD 23X23 |         |
| 4    | ICJFUSION2324P | INN COV FUSION POS FLD 23X24 |         |
| 4    | ICJFUSION2419P | INN COV FUSION POS FLD 24X19 |         |
| 4    | ICJFUSION2420P | INN COV FUSION POS FLD 24X20 |         |
| 4    | ICJFUSION2421P | INN COV FUSION POS FLD 24X21 |         |
| 4    | ICJFUSION2422P | INN COV FUSION POS FLD 24X22 |         |
| 4    | ICJFUSION2423P | INN COV FUSION POS FLD 24X23 |         |
| 4    | ICJFUSION2424P | INN COV FUSION POS FLD 24X24 |         |
| 5    | BJFUSION1414   | BASE JFUSION1414             |         |
| 5    | BJFUSION1415   | BASE JFUSION1415             |         |
| 5    | BJFUSION1416   | BASE JFUSION1416             |         |
| 5    | BJFUSION1417   | BASE JFUSION1417             |         |
| 5    | BJFUSION1418   | BASE JFUSION1418             |         |
| 5    | BJFUSION1514   | BASE JFUSION1514             |         |
| 5    | BJFUSION1515   | BASE JFUSION1515             |         |
| 5    | BJFUSION1516   | BASE JFUSION1516             |         |
| 5    | BJFUSION1517   | BASE JFUSION1517             |         |
| 5    | BJFUSION1518   | BASE JFUSION1518             |         |
| 5    | BJFUSION1519   | BASE JFUSION1519             |         |

| Pos. | Item Number  | Description      | Remarks |
|------|--------------|------------------|---------|
| 5    | BJFUSION1520 | BASE JFUSION1520 |         |
| 5    | BJFUSION1521 | BASE JFUSION1521 |         |
| 5    | BJFUSION1614 | BASE JFUSION1614 |         |
| 5    | BJFUSION1615 | BASE JFUSION1615 |         |
| 5    | BJFUSION1616 | BASE JFUSION1616 |         |
| 5    | BJFUSION1617 | BASE JFUSION1617 |         |
| 5    | BJFUSION1618 | BASE JFUSION1618 |         |
| 5    | BJFUSION1619 | BASE JFUSION1619 |         |
| 5    | BJFUSION1620 | BASE JFUSION1620 |         |
| 5    | BJFUSION1621 | BASE JFUSION1621 |         |
| 5    | BJFUSION1622 | BASE JFUSION1622 |         |
| 5    | BJFUSION1714 | BASE JFUSION1714 |         |
| 5    | BJFUSION1715 | BASE JFUSION1715 |         |
| 5    | BJFUSION1716 | BASE JFUSION1716 |         |
| 5    | BJFUSION1717 | BASE JFUSION1717 |         |
| 5    | BJFUSION1718 | BASE JFUSION1718 |         |
| 5    | BJFUSION1719 | BASE JFUSION1719 |         |
| 5    | BJFUSION1720 | BASE JFUSION1720 |         |
| 5    | BJFUSION1721 | BASE JFUSION1721 |         |
| 5    | BJFUSION1722 | BASE JFUSION1722 |         |
| 5    | BJFUSION1723 | BASE JFUSION1723 |         |
| 5    | BJFUSION1814 | BASE JFUSION1814 |         |
| 5    | BJFUSION1815 | BASE JFUSION1815 |         |
| 5    | BJFUSION1816 | BASE JFUSION1816 |         |
| 5    | BJFUSION1817 | BASE JFUSION1817 |         |
| 5    | BJFUSION1818 | BASE JFUSION1818 |         |
| 5    | BJFUSION1819 | BASE JFUSION1819 |         |
| 5    | BJFUSION1820 | BASE JFUSION1820 |         |
| 5    | BJFUSION1821 | BASE JFUSION1821 |         |
| 5    | BJFUSION1822 | BASE JFUSION1822 |         |
| 5    | BJFUSION1823 | BASE JFUSION1823 |         |
| 5    | BJFUSION1824 | BASE JFUSION1824 |         |
| 5    | BJFUSION1915 | BASE JFUSION1915 |         |
| 5    | BJFUSION1916 | BASE JFUSION1916 |         |
| 5    | BJFUSION1917 | BASE JFUSION1917 |         |
| 5    | BJFUSION1918 | BASE JFUSION1918 |         |
| 5    | BJFUSION1919 | BASE JFUSION1919 |         |
| 5    | BJFUSION1920 | BASE JFUSION1920 |         |
| 5    | BJFUSION1921 | BASE JFUSION1921 |         |
| 5    | BJFUSION1922 | BASE JFUSION1922 |         |
| 5    | BJFUSION1923 | BASE JFUSION1923 |         |
| 5    | BJFUSION1924 | BASE JFUSION1924 |         |
| 5    | BJFUSION2015 | BASE JFUSION2015 |         |
| 5    | BJFUSION2016 | BASE JFUSION2016 |         |
| 5    | BJFUSION2017 | BASE JFUSION2017 |         |
| 5    | BJFUSION2018 | BASE JFUSION2018 |         |
| 5    | BJFUSION2019 | BASE JFUSION2019 |         |
| 5    | BJFUSION2020 | BASE JFUSION2020 |         |
| 5    | BJFUSION2021 | BASE JFUSION2021 |         |
| 5    | BJFUSION2022 | BASE JFUSION2022 |         |
| 5    | BJFUSION2023 | BASE JFUSION2023 |         |
| 5    | BJFUSION2024 | BASE JFUSION2024 |         |
| 5    | BJFUSION2117 | BASE JFUSION2117 |         |
| 5    | BJFUSION2118 | BASE JFUSION2118 |         |
| 5    | BJFUSION2119 | BASE JFUSION2119 |         |
| 5    | BJFUSION2120 | BASE JFUSION2120 |         |
| 5    | BJFUSION2121 | BASE JFUSION2121 |         |

| Pos. | Item Number  | Description      | Remarks  |
|------|--------------|------------------|----------|
| 5    | BJFUSION2122 | BASE JFUSION2122 |          |
| 5    | BJFUSION2123 | BASE JFUSION2123 |          |
| 5    | BJFUSION2124 | BASE JFUSION2124 |          |
| 5    | BJFUSION2217 | BASE JFUSION2217 |          |
| 5    | BJFUSION2218 | BASE JFUSION2218 |          |
| 5    | BJFUSION2219 | BASE JFUSION2219 |          |
| 5    | BJFUSION2220 | BASE JFUSION2220 |          |
| 5    | BJFUSION2221 | BASE JFUSION2221 |          |
| 5    | BJFUSION2222 | BASE JFUSION2222 |          |
| 5    | BJFUSION2223 | BASE JFUSION2223 |          |
| 5    | BJFUSION2224 | BASE JFUSION2224 |          |
| 5    | BJFUSION2319 | BASE JFUSION2319 |          |
| 5    | BJFUSION2320 | BASE JFUSION2320 |          |
| 5    | BJFUSION2321 | BASE JFUSION2321 |          |
| 5    | BJFUSION2322 | BASE JFUSION2322 |          |
| 5    | BJFUSION2323 | BASE JFUSION2323 |          |
| 5    | BJFUSION2324 | BASE JFUSION2324 |          |
| 5    | BJFUSION2419 | BASE JFUSION2419 |          |
| 5    | BJFUSION2420 | BASE JFUSION2420 |          |
| 5    | BJFUSION2421 | BASE JFUSION2421 |          |
| 5    | BJFUSION2422 | BASE JFUSION2422 |          |
| 5    | BJFUSION2423 | BASE JFUSION2423 |          |
| 5    | BJFUSION2424 | BASE JFUSION2424 |          |
| 3    |              | Part             | 5 pieces |

Note: Cushion widths 14" & 15" use A size Fluid Bladders, 16" - 20" use B size bladders and 21"-24" use C size bladders.




| Pos. | Item Number   | Description                 | Remarks |
|------|---------------|-----------------------------|---------|
| 1    | CJFUSION1414C | COVER OUT FUSION CORE 14X14 |         |
| 1    | CJFUSION1415C | COVER OUT FUSION CORE 14X15 |         |
| 1    | CJFUSION1416C | COVER OUT FUSION CORE 14X16 |         |
| 1    | CJFUSION1417C | COVER OUT FUSION CORE 14X17 |         |
| 1    | CJFUSION1418C | COVER OUT FUSION CORE 14X18 |         |
| 1    | CJFUSION1514C | COVER OUT FUSION CORE 15X14 |         |
| 1    | CJFUSION1515C | COVER OUT FUSION CORE 15X15 |         |
| 1    | CJFUSION1516C | COVER OUT FUSION CORE 15X16 |         |
| 1    | CJFUSION1517C | COVER OUT FUSION CORE 15X17 |         |
| 1    | CJFUSION1518C | COVER OUT FUSION CORE 15X18 |         |
| 1    | CJFUSION1519C | COVER OUT FUSION CORE 15X19 |         |
| 1    | CJFUSION1520C | COVER OUT FUSION CORE 15X20 |         |
| 1    | CJFUSION1521C | COVER OUT FUSION CORE 15X21 |         |
| 1    | CJFUSION1614C | COVER OUT FUSION CORE 16X14 |         |
| 1    | CJFUSION1615C | COVER OUT FUSION CORE 16X15 |         |
| 1    | CJFUSION1616C | COVER OUT FUSION CORE 16X16 |         |
| 1    | CJFUSION1617C | COVER OUT FUSION CORE 16X17 |         |
| 1    | CJFUSION1618C | COVER OUT FUSION CORE 16X18 |         |

| Pos. | Item Number   | Description                 | Remarks |
|------|---------------|-----------------------------|---------|
| 1    | CJFUSION1619C | COVER OUT FUSION CORE 16X19 |         |
| 1    | CJFUSION1620C | COVER OUT FUSION CORE 16X20 |         |
| 1    | CJFUSION1621C | COVER OUT FUSION CORE 16X21 |         |
| 1    | CJFUSION1622C | COVER OUT FUSION CORE 16X22 |         |
| 1    | CJFUSION1714C | COVER OUT FUSION CORE 17X14 |         |
| 1    | CJFUSION1715C | COVER OUT FUSION CORE 17X15 |         |
| 1    | CJFUSION1716C | COVER OUT FUSION CORE 17X16 |         |
| 1    | CJFUSION1717C | COVER OUT FUSION CORE 17X17 |         |
| 1    | CJFUSION1718C | COVER OUT FUSION CORE 17X18 |         |
| 1    | CJFUSION1719C | COVER OUT FUSION CORE 17X19 |         |
| 1    | CJFUSION1720C | COVER OUT FUSION CORE 17X20 |         |
| 1    | CJFUSION1721C | COVER OUT FUSION CORE 17X21 |         |
| 1    | CJFUSION1722C | COVER OUT FUSION CORE 17X22 |         |
| 1    | CJFUSION1723C | COVER OUT FUSION CORE 17X23 |         |
| 1    | CJFUSION1814C | COVER OUT FUSION CORE 18X14 |         |
| 1    | CJFUSION1815C | COVER OUT FUSION CORE 18X15 |         |
| 1    | CJFUSION1816C | COVER OUT FUSION CORE 18X16 |         |
| 1    | CJFUSION1817C | COVER OUT FUSION CORE 18X17 |         |
| 1    | CJFUSION1818C | COVER OUT FUSION CORE 18X18 |         |
| 1    | CJFUSION1819C | COVER OUT FUSION CORE 18X19 |         |
| 1    | CJFUSION1820C | COVER OUT FUSION CORE 18X20 |         |
| 1    | CJFUSION1821C | COVER OUT FUSION CORE 18X21 |         |
| 1    | CJFUSION1822C | COVER OUT FUSION CORE 18X22 |         |
| 1    | CJFUSION1823C | COVER OUT FUSION CORE 18X23 |         |
| 1    | CJFUSION1824C | COVER OUT FUSION CORE 18X24 |         |
| 1    | CJFUSION1915C | COVER OUT FUSION CORE 19X15 |         |
| 1    | CJFUSION1916C | COVER OUT FUSION CORE 19X16 |         |
| 1    | CJFUSION1917C | COVER OUT FUSION CORE 19X17 |         |
| 1    | CJFUSION1918C | COVER OUT FUSION CORE 19X18 |         |
| 1    | CJFUSION1919C | COVER OUT FUSION CORE 19X19 |         |
| 1    | CJFUSION1920C | COVER OUT FUSION CORE 19X20 |         |
| 1    | CJFUSION1921C | COVER OUT FUSION CORE 19X21 |         |
| 1    | CJFUSION1922C | COVER OUT FUSION CORE 19X22 |         |
| 1    | CJFUSION1923C | COVER OUT FUSION CORE 19X23 |         |
| 1    | CJFUSION1924C | COVER OUT FUSION CORE 19X24 |         |
| 1    | CJFUSION2015C | COVER OUT FUSION CORE 20X15 |         |
| 1    | CJFUSION2016C | COVER OUT FUSION CORE 20X16 |         |
| 1    | CJFUSION2017C | COVER OUT FUSION CORE 20X17 |         |
| 1    | CJFUSION2018C | COVER OUT FUSION CORE 20X18 |         |
| 1    | CJFUSION2019C | COVER OUT FUSION CORE 20X19 |         |
| 1    | CJFUSION2020C | COVER OUT FUSION CORE 20X20 |         |
| 1    | CJFUSION2021C | COVER OUT FUSION CORE 20X21 |         |
| 1    | CJFUSION2022C | COVER OUT FUSION CORE 20X22 |         |
| 1    | CJFUSION2023C | COVER OUT FUSION CORE 20X23 |         |
| 1    | CJFUSION2024C | COVER OUT FUSION CORE 20X24 |         |
| 1    | CJFUSION2117C | COVER OUT FUSION CORE 21X17 |         |
| 1    | CJFUSION2118C | COVER OUT FUSION CORE 21X18 |         |
| 1    | CJFUSION2119C | COVER OUT FUSION CORE 21X19 |         |
| 1    | CJFUSION2120C | COVER OUT FUSION CORE 21X20 |         |
| 1    | CJFUSION2121C | COVER OUT FUSION CORE 21X21 |         |
| 1    | CJFUSION2122C | COVER OUT FUSION CORE 21X22 |         |
| 1    | CJFUSION2123C | COVER OUT FUSION CORE 21X23 |         |
| 1    | CJFUSION2124C | COVER OUT FUSION CORE 21X24 |         |
| 1    | CJFUSION2217C | COVER OUT FUSION CORE 22X17 |         |
| 1    | CJFUSION2218C | COVER OUT FUSION CORE 22X18 |         |
| 1    | CJFUSION2219C | COVER OUT FUSION CORE 22X19 |         |
| 1    | CJFUSION2220C | COVER OUT FUSION CORE 22X20 |         |

| Pos. | Item Number   | Description                    | Remarks             |
|------|---------------|--------------------------------|---------------------|
| 1    | CJFUSION2221C | COVER OUT FUSION CORE 22X21    |                     |
| 1    | CJFUSION2222C | COVER OUT FUSION CORE 22X22    |                     |
| 1    | CJFUSION2223C | COVER OUT FUSION CORE 22X23    |                     |
| 1    | CJFUSION2224C | COVER OUT FUSION CORE 22X24    |                     |
| 1    | CJFUSION2319C | COVER OUT FUSION CORE 23X19    |                     |
| 1    | CJFUSION2320C | COVER OUT FUSION CORE 23X20    |                     |
| 1    | CJFUSION2321C | COVER OUT FUSION CORE 23X21    |                     |
| 1    | CJFUSION2322C | COVER OUT FUSION CORE 23X22    |                     |
| 1    | CJFUSION2323C | COVER OUT FUSION CORE 23X23    |                     |
| 1    | CJFUSION2324C | COVER OUT FUSION CORE 23X24    |                     |
| 1    | CJFUSION2419C | COVER OUT FUSION CORE 24X19    |                     |
| 1    | CJFUSION2420C | COVER OUT FUSION CORE 24X20    |                     |
| 1    | CJFUSION2421C | COVER OUT FUSION CORE 24X21    |                     |
| 1    | CJFUSION2422C | COVER OUT FUSION CORE 24X22    |                     |
| 1    | CJFUSION2423C | COVER OUT FUSION CORE 24X23    |                     |
| 1    | CJFUSION2424C | COVER OUT FUSION CORE 24X24    |                     |
| 2    | FJX2-1        | FL PD JAY X2 14X14 14X16       | all 14 & 15 wide    |
| 2    | FJX2-5        | FL PD JAY X2 17X17 18X18 18X20 | all 16 thru 21 wide |
| 2    | FJX2-6        | FL PD JAY X2 20X20             | all 22 thru 24 wide |
| 3    | BJFUSION1414  | BASE JFUSION1414               |                     |
| 3    | BJFUSION1415  | BASE JFUSION1415               |                     |
| 3    | BJFUSION1416  | BASE JFUSION1416               |                     |
| 3    | BJFUSION1417  | BASE JFUSION1417               |                     |
| 3    | BJFUSION1418  | BASE JFUSION1418               |                     |
| 3    | BJFUSION1514  | BASE JFUSION1514               |                     |
| 3    | BJFUSION1515  | BASE JFUSION1515               |                     |
| 3    | BJFUSION1516  | BASE JFUSION1516               |                     |
| 3    | BJFUSION1517  | BASE JFUSION1517               |                     |
| 3    | BJFUSION1518  | BASE JFUSION1518               |                     |
| 3    | BJFUSION1519  | BASE JFUSION1519               |                     |
| 3    | BJFUSION1520  | BASE JFUSION1520               |                     |
| 3    | BJFUSION1521  | BASE JFUSION1521               |                     |
| 3    | BJFUSION1614  | BASE JFUSION1614               |                     |
| 3    | BJFUSION1615  | BASE JFUSION1615               |                     |
| 3    | BJFUSION1616  | BASE JFUSION1616               |                     |
| 3    | BJFUSION1617  | BASE JFUSION1617               |                     |
| 3    | BJFUSION1618  | BASE JFUSION1618               |                     |
| 3    | BJFUSION1619  | BASE JFUSION1619               |                     |
| 3    | BJFUSION1620  | BASE JFUSION1620               |                     |
| 3    | BJFUSION1621  | BASE JFUSION1621               |                     |
| 3    | BJFUSION1622  | BASE JFUSION1622               |                     |
| 3    | BJFUSION1714  | BASE JFUSION1714               |                     |
| 3    | BJFUSION1715  | BASE JFUSION1715               |                     |
| 3    | BJFUSION1716  | BASE JFUSION1716               |                     |
| 3    | BJFUSION1717  | BASE JFUSION1717               |                     |
| 3    | BJFUSION1718  | BASE JFUSION1718               |                     |
| 3    | BJFUSION1719  | BASE JFUSION1719               |                     |
| 3    | BJFUSION1720  | BASE JFUSION1720               |                     |
| 3    | BJFUSION1721  | BASE JFUSION1721               |                     |
| 3    | BJFUSION1722  | BASE JFUSION1722               |                     |
| 3    | BJFUSION1723  | BASE JFUSION1723               |                     |
| 3    | BJFUSION1814  | BASE JFUSION1814               |                     |
| 3    | BJFUSION1815  | BASE JFUSION1815               |                     |
| 3    | BJFUSION1816  | BASE JFUSION1816               |                     |
| 3    | BJFUSION1817  | BASE JFUSION1817               |                     |
| 3    | BJFUSION1818  | BASE JFUSION1818               |                     |
| 3    | BJFUSION1819  | BASE JFUSION1819               |                     |

| Pos. | Item Number  | Description      | Remarks |
|------|--------------|------------------|---------|
| 3    | BJFUSION1820 | BASE JFUSION1820 |         |
| 3    | BJFUSION1821 | BASE JFUSION1821 |         |
| 3    | BJFUSION1822 | BASE JFUSION1822 |         |
| 3    | BJFUSION1823 | BASE JFUSION1823 |         |
| 3    | BJFUSION1824 | BASE JFUSION1824 |         |
| 3    | BJFUSION1915 | BASE JFUSION1915 |         |
| 3    | BJFUSION1916 | BASE JFUSION1916 |         |
| 3    | BJFUSION1917 | BASE JFUSION1917 |         |
| 3    | BJFUSION1918 | BASE JFUSION1918 |         |
| 3    | BJFUSION1919 | BASE JFUSION1919 |         |
| 3    | BJFUSION1920 | BASE JFUSION1920 |         |
| 3    | BJFUSION1921 | BASE JFUSION1921 |         |
| 3    | BJFUSION1922 | BASE JFUSION1922 |         |
| 3    | BJFUSION1923 | BASE JFUSION1923 |         |
| 3    | BJFUSION1924 | BASE JFUSION1924 |         |
| 3    | BJFUSION2015 | BASE JFUSION2015 |         |
| 3    | BJFUSION2016 | BASE JFUSION2016 |         |
| 3    | BJFUSION2017 | BASE JFUSION2017 |         |
| 3    | BJFUSION2018 | BASE JFUSION2018 |         |
| 3    | BJFUSION2019 | BASE JFUSION2019 |         |
| 3    | BJFUSION2020 | BASE JFUSION2020 |         |
| 3    | BJFUSION2021 | BASE JFUSION2021 |         |
| 3    | BJFUSION2022 | BASE JFUSION2022 |         |
| 3    | BJFUSION2023 | BASE JFUSION2023 |         |
| 3    | BJFUSION2024 | BASE JFUSION2024 |         |
| 3    | BJFUSION2117 | BASE JFUSION2117 |         |
| 3    | BJFUSION2118 | BASE JFUSION2118 |         |
| 3    | BJFUSION2119 | BASE JFUSION2119 |         |
| 3    | BJFUSION2120 | BASE JFUSION2120 |         |
| 3    | BJFUSION2121 | BASE JFUSION2121 |         |
| 3    | BJFUSION2122 | BASE JFUSION2122 |         |
| 3    | BJFUSION2123 | BASE JFUSION2123 |         |
| 3    | BJFUSION2124 | BASE JFUSION2124 |         |
| 3    | BJFUSION2217 | BASE JFUSION2217 |         |
| 3    | BJFUSION2218 | BASE JFUSION2218 |         |
| 3    | BJFUSION2219 | BASE JFUSION2219 |         |
| 3    | BJFUSION2220 | BASE JFUSION2220 |         |
| 3    | BJFUSION2221 | BASE JFUSION2221 |         |
| 3    | BJFUSION2222 | BASE JFUSION2222 |         |
| 3    | BJFUSION2223 | BASE JFUSION2223 |         |
| 3    | BJFUSION2224 | BASE JFUSION2224 |         |
| 3    | BJFUSION2319 | BASE JFUSION2319 |         |
| 3    | BJFUSION2320 | BASE JFUSION2320 |         |
| 3    | BJFUSION2321 | BASE JFUSION2321 |         |
| 3    | BJFUSION2322 | BASE JFUSION2322 |         |
| 3    | BJFUSION2323 | BASE JFUSION2323 |         |
| 3    | BJFUSION2324 | BASE JFUSION2324 |         |
| 3    | BJFUSION2419 | BASE JFUSION2419 |         |
| 3    | BJFUSION2420 | BASE JFUSION2420 |         |
| 3    | BJFUSION2421 | BASE JFUSION2421 |         |
| 3    | BJFUSION2422 | BASE JFUSION2422 |         |
| 3    | BJFUSION2423 | BASE JFUSION2423 |         |
| 3    | BJFUSION2424 | BASE JFUSION2424 |         |




| Pos. | Item Number | Description        | Pos. | Item Number | Description        |
|------|-------------|--------------------|------|-------------|--------------------|
| 1    | CJGO1414    | COVER JAY GO 14X14 | 1    | CJGO1620    | COVER JAY GO 16X20 |
| 1    | CJGO1415    | COVER JAY GO 14X15 | 1    | CJGO1621    | COVER JAY GO 16X21 |
| 1    | CJGO1416    | COVER JAY GO 14X16 | 1    | CJGO1622    | COVER JAY GO 16X22 |
| 1    | CJGO1417    | COVER JAY GO 14X17 | 1    | CJGO1714    | COVER JAY GO 17X14 |
| 1    | CJGO1418    | COVER JAY GO 14X18 | 1    | CJGO1715    | COVER JAY GO 17X15 |
| 1    | CJGO1419    | COVER JAY GO 14X19 | 1    | CJGO1716    | COVER JAY GO 17X16 |
| 1    | CJGO1420    | COVER JAY GO 14X20 | 1    | CJGO1717    | COVER JAY GO 17X17 |
| 1    | CJGO1514    | COVER JAY GO 15X14 | 1    | CJGO1718    | COVER JAY GO 17X18 |
| 1    | CJGO1515    | COVER JAY GO 15X15 | 1    | CJGO1719    | COVER JAY GO 17X19 |
| 1    | CJGO1516    | COVER JAY GO 15X16 | 1    | CJGO1720    | COVER JAY GO 17X20 |
| 1    | CJGO1517    | COVER JAY GO 15X17 | 1    | CJGO1721    | COVER JAY GO 17X21 |
| 1    | CJGO1518    | COVER JAY GO 15X18 | 1    | CJGO1722    | COVER JAY GO 17X22 |
| 1    | CJGO1519    | COVER JAY GO 15X19 | 1    | CJGO1723    | COVER JAY GO 17X23 |
| 1    | CJGO1520    | COVER JAY GO 15X20 | 1    | CJGO1814    | COVER JAY GO 18X14 |
| 1    | CJGO1521    | COVER JAY GO 15X21 | 1    | CJGO1815    | COVER JAY GO 18X15 |
| 1    | CJGO1614    | COVER JAY GO 16X14 | 1    | CJGO1816    | COVER JAY GO 18X16 |
| 1    | CJGO1615    | COVER JAY GO 16X15 | 1    | CJGO1817    | COVER JAY GO 18X17 |
| 1    | CJGO1616    | COVER JAY GO 16X16 | 1    | CJGO1818    | COVER JAY GO 18X18 |
| 1    | CJGO1617    | COVER JAY GO 16X17 | 1    | CJGO1819    | COVER JAY GO 18X19 |
| 1    | CJGO1618    | COVER JAY GO 16X18 | 1    | CJGO1820    | COVER JAY GO 18X20 |
| 1    | CJGO1619    | COVER JAY GO 16X19 | 1    | CJGO1821    | COVER JAY GO 18X21 |

| Pos. | Item Number | Description            | Pos. | Item Number | Description            |
|------|-------------|------------------------|------|-------------|------------------------|
| 1    | CJGO1822    | COVER JAY GO 18X22     | 2    | BJGO1420    | BASE FOAM JAY GO 14X20 |
| 1    | CJGO1823    | COVER JAY GO 18X23     | 2    | BJGO1514    | BASE FOAM JAY GO 15X14 |
| 1    | CJGO1824    | COVER JAY GO 18X24     | 2    | BJGO1515    | BASE FOAM JAY GO 15X15 |
| 1    | CJGO1915    | COVER JAY GO 19X15     | 2    | BJGO1516    | BASE FOAM JAY GO 15X16 |
| 1    | CJGO1916    | COVER JAY GO 19X16     | 2    | BJGO1517    | BASE FOAM JAY GO 15X17 |
| 1    | CJGO1917    | COVER JAY GO 19X17     | 2    | BJGO1518    | BASE FOAM JAY GO 15X18 |
| 1    | CJGO1918    | COVER JAY GO 19X18     | 2    | BJGO1519    | BASE FOAM JAY GO 15X19 |
| 1    | CJGO1919    | COVER JAY GO 19X19     | 2    | BJGO1520    | BASE FOAM JAY GO 15X20 |
| 1    | CJGO1920    | COVER JAY GO 19X20     | 2    | BJGO1521    | BASE FOAM JAY GO 15X21 |
| 1    | CJGO1921    | COVER JAY GO 19X21     | 2    | BJGO1614    | BASE FOAM JAY GO 16X14 |
| 1    | CJGO1922    | COVER JAY GO 19X22     | 2    | BJGO1615    | BASE FOAM JAY GO 16X15 |
| 1    | CJGO1923    | COVER JAY GO 19X23     | 2    | BJGO1616    | BASE FOAM JAY GO 16X16 |
| 1    | CJGO1924    | COVER JAY GO 19X24     | 2    | BJGO1617    | BASE FOAM JAY GO 16X17 |
| 1    | CJGO2015    | COVER JAY GO 20X15     | 2    | BJGO1618    | BASE FOAM JAY GO 16X18 |
| 1    | CJGO2016    | COVER JAY GO 20X16     | 2    | BJGO1619    | BASE FOAM JAY GO 16X19 |
| 1    | CJGO2017    | COVER JAY GO 20X17     | 2    | BJGO1620    | BASE FOAM JAY GO 16X20 |
| 1    | CJGO2018    | COVER JAY GO 20X18     | 2    | BJGO1621    | BASE FOAM JAY GO 16X21 |
| 1    | CJGO2019    | COVER JAY GO 20X19     | 2    | BJGO1622    | BASE FOAM JAY GO 16X22 |
| 1    | CJGO2020    | COVER JAY GO 20X20     | 2    | BJGO1714    | BASE FOAM JAY GO 17X14 |
| 1    | CJGO2021    | COVER JAY GO 20X21     | 2    | BJGO1715    | BASE FOAM JAY GO 17X15 |
| 1    | CJGO2022    | COVER JAY GO 20X22     | 2    | BJGO1716    | BASE FOAM JAY GO 17X16 |
| 1    | CJGO2023    | COVER JAY GO 20X23     | 2    | BJGO1717    | BASE FOAM JAY GO 17X17 |
| 1    | CJGO2024    | COVER JAY GO 20X24     | 2    | BJGO1718    | BASE FOAM JAY GO 17X18 |
| 1    | CJGO2117    | COVER JAY GO 21X17     | 2    | BJGO1719    | BASE FOAM JAY GO 17X19 |
| 1    | CJGO2118    | COVER JAY GO 21X18     | 2    | BJGO1720    | BASE FOAM JAY GO 17X20 |
| 1    | CJGO2119    | COVER JAY GO 21X19     | 2    | BJGO1721    | BASE FOAM JAY GO 17X21 |
| 1    | CJGO2120    | COVER JAY GO 21X20     | 2    | BJGO1722    | BASE FOAM JAY GO 17X22 |
| 1    | CJGO2121    | COVER JAY GO 21X21     | 2    | BJGO1723    | BASE FOAM JAY GO 17X23 |
| 1    | CJGO2122    | COVER JAY GO 21X22     | 2    | BJGO1814    | BASE FOAM JAY GO 18X14 |
| 1    | CJGO2123    | COVER JAY GO 21X23     | 2    | BJGO1815    | BASE FOAM JAY GO 18X15 |
| 1    | CJGO2124    | COVER JAY GO 21X24     | 2    | BJGO1816    | BASE FOAM JAY GO 18X16 |
| 1    | CJGO2217    | COVER JAY GO 22X17     | 2    | BJGO1817    | BASE FOAM JAY GO 18X17 |
| 1    | CJGO2218    | COVER JAY GO 22X18     | 2    | BJGO1818    | BASE FOAM JAY GO 18X18 |
| 1    | CJGO2219    | COVER JAY GO 22X19     | 2    | BJGO1819    | BASE FOAM JAY GO 18X19 |
| 1    | CJGO2220    | COVER JAY GO 22X20     | 2    | BJGO1820    | BASE FOAM JAY GO 18X20 |
| 1    | CJGO2221    | COVER JAY GO 22X21     | 2    | BJGO1821    | BASE FOAM JAY GO 18X21 |
| 1    | CJGO2222    | COVER JAY GO 22X22     | 2    | BJGO1822    | BASE FOAM JAY GO 18X22 |
| 1    | CJGO2223    | COVER JAY GO 22X23     | 2    | BJGO1823    | BASE FOAM JAY GO 18X23 |
| 1    | CJGO2224    | COVER JAY GO 22X24     | 2    | BJGO1824    | BASE FOAM JAY GO 18X24 |
| 1    | CJGO2319    | COVER JAY GO 23X19     | 2    | BJGO1915    | BASE FOAM JAY GO 19X15 |
| 1    | CJGO2320    | COVER JAY GO 23X20     | 2    | BJGO1916    | BASE FOAM JAY GO 19X16 |
| 1    | CJGO2321    | COVER JAY GO 23X21     | 2    | BJGO1917    | BASE FOAM JAY GO 19X17 |
| 1    | CJGO2322    | COVER JAY GO 23X22     | 2    | BJGO1918    | BASE FOAM JAY GO 19X18 |
| 1    | CJGO2323    | COVER JAY GO 23X23     | 2    | BJGO1919    | BASE FOAM JAY GO 19X19 |
| 1    | CJGO2324    | COVER JAY GO 23X24     | 2    | BJGO1920    | BASE FOAM JAY GO 19X20 |
| 1    | CJGO2419    | COVER JAY GO 24X19     | 2    | BJGO1921    | BASE FOAM JAY GO 19X21 |
| 1    | CJGO2420    | COVER JAY GO 24X20     | 2    | BJGO1922    | BASE FOAM JAY GO 19X22 |
| 1    | CJGO2421    | COVER JAY GO 24X21     | 2    | BJGO1923    | BASE FOAM JAY GO 19X23 |
| 1    | CJGO2422    | COVER JAY GO 24X22     | 2    | BJGO1924    | BASE FOAM JAY GO 19X24 |
| 1    | CJGO2423    | COVER JAY GO 24X23     | 2    | BJGO2015    | BASE FOAM JAY GO 20X15 |
| 1    | CJGO2424    | COVER JAY GO 24X24     | 2    | BJGO2016    | BASE FOAM JAY GO 20X16 |
| 2    | BJGO1414    | BASE FOAM JAY GO 14X14 | 2    | BJGO2017    | BASE FOAM JAY GO 20X17 |
| 2    | BJGO1415    | BASE FOAM JAY GO 14X15 | 2    | BJGO2018    | BASE FOAM JAY GO 20X18 |
| 2    | BJGO1416    | BASE FOAM JAY GO 14X16 | 2    | BJGO2019    | BASE FOAM JAY GO 20X19 |
| 2    | BJGO1417    | BASE FOAM JAY GO 14X17 | 2    | BJGO2020    | BASE FOAM JAY GO 20X20 |
| 2    | BJGO1418    | BASE FOAM JAY GO 14X18 | 2    | BJGO2021    | BASE FOAM JAY GO 20X21 |
| 2    | BJGO1419    | BASE FOAM JAY GO 14X19 | 2    | BJGO2022    | BASE FOAM JAY GO 20X22 |

| Pos. | Item Number | Description            | Pos. | Item Number | Description            |
|------|-------------|------------------------|------|-------------|------------------------|
| 2    | BJGO2023    | BASE FOAM JAY GO 20X23 | 2    | BJGO2222    | BASE FOAM JAY GO 22X22 |
| 2    | BJGO2024    | BASE FOAM JAY GO 20X24 | 2    | BJGO2223    | BASE FOAM JAY GO 22X23 |
| 2    | BJGO2117    | BASE FOAM JAY GO 21X17 | 2    | BJGO2224    | BASE FOAM JAY GO 22X24 |
| 2    | BJGO2118    | BASE FOAM JAY GO 21X18 | 2    | BJGO2319    | BASE FOAM JAY GO 23X19 |
| 2    | BJGO2119    | BASE FOAM JAY GO 21X19 | 2    | BJGO2320    | BASE FOAM JAY GO 23X20 |
| 2    | BJGO2120    | BASE FOAM JAY GO 21X20 | 2    | BJGO2321    | BASE FOAM JAY GO 23X21 |
| 2    | BJGO2121    | BASE FOAM JAY GO 21X21 | 2    | BJGO2322    | BASE FOAM JAY GO 23X22 |
| 2    | BJGO2122    | BASE FOAM JAY GO 21X22 | 2    | BJGO2323    | BASE FOAM JAY GO 23X23 |
| 2    | BJGO2123    | BASE FOAM JAY GO 21X23 | 2    | BJGO2324    | BASE FOAM JAY GO 23X24 |
| 2    | BJGO2124    | BASE FOAM JAY GO 21X24 | 2    | BJGO2419    | BASE FOAM JAY GO 24X19 |
| 2    | BJGO2217    | BASE FOAM JAY GO 22X17 | 2    | BJGO2420    | BASE FOAM JAY GO 24X20 |
| 2    | BJGO2218    | BASE FOAM JAY GO 22X18 | 2    | BJGO2421    | BASE FOAM JAY GO 24X21 |
| 2    | BJGO2219    | BASE FOAM JAY GO 22X19 | 2    | BJGO2422    | BASE FOAM JAY GO 24X22 |
| 2    | BJGO2220    | BASE FOAM JAY GO 22X20 | 2    | BJGO2423    | BASE FOAM JAY GO 24X23 |
| 2    | BJGO2221    | BASE FOAM JAY GO 22X21 | 2    | BJGO2424    | BASE FOAM JAY GO 24X24 |


Note: For cushions purchased after 3/7/16 adapter Velcro will be shipped to allow the replacement parts to attach to existing base, fluid bladder, or positioning pieces.

Note: For cushions purchased after 3/7/16 that are needing a replacement cover a new style replacement cover will need to be ordered. If positioning pieces are utilized in the existing cushion a P2 or P3 cushion is required, and if no positioning a non P cover will be adequate.



| Pos. | Item Number | Description                       | Remarks |
|------|-------------|-----------------------------------|---------|
| 1A   | CJGS10      | CVR, JAY GS, 10IN X 10-13IN       |         |
| 1A   | CJGS12      | CVR, JAY GS, 12IN X 12-15IN       |         |
| 1A   | CJGS14      | CVR, JAY GS, 14IN X 14-17IN       |         |
| 1A   | CJGS16      | CVR, JAY GS, 16IN X 16-19IN       |         |
| 1A   | CJGS18      | CVR, JAY GS, 18IN X 17-20IN       |         |
| 1B   | CJGS10P2    | CVR, JAY GS, 2POS, 10IN X 10-13IN |         |
| 1B   | CJGS12P2    | CVR, JAY GS, 2POS, 12IN X 12-15IN |         |
| 1B   | CJGS14P2    | CVR, JAY GS, 2POS, 14IN X 14-17IN |         |
| 1B   | CJGS16P2    | CVR, JAY GS, 2POS, 16IN X 16-19IN |         |
| 1B   | CJGS18P2    | CVR, JAY GS, 2POS, 18IN X 17-20IN |         |

| Pos. | Item Number | Description                                | Remarks |
|------|-------------|--------------------------------------------|---------|
| 1C   | CJGS14P3    | CVR,JAY GS,3POS,14INx14-17IN               |         |
| 1C   | CJGS16P3    | CVR,JAY GS,3POS,16INx16-19IN               |         |
| 1C   | CJGS18P3    | CVR,JAY GS,3POS,18INx17-20IN               |         |
| 2A   | F1100       | FLUID PAD JAY GS JUNIOR 18 X 20            |         |
| 2A   | F1105       | FLUID PAD JAY GS JUNIOR 16 X 19            |         |
| 2A   | F1110       | FLUID PAD JAY GS JUNIOR 14 X 17            |         |
| 2A   | F1111       | FLUID PAD JAY GS JUNIOR 10 X 13            |         |
| 2A   | F1115       | FLUID PAD JAY GS JUNIOR 12 X 15            |         |
| 2B   | F1100.1     | PRESSURE RELIEVING PAD 18 X 20             |         |
| 2B   | F1105.1     | PRESSURE RELIEVING PAD 16 X 19             |         |
| 2B   | F1110.1     | PRESSURE RELIEVING PAD 14 X 17             |         |
| 2B   | F1111.1     | PRESSURE RELIEVING PAD 10 X 13             |         |
| 2B   | F1115.1     | PRESSURE RELIEVING PAD 12 X 15             |         |
| 3    | B1171       | PRE ISCHIAL SHELF EXTENSION GS (10" x 13") |         |
| 3    | B1173       | PRE ISCHIAL SHELF EXTENSION GS (12" x 15") |         |
| 3    | B1174       | PRE ISCHIAL SHELF EXTENSION GS (14" x 17") |         |
| 3    | B1175       | PRE ISCHIAL SHELF EXTENSION GS (16" x 19") |         |
| 3    | B1176       | PRE ISCHIAL SHELF EXTENSION GS (18" x 20") |         |
| 4    | B1100       | JAY GS BASE 18"                            |         |
| 4    | B1105       | JAY GS BASE 16"                            |         |
| 4    | B1110       | JAY GS BASE JUNIOR 14"                     |         |
| 4    | B1111       | JAY GS BASE SMALL CHILD 10"                |         |
| 4    | B1115       | JAY GS BASE CHILD 12"                      |         |
| 5A   | B163        | 2" JUNIOR ADDUCTOR WEDGES                  | Pair*   |
| 5B   | B1113       | WEDGES,3" ADDUCTOR                         | Pair*   |
| 6A   | B1118       | 3" ABDUCTOR BUILD-UP                       | *       |
| 6B   | B168        | BUILD-UP 2" JUNIOR ABDUCTOR                | *       |
| 7    | B1167       | HIP GUIDES,2",JAY & GS                     | Pair*   |



| Pos. | Item Number | Description           | Pos. | Item Number | Description           |
|------|-------------|-----------------------|------|-------------|-----------------------|
| 1    | CJI1414     | COVER OUTER ION 14X14 | 1    | CJI1620     | COVER OUTER ION 16X20 |
| 1    | CJI1415     | COVER OUTER ION 14X15 | 1    | CJI1621     | COVER OUTER ION 16X21 |
| 1    | CJI1416     | COVER OUTER ION 14X16 | 1    | CJI1622     | COVER OUTER ION 16X22 |
| 1    | CJI1417     | COVER OUTER ION 14X17 | 1    | CJI1714     | COVER OUTER ION 17X14 |
| 1    | CJI1418     | COVER OUTER ION 14X18 | 1    | CJI1715     | COVER OUTER ION 17X15 |
| 1    | CJI1514     | COVER OUTER ION 15X14 | 1    | CJI1716     | COVER OUTER ION 17X16 |
| 1    | CJI1515     | COVER OUTER ION 15X15 | 1    | CJI1717     | COVER OUTER ION 17X17 |
| 1    | CJI1516     | COVER OUTER ION 15X16 | 1    | CJI1718     | COVER OUTER ION 17X18 |
| 1    | CJI1517     | COVER OUTER ION 15X17 | 1    | CJI1719     | COVER OUTER ION 17X19 |
| 1    | CJI1518     | COVER OUTER ION 15X18 | 1    | CJI1720     | COVER OUTER ION 17X20 |
| 1    | CJI1519     | COVER OUTER ION 15X19 | 1    | CJI1721     | COVER OUTER ION 17X21 |
| 1    | CJI1520     | COVER OUTER ION 15X20 | 1    | CJI1722     | COVER OUTER ION 17X22 |
| 1    | CJI1521     | COVER OUTER ION 15X21 | 1    | CJI1723     | COVER OUTER ION 17X23 |
| 1    | CJI1614     | COVER OUTER ION 16X14 | 1    | CJI1814     | COVER OUTER ION 18X14 |
| 1    | CJI1615     | COVER OUTER ION 16X15 | 1    | CJI1815     | COVER OUTER ION 18X15 |
| 1    | CJI1616     | COVER OUTER ION 16X16 | 1    | CJI1816     | COVER OUTER ION 18X16 |
| 1    | CJI1617     | COVER OUTER ION 16X17 | 1    | CJI1817     | COVER OUTER ION 18X17 |
| 1    | CJI1618     | COVER OUTER ION 16X18 | 1    | CJI1818     | COVER OUTER ION 18X18 |
| 1    | CJI1619     | COVER OUTER ION 16X19 | 1    | CJI1819     | COVER OUTER ION 18X19 |

| Pos. | Item Number | Description           | Pos. | Item Number | Description           |
|------|-------------|-----------------------|------|-------------|-----------------------|
| 1    | CJI1820     | COVER OUTER ION 18X20 | 2    | ICJI1418    | COVER INNER ION 14X18 |
| 1    | CJI1821     | COVER OUTER ION 18X21 | 2    | ICJI1514    | COVER INNER ION 15X14 |
| 1    | CJI1822     | COVER OUTER ION 18X22 | 2    | ICJI1515    | COVER INNER ION 15X15 |
| 1    | CJI1823     | COVER OUTER ION 18X23 | 2    | ICJI1516    | COVER INNER ION 15X16 |
| 1    | CJI1824     | COVER OUTER ION 18X24 | 2    | ICJI1517    | COVER INNER ION 15X17 |
| 1    | CJI1915     | COVER OUTER ION 19X15 | 2    | ICJI1518    | COVER INNER ION 15X18 |
| 1    | CJI1916     | COVER OUTER ION 19X16 | 2    | ICJI1519    | COVER INNER ION 15X19 |
| 1    | CJI1917     | COVER OUTER ION 19X17 | 2    | ICJI1520    | COVER INNER ION 15X20 |
| 1    | CJI1918     | COVER OUTER ION 19X18 | 2    | ICJI1521    | COVER INNER ION 15X21 |
| 1    | CJI1919     | COVER OUTER ION 19X19 | 2    | ICJI1614    | COVER INNER ION 16X14 |
| 1    | CJI1920     | COVER OUTER ION 19X20 | 2    | ICJI1615    | COVER INNER ION 16X15 |
| 1    | CJI1921     | COVER OUTER ION 19X21 | 2    | ICJI1616    | COVER INNER ION 16X16 |
| 1    | CJI1922     | COVER OUTER ION 19X22 | 2    | ICJI1617    | COVER INNER ION 16X17 |
| 1    | CJI1923     | COVER OUTER ION 19X23 | 2    | ICJI1618    | COVER INNER ION 16X18 |
| 1    | CJI1924     | COVER OUTER ION 19X24 | 2    | ICJI1619    | COVER INNER ION 16X19 |
| 1    | CJI2015     | COVER OUTER ION 20X15 | 2    | ICJI1620    | COVER INNER ION 16X20 |
| 1    | CJI2016     | COVER OUTER ION 20X16 | 2    | ICJI1621    | COVER INNER ION 16X21 |
| 1    | CJI2017     | COVER OUTER ION 20X17 | 2    | ICJI1622    | COVER INNER ION 16X22 |
| 1    | CJI2018     | COVER OUTER ION 20X18 | 2    | ICJI1714    | COVER INNER ION 17X14 |
| 1    | CJI2019     | COVER OUTER ION 20X19 | 2    | ICJI1715    | COVER INNER ION 17X15 |
| 1    | CJI2020     | COVER OUTER ION 20X20 | 2    | ICJI1716    | COVER INNER ION 17X16 |
| 1    | CJI2021     | COVER OUTER ION 20X21 | 2    | ICJI1717    | COVER INNER ION 17X17 |
| 1    | CJI2022     | COVER OUTER ION 20X22 | 2    | ICJI1718    | COVER INNER ION 17X18 |
| 1    | CJI2023     | COVER OUTER ION 20X23 | 2    | ICJI1719    | COVER INNER ION 17X19 |
| 1    | CJI2024     | COVER OUTER ION 20X24 | 2    | ICJI1720    | COVER INNER ION 17X20 |
| 1    | CJI2117     | COVER OUTER ION 21X17 | 2    | ICJI1721    | COVER INNER ION 17X21 |
| 1    | CJI2118     | COVER OUTER ION 21X18 | 2    | ICJI1722    | COVER INNER ION 17X22 |
| 1    | CJI2119     | COVER OUTER ION 21X19 | 2    | ICJI1723    | COVER INNER ION 17X23 |
| 1    | CJI2120     | COVER OUTER ION 21X20 | 2    | ICJI1814    | COVER INNER ION 18X14 |
| 1    | CJI2121     | COVER OUTER ION 21X21 | 2    | ICJI1815    | COVER INNER ION 18X15 |
| 1    | CJI2122     | COVER OUTER ION 21X22 | 2    | ICJI1816    | COVER INNER ION 18X16 |
| 1    | CJI2123     | COVER OUTER ION 21X23 | 2    | ICJI1817    | COVER INNER ION 18X17 |
| 1    | CJI2124     | COVER OUTER ION 21X24 | 2    | ICJI1818    | COVER INNER ION 18X18 |
| 1    | CJI2217     | COVER OUTER ION 22X17 | 2    | ICJI1819    | COVER INNER ION 18X19 |
| 1    | CJI2218     | COVER OUTER ION 22X18 | 2    | ICJI1820    | COVER INNER ION 18X20 |
| 1    | CJI2219     | COVER OUTER ION 22X19 | 2    | ICJI1821    | COVER INNER ION 18X21 |
| 1    | CJI2220     | COVER OUTER ION 22X20 | 2    | ICJI1822    | COVER INNER ION 18X22 |
| 1    | CJI2221     | COVER OUTER ION 22X21 | 2    | ICJI1823    | COVER INNER ION 18X23 |
| 1    | CJI2222     | COVER OUTER ION 22X22 | 2    | ICJI1824    | COVER INNER ION 18X24 |
| 1    | CJI2223     | COVER OUTER ION 22X23 | 2    | ICJI1915    | COVER INNER ION 19X15 |
| 1    | CJI2224     | COVER OUTER ION 22X24 | 2    | ICJI1916    | COVER INNER ION 19X16 |
| 1    | CJI2319     | COVER OUTER ION 23X19 | 2    | ICJI1917    | COVER INNER ION 19X17 |
| 1    | CJI2320     | COVER OUTER ION 23X20 | 2    | ICJI1918    | COVER INNER ION 19X18 |
| 1    | CJI2321     | COVER OUTER ION 23X21 | 2    | ICJI1919    | COVER INNER ION 19X19 |
| 1    | CJI2322     | COVER OUTER ION 23X22 | 2    | ICJI1920    | COVER INNER ION 19X20 |
| 1    | CJI2323     | COVER OUTER ION 23X23 | 2    | ICJI1921    | COVER INNER ION 19X21 |
| 1    | CJI2324     | COVER OUTER ION 23X24 | 2    | ICJI1922    | COVER INNER ION 19X22 |
| 1    | CJI2419     | COVER OUTER ION 24X19 | 2    | ICJI1923    | COVER INNER ION 19X23 |
| 1    | CJI2420     | COVER OUTER ION 24X20 | 2    | ICJI1924    | COVER INNER ION 19X24 |
| 1    | CJI2421     | COVER OUTER ION 24X21 | 2    | ICJI2015    | COVER INNER ION 20X15 |
| 1    | CJI2422     | COVER OUTER ION 24X22 | 2    | ICJI2016    | COVER INNER ION 20X16 |
| 1    | CJI2423     | COVER OUTER ION 24X23 | 2    | ICJI2017    | COVER INNER ION 20X17 |
| 1    | CJI2424     | COVER OUTER ION 24X24 | 2    | ICJI2018    | COVER INNER ION 20X18 |
| 2    | ICJI1414    | COVER INNER ION 14X14 | 2    | ICJI2019    | COVER INNER ION 20X19 |
| 2    | ICJI1415    | COVER INNER ION 14X15 | 2    | ICJI2020    | COVER INNER ION 20X20 |
| 2    | ICJI1416    | COVER INNER ION 14X16 | 2    | ICJI2021    | COVER INNER ION 20X21 |
| 2    | ICJI1417    | COVER INNER ION 14X17 | 2    | ICJI2022    | COVER INNER ION 20X22 |

| Pos. | Item Number | Description             | Pos. | Item Number | Description             |
|------|-------------|-------------------------|------|-------------|-------------------------|
| 2    | ICJI2023    | COVER INNER ION 20X23   | 3    | BJI1719     | BASE FOAM JAY ION 17X19 |
| 2    | ICJI2024    | COVER INNER ION 20X24   | 3    | BJI1720     | BASE FOAM JAY ION 17X20 |
| 2    | ICJI2117    | COVER INNER ION 21X17   | 3    | BJI1721     | BASE FOAM JAY ION 17X21 |
| 2    | ICJI2118    | COVER INNER ION 21X18   | 3    | BJI1722     | BASE FOAM JAY ION 17X22 |
| 2    | ICJI2119    | COVER INNER ION 21X19   | 3    | BJI1723     | BASE FOAM JAY ION 17X23 |
| 2    | ICJI2120    | COVER INNER ION 21X20   | 3    | BJI1814     | BASE FOAM JAY ION 18X14 |
| 2    | ICJI2121    | COVER INNER ION 21X21   | 3    | BJI1815     | BASE FOAM JAY ION 18X15 |
| 2    | ICJI2122    | COVER INNER ION 21X22   | 3    | BJI1816     | BASE FOAM JAY ION 18X16 |
| 2    | ICJI2123    | COVER INNER ION 21X23   | 3    | BJI1817     | BASE FOAM JAY ION 18X17 |
| 2    | ICJI2124    | COVER INNER ION 21X24   | 3    | BJI1818     | BASE FOAM JAY ION 18X18 |
| 2    | ICJI2217    | COVER INNER ION 22X17   | 3    | BJI1819     | BASE FOAM JAY ION 18X19 |
| 2    | ICJI2218    | COVER INNER ION 22X18   | 3    | BJI1820     | BASE FOAM JAY ION 18X20 |
| 2    | ICJI2219    | COVER INNER ION 22X19   | 3    | BJI1821     | BASE FOAM JAY ION 18X21 |
| 2    | ICJI2220    | COVER INNER ION 22X20   | 3    | BJI1822     | BASE FOAM JAY ION 18X22 |
| 2    | ICJI2221    | COVER INNER ION 22X21   | 3    | BJI1823     | BASE FOAM JAY ION 18X23 |
| 2    | ICJI2222    | COVER INNER ION 22X22   | 3    | BJI1824     | BASE FOAM JAY ION 18X24 |
| 2    | ICJI2223    | COVER INNER ION 22X23   | 3    | BJI1915     | BASE FOAM JAY ION 19X15 |
| 2    | ICJI2224    | COVER INNER ION 22X24   | 3    | BJI1916     | BASE FOAM JAY ION 19X16 |
| 2    | ICJI2319    | COVER INNER ION 23X19   | 3    | BJI1917     | BASE FOAM JAY ION 19X17 |
| 2    | ICJI2320    | COVER INNER ION 23X20   | 3    | BJI1918     | BASE FOAM JAY ION 19X18 |
| 2    | ICJI2321    | COVER INNER ION 23X21   | 3    | BJI1919     | BASE FOAM JAY ION 19X19 |
| 2    | ICJI2322    | COVER INNER ION 23X22   | 3    | BJI1920     | BASE FOAM JAY ION 19X20 |
| 2    | ICJI2323    | COVER INNER ION 23X23   | 3    | BJI1921     | BASE FOAM JAY ION 19X21 |
| 2    | ICJI2324    | COVER INNER ION 23X24   | 3    | BJI1922     | BASE FOAM JAY ION 19X22 |
| 2    | ICJI2419    | COVER INNER ION 24X19   | 3    | BJI1923     | BASE FOAM JAY ION 19X23 |
| 2    | ICJI2420    | COVER INNER ION 24X20   | 3    | BJI1924     | BASE FOAM JAY ION 19X24 |
| 2    | ICJI2421    | COVER INNER ION 24X21   | 3    | BJI2015     | BASE FOAM JAY ION 20X15 |
| 2    | ICJI2422    | COVER INNER ION 24X22   | 3    | BJI2016     | BASE FOAM JAY ION 20X16 |
| 2    | ICJI2423    | COVER INNER ION 24X23   | 3    | BJI2017     | BASE FOAM JAY ION 20X17 |
| 2    | ICJI2424    | COVER INNER ION 24X24   | 3    | BJI2018     | BASE FOAM JAY ION 20X18 |
| 3    | BJI1414     | BASE FOAM JAY ION 14X14 | 3    | BJI2019     | BASE FOAM JAY ION 20X19 |
| 3    | BJI1415     | BASE FOAM JAY ION 14X15 | 3    | BJI2020     | BASE FOAM JAY ION 20X20 |
| 3    | BJI1416     | BASE FOAM JAY ION 14X16 | 3    | BJI2021     | BASE FOAM JAY ION 20X21 |
| 3    | BJI1417     | BASE FOAM JAY ION 14X17 | 3    | BJI2022     | BASE FOAM JAY ION 20X22 |
| 3    | BJI1418     | BASE FOAM JAY ION 14X18 | 3    | BJI2023     | BASE FOAM JAY ION 20X23 |
| 3    | BJI1514     | BASE FOAM JAY ION 15X14 | 3    | BJI2024     | BASE FOAM JAY ION 20X24 |
| 3    | BJI1515     | BASE FOAM JAY ION 15X15 | 3    | BJI2117     | BASE FOAM JAY ION 21X17 |
| 3    | BJI1516     | BASE FOAM JAY ION 15X16 | 3    | BJI2118     | BASE FOAM JAY ION 21X18 |
| 3    | BJI1517     | BASE FOAM JAY ION 15X17 | 3    | BJI2119     | BASE FOAM JAY ION 21X19 |
| 3    | BJI1518     | BASE FOAM JAY ION 15X18 | 3    | BJI2120     | BASE FOAM JAY ION 21X20 |
| 3    | BJI1519     | BASE FOAM JAY ION 15X19 | 3    | BJI2121     | BASE FOAM JAY ION 21X21 |
| 3    | BJI1520     | BASE FOAM JAY ION 15X20 | 3    | BJI2122     | BASE FOAM JAY ION 21X22 |
| 3    | BJI1521     | BASE FOAM JAY ION 15X21 | 3    | BJI2123     | BASE FOAM JAY ION 21X23 |
| 3    | BJI1614     | BASE FOAM JAY ION 16X14 | 3    | BJI2124     | BASE FOAM JAY ION 21X24 |
| 3    | BJI1615     | BASE FOAM JAY ION 16X15 | 3    | BJI2217     | BASE FOAM JAY ION 22X17 |
| 3    | BJI1616     | BASE FOAM JAY ION 16X16 | 3    | BJI2218     | BASE FOAM JAY ION 22X18 |
| 3    | BJI1617     | BASE FOAM JAY ION 16X17 | 3    | BJI2219     | BASE FOAM JAY ION 22X19 |
| 3    | BJI1618     | BASE FOAM JAY ION 16X18 | 3    | BJI2220     | BASE FOAM JAY ION 22X20 |
| 3    | BJI1619     | BASE FOAM JAY ION 16X19 | 3    | BJI2221     | BASE FOAM JAY ION 22X21 |
| 3    | BJI1620     | BASE FOAM JAY ION 16X20 | 3    | BJI2222     | BASE FOAM JAY ION 22X22 |
| 3    | BJI1621     | BASE FOAM JAY ION 16X21 | 3    | BJI2223     | BASE FOAM JAY ION 22X23 |
| 3    | BJI1622     | BASE FOAM JAY ION 16X22 | 3    | BJI2224     | BASE FOAM JAY ION 22X24 |
| 3    | BJI1714     | BASE FOAM JAY ION 17X14 | 3    | BJI2319     | BASE FOAM JAY ION 23X19 |
| 3    | BJI1715     | BASE FOAM JAY ION 17X15 | 3    | BJI2320     | BASE FOAM JAY ION 23X20 |
| 3    | BJI1716     | BASE FOAM JAY ION 17X16 | 3    | BJI2321     | BASE FOAM JAY ION 23X21 |
| 3    | BJI1717     | BASE FOAM JAY ION 17X17 | 3    | BJI2322     | BASE FOAM JAY ION 23X22 |
| 3    | BJI1718     | BASE FOAM JAY ION 17X18 | 3    | BJI2323     | BASE FOAM JAY ION 23X23 |


| Pos. | Item Number | Description             | Pos. | Item Number | Description             |
|------|-------------|-------------------------|------|-------------|-------------------------|
| 3    | BJI2324     | BASE FOAM JAY ION 23X24 | 3    | BJI2422     | BASE FOAM JAY ION 24X22 |
| 3    | BJI2419     | BASE FOAM JAY ION 24X19 | 3    | BJI2423     | BASE FOAM JAY ION 24X23 |
| 3    | BJI2420     | BASE FOAM JAY ION 24X20 | 3    | BJI2424     | BASE FOAM JAY ION 24X24 |
| 3    | BJI2421     | BASE FOAM JAY ION 24X21 |      |             |                         |

NOTE: SACRAL SEAL HELPS PREVENT FLUID FROM MIGRATING TO THE REAR OF THE FLUID PAD WHEN INDIVIDUALS ARE IN RECLINED OR POSTERIOR PELVIC TILT POSITION.



| Pos. | Item Number | Description                   | Remarks | Pos. | Item Number | Description                     | Remarks |
|------|-------------|-------------------------------|---------|------|-------------|---------------------------------|---------|
| 1    | C2100-MJ    | COVER J2 18 x 16              |         | 1A   | C2104A      | COVER J2 AIR EXCHANGE 24 x 20   |         |
| 1    | C2101       | COVER J2 18 x 18              |         | 1A   | C2105A      | COVER J2 AIR EXCHANGE 15.5 x 16 |         |
| 1    | C2102       | COVER J2 20 x 20              |         | 1A   | C2106A      | COVER J2 AIR EXCHANGE 15.5 x 18 |         |
| 1    | C2103       | COVER J2 18 x 20              |         | 1A   | C2107A      | COVER J2 AIR EXCHANGE 15.5 x 20 |         |
| 1    | C2104       | COVER J2 24 x 20              |         | 1A   | C2108A      | COVER J2 AIR EXCHANGE 14 x 16   |         |
| 1    | C2105       | COVER J2 15.5 x 16            |         | 1A   | C2116A      | COVER J2 AIR EXCHANGE 20 x 16   |         |
| 1    | C2106       | COVER J2 15.5 x 18            |         | 1A   | C2117A      | COVER J2 AIR EXCHANGE 17 x 17   |         |
| 1    | C2107       | COVER J2 15.5 x 20            |         | 1A   | C2118A      | COVER J2 AIR EXCHANGE 20 x 18   |         |
| 1    | C2108       | COVER J2 14 x 16              |         | 1B   | C2100N      | COVER J2 INCONTINENT 18 x 16    |         |
| 1    | C2116       | COVER J2 20 x 16              |         | 1B   | C2101N      | COVER J2 INCONTINENT 18 x 18    |         |
| 1    | C2117       | COVER J2 17 x 17              |         | 1B   | C2102N      | COVER J2 INCONTINENT 20 x 20    |         |
| 1    | C2118       | COVER J2 20 x 18              |         | 1B   | C2103N      | COVER J2 INCONTINENT 18 x 20    |         |
| 1A   | C2100A      | COVER J2 AIR EXCHANGE 18 x 16 |         | 1B   | C2104N      | COVER J2 INCONTINENT 24 x 20    |         |
| 1A   | C2101A      | COVER J2 AIR EXCHANGE 18 x 18 |         | 1B   | C2105N      | COVER J2 INCONTINENT 15.5 x 16  |         |
| 1A   | C2102A      | COVER J2 AIR EXCHANGE 20 x 20 |         | 1B   | C2106N      | COVER J2 INCONTINENT 15.5 x 18  |         |
| 1A   | C2103A      | COVER J2 AIR EXCHANGE 18 x 20 |         | 1B   | C2107N      | COVER J2 INCONTINENT 15.5 x 20  |         |

| Pos. | Item Number | Description                    | Remarks | Pos. | Item Number | Description                      | Remarks |
|------|-------------|--------------------------------|---------|------|-------------|----------------------------------|---------|
| 1B   | C2108N      | COVER J2 INCONTINEN T14 x 16   |         | 2A   | F2117.1     | FLUID PAD J2 RECLINE 17 x 17     | *       |
| 1B   | C2116N      | COVER J2 INCONTINENT 20 x 16   |         | 2A   | F2118.1     | FLUID PAD J2 RECLINE 20 x 18     | *       |
| 1B   | C2117N      | COVER J2 INCONTINENT 17 x 17   |         | 3    | B2100       | BASE J2 18 x 16                  |         |
| 1B   | C2118N      | COVER J2 INCONTINENT 20 x 18   |         | 3    | B2101       | BASE J2 18 x 18                  |         |
| 2    | F2100       | FLUID PAD J2 18 x 16           |         | 3    | B2102       | BASE J2 20 x20                   |         |
| 2    | F2101       | FLUID PAD J2 18 x 18           |         | 3    | B2103       | BASE J2 18 x 20                  |         |
| 2    | F2102       | FLUID PAD J2 20 x20            |         | 3    | B2104       | BASE J2 24 x 20                  |         |
| 2    | F2103       | FLUID PAD J2 18 x 20           |         | 3    | B2105       | BASE J2 15.5 x 16                |         |
| 2    | F2104       | FLUID PAD J2 24 x 20           |         | 3    | B2106       | BASE J2 15.5 x 18                |         |
| 2    | F2105       | FLUID PAD J2 15.5 x 16         |         | 3    | B2107       | BASE J2 15.5 x 20                |         |
| 2    | F2106       | FLUID PAD J2 15.5 x 18         |         | 3    | B2108       | BASE J2 14 x 16                  |         |
| 2    | F2107       | FLUID PAD J2 15.5 x 20         |         | 3    | B2116       | BASE J2 20 x 16                  |         |
| 2    | F2108       | FLUID PAD J2 14 x 16           |         | 3    | B2117       | BASE J2 17 x 17                  |         |
| 2    | F2116       | FLUID PAD J2 20 x 16           |         | 3    | B2118       | BASE J2 20 x 18                  |         |
| 2    | F2117       | FLUID PAD J2 17 x 17           |         | 4    | S2100       | SOLID SEAT J2 NON SKID 18 x 16   |         |
| 2    | F2118       | FLUID PAD J2 20 x 18           |         | 4    | S2101       | SOLID SEAT J2 NON SKID 18 x 18   |         |
| 2A   | F2100.1     | FLUID PAD J2 RECLINE 18 x 16   | *       | 4    | S2102       | SOLID SEAT J2 NON SKID 20 x20    |         |
| 2A   | F2101.1     | FLUID PAD J2 RECLINE 18 x 18   | *       | 4    | S2103       | SOLID SEAT J2 NON SKID 18 x 20   |         |
| 2A   | F2102.1     | FLUID PAD J2 RECLINE 20 x20    | *       | 4    | S2104       | SOLID SEAT J2 NON SKID 24 x 20   |         |
| 2A   | F2103.1     | FLUID PAD J2 RECLINE 18 x 20   | *       | 4    | S2105       | SOLID SEAT J2 NON SKID 15.5 x 16 |         |
| 2A   | F2104.1     | FLUID PAD J2 RECLINE 24 x 20   | *       | 4    | S2106       | SOLID SEAT J2 NON SKID 15.5 x 18 |         |
| 2A   | F2105.1     | FLUID PAD J2 RECLINE 15.5 x 16 | *       | 4    | S2107       | SOLID SEAT J2 NON SKID 15.5 x 20 |         |
| 2A   | F2106.1     | FLUID PAD J2 RECLINE 15.5 x 18 | *       | 4    | S2108       | SOLID SEAT J2 NON SKID 14 x 16   |         |
| 2A   | F2107.1     | FLUID PAD J2 RECLINE 15.5 x 20 | *       | 4    | S2116       | SOLID SEAT J2 NON SKID 20 x 16   |         |
| 2A   | F2108.1     | FLUID PAD J2 RECLINE 14 x 16   | *       | 4    | S2117       | SOLID SEAT J2 NON SKID 17 x 17   |         |
| 2A   | F2116.1     | FLUID PAD J2 RECLINE 20 x 16   | *       | 4    | S2118       | SOLID SEAT J2 NON SKID 20 x 18   |         |



| Pos. | Item Number | Description                                | Remarks        |
|------|-------------|--------------------------------------------|----------------|
| 1A   | C2444       | COVER J2 DEEP CONTUR 14 x 14               | ZIPPERED COVER |
| 1A   | C2446       | COVER J2 DEEP CONTUR 14 x 16               | ZIPPERED COVER |
| 1A   | C2455       | COVER J2 DEEP CONTUR 15 x 15               | ZIPPERED COVER |
| 1A   | C2457       | COVER J2 DEEP CONTUR 15 x 17               | ZIPPERED COVER |
| 1A   | C2466       | COVER J2 DEEP CONTUR 16 x 16               | ZIPPERED COVER |
| 1A   | C2468       | COVER J2 DEEP CONTUR 16 x 18               | ZIPPERED COVER |
| 1A   | C2475       | COVER J2 DEEP CONTUR 17 x 15               | ZIPPERED COVER |
| 1A   | C2477       | COVER J2 DEEP CONTUR 17 x 17               | ZIPPERED COVER |
| 1A   | C2486       | COVER J2 DEEP CONTUR 18 x 16               | ZIPPERED COVER |
| 1A   | C2488       | COVER J2 DEEP CONTUR 18 x 18               | ZIPPERED COVER |
| 1A   | C2406       | COVER J2 DEEP CONTUR 20 x 16               | ZIPPERED COVER |
| 1A   | C2408       | COVER J2 DEEP CONTUR 20 x 18               | ZIPPERED COVER |
| 1A   | C2426       | COVER J2 DEEP CONTUR 22 x 16               | ZIPPERED COVER |
| 1A   | C2428       | COVER J2 DEEP CONTUR 22 x 18               | ZIPPERED COVER |
| 1A   | C2460       | COVER J2 DEEP CONTUR 16 x 20               | ZIPPERED COVER |
| 1A   | C2480       | COVER J2 DEEP CONTUR 18 x 20               | ZIPPERED COVER |
| 1A   | C2400       | COVER J2 DEEP CONTUR 20 x 20               | ZIPPERED COVER |
| 1A   | C2448       | COVER J2 DEEP CONTUR 24 x 18               | ZIPPERED COVER |
| 1B   | C2444A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 14 x 14 | ZIPPERED COVER |
| 1B   | C2446A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 14 x 16 | ZIPPERED COVER |

| Pos. | Item Number | Description                                               | Remarks                               |
|------|-------------|-----------------------------------------------------------|---------------------------------------|
| 1B   | C2455A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 15 x 15                | ZIPPERED COVER                        |
| 1B   | C2457A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 15 x 17                | ZIPPERED COVER                        |
| 1B   | C2466A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 16 x 16                | ZIPPERED COVER                        |
| 1B   | C2468A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 16 x 18                | ZIPPERED COVER                        |
| 1B   | C2475A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 17 x 15                | ZIPPERED COVER                        |
| 1B   | C2477A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 17 x 17                | ZIPPERED COVER                        |
| 1B   | C2486A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 18 x 16                | ZIPPERED COVER                        |
| 1B   | C2488A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 18 x 18                | ZIPPERED COVER                        |
| 1B   | C2406A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 20 x 16                | ZIPPERED COVER                        |
| 1B   | C2408A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 20 x 18                | ZIPPERED COVER                        |
| 1B   | C2426A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 22 x 16                | ZIPPERED COVER                        |
| 1B   | C2428A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 22 x 18                | ZIPPERED COVER                        |
| 1B   | C2460A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 16 x 20                | ZIPPERED COVER                        |
| 1B   | C2480A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 18 x 20                | ZIPPERED COVER                        |
| 1B   | C2400A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 20 x 20                | ZIPPERED COVER                        |
| 1B   | C2448A      | COVER J2 DEEP CONTOUR AIR EXCHANGE 24 x 18                | ZIPPERED COVER                        |
| 1C   | C2444N      | COVER J2 DEEP CONTOUR INCONTINENT 14 x 14                 | ZIPPERED COVER                        |
| 1C   | C2446N      | COVER J2 DEEP CONTOUR INCONTINENT 14 x 16                 | ZIPPERED COVER                        |
| 1C   | C2455N      | COVER J2 DEEP CONTOUR INCONTINENT 15 x 15                 | ZIPPERED COVER                        |
| 1C   | C2457N      | COVER J2 DEEP CONTOUR INCONTINENT 15 x 17                 | ZIPPERED COVER                        |
| 1C   | C2466N      | COVER J2 DEEP CONTOUR INCONTINENT 16 x 16                 | ZIPPERED COVER                        |
| 1C   | C2468N      | COVER J2 DEEP CONTOUR INCONTINENT 16 x 18                 | ZIPPERED COVER                        |
| 1C   | C2475N      | COVER J2 DEEP CONTOUR INCONTINENT 17 x 15                 | ZIPPERED COVER                        |
| 1C   | C2477N      | COVER J2 DEEP CONTOUR INCONTINENT 17 x 17                 | ZIPPERED COVER                        |
| 1C   | C2486N      | COVER J2 DEEP CONTOUR INCONTINENT 18 x 16                 | ZIPPERED COVER                        |
| 1C   | C2488N      | COVER J2 DEEP CONTOUR INCONTINENT 18 x 18                 | ZIPPERED COVER                        |
| 1C   | C2406N      | COVER J2 DEEP CONTOUR INCONTINENT 20 x 16                 | ZIPPERED COVER                        |
| 1C   | C2408N      | COVER J2 DEEP CONTOUR INCONTINENT 20 x 18                 | ZIPPERED COVER                        |
| 1C   | C2426N      | COVER J2 DEEP CONTOUR INCONTINENT 22 x 16                 | ZIPPERED COVER                        |
| 1C   | C2428N      | COVER J2 DEEP CONTOUR INCONTINENT 22 x 18                 | ZIPPERED COVER                        |
| 1C   | C2460N      | COVER J2 DEEP CONTOUR INCONTINENT 16 x 20                 | ZIPPERED COVER                        |
| 1C   | C2480N      | COVER J2 DEEP CONTOUR INCONTINENT 18 x 20                 | ZIPPERED COVER                        |
| 1C   | C2400N      | COVER J2 DEEP CONTOUR INCONTINENT 20 x 20                 | ZIPPERED COVER                        |
| 1C   | C2448N      | COVER J2 DEEP CONTOUR INCONTINENT 24 x 18                 | ZIPPERED COVER                        |
| 1D   | C2444D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 14 x 14              | Non-zippered Cover - Standard Stretch |
| 1D   | C2446D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 14 x 16              | Non-zippered Cover - Standard Stretch |
| 1D   | C2455D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 15 x 15              | Non-zippered Cover - Standard Stretch |
| 1D   | C2457D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 15 x 17              | Non-zippered Cover - Standard Stretch |
| 1D   | C2466D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 16 x 16              | Non-zippered Cover - Standard Stretch |
| 1D   | C2468D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 16 x 18              | Non-zippered Cover - Standard Stretch |
| 1D   | C2475D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 17 x 15              | Non-zippered Cover - Standard Stretch |
| 1D   | C2477D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 17 x 17              | Non-zippered Cover - Standard Stretch |
| 1D   | C2486D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 18 x 16              | Non-zippered Cover - Standard Stretch |
| 1D   | C2488D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 18 x 18              | Non-zippered Cover - Standard Stretch |
| 1D   | C2406D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 20 x 16              | Non-zippered Cover - Standard Stretch |
| 1D   | C2408D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 20 x 18              | Non-zippered Cover - Standard Stretch |
| 1D   | C2426D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 22 x 16              | Non-zippered Cover - Standard Stretch |
| 1D   | C2428D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 22 x 18              | Non-zippered Cover - Standard Stretch |
| 1D   | C2460D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 16 x 20              | Non-zippered Cover - Standard Stretch |
| 1D   | C2480D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 18 x 20              | Non-zippered Cover - Standard Stretch |
| 1D   | C2400D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 20 x 20              | Non-zippered Cover - Standard Stretch |
| 1D   | C2448D      | COVER J2 DEEP CONTOUR WITH DROP SEAT 24 x 18              | Non-zippered Cover - Standard Stretch |
| 1E   | C2444AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 14 x 14 | NON-ZIPPERD COVER                     |
| 1E   | C2446AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 14 x 16 | NON-ZIPPERD COVER                     |
| 1E   | C2455AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 15 x 15 | NON-ZIPPERD COVER                     |
| 1E   | C2457AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 15 x 17 | NON-ZIPPERD COVER                     |
| 1E   | C2466AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 16 x 16 | NON-ZIPPERD COVER                     |


| Pos. | Item Number | Description                                               | Remarks           |
|------|-------------|-----------------------------------------------------------|-------------------|
| 1E   | C2468AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 16 x 18 | NON-ZIPPERD COVER |
| 1E   | C2475AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 17 x 15 | NON-ZIPPERD COVER |
| 1E   | C2477AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 17 x 17 | NON-ZIPPERD COVER |
| 1E   | C2486AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 18 x 16 | NON-ZIPPERD COVER |
| 1E   | C2488AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 18 x 18 | NON-ZIPPERD COVER |
| 1E   | C2406AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 20 x 16 | NON-ZIPPERD COVER |
| 1E   | C2408AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 20 x 18 | NON-ZIPPERD COVER |
| 1E   | C2426AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 22 x 16 | NON-ZIPPERD COVER |
| 1E   | C2428AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 22 x 18 | NON-ZIPPERD COVER |
| 1E   | C2460AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 16 x 20 | NON-ZIPPERD COVER |
| 1E   | C2480AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 18 x 20 | NON-ZIPPERD COVER |
| 1E   | C2400AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 20 x 20 | NON-ZIPPERD COVER |
| 1E   | C2448AD     | COVER J2 DEEP CONTOUR AIR EXCHANGE WITH DROP SEAT 24 x 18 | NON-ZIPPERD COVER |
| 1F   | C2444ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 14 x 14  | NON-ZIPPERD COVER |
| 1F   | C2446ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 14 x 16  | NON-ZIPPERD COVER |
| 1F   | C2455ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 15 x 15  | NON-ZIPPERD COVER |
| 1F   | C2457ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 15 x 17  | NON-ZIPPERD COVER |
| 1F   | C2466ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 16 x 16  | NON-ZIPPERD COVER |
| 1F   | C2468ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 16 x 18  | NON-ZIPPERD COVER |
| 1F   | C2475ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 17 x 15  | NON-ZIPPERD COVER |
| 1F   | C2477ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 17 x 17  | NON-ZIPPERD COVER |
| 1F   | C2486ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 18 x 16  | NON-ZIPPERD COVER |
| 1F   | C2488ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 18 x 18  | NON-ZIPPERD COVER |
| 1F   | C2406ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 20 x 16  | NON-ZIPPERD COVER |
| 1F   | C2408ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 20 x 18  | NON-ZIPPERD COVER |
| 1F   | C2426ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 22 x 16  | NON-ZIPPERD COVER |
| 1F   | C2428ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 22 x 18  | NON-ZIPPERD COVER |
| 1F   | C2460ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 16 x 20  | NON-ZIPPERD COVER |
| 1F   | C2480ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 18 x 20  | NON-ZIPPERD COVER |
| 1F   | C2400ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 20 x 20  | NON-ZIPPERD COVER |
| 1F   | C2448ND     | COVER J2 DEEP CONTOUR INCONTINENT WITH DROP SEAT 24 x 18  | NON-ZIPPERD COVER |
| 2    | F2444       | FLUID PAD J2 DEEP CONTOUR 14 x 14                         |                   |
| 2    | F2446       | FLUID PAD J2 DEEP CONTOUR 14 x 16                         |                   |
| 2    | F2455       | FLUID PAD J2 DEEP CONTOUR 15 x 15                         |                   |
| 2    | F2457       | FLUID PAD J2 DEEP CONTOUR 15 x 17                         |                   |
| 2    | F2466       | FLUID PAD J2 DEEP CONTOUR 16 x 16                         |                   |
| 2    | F2468       | FLUID PAD J2 DEEP CONTOUR 16 x 18                         |                   |
| 2    | F2475       | FLUID PAD J2 DEEP CONTOUR 17 x 15                         |                   |
| 2    | F2477       | FLUID PAD J2 DEEP CONTOUR 17 x 17                         |                   |
| 2    | F2486       | FLUID PAD J2 DEEP CONTOUR 18 x 16                         |                   |
| 2    | F2488       | FLUID PAD J2 DEEP CONTOUR 18 x 18                         |                   |
| 2    | F2406       | FLUID PAD J2 DEEP CONTOUR 20 x 16                         |                   |
| 2    | F2408       | FLUID PAD J2 DEEP CONTOUR 20 x 18                         |                   |
| 2    | F2426       | FLUID PAD J2 DEEP CONTOUR 22 x 16                         |                   |
| 2    | F2428       | FLUID PAD J2 DEEP CONTOUR 22 x 18                         |                   |
| 2    | F2460       | FLUID PAD J2 DEEP CONTOUR 16 x 20                         |                   |
| 2    | F2480       | FLUID PAD J2 DEEP CONTOUR 18 x 20                         |                   |
| 2    | F2400       | FLUID PAD J2 DEEP CONTOUR 20 x 20                         |                   |
| 2    | F2448       | FLUID PAD J2 DEEP CONTOUR 24 x 18                         |                   |
| 3    | B2444       | BASE J2 DEEP CONTOUR 14 x 14                              |                   |
| 3    | B2446       | BASE J2 DEEP CONTOUR 14 x 16                              |                   |
| 3    | B2455       | BASE J2 DEEP CONTOUR 15 x 15                              |                   |
| 3    | B2457       | BASE J2 DEEP CONTOUR 15 x 17                              |                   |
| 3    | B2466       | BASE J2 DEEP CONTOUR 16 x 16                              |                   |
| 3    | B2468       | BASE J2 DEEP CONTOUR 16 x 18                              |                   |
| 3    | B2475       | BASE J2 DEEP CONTOUR 17 x 15                              |                   |
| 3    | B2477       | BASE J2 DEEP CONTOUR 17 x 17                              |                   |

| Pos. | Item Number | Description                               | Remarks                             |
|------|-------------|-------------------------------------------|-------------------------------------|
| 3    | B2486       | BASE J2 DEEP CONTOUR 18 x 16              |                                     |
| 3    | B2488       | BASE J2 DEEP CONTOUR 18 x 18              |                                     |
| 3    | B2406       | BASE J2 DEEP CONTOUR 20 x 16              |                                     |
| 3    | B2408       | BASE J2 DEEP CONTOUR 20 x 18              |                                     |
| 3    | B2426       | BASE J2 DEEP CONTOUR 22 x 16              |                                     |
| 3    | B2428       | BASE J2 DEEP CONTOUR 22 x 18              |                                     |
| 3    | B2460       | BASE J2 DEEP CONTOUR 16 x 20              |                                     |
| 3    | B2480       | BASE J2 DEEP CONTOUR 18 x 20              |                                     |
| 3    | B2400       | BASE J2 DEEP CONTOUR 20 x 20              |                                     |
| 3    | B2448       | BASE J2 DEEP CONTOUR 24 x 18              |                                     |
| 4A   | S2444       | SOLID SEAT INSERT J2 DEEP CONTOUR 14 x 14 |                                     |
| 4A   | S2446       | SOLID SEAT INSERT J2 DEEP CONTOUR 14 x 16 |                                     |
| 4A   | S2455       | SOLID SEAT INSERT J2 DEEP CONTOUR 15 x 15 |                                     |
| 4A   | S2457       | SOLID SEAT INSERT J2 DEEP CONTOUR 15 x 17 |                                     |
| 4A   | S2466       | SOLID SEAT INSERT J2 DEEP CONTOUR 16 x 16 |                                     |
| 4A   | S2468       | SOLID SEAT INSERT J2 DEEP CONTOUR 16 x 18 |                                     |
| 4A   | S2475       | SOLID SEAT INSERT J2 DEEP CONTOUR 17 x 15 |                                     |
| 4A   | S2477       | SOLID SEAT INSERT J2 DEEP CONTOUR 17 x 17 |                                     |
| 4A   | S2486       | SOLID SEAT INSERT J2 DEEP CONTOUR 18 x 16 |                                     |
| 4A   | S2488       | SOLID SEAT INSERT J2 DEEP CONTOUR 18 x 18 |                                     |
| 4A   | S2406       | SOLID SEAT INSERT J2 DEEP CONTOUR 20 x 16 |                                     |
| 4A   | S2408       | SOLID SEAT INSERT J2 DEEP CONTOUR 20 x 18 |                                     |
| 4A   | S2426       | SOLID SEAT INSERT J2 DEEP CONTOUR 22 x 16 |                                     |
| 4A   | S2428       | SOLID SEAT INSERT J2 DEEP CONTOUR 22 x 18 |                                     |
| 4A   | S2460       | SOLID SEAT INSERT J2 DEEP CONTOUR 16 x 20 |                                     |
| 4A   | S2480       | SOLID SEAT INSERT J2 DEEP CONTOUR 18 x 20 |                                     |
| 4A   | S2400       | SOLID SEAT INSERT J2 DEEP CONTOUR 20 x 20 |                                     |
| 4A   | S2448       | SOLID SEAT INSERT J2 DEEP CONTOUR 24 x 18 |                                     |
| 5    | M-24S       | SHORT 1/2" DROP HOOKS                     | Pack of 4                           |
| 5    | M-24M       | MEDIUM 1" DROP HOOKS                      | Pack of 4                           |
| 5    | M-24L       | LONG 1 1/2" DROP HOOKS                    | Pack of 4                           |
| 4-6  | SD2444      | SOLID DROP SEAT WITH HARDWARE 14 x 14     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2446      | SOLID DROP SEAT WITH HARDWARE 14 x 16     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2455      | SOLID DROP SEAT WITH HARDWARE 15 x 15     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2457      | SOLID DROP SEAT WITH HARDWARE 15 x 17     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2466      | SOLID DROP SEAT WITH HARDWARE 16 x 16     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2468      | SOLID DROP SEAT WITH HARDWARE 16 x 18     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2475      | SOLID DROP SEAT WITH HARDWARE 17 x 15     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2477      | SOLID DROP SEAT WITH HARDWARE 17 x 17     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2486      | SOLID DROP SEAT WITH HARDWARE 18 x 16     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2488      | SOLID DROP SEAT WITH HARDWARE 18 x 18     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2406      | SOLID DROP SEAT WITH HARDWARE 20 x 16     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2408      | SOLID DROP SEAT WITH HARDWARE 20 x 18     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2426      | SOLID DROP SEAT WITH HARDWARE 22 x 16     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2428      | SOLID DROP SEAT WITH HARDWARE 22 x 18     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2460      | SOLID DROP SEAT WITH HARDWARE 16 x 20     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2480      | SOLID DROP SEAT WITH HARDWARE 18 x 20     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2400      | SOLID DROP SEAT WITH HARDWARE 20 x 20     | Ships standard w/ 1 1/2' drop hooks |
| 4-6  | SD2448      | SOLID DROP SEAT WITH HARDWARE 24 x 18     | Ships standard w/ 1 1/2' drop hooks |
| 6    | M-2404H     | HANDLE KIT J2 DEEP CONTOUR CUSHION 14"    |                                     |
| 6    | M-2405H     | HANDLE KIT J2 DEEP CONTOUR CUSHION 15"    |                                     |
| 6    | M-2406H     | HANDLE KIT J2 DEEP CONTOUR CUSHION 16"    |                                     |
| 6    | M-2407H     | HANDLE KIT J2 DEEP CONTOUR CUSHION 17"    |                                     |
| 6    | M-2408H     | HANDLE KIT J2 DEEP CONTOUR CUSHION 18"    |                                     |
| 6    | M-2409H     | HANDLE KIT J2 DEEP CONTOUR CUSHION 20"    |                                     |
| 6    | M-2410H     | HANDLE KIT J2 DEEP CONTOUR CUSHION 22"    |                                     |
| 6    | M-2411H     | HANDLE KIT J2 DEEP CONTOUR CUSHION 24"    |                                     |

| Pos.     | Item Number | Description                                 | Remarks                             |
|----------|-------------|---------------------------------------------|-------------------------------------|
| 3,4,5, 6 | B2444       | BASE J2 DEEP CONTOUR WITH DROP SEAT 14 x 14 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2446       | BASE J2 DEEP CONTOUR WITH DROP SEAT 14 x 16 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2455       | BASE J2 DEEP CONTOUR WITH DROP SEAT 15 x 15 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2457       | BASE J2 DEEP CONTOUR WITH DROP SEAT 15 x 17 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2466       | BASE J2 DEEP CONTOUR WITH DROP SEAT 16 x 16 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2468       | BASE J2 DEEP CONTOUR WITH DROP SEAT 16 x 18 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2475       | BASE J2 DEEP CONTOUR WITH DROP SEAT 17 x 15 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2477       | BASE J2 DEEP CONTOUR WITH DROP SEAT 17 x 17 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2486       | BASE J2 DEEP CONTOUR WITH DROP SEAT 18 x 16 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2488       | BASE J2 DEEP CONTOUR WITH DROP SEAT 18 x 18 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2406       | BASE J2 DEEP CONTOUR WITH DROP SEAT 20 x 16 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2408       | BASE J2 DEEP CONTOUR WITH DROP SEAT 20 x 18 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2426       | BASE J2 DEEP CONTOUR WITH DROP SEAT 22 x 16 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2428       | BASE J2 DEEP CONTOUR WITH DROP SEAT 22 x 18 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2460       | BASE J2 DEEP CONTOUR WITH DROP SEAT 16 x 20 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2480       | BASE J2 DEEP CONTOUR WITH DROP SEAT 18 x 20 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2400       | BASE J2 DEEP CONTOUR WITH DROP SEAT 20 x 20 | Ships standard w/ 1 1/2' drop hooks |
| 3,4,5, 6 | B2448       | BASE J2 DEEP CONTOUR WITH DROP SEAT 24 x 18 | Ships standard w/ 1 1/2' drop hooks |

**NOTE:** SHIPS WITH FULL SET OF 6 FOAM PIECES.

**NOTE:** EFFECTIVE MAY 1, 2016, ODD SIZES AND CUSTOM FOAM KIT COMPLETE CUSHIONS WILL NO LONGER BE AVAILABLE. PARTS WILL CONTINUE TO BE OFFERED.



**WELL FOAM KIT DENSITY GUIDE**

| ITEM | LAYER  | FOAM DENSITY |
|------|--------|--------------|
| 5A   | BOTTOM | MEDIUM       |
| 5B   | BOTTOM | FIRM         |
| 5A   | MIDDLE | MED/SOFT     |
| 5B   | MIDDLE | MEDIUM       |
| 5A   | TOP    | MED/SOFT     |
| 5B   | TOP    | SOFT         |

**FLUID PAD/WELL FOAM INSERT GUIDE**

| CUSHION DEPTH | FLUID PAD | STANDARD WELL FOAM | CUSTOM WELL FOAM |
|---------------|-----------|--------------------|------------------|
| 18"           | F300018   | 300018-MJ          | 300018-MJ        |
| 19"           | F300020   | 300020-MJ          | 300020-MJ        |
| 20"           | F300020   | 300020-MJ          | 300020-MJ        |
| 21"           | F300022   | 300022-MJ          | 300022-MJ        |
| 22"           | F300022   | 300022-MJ          | 300022-MJ        |

| Pos. | Item Number | Description                                | Remarks |
|------|-------------|--------------------------------------------|---------|
| 1A   | C302018A    | J2 PLUS CUSHION COVER AIR EXCHANGE 20 x 18 |         |
| 1A   | C302019A    | J2 PLUS CUSHION COVER AIR EXCHANGE 20 x 19 |         |
| 1A   | C302020A    | J2 PLUS CUSHION COVER AIR EXCHANGE 20 x 20 |         |
| 1A   | C302021A    | J2 PLUS CUSHION COVER AIR EXCHANGE 20 x 21 |         |
| 1A   | C302022A    | J2 PLUS CUSHION COVER AIR EXCHANGE 20 x 22 |         |
| 1A   | C302118A    | J2 PLUS CUSHION COVER AIR EXCHANGE 21 x 18 |         |
| 1A   | C302119A    | J2 PLUS CUSHION COVER AIR EXCHANGE 21 x 19 |         |
| 1A   | C302120A    | J2 PLUS CUSHION COVER AIR EXCHANGE 21 x 20 |         |
| 1A   | C302121A    | J2 PLUS CUSHION COVER AIR EXCHANGE 21 x 21 |         |
| 1A   | C302122A    | J2 PLUS CUSHION COVER AIR EXCHANGE 21 x 22 |         |
| 1A   | C302218A    | J2 PLUS CUSHION COVER AIR EXCHANGE 22 x 18 |         |

| Pos. | Item Number | Description                                | Remarks |
|------|-------------|--------------------------------------------|---------|
| 1A   | C302219A    | J2 PLUS CUSHION COVER AIR EXCHANGE 22 x 19 |         |
| 1A   | C302220A    | J2 PLUS CUSHION COVER AIR EXCHANGE 22 x 20 |         |
| 1A   | C302221A    | J2 PLUS CUSHION COVER AIR EXCHANGE 22 x 21 |         |
| 1A   | C302222A    | J2 PLUS CUSHION COVER AIR EXCHANGE 22 x 22 |         |
| 1A   | C302318A    | J2 PLUS CUSHION COVER AIR EXCHANGE 23 x 18 |         |
| 1A   | C302319A    | J2 PLUS CUSHION COVER AIR EXCHANGE 23 x 19 |         |
| 1A   | C302320A    | J2 PLUS CUSHION COVER AIR EXCHANGE 23 x 20 |         |
| 1A   | C302321A    | J2 PLUS CUSHION COVER AIR EXCHANGE 23 x 21 |         |
| 1A   | C302322A    | J2 PLUS CUSHION COVER AIR EXCHANGE 23 x 22 |         |
| 1A   | C302418A    | J2 PLUS CUSHION COVER AIR EXCHANGE 24 x 18 |         |
| 1A   | C302419A    | J2 PLUS CUSHION COVER AIR EXCHANGE 24 x 19 |         |
| 1A   | C302420A    | J2 PLUS CUSHION COVER AIR EXCHANGE 24 x 20 |         |
| 1A   | C302421A    | J2 PLUS CUSHION COVER AIR EXCHANGE 24 x 21 |         |
| 1A   | C302422A    | J2 PLUS CUSHION COVER AIR EXCHANGE 24 x 22 |         |
| 1A   | C302518A    | J2 PLUS CUSHION COVER AIR EXCHANGE 25 x 18 |         |
| 1A   | C302519A    | J2 PLUS CUSHION COVER AIR EXCHANGE 25 x 19 |         |
| 1A   | C302520A    | J2 PLUS CUSHION COVER AIR EXCHANGE 25 x 20 |         |
| 1A   | C302521A    | J2 PLUS CUSHION COVER AIR EXCHANGE 25 x 21 |         |
| 1A   | C302522A    | J2 PLUS CUSHION COVER AIR EXCHANGE 25 x 22 |         |
| 1A   | C302618A    | J2 PLUS CUSHION COVER AIR EXCHANGE 26 x 18 |         |
| 1A   | C302619A    | J2 PLUS CUSHION COVER AIR EXCHANGE 26 x 19 |         |
| 1A   | C302620A    | J2 PLUS CUSHION COVER AIR EXCHANGE 26 x 20 |         |
| 1A   | C302621A    | J2 PLUS CUSHION COVER AIR EXCHANGE 26 x 21 |         |
| 1A   | C302622A    | J2 PLUS CUSHION COVER AIR EXCHANGE 26 x 22 |         |
| 1B   | C302018LS   | J2 PLUS CUSHION COVER LOW SHEAR 20 x 18    |         |
| 1B   | C302019LS   | J2 PLUS CUSHION COVER LOW SHEAR 20 x 19    |         |
| 1B   | C302020LS   | J2 PLUS CUSHION COVER LOW SHEAR 20 x 20    |         |
| 1B   | C302021LS   | J2 PLUS CUSHION COVER LOW SHEAR 20 x 21    |         |
| 1B   | C302022LS   | J2 PLUS CUSHION COVER LOW SHEAR 20 x 22    |         |
| 1B   | C302118LS   | J2 PLUS CUSHION COVER LOW SHEAR 21 x 18    |         |
| 1B   | C302119LS   | J2 PLUS CUSHION COVER LOW SHEAR 21 x 19    |         |
| 1B   | C302120LS   | J2 PLUS CUSHION COVER LOW SHEAR 21 x 20    |         |
| 1B   | C302121LS   | J2 PLUS CUSHION COVER LOW SHEAR 21 x 21    |         |
| 1B   | C302122LS   | J2 PLUS CUSHION COVER LOW SHEAR 21 x 22    |         |
| 1B   | C302218LS   | J2 PLUS CUSHION COVER LOW SHEAR 22 x 18    |         |
| 1B   | C302219LS   | J2 PLUS CUSHION COVER LOW SHEAR 22 x 19    |         |
| 1B   | C302220LS   | J2 PLUS CUSHION COVER LOW SHEAR 22 x 20    |         |
| 1B   | C302221LS   | J2 PLUS CUSHION COVER LOW SHEAR 22 x 21    |         |
| 1B   | C302222LS   | J2 PLUS CUSHION COVER LOW SHEAR 22 x 22    |         |
| 1B   | C302318LS   | J2 PLUS CUSHION COVER LOW SHEAR 23 x 18    |         |
| 1B   | C302319LS   | J2 PLUS CUSHION COVER LOW SHEAR 23 x 19    |         |
| 1B   | C302320LS   | J2 PLUS CUSHION COVER LOW SHEAR 23 x 20    |         |
| 1B   | C302321LS   | J2 PLUS CUSHION COVER LOW SHEAR 23 x 21    |         |
| 1B   | C302322LS   | J2 PLUS CUSHION COVER LOW SHEAR 23 x 22    |         |
| 1B   | C302418LS   | J2 PLUS CUSHION COVER LOW SHEAR 24 x 18    |         |
| 1B   | C302419LS   | J2 PLUS CUSHION COVER LOW SHEAR 24 x 19    |         |
| 1B   | C302420LS   | J2 PLUS CUSHION COVER LOW SHEAR 24 x 20    |         |
| 1B   | C302421LS   | J2 PLUS CUSHION COVER LOW SHEAR 24 x 21    |         |
| 1B   | C302422LS   | J2 PLUS CUSHION COVER LOW SHEAR 24 x 22    |         |
| 1B   | C302518LS   | J2 PLUS CUSHION COVER LOW SHEAR 25 x 18    |         |
| 1B   | C302519LS   | J2 PLUS CUSHION COVER LOW SHEAR 25 x 19    |         |
| 1B   | C302520LS   | J2 PLUS CUSHION COVER LOW SHEAR 25 x 20    |         |
| 1B   | C302521LS   | J2 PLUS CUSHION COVER LOW SHEAR 25 x 21    |         |
| 1B   | C302522LS   | J2 PLUS CUSHION COVER LOW SHEAR 25 x 22    |         |
| 1B   | C302618LS   | J2 PLUS CUSHION COVER LOW SHEAR 26 x 18    |         |
| 1B   | C302619LS   | J2 PLUS CUSHION COVER LOW SHEAR 26 x 19    |         |
| 1B   | C302620LS   | J2 PLUS CUSHION COVER LOW SHEAR 26 x 20    |         |

| Pos. | Item Number | Description                             | Remarks |
|------|-------------|-----------------------------------------|---------|
| 1B   | C302621LS   | J2 PLUS CUSHION COVER LOW SHEAR 26 x 21 |         |
| 1B   | C302622LS   | J2 PLUS CUSHION COVER LOW SHEAR 26 x 22 |         |
| 2    | E302018     | ENVELOPE J2 PLUS 20 x 18                |         |
| 2    | E302019     | ENVELOPE J2 PLUS 20 x 19                |         |
| 2    | E302020     | ENVELOPE J2 PLUS 20 x 20                |         |
| 2    | E302021     | ENVELOPE J2 PLUS 20 x 21                |         |
| 2    | E302022     | ENVELOPE J2 PLUS 20 x 22                |         |
| 2    | E302118     | ENVELOPE J2 PLUS 21 x 18                |         |
| 2    | E302119     | ENVELOPE J2 PLUS 21 x 19                |         |
| 2    | E302120     | ENVELOPE J2 PLUS 21 x 20                |         |
| 2    | E302121     | ENVELOPE J2 PLUS 21 x 21                |         |
| 2    | E302122     | ENVELOPE J2 PLUS 21 x 22                |         |
| 2    | E302218     | ENVELOPE J2 PLUS 22 x 18                |         |
| 2    | E302219     | ENVELOPE J2 PLUS 22 x 19                |         |
| 2    | E302220     | ENVELOPE J2 PLUS 22 x 20                |         |
| 2    | E302221     | ENVELOPE J2 PLUS 22 x 21                |         |
| 2    | E302222     | ENVELOPE J2 PLUS 22 x 22                |         |
| 2    | E302318     | ENVELOPE J2 PLUS 23 x 18                |         |
| 2    | E302319     | ENVELOPE J2 PLUS 23 x 19                |         |
| 2    | E302320     | ENVELOPE J2 PLUS 23 x 20                |         |
| 2    | E302321     | ENVELOPE J2 PLUS 23 x 21                |         |
| 2    | E302322     | ENVELOPE J2 PLUS 23 x 22                |         |
| 2    | E302418     | ENVELOPE J2 PLUS 24 x 18                |         |
| 2    | E302419     | ENVELOPE J2 PLUS 24 x 19                |         |
| 2    | E302420     | ENVELOPE J2 PLUS 24 x 20                |         |
| 2    | E302421     | ENVELOPE J2 PLUS 24 x 21                |         |
| 2    | E302422     | ENVELOPE J2 PLUS 24 x 22                |         |
| 2    | E302518     | ENVELOPE J2 PLUS 25 x 18                |         |
| 2    | E302519     | ENVELOPE J2 PLUS 25 x 19                |         |
| 2    | E302520     | ENVELOPE J2 PLUS 25 x 20                |         |
| 2    | E302521     | ENVELOPE J2 PLUS 25 x 21                |         |
| 2    | E302522     | ENVELOPE J2 PLUS 25 x 22                |         |
| 2    | E302618     | ENVELOPE J2 PLUS 26 x 18                |         |
| 2    | E302619     | ENVELOPE J2 PLUS 26 x 19                |         |
| 2    | E302620     | ENVELOPE J2 PLUS 26 x 20                |         |
| 2    | E302621     | ENVELOPE J2 PLUS 26 x 21                |         |
| 2    | E302622     | ENVELOPE J2 PLUS 26 x 22                |         |
| 3    | B302018     | BASE FOAM J2 PLUS 20 x 18               |         |
| 3    | B302019     | BASE FOAM J2 PLUS 20 x 19               |         |
| 3    | B302020     | BASE FOAM J2 PLUS 20 x 20               |         |
| 3    | B302021     | BASE FOAM J2 PLUS 20 x 21               |         |
| 3    | B302022     | BASE FOAM J2 PLUS 20 x 22               |         |
| 3    | B302118     | BASE FOAM J2 PLUS 21 x 18               |         |
| 3    | B302119     | BASE FOAM J2 PLUS 21 x 19               |         |
| 3    | B302120     | BASE FOAM J2 PLUS 21 x 20               |         |
| 3    | B302121     | BASE FOAM J2 PLUS 21 x 21               |         |
| 3    | B302122     | BASE FOAM J2 PLUS 21 x 22               |         |
| 3    | B302218     | BASE FOAM J2 PLUS 22 x 18               |         |
| 3    | B302219     | BASE FOAM J2 PLUS 22 x 19               |         |
| 3    | B302220     | BASE FOAM J2 PLUS 22 x 20               |         |
| 3    | B302221     | BASE FOAM J2 PLUS 22 x 21               |         |
| 3    | B302222     | BASE FOAM J2 PLUS 22 x 22               |         |
| 3    | B302318     | BASE FOAM J2 PLUS 23 x 18               |         |
| 3    | B302319     | BASE FOAM J2 PLUS 23 x 19               |         |
| 3    | B302320     | BASE FOAM J2 PLUS 23 x 20               |         |
| 3    | B302321     | BASE FOAM J2 PLUS 23 x 21               |         |
| 3    | B302322     | BASE FOAM J2 PLUS 23 x 22               |         |

| Pos. | Item Number | Description                     | Remarks    |
|------|-------------|---------------------------------|------------|
| 3    | B302418     | BASE FOAM J2 PLUS 24 x 18       |            |
| 3    | B302419     | BASE FOAM J2 PLUS 24 x 19       |            |
| 3    | B302420     | BASE FOAM J2 PLUS 24 x 20       |            |
| 3    | B302421     | BASE FOAM J2 PLUS 24 x 21       |            |
| 3    | B302422     | BASE FOAM J2 PLUS 24 x 22       |            |
| 3    | B302518     | BASE FOAM J2 PLUS 25 x 18       |            |
| 3    | B302519     | BASE FOAM J2 PLUS 25 x 19       |            |
| 3    | B302520     | BASE FOAM J2 PLUS 25 x 20       |            |
| 3    | B302521     | BASE FOAM J2 PLUS 25 x 21       |            |
| 3    | B302522     | BASE FOAM J2 PLUS 25 x 22       |            |
| 3    | B302618     | BASE FOAM J2 PLUS 26 x 18       |            |
| 3    | B302619     | BASE FOAM J2 PLUS 26 x 19       |            |
| 3    | B302620     | BASE FOAM J2 PLUS 26 x 20       |            |
| 3    | B302621     | BASE FOAM J2 PLUS 26 x 21       |            |
| 3    | B302622     | BASE FOAM J2 PLUS 26 x 22       |            |
| 4    | F300018     | FLUID PAD INSERT J2 PLUS 18"    |            |
| 4    | F300020     | FLUID PAD INSERT J2 PLUS 20"    |            |
| 4    | F300022     | FLUID PAD INSERT J2 PLUS 22"    |            |
| 5    | 300018-MJ   | WELL FOAM KIT J2 PLUS 18"       |            |
| 5    | 300020-MJ   | WELL FOAM KIT J2 PLUS 20"       |            |
| 5    | 300022-MJ   | WELL FOAM KIT J2 PLUS 22"       |            |
| 5A   | 300118-MJ   | WELL FOAM KIT J2 PLUS 18" C-KIT |            |
| 5A   | 300120-MJ   | WELL FOAM KIT J2 PLUS 20" C-KIT |            |
| 5A   | 300122-MJ   | WELL FOAM KIT J2 PLUS 22" C-KIT |            |
| 6    | S302018     | SOLID SEAT J2 PLUS 20 x 18      | 1/2" THICK |
| 6    | S302019     | SOLID SEAT J2 PLUS 20 x 19      | 1/2" THICK |
| 6    | S302020     | SOLID SEAT J2 PLUS 20 x 20      | 1/2" THICK |
| 6    | S302021     | SOLID SEAT J2 PLUS 20 x 21      | 1/2" THICK |
| 6    | S302022     | SOLID SEAT J2 PLUS 20 x 22      | 1/2" THICK |
| 6    | S302118     | SOLID SEAT J2 PLUS 21 x 18      | 1/2" THICK |
| 6    | S302119     | SOLID SEAT J2 PLUS 21 x 19      | 1/2" THICK |
| 6    | S302120     | SOLID SEAT J2 PLUS 21 x 20      | 1/2" THICK |
| 6    | S302121     | SOLID SEAT J2 PLUS 21 x 21      | 1/2" THICK |
| 6    | S302122     | SOLID SEAT J2 PLUS 21 x 22      | 1/2" THICK |
| 6    | S302218     | SOLID SEAT J2 PLUS 22 x 18      | 1/2" THICK |
| 6    | S302219     | SOLID SEAT J2 PLUS 22 x 19      | 1/2" THICK |
| 6    | S302220     | SOLID SEAT J2 PLUS 22 x 20      | 1/2" THICK |
| 6    | S302221     | SOLID SEAT J2 PLUS 22 x 21      | 1/2" THICK |
| 6    | S302222     | SOLID SEAT J2 PLUS 22 x 22      | 1/2" THICK |
| 6    | S302318     | SOLID SEAT J2 PLUS 23 x 18      | 1/2" THICK |
| 6    | S302319     | SOLID SEAT J2 PLUS 23 x 19      | 1/2" THICK |
| 6    | S302320     | SOLID SEAT J2 PLUS 23 x 20      | 1/2" THICK |
| 6    | S302321     | SOLID SEAT J2 PLUS 23 x 21      | 1/2" THICK |
| 6    | S302322     | SOLID SEAT J2 PLUS 23 x 22      | 1/2" THICK |
| 6    | S302418     | SOLID SEAT J2 PLUS 24 x 18      | 1/2" THICK |
| 6    | S302419     | SOLID SEAT J2 PLUS 24 x 19      | 1/2" THICK |
| 6    | S302420     | SOLID SEAT J2 PLUS 24 x 20      | 1/2" THICK |
| 6    | S302421     | SOLID SEAT J2 PLUS 24 x 21      | 1/2" THICK |
| 6    | S302422     | SOLID SEAT J2 PLUS 24 x 22      | 1/2" THICK |
| 6    | S302518     | SOLID SEAT J2 PLUS 25 x 18      | 1/2" THICK |
| 6    | S302519     | SOLID SEAT J2 PLUS 25 x 19      | 1/2" THICK |
| 6    | S302520     | SOLID SEAT J2 PLUS 25 x 20      | 1/2" THICK |
| 6    | S302521     | SOLID SEAT J2 PLUS 25 x 21      | 1/2" THICK |
| 6    | S302522     | SOLID SEAT J2 PLUS 25 x 22      | 1/2" THICK |
| 6    | S302618     | SOLID SEAT J2 PLUS 26 x 18      | 1/2" THICK |
| 6    | S302619     | SOLID SEAT J2 PLUS 26 x 19      | 1/2" THICK |
| 6    | S302620     | SOLID SEAT J2 PLUS 26 x 20      | 1/2" THICK |

| Pos. | Item Number | Description                | Remarks    |
|------|-------------|----------------------------|------------|
| 6    | S302621     | SOLID SEAT J2 PLUS 26 x 21 | 1/2" THICK |
| 6    | S302622     | SOLID SEAT J2 PLUS 26 x 22 | 1/2" THICK |

NOTE: 1 WHEN REDUCING THE P.L.A. SIZE, A NEW P.L.A. INSERT MUST BE ORDERED SEPARATELY IN THE MATCHING P.L.A. SIZE.

NOTE: 2. 12 THRU 15 WIDTH CUSHIONS ARE STANDARD WITH SIZE "A" P.L.A. RELIEF AND CAN NOT BE REDUCED.

NOTE: 3. 16 THRU 19 WIDTH CUSHIONS ARE STANDARD WITH SIZE "B" P.L.A. RELIEF AND CAN BE REDUCED TO SIZE "A" ONLY.

NOTE: 4. 20 THRU 24 WIDTH CUSHIONS ARE STANDARD WITH SIZE "C" P.L.A. RELIEF AND CAN BE REDUCED TO SIZE "B" OR SIZE "A".

NOTE: 5. WHEN CHANGING THE P.L.A. SIZE FROM "A" TO "C" OR FROM "C" TO "A" A NEW COVER MUST BE ORDERED SEPERATELY.

NOTE: 6. WHEN CHANGING THE P.L.A. SIZE FROM "A" TO "B", FROM "B" TO "A", FROM "B" TO "C", OR FROM "C" TO "B" A NEW COVER IS NOT REQUIRED.

NOTE: 7. SIZE "A" FIELD VARIABLE P.L.A. INSERT COMES WITH QTY (2) MEDIUM SUPPLEMENT FLUID PADS.

NOTE: 8. SIZE "B" FIELD VARIABLE P.L.A. INSERT COMES WITH QTY (2) LARGE SUPPLEMENT FLUID PADS.

NOTE: 9. SIZE "C" FIELD VARIABLE P.L.A. INSERT COMES WITH QTY (2) XLARGE SUPPLEMENT FLUID PADS.

NOTE: 10. 12 THRU 15 WIDTH REPLACEMENT BASE'S WILL ALWAYS SHIP AS SIZE "A" P.L.A. RELIEF.

NOTE: 11. 16 THRU 19 WIDTH REPLACEMENT BASE'S WILL SHIP STANDARD WITH SIZE "B" P.L.A. RELIEF AND WILL INCLUDE THE SIZE "A" P.L.A. REDUCER RING IF SIZE "A" IS ORDERED.

NOTE: 12. 20 THRU 24 WIDTH REPLACEMENT BASE'S WILL SHIP STANDARD WITH SIZE "C" P.L.A. RELIEF AND WILL INCLUDE THE SIZE "B" OR "A" P.L.A. REDUCER RING IF SIZE "B" OR "A" IS ORDERED.



| Pos. | Item Number | Description               | Remarks                                                                                 |
|------|-------------|---------------------------|-----------------------------------------------------------------------------------------|
| 1A   | CJ3C        | J3 CUSHION COVER          | You will be prompted for Height, Width, Depth, P.L.A. Size & Cover Type                 |
| 1B   | MCJ3C       | MODIFIED J3 CUSHION COVER | Order Modified Cushion Cover's Thru Customer Service ONLY.                              |
| 2A   | BJ3C        | J3 CUSHION BASE           | SEE NOTES 10-12 ABOVE -<br>You will be prompted for Height, Width, Depth, & P.L.A. Size |
| 2B   | MBJ3C       | MODIFIED J3 CUSHION BASE  | Order Modified Cushion Base's Thru Customer Service ONLY.                               |
| 3    | J3CARSCA    | STD P.L.A REDUCER C-A     | SEE NOTES 1-6,10-12 ABOVE                                                               |
| 3    | J3CARSCB    | STD P.L.A REDUCER C-B     | SEE NOTES 1-6,10-12 ABOVE                                                               |
| 3    | J3CARSBA    | STD P.L.A REDUCER B-A     | SEE NOTES 1-6,10-12 ABOVE                                                               |
| 3    | J3CARDCA    | DEEP P.L.A REDUCER C-A    | SEE NOTES 1-6,10-12 ABOVE                                                               |
| 3    | J3CARDCB    | DEEP P.L.A REDUCER C-B    | SEE NOTES 1-6,10-12 ABOVE                                                               |
| 3    | J3CARDBA    | DEEP P.L.A REDUCER B-A    | SEE NOTES 1-6,10-12 ABOVE                                                               |

| Pos. | Item Number | Description              | Remarks                                                                                      |
|------|-------------|--------------------------|----------------------------------------------------------------------------------------------|
| 4A   | FJ3C        | J3 P.L.A. INSERT BLADDER | SEE NOTES 7-9 ABOVE -<br>You will be prompted for Height, P.L.A. size, and P.L.A. Fill Type. |
| 4B   | AJ3C        | J3 AIR P.L.A. INSERT     | You will be prompted for height,PLA size and single or dual valve                            |
| 5    | J3CAFSM     | FIELD ADJ SUPPL SM       | SMALL SUPPLEMENT FLUID PAD                                                                   |
| 5    | J3CAFMD     | FIELD ADJ SUPPL MED      | MEDIUM SUPPLEMENT FLUID PAD                                                                  |
| 5    | J3CAFLG     | FIELD ADJ SUPPL LG       | LARGE SUPPLEMENT FLUID PAD                                                                   |
| 5    | J3CAFXL     | FIELD ADJ SUPPL XL       | X-LARGE SUPPLEMENT FLUID PAD                                                                 |
| 6    | J3CAR1212   | J3 SEAT RIGIDIZER 12X12  | 12" WIDTH X 12" DEPTH                                                                        |
| 6    | J3CAR1214   | J3 SEAT RIGIDIZER 12X14  | 12" WIDTH X 14" DEPTH                                                                        |
| 6    | J3CAR1412   | J3 SEAT RIGIDIZER 14X12  | 14" WIDTH X 12" DEPTH                                                                        |
| 6    | J3CAR1414   | J3 SEAT RIGIDIZER 14X14  | 14" WIDTH X 14" DEPTH                                                                        |
| 6    | J3CAR1416   | J3 SEAT RIGIDIZER 14X16  | 14" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR1418   | J3 SEAT RIGIDIZER 14X18  | 14" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR1515   | J3 SEAT RIGIDIZER 15X15  | 15" WIDTH X 15" DEPTH                                                                        |
| 6    | J3CAR1516   | J3 SEAT RIGIDIZER 15X16  | 15" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR1517   | J3 SEAT RIGIDIZER 15X17  | 15" WIDTH X 17" DEPTH                                                                        |
| 6    | J3CAR1615   | J3 SEAT RIGIDIZER 16X15  | 16" WIDTH X 15" DEPTH                                                                        |
| 6    | J3CAR1616   | J3 SEAT RIGIDIZER 16X16  | 16" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR1617   | J3 SEAT RIGIDIZER 16X17  | 16" WIDTH X 17" DEPTH                                                                        |
| 6    | J3CAR1618   | J3 SEAT RIGIDIZER 16X18  | 16" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR1620   | J3 SEAT RIGIDIZER 16X20  | 16" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR1715   | J3 SEAT RIGIDIZER 17X15  | 17" WIDTH X 15" DEPTH                                                                        |
| 6    | J3CAR1716   | J3 SEAT RIGIDIZER 17X16  | 17" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR1717   | J3 SEAT RIGIDIZER 17X17  | 17" WIDTH X 17" DEPTH                                                                        |
| 6    | J3CAR1718   | J3 SEAT RIGIDIZER 17X18  | 17" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR1816   | J3 SEAT RIGIDIZER 18X16  | 18" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR1817   | J3 SEAT RIGIDIZER 18X17  | 18" WIDTH X 17" DEPTH                                                                        |
| 6    | J3CAR1818   | J3 SEAT RIGIDIZER 18X18  | 18" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR1820   | J3 SEAT RIGIDIZER 18X20  | 18" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR1916   | J3 SEAT RIGIDIZER 19X16  | 19" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR1917   | J3 SEAT RIGIDIZER 19X17  | 19" WIDTH X 17" DEPTH                                                                        |
| 6    | J3CAR1918   | J3 SEAT RIGIDIZER 19X18  | 19" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR1920   | J3 SEAT RIGIDIZER 19X20  | 19" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR2016   | J3 SEAT RIGIDIZER 20X16  | 20" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR2018   | J3 SEAT RIGIDIZER 20X18  | 20" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR2020   | J3 SEAT RIGIDIZER 20X20  | 20" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR2118   | J3 SEAT RIGIDIZER 21X18  | 21" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR2119   | J3 SEAT RIGIDIZER 21X19  | 21" WIDTH X 19" DEPTH                                                                        |
| 6    | J3CAR2120   | J3 SEAT RIGIDIZER 21X20  | 21" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR2121   | J3 SEAT RIGIDIZER 21X21  | 21" WIDTH X 21" DEPTH                                                                        |
| 6    | J3CAR2122   | J3 SEAT RIGIDIZER 21X22  | 21" WIDTH X 22" DEPTH                                                                        |
| 6    | J3CAR2218   | J3 SEAT RIGIDIZER 22X18  | 22" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR2219   | J3 SEAT RIGIDIZER 22X19  | 22" WIDTH X 19" DEPTH                                                                        |
| 6    | J3CAR2220   | J3 SEAT RIGIDIZER 22X20  | 22" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR2221   | J3 SEAT RIGIDIZER 22X21  | 22" WIDTH X 21" DEPTH                                                                        |
| 6    | J3CAR2222   | J3 SEAT RIGIDIZER 22X22  | 22" WIDTH X 22" DEPTH                                                                        |
| 6    | J3CAR2318   | J3 SEAT RIGIDIZER 23X18  | 23" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR2319   | J3 SEAT RIGIDIZER 23X19  | 23" WIDTH X 19" DEPTH                                                                        |
| 6    | J3CAR2320   | J3 SEAT RIGIDIZER 23X20  | 23" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR2321   | J3 SEAT RIGIDIZER 23X21  | 23" WIDTH X 21" DEPTH                                                                        |
| 6    | J3CAR2322   | J3 SEAT RIGIDIZER 23X22  | 23" WIDTH X 22" DEPTH                                                                        |
| 6    | J3CAR2418   | J3 SEAT RIGIDIZER 24X18  | 24" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR2419   | J3 SEAT RIGIDIZER 24X19  | 24" WIDTH X 19" DEPTH                                                                        |
| 6    | J3CAR2420   | J3 SEAT RIGIDIZER 24X20  | 24" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR2421   | J3 SEAT RIGIDIZER 24X21  | 24" WIDTH X 21" DEPTH                                                                        |
| 6    | J3CAR2422   | J3 SEAT RIGIDIZER 24X22  | 24" WIDTH X 22" DEPTH                                                                        |

| Pos. | Item Number | Description             | Remarks               |
|------|-------------|-------------------------|-----------------------|
| 6    | J3CAR2424   | J3 SEAT RIGIDIZER 24X24 | 24" WIDTH X 24" DEPTH |
| 7    | PATCHKIT    | JAY AIR PATCH KIT       |                       |
| 8    | AIRPUMP     | JAY AIR PUMP            |                       |



| Pos. | Item Number | Description               | Remarks                                                                                 |
|------|-------------|---------------------------|-----------------------------------------------------------------------------------------|
| 1A   | CJ3CP       | J3 CUSHION COVER          | You will be prompted for Height, Width, Depth, P.L.A. Size & Cover Type                 |
| 1B   | MCJ3CP      | MODIFIED J3 CUSHION COVER | Order Modified Cushion Cover's Thru Customer Service ONLY.                              |
| 2A   | BJ3CP       | J3 CUSHION BASE           | SEE NOTES 10-12 ABOVE -<br>You will be prompted for Height, Width, Depth, & P.L.A. Size |
| 2B   | MBJ3CP      | MODIFIED J3 CUSHION BASE  | Order Modified Cushion Base's Thru Customer Service ONLY.                               |
| 3    | J3CARSCA    | STD P.L.A REDUCER C-A     | SEE NOTES 1-6,10-12 ABOVE                                                               |
| 3    | J3CARSCB    | STD P.L.A REDUCER C-B     | SEE NOTES 1-6,10-12 ABOVE                                                               |
| 3    | J3CARSBA    | STD P.L.A REDUCER B-A     | SEE NOTES 1-6,10-12 ABOVE                                                               |
| 3    | J3CARDCA    | DEEP P.L.A REDUCER C-A    | SEE NOTES 1-6,10-12 ABOVE                                                               |
| 3    | J3CARD CB   | DEEP P.L.A REDUCER C-B    | SEE NOTES 1-6,10-12 ABOVE                                                               |

| Pos. | Item Number | Description              | Remarks                                                                                      |
|------|-------------|--------------------------|----------------------------------------------------------------------------------------------|
| 3    | J3CARDBA    | DEEP P.L.A REDUCER B-A   | SEE NOTES 1-6,10-12 ABOVE                                                                    |
| 4A   | FJ3C        | J3 P.L.A. INSERT BLADDER | SEE NOTES 7-9 ABOVE -<br>You will be prompted for Height, P.L.A. size, and P.L.A. Fill Type. |
| 4B   | AJ3C        | J3 AIR P.L.A. INSERT     | You will be prompted for height,PLA size and single or dual valve                            |
| 5    | J3CAFSM     | FIELD ADJ SUPPL SM       | SMALL SUPPLEMENT FLUID PAD                                                                   |
| 5    | J3CAFMD     | FIELD ADJ SUPPL MED      | MEDIUM SUPPLEMENT FLUID PAD                                                                  |
| 5    | J3CAFLG     | FIELD ADJ SUPPL LG       | LARGE SUPPLEMENT FLUID PAD                                                                   |
| 5    | J3CAFXL     | FIELD ADJ SUPPL XL       | X-LARGE SUPPLEMENT FLUID PAD                                                                 |
| 6    | J3CAR1212   | J3 SEAT RIGIDIZER 12X12  | 12" WIDTH X 12" DEPTH                                                                        |
| 6    | J3CAR1214   | J3 SEAT RIGIDIZER 12X14  | 12" WIDTH X 14" DEPTH                                                                        |
| 6    | J3CAR1412   | J3 SEAT RIGIDIZER 14X12  | 14" WIDTH X 12" DEPTH                                                                        |
| 6    | J3CAR1414   | J3 SEAT RIGIDIZER 14X14  | 14" WIDTH X 14" DEPTH                                                                        |
| 6    | J3CAR1416   | J3 SEAT RIGIDIZER 14X16  | 14" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR1418   | J3 SEAT RIGIDIZER 14X18  | 14" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR1515   | J3 SEAT RIGIDIZER 15X15  | 15" WIDTH X 15" DEPTH                                                                        |
| 6    | J3CAR1516   | J3 SEAT RIGIDIZER 15X16  | 15" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR1517   | J3 SEAT RIGIDIZER 15X17  | 15" WIDTH X 17" DEPTH                                                                        |
| 6    | J3CAR1615   | J3 SEAT RIGIDIZER 16X15  | 16" WIDTH X 15" DEPTH                                                                        |
| 6    | J3CAR1616   | J3 SEAT RIGIDIZER 16X16  | 16" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR1617   | J3 SEAT RIGIDIZER 16X17  | 16" WIDTH X 17" DEPTH                                                                        |
| 6    | J3CAR1618   | J3 SEAT RIGIDIZER 16X18  | 16" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR1620   | J3 SEAT RIGIDIZER 16X20  | 16" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR1715   | J3 SEAT RIGIDIZER 17X15  | 17" WIDTH X 15" DEPTH                                                                        |
| 6    | J3CAR1716   | J3 SEAT RIGIDIZER 17X16  | 17" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR1717   | J3 SEAT RIGIDIZER 17X17  | 17" WIDTH X 17" DEPTH                                                                        |
| 6    | J3CAR1718   | J3 SEAT RIGIDIZER 17X18  | 17" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR1816   | J3 SEAT RIGIDIZER 18X16  | 18" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR1817   | J3 SEAT RIGIDIZER 18X17  | 18" WIDTH X 17" DEPTH                                                                        |
| 6    | J3CAR1818   | J3 SEAT RIGIDIZER 18X18  | 18" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR1820   | J3 SEAT RIGIDIZER 18X20  | 18" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR1916   | J3 SEAT RIGIDIZER 19X16  | 19" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR1917   | J3 SEAT RIGIDIZER 19X17  | 19" WIDTH X 17" DEPTH                                                                        |
| 6    | J3CAR1918   | J3 SEAT RIGIDIZER 19X18  | 19" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR1920   | J3 SEAT RIGIDIZER 19X20  | 19" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR2016   | J3 SEAT RIGIDIZER 20X16  | 20" WIDTH X 16" DEPTH                                                                        |
| 6    | J3CAR2018   | J3 SEAT RIGIDIZER 20X18  | 20" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR2020   | J3 SEAT RIGIDIZER 20X20  | 20" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR2118   | J3 SEAT RIGIDIZER 21X18  | 21" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR2119   | J3 SEAT RIGIDIZER 21X19  | 21" WIDTH X 19" DEPTH                                                                        |
| 6    | J3CAR2120   | J3 SEAT RIGIDIZER 21X20  | 21" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR2121   | J3 SEAT RIGIDIZER 21X21  | 21" WIDTH X 21" DEPTH                                                                        |
| 6    | J3CAR2122   | J3 SEAT RIGIDIZER 21X22  | 21" WIDTH X 22" DEPTH                                                                        |
| 6    | J3CAR2218   | J3 SEAT RIGIDIZER 22X18  | 22" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR2219   | J3 SEAT RIGIDIZER 22X19  | 22" WIDTH X 19" DEPTH                                                                        |
| 6    | J3CAR2220   | J3 SEAT RIGIDIZER 22X20  | 22" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR2221   | J3 SEAT RIGIDIZER 22X21  | 22" WIDTH X 21" DEPTH                                                                        |
| 6    | J3CAR2222   | J3 SEAT RIGIDIZER 22X22  | 22" WIDTH X 22" DEPTH                                                                        |
| 6    | J3CAR2318   | J3 SEAT RIGIDIZER 23X18  | 23" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR2319   | J3 SEAT RIGIDIZER 23X19  | 23" WIDTH X 19" DEPTH                                                                        |
| 6    | J3CAR2320   | J3 SEAT RIGIDIZER 23X20  | 23" WIDTH X 20" DEPTH                                                                        |
| 6    | J3CAR2321   | J3 SEAT RIGIDIZER 23X21  | 23" WIDTH X 21" DEPTH                                                                        |
| 6    | J3CAR2322   | J3 SEAT RIGIDIZER 23X22  | 23" WIDTH X 22" DEPTH                                                                        |
| 6    | J3CAR2418   | J3 SEAT RIGIDIZER 24X18  | 24" WIDTH X 18" DEPTH                                                                        |
| 6    | J3CAR2419   | J3 SEAT RIGIDIZER 24X19  | 24" WIDTH X 19" DEPTH                                                                        |

| Pos. | Item Number | Description                  | Remarks                           |
|------|-------------|------------------------------|-----------------------------------|
| 6    | J3CAR2420   | J3 SEAT RIGIDIZER 24X20      | 24" WIDTH X 20" DEPTH             |
| 6    | J3CAR2421   | J3 SEAT RIGIDIZER 24X21      | 24" WIDTH X 21" DEPTH             |
| 6    | J3CAR2422   | J3 SEAT RIGIDIZER 24X22      | 24" WIDTH X 22" DEPTH             |
| 6    | J3CAR2424   | J3 SEAT RIGIDIZER 24X24      | 24" WIDTH X 24" DEPTH             |
| 7    | J3CAPS2     | LATPEL 2"H (SMALL)           | Pair- Used on 12 deep to 18 deep  |
| 7    | J3CAPM2     | LATPEL 2"H (MEDIUM)          | Pair - Used on 19 deep to 24 deep |
| 8    | J3CAMS2     | MEDTHI 2"H (CUSHIONS12W-14W) | Used on 12 deep to 16 deep        |
| 8    | J3CMM2      | MEDTHI 2"H (CUSHIONS16W-19W) | Used on 17 deep to 19 deep        |
| 8    | J3CAML2     | MEDTHI 2"H (CUSHIONS20W-24W) | Used on 20 deep to 24 deep        |
| 9    | J3CALS2     | LATTHI 2"H (CUSHIONS12W-14W) | Pair- Used on 12 deep to 16 deep  |
| 9    | J3CALM2     | LATTHI 2"H (CUSHIONS16W-19W) | Pair - Used on 17 deep to 19 deep |
| 9    | J3CALL2     | LATTHI 2"H (CUSHIONS20W-24W) | Pair - Used on 20 deep to 24 deep |
| 10   | PATCHKIT    | JAY AIR PATCH KIT            |                                   |
| 11   | AIRPUMP     | JAY AIR PUMP                 |                                   |

NOTE: 1 WHEN REDUCING THE P.L.A. SIZE, A NEW P.L.A. INSERT MUST BE ORDERED SEPARATELY IN THE MATCHING P.L.A. SIZE.

NOTE: 2. 12 THRU 15 WIDTH CUSHIONS ARE STANDARD WITH SIZE "A" P.L.A. RELIEF AND CAN NOT BE REDUCED.

NOTE: 3. 16 THRU 19 WIDTH CUSHIONS ARE STANDARD WITH SIZE "B" P.L.A. RELIEF AND CAN BE REDUCED TO SIZE "A" ONLY.

NOTE: 4. 20 THRU 24 WIDTH CUSHIONS ARE STANDARD WITH SIZE "C" P.L.A. RELIEF AND CAN BE REDUCED TO SIZE "B" OR SIZE "A".

NOTE: 5. WHEN CHANGING THE P.L.A. SIZE FROM "A" TO "C" OR FROM "C" TO "A" A NEW COVER MUST BE ORDERED SEPERATELY.

NOTE: 6. WHEN CHANGING THE P.L.A. SIZE FROM "A" TO "B", FROM "B" TO "A", FROM "B" TO "C", OR FROM "C" TO "B" A NEW COVER IS NOT REQUIRED.

NOTE: 7. SIZE "A" FIELD VARIABLE P.L.A. INSERT COMES WITH QTY (2) MEDIUM SUPPLEMENT FLUID PADS.

NOTE: 8. SIZE "B" FIELD VARIABLE P.L.A. INSERT COMES WITH QTY (2) LARGE SUPPLEMENT FLUID PADS.


NOTE: 9. SIZE "C" FIELD VARIABLE P.L.A. INSERT COMES WITH QTY (2) XLARGE SUPPLEMENT FLUID PADS.

[09/2010]


NOTE: 10. 12 THRU 15 WIDTH REPLACEMENT BASE'S WILL ALWAYS SHIP AS SIZE "A" P.L.A. RELIEF.

NOTE: 11. 16 THRU 19 WIDTH REPLACEMENT BASE'S WILL SHIP STANDARD WITH SIZE "B" P.L.A. RELIEF AND WILL INCLUDE THE SIZE "A" P.L.A.REDUCER RING IF SIZE "A" IS ORDERED.

NOTE: 12. 20 THRU 24 WIDTH REPLACEMENT BASE'S WILL SHIP STANDARD WITH SIZE "C" P.L.A. RELIEF AND WILL INCLUDE THE SIZE "B" OR "A" P.L.A. REDUCER RING IF SIZE "B" OR "A" IS ORDERED.



| Pos.        | Item Number | Description              |
|-------------|-------------|--------------------------|
| 1,2A,3A,4-7 | J3CEKS      | J3 CUSH STD PLA EVAL KIT |
| 1,2B,3B,4-7 | J3CEKD      | J3 CUSH DP PLA EVAL KIT  |



| Pos. | Item Number | Description                    | Remarks        |
|------|-------------|--------------------------------|----------------|
| 2    | J3CALM1     | LATERAL THIGH 1" HEIGHT MEDIUM | 8.25" LONG     |
| 2    | J3CALL1     | LATERAL THIGH 1" HEIGHT LARGE  | 11" LONG       |
| 2    | J3CALS2     | LATERAL THIGH 2" HEIGHT SMALL  | 5.5" LONG      |
| 2    | J3CALM2     | LATERAL THIGH 2" HEIGHT MEDIUM | 8.25" LONG     |
| 2    | J3CALL2     | LATERAL THIGH 2" HEIGHT LARGE  | 11" LONG       |
| 3    | J3CAMP1     | MED THIGH 1" HEIGHT MEDIUM     | 6" LONG        |
| 3    | J3CAML1     | MED THIGH 1" HEIGHT LARGE      | 8.5" LONG      |
| 3    | J3CAMS2     | MED THIGH 2" HEIGHT SMALL      | 3.5" LONG      |
| 3    | J3CAMP2     | MED THIGH 2" HEIGHT MEDIUM     | 6" LONG        |
| 3    | J3CAML2     | MED THIGH 2" HEIGHT LARGE      | 8.5" LONG      |
| 4    | J3CAPS2     | LAT PEL 2" HEIGHT SMALL        | 4.75" LONG     |
| 4    | J3CAPM2     | LAT PEL 2" HEIGHT MEDIUM       | 6.25" LONG     |
| 5    | J3CSAOR     | OBLIQUITY PAD STD A RIGHT      | 1/2" in height |
| 5    | J3CSAOL     | OBLIQUITY PAD STD A LEFT       | 1/2" in height |
| 5    | J3CSBOR     | OBLIQUITY PAD STD B RIGHT      | 1/2" in height |
| 5    | J3CSBOL     | OBLIQUITY PAD STD B LEFT       | 1/2" in height |
| 5    | J3CSCOR     | OBLIQUITY PAD STD C RIGHT      | 1/2" in height |
| 5    | J3CSCOL     | OBLIQUITY PAD STD C LEFT       | 1/2" in height |
| 5    | J3CDAOR     | OBLIQUITY PAD DP A RIGHT       | 1/2" in height |
| 5    | J3CDAOL     | OBLIQUITY PAD DP A LEFT        | 1/2" in height |
| 5    | J3CDBOR     | OBLIQUITY PAD DP B RIGHT       | 1/2" in height |

| Pos. | Item Number | Description                     | Remarks        |
|------|-------------|---------------------------------|----------------|
| 5    | J3CDBOL     | OBLIQUITY PAD DP B LEFT         | 1/2" in height |
| 5    | J3CDCOR     | OBLIQUITY PAD DP C RIGHT        | 1/2" in height |
| 5    | J3CDCOL     | OBLIQUITY PAD DP C LEFT         | 1/2" in height |
| 6    | J3CPRT      | REPLACEMENT TAPE 30' ROLL       |                |
| 1-6  | J3CPEK      | J3 CUSHION POSITIONING EVAL KIT |                |



| Pos. | Item Number | Description      | Remarks                             |
|------|-------------|------------------|-------------------------------------|
| >>>> |             | CUSHION COVERS   |                                     |
| N/A  | 109422-xxx  | J3 CUSHION COVER | PLEASE SELL<br>"CJ3C" OR<br>"CJ3CP" |
| N/A  | 109423-xxx  | J3 CUSHION COVER | PLEASE SELL<br>"CJ3C" OR<br>"CJ3CP" |
| N/A  | 109424-xxx  | J3 CUSHION COVER | PLEASE SELL<br>"CJ3C" OR<br>"CJ3CP" |
| N/A  | 109425-xxx  | J3 CUSHION COVER | PLEASE SELL<br>"CJ3C" OR<br>"CJ3CP" |


| Pos. | Item Number | Description                    | Remarks                             |
|------|-------------|--------------------------------|-------------------------------------|
| N/A  | 109426-xxx  | J3 CUSHION COVER               | PLEASE SELL<br>"CJ3C" OR<br>"CJ3CP" |
| N/A  | 109427-xxx  | J3 CUSHION COVER               | PLEASE SELL<br>"CJ3C" OR<br>"CJ3CP" |
| >>>> |             | FLUID P.L.A. INSERT            |                                     |
| N/A  | J3CSAFF     | J3 P.L.A. INSERT FLUID BLADDER | PLEASE SELL<br>"FJ3C"               |
| N/A  | J3CSBFF     | J3 P.L.A. INSERT FLUID BLADDER | PLEASE SELL<br>"FJ3C"               |
| N/A  | J3CSCFF     | J3 P.L.A. INSERT FLUID BLADDER | PLEASE SELL<br>"FJ3C"               |
| N/A  | J3CSAFV     | J3 P.L.A. INSERT FLUID BLADDER | PLEASE SELL<br>"FJ3C"               |
| N/A  | J3CSBFV     | J3 P.L.A. INSERT FLUID BLADDER | PLEASE SELL<br>"FJ3C"               |
| N/A  | J3CSCFV     | J3 P.L.A. INSERT FLUID BLADDER | PLEASE SELL<br>"FJ3C"               |
| N/A  | J3CDAFF     | J3 P.L.A. INSERT FLUID BLADDER | PLEASE SELL<br>"FJ3C"               |
| N/A  | J3CDBFF     | J3 P.L.A. INSERT FLUID BLADDER | PLEASE SELL<br>"FJ3C"               |
| N/A  | J3CDCFF     | J3 P.L.A. INSERT FLUID BLADDER | PLEASE SELL<br>"FJ3C"               |
| N/A  | J3CDAFV     | J3 P.L.A. INSERT FLUID BLADDER | PLEASE SELL<br>"FJ3C"               |
| N/A  | J3CDBFV     | J3 P.L.A. INSERT FLUID BLADDER | PLEASE SELL<br>"FJ3C"               |
| >>>> |             | AIR P.L.A. INSERT              |                                     |
| N/A  | J3CDCFV     | J3 P.L.A. INSERT AIR BLADDER   | PLEASE SELL<br>"FJ3C"               |
| N/A  | J3CSAAS     | J3 P.L.A. INSERT AIR BLADDER   | PLEASE SELL<br>"AJ3C"               |
| N/A  | J3CSBAS     | J3 P.L.A. INSERT AIR BLADDER   | PLEASE SELL<br>"AJ3C"               |
| N/A  | J3CSCAS     | J3 P.L.A. INSERT AIR BLADDER   | PLEASE SELL<br>"AJ3C"               |
| N/A  | J3CSAAD     | J3 P.L.A. INSERT AIR BLADDER   | PLEASE SELL<br>"AJ3C"               |
| N/A  | J3CSBAD     | J3 P.L.A. INSERT AIR BLADDER   | PLEASE SELL<br>"AJ3C"               |
| N/A  | J3CSCAD     | J3 P.L.A. INSERT AIR BLADDER   | PLEASE SELL<br>"AJ3C"               |
| N/A  | J3CDAAS     | J3 P.L.A. INSERT AIR BLADDER   | PLEASE SELL<br>"AJ3C"               |
| N/A  | J3CDBAS     | J3 P.L.A. INSERT AIR BLADDER   | PLEASE SELL<br>"AJ3C"               |

| Pos. | Item Number | Description                                            | Remarks                                       |
|------|-------------|--------------------------------------------------------|-----------------------------------------------|
| N/A  | J3CDCAS     | J3 P.L.A. INSERT AIR BLADDER                           | PLEASE SELL "AJ3C"                            |
| N/A  | J3CDAAD     | J3 P.L.A. INSERT AIR BLADDER                           | PLEASE SELL "AJ3C"                            |
| N/A  | J3CDBAD     | J3 P.L.A. INSERT AIR BLADDER                           | PLEASE SELL "AJ3C"                            |
| N/A  | J3CDCAD     | J3 P.L.A. INSERT AIR BLADDER                           | PLEASE SELL "AJ3C"                            |
| >>>> |             | SUPPLEMENTAL FLUID PADS                                |                                               |
| N/A  | 090509      | SUPPLEMENTAL FLUID PAD SMALL                           | PLEASE SELL "J3CAFSM"                         |
| N/A  | 012010      | SUPPLEMENTAL FLUID PAD MEDIUM                          | PLEASE SELL "J3CAFMD"                         |
| N/A  | 010510      | SUPPLEMENTAL FLUID PAD LARGE                           | PLEASE SELL "J3CAFLG"                         |
| N/A  | 091409      | SUPPLEMENTAL FLUID PAD X-LARGE                         | PLEASE SELL "J3CAFXL"                         |
| >>>> |             | LATERAL THIGH SUPPORTS (PAIR)                          |                                               |
| N/A  | 108821-001  | LATERAL THIGH SUPPORTS 1" HIGH (CUSHIONS12W-14W) LEFT  | PLEASE SELL "J3CALS1" (INCLUDES LEFT & RIGHT) |
| N/A  | 108822-001  | LATERAL THIGH SUPPORTS 1" HIGH (CUSHIONS12W-14W) RIGHT | PLEASE SELL "J3CALS1" (INCLUDES LEFT & RIGHT) |
| N/A  | 108821-002  | LATERAL THIGH SUPPORTS 1" HIGH (CUSHIONS16W-19W) LEFT  | PLEASE SELL "J3CALM1" (INCLUDES LEFT & RIGHT) |
| N/A  | 108822-002  | LATERAL THIGH SUPPORTS 1" HIGH (CUSHIONS16W-19W) RIGHT | PLEASE SELL "J3CALM1" (INCLUDES LEFT & RIGHT) |
| N/A  | 108821-003  | LATERAL THIGH SUPPORTS 1" HIGH (CUSHIONS20W-24W) LEFT  | PLEASE SELL "J3CALL1" (INCLUDES LEFT & RIGHT) |
| N/A  | 108822-003  | LATERAL THIGH SUPPORTS 1" HIGH (CUSHIONS20W-24W) RIGHT | PLEASE SELL "J3CALL1" (INCLUDES LEFT & RIGHT) |
| N/A  | 108823-001  | LATERAL THIGH SUPPORTS 2" HIGH (CUSHIONS12W-14W) LEFT  | PLEASE SELL "J3CALS2" (INCLUDES LEFT & RIGHT) |
| N/A  | 108824-001  | LATERAL THIGH SUPPORTS 2" HIGH (CUSHIONS12W-14W) RIGHT | PLEASE SELL "J3CALS2" (INCLUDES LEFT & RIGHT) |

| Pos. | Item Number | Description                                            | Remarks                                                |
|------|-------------|--------------------------------------------------------|--------------------------------------------------------|
| N/A  | 108823-002  | LATERAL THIGH SUPPORTS 2" HIGH (CUSHIONS16W-19W) LEFT  | PLEASE SELL<br>"J3CALM2"<br>(INCLUDES LEFT<br>& RIGHT) |
| N/A  | 108824-002  | LATERAL THIGH SUPPORTS 2" HIGH (CUSHIONS16W-19W) RIGHT | PLEASE SELL<br>"J3CALM2"<br>(INCLUDES LEFT<br>& RIGHT) |
| N/A  | 108823-003  | LATERAL THIGH SUPPORTS 2" HIGH (CUSHIONS20W-24W) LEFT  | PLEASE SELL<br>"J3CALL2"<br>(INCLUDES LEFT<br>& RIGHT) |
| N/A  | 108824-003  | LATERAL THIGH SUPPORTS 2" HIGH (CUSHIONS20W-24W) RIGHT | PLEASE SELL<br>"J3CALL2"<br>(INCLUDES LEFT<br>& RIGHT) |
| >>>> |             | MEDIAL THIGH SUPPORTS (EACH)                           |                                                        |
| N/A  | 108825-001  | MEDIAL THIGH SUPPORT 1" HIGH (CUSHIONS 12W-14W)        | PLEASE SELL<br>"J3CAMS1"<br>(EACH)                     |
| N/A  | 108825-002  | MEDIAL THIGH SUPPORT 1" HIGH (CUSHIONS 16W-19W)        | PLEASE SELL<br>"J3CMM1"<br>(EACH)                      |
| N/A  | 108825-003  | MEDIAL THIGH SUPPORT 1" HIGH (CUSHIONS 20W-24W)        | PLEASE SELL<br>"J3CAML1" (EACH)                        |
| N/A  | 108826-001  | MEDIAL THIGH SUPPORT 2" HIGH (CUSHIONS 12W-14W)        | PLEASE SELL<br>"J3CAMS2"<br>(EACH)                     |
| N/A  | 108826-002  | MEDIAL THIGH SUPPORT 2" HIGH (CUSHIONS 16W-19W)        | PLEASE SELL<br>"J3CMM2"<br>(EACH)                      |
| N/A  | 108826-003  | MEDIAL THIGH SUPPORT 2" HIGH (CUSHIONS 20W-24W)        | PLEASE SELL<br>"J3CAML2" (EACH)                        |
| >>>> |             | LATERAL PELVIC SUPPORTS (PAIR)                         |                                                        |
| N/A  | 108227-001  | LATERAL PELVIC SUPPORTS 2" HIGH (CUSHIONS 12W-14W)     | PLEASE SELL<br>"J3CAPS2"<br>(INCLUDES LEFT<br>& RIGHT) |
| N/A  | 108827-002  | LATERAL PELVIC SUPPORTS 2" HIGH (CUSHIONS 16W-19W)     | PLEASE SELL<br>"J3CAPM2"<br>(INCLUDES LEFT<br>& RIGHT) |
| >>>> |             | P.L.A. REDUCER                                         |                                                        |
| N/A  | 109368      | STANDARD P.L.A REDUCER - SIZE B TO A                   | PLEASE SELL<br>"J3CARsBA"                              |
| N/A  | 109440      | STANDARD P.L.A REDUCER - SIZE C TO A                   | PLEASE SELL<br>"J3CARsCA"                              |
| N/A  | 109369      | STANDARD P.L.A REDUCER - SIZE C TO B                   | PLEASE SELL<br>"J3CARsCB"                              |
| N/A  | 109370      | DEEP P.L.A REDUCER - SIZE B TO A                       | PLEASE SELL<br>"J3CARDBA"                              |

| Pos. | Item Number | Description                      | Remarks                    |
|------|-------------|----------------------------------|----------------------------|
| N/A  | 109441      | DEEP P.L.A REDUCER - SIZE C TO A | PLEASE SELL<br>"J3CARDCA"  |
| N/A  | 109371      | DEEP P.L.A REDUCER - SIZE C TO B | PLEASE SELL<br>"J3CARD CB" |
| >>>> |             | OBLIQUITY PADS                   |                            |
| N/A  | 108830-001  | OBLIQUITY PAD STANDARD A RIGHT   | PLEASE SELL<br>"J3CSAOR"   |
| N/A  | 108831-001  | OBLIQUITY PAD STANDARD A LEFT    | PLEASE SELL<br>"J3CSAOL"   |
| N/A  | 108830-002  | OBLIQUITY PAD STANDARD B RIGHT   | PLEASE SELL<br>"J3CSBOR"   |
| N/A  | 108831-002  | OBLIQUITY PAD STANDARD B LEFT    | PLEASE SELL<br>"J3CSBOL"   |
| N/A  | 108830-003  | OBLIQUITY PAD STANDARD C RIGHT   | PLEASE SELL<br>"J3CSCOR"   |
| N/A  | 108831-003  | OBLIQUITY PAD STANDARD C LEFT    | PLEASE SELL<br>"J3CSCOL"   |
| N/A  | 108846-001  | OBLIQUITY PAD DEEP A RIGHT       | PLEASE SELL<br>"J3CDAOR"   |
| N/A  | 108845-001  | OBLIQUITY PAD DEEP A LEFT        | PLEASE SELL<br>"J3CDAOL"   |
| N/A  | 108846-002  | OBLIQUITY PAD DEEP B RIGHT       | PLEASE SELL<br>"J3CDBOR"   |
| N/A  | 108845-002  | OBLIQUITY PAD DEEP B LEFT        | PLEASE SELL<br>"J3CDBOL"   |
| N/A  | 108846-003  | OBLIQUITY PAD DEEP C RIGHT       | PLEASE SELL<br>"J3CDCOR"   |
| N/A  | 108845-003  | OBLIQUITY PAD DEEP C LEFT        | PLEASE SELL<br>"J3CDCOL"   |

Note: Positioning attachments are 1 inch tall and are not replaceable.





| Pos. | Item Number | Description                  | Remarks |
|------|-------------|------------------------------|---------|
| 1A   | CJL1414     | COVER JAY LITE 1414IN 3535CM |         |
| 1A   | CJL1416     | COVER JAY LITE 1416IN 3535CM |         |
| 1A   | CJL1418     | COVER JAY LITE 1418IN 3535CM |         |
| 1A   | CJL1515     | COVER JAY LITE 1515IN        |         |
| 1A   | CJL1516     | COVER JAY LITE 1516IN 3840CM |         |
| 1A   | CJL1517     | COVER JAY LITE 1517IN 3844CM |         |
| 1A   | CJL1518     | COVER JAY LITE 1518IN        |         |
| 1A   | CJL1615     | COVER JAY LITE 1615IN 4038CM |         |
| 1A   | CJL1616     | COVER JAY LITE 1616IN 4040CM |         |
| 1A   | CJL1617     | COVER JAY LITE 1617IN 4043CM |         |
| 1A   | CJL1618     | COVER JAY LITE 1618IN 4045CM |         |
| 1A   | CJL1620     | COVER JAY LITE 1620IN 4050CM |         |
| 1A   | CJL1715     | COVER JAY LITE 1715IN 4338CM |         |
| 1A   | CJL1716     | COVER JAY LITE 1716IN 4340   |         |

| Pos. | Item Number | Description                        | Remarks |
|------|-------------|------------------------------------|---------|
| 1A   | CJL1717     | COVER JAY LITE 1717IN 4343CM       |         |
| 1A   | CJL1718     | COVER JAY LITE 1718IN 4345CM       |         |
| 1A   | CJL1816     | COVER JAY LITE 1816IN 4540CM       |         |
| 1A   | CJL1817     | COVER JAY LITE 1817IN 4543CM       |         |
| 1A   | CJL1818     | COVER JAY LITE 1818IN 4545CM       |         |
| 1A   | CJL1820     | COVER JAY LITE 1820IN 4550CM       |         |
| 1A   | CJL1916     | COVER JAY LITE 1916IN 4840CM       |         |
| 1A   | CJL1917     | COVER JAY LITE 1917IN 4843CM       |         |
| 1A   | CJL1918     | COVER JAY LITE 1918IN 4845CM       |         |
| 1A   | CJL1920     | COVER JAY LITE 1920IN 4850CM       |         |
| 1A   | CJL2016     | COVER JAY LITE 2016IN 5040CM       |         |
| 1A   | CJL2018     | COVER JAY LITE 2018IN 5045CM       |         |
| 1A   | CJL2020     | COVER JAY LITE 2020IN 5050CM       |         |
| 1A   | CJL2218     | COVER JAY LITE 2218IN 5545CM       |         |
| 1A   | CJL2220     | COVER JAY LITE 2220IN 5550CM       |         |
| 1A   | CJL2222     | COVER JAY LITE 2222IN 5555CM       |         |
| 1A   | CJL2418     | COVER JAY LITE 2418IN 6045CM       |         |
| 1A   | CJL2420     | COVER JAY LITE 2420IN 6050CM       |         |
| 1A   | CJL2422     | COVER JAY LITE 2422IN 6055CM       |         |
| 1A   | CJL2424     | COVER JAY LITE 2424IN 6060CM       |         |
| 1B   | CJL1414P    | COVER, POS, JAY LITE 1414IN 3535CM |         |
| 1B   | CJL1416P    | COVER, POS, JAY LITE 1416IN 3535CM |         |
| 1B   | CJL1418P    | COVER, POS, JAY LITE 1418IN 3535CM |         |
| 1B   | CJL1515P    | COVER, POS, JAY LITE 1515IN        |         |
| 1B   | CJL1516P    | COVER, POS, JAY LITE 1516IN 3840CM |         |
| 1B   | CJL1517P    | COVER, POS, JAY LITE 1517IN 3844CM |         |
| 1B   | CJL1518P    | COVER, POS, JAY LITE 1518IN        |         |
| 1B   | CJL1615P    | COVER, POS, JAY LITE 1615IN 4038CM |         |
| 1B   | CJL1616P    | COVER, POS, JAY LITE 1616IN 4040CM |         |
| 1B   | CJL1617P    | COVER, POS, JAY LITE 1617IN 4043CM |         |
| 1B   | CJL1618P    | COVER, POS, JAY LITE 1618IN 4045CM |         |
| 1B   | CJL1620P    | COVER, POS, JAY LITE 1620IN 4050CM |         |
| 1B   | CJL1715P    | COVER, POS, JAY LITE 1715IN 4338CM |         |
| 1B   | CJL1716P    | COVER, POS, JAY LITE 1716IN 4340   |         |
| 1B   | CJL1717P    | COVER, POS, JAY LITE 1717IN 4343CM |         |
| 1B   | CJL1718P    | COVER, POS, JAY LITE 1718IN 4345CM |         |
| 1B   | CJL1816P    | COVER, POS, JAY LITE 1816IN 4540CM |         |
| 1B   | CJL1817P    | COVER, POS, JAY LITE 1817IN 4543CM |         |
| 1B   | CJL1818P    | COVER, POS, JAY LITE 1818IN 4545CM |         |
| 1B   | CJL1820P    | COVER, POS, JAY LITE 1820IN 4550CM |         |
| 1B   | CJL1916P    | COVER, POS, JAY LITE 1916IN 4840CM |         |
| 1B   | CJL1917P    | COVER, POS, JAY LITE 1917IN 4843CM |         |
| 1B   | CJL1918P    | COVER, POS, JAY LITE 1918IN 4845CM |         |
| 1B   | CJL1920P    | COVER, POS, JAY LITE 1920IN 4850CM |         |
| 1B   | CJL2016P    | COVER, POS, JAY LITE 2016IN 5040CM |         |
| 1B   | CJL2018P    | COVER, POS, JAY LITE 2018IN 5045CM |         |
| 1B   | CJL2020P    | COVER, POS, JAY LITE 2020IN 5050CM |         |
| 1B   | CJL2218P    | COVER, POS, JAY LITE 2218IN 5545CM |         |
| 1B   | CJL2220P    | COVER, POS, JAY LITE 2220IN 5550CM |         |
| 1B   | CJL2222P    | COVER, POS, JAY LITE 2222IN 5555CM |         |
| 1B   | CJL2418P    | COVER, POS, JAY LITE 2418IN 6045CM |         |
| 1B   | CJL2420P    | COVER, POS, JAY LITE 2420IN 6050CM |         |
| 1B   | CJL2422P    | COVER, POS, JAY LITE 2422IN 6055CM |         |
| 1B   | CJL2424P    | COVER, POS, JAY LITE 2424IN 6060CM |         |
| 2A   | BJL1414     | BASE JAY LITE 1414IN 3535CM        |         |
| 2A   | BJL1416     | BASE JAY LITE 1416IN 3535CM        |         |
| 2A   | BJL1418     | BASE JAY LITE 1418IN 3535CM        |         |

| Pos. | Item Number | Description                       | Remarks |
|------|-------------|-----------------------------------|---------|
| 2A   | BJL1515     | BASE JAY LITE 1515IN              |         |
| 2A   | BJL1516     | BASE JAY LITE 1516IN 3840CM       |         |
| 2A   | BJL1517     | BASE JAY LITE 1517IN 3844CM       |         |
| 2A   | BJL1518     | BASE JAY LITE 1518IN              |         |
| 2A   | BJL1615     | BASE JAY LITE 1615IN 4038CM       |         |
| 2A   | BJL1616     | BASE JAY LITE 1616IN 4040CM       |         |
| 2A   | BJL1617     | BASE JAY LITE 1617IN 4043CM       |         |
| 2A   | BJL1618     | BASE JAY LITE 1618IN 4045CM       |         |
| 2A   | BJL1620     | BASE JAY LITE 1620IN 4050CM       |         |
| 2A   | BJL1715     | BASE JAY LITE 1715IN 4338CM       |         |
| 2A   | BJL1716     | BASE JAY LITE 1716IN 4340         |         |
| 2A   | BJL1717     | BASE JAY LITE 1717IN 4343CM       |         |
| 2A   | BJL1718     | BASE JAY LITE 1718IN 4345CM       |         |
| 2A   | BJL1816     | BASE JAY LITE 1816IN 4540CM       |         |
| 2A   | BJL1817     | BASE JAY LITE 1817IN 4543CM       |         |
| 2A   | BJL1818     | BASE JAY LITE 1818IN 4545CM       |         |
| 2A   | BJL1820     | BASE JAY LITE 1820IN 4550CM       |         |
| 2A   | BJL1916     | BASE JAY LITE 1916IN 4840CM       |         |
| 2A   | BJL1917     | BASE JAY LITE 1917IN 4843CM       |         |
| 2A   | BJL1918     | BASE JAY LITE 1918IN 4845CM       |         |
| 2A   | BJL1920     | BASE JAY LITE 1920IN 4850CM       |         |
| 2A   | BJL2016     | BASE JAY LITE 2016IN 5040CM       |         |
| 2A   | BJL2018     | BASE JAY LITE 2018IN 5045CM       |         |
| 2A   | BJL2020     | BASE JAY LITE 2020IN 5050CM       |         |
| 2A   | BJL2218     | BASE JAY LITE 2218IN 5545CM       |         |
| 2A   | BJL2220     | BASE JAY LITE 2220IN 5550CM       |         |
| 2A   | BJL2222     | BASE JAY LITE 2222IN 5555CM       |         |
| 2A   | BJL2418     | BASE JAY LITE 2418IN 6045CM       |         |
| 2A   | BJL2420     | BASE JAY LITE 2420IN 6050CM       |         |
| 2A   | BJL2422     | BASE JAY LITE 2422IN 6055CM       |         |
| 2A   | BJL2424     | BASE JAY LITE 2424IN 6060CM       |         |
| 2B   | BJL1414P    | BASE, POS, JAY LITE 1414IN 3535CM |         |
| 2B   | BJL1416P    | BASE, POS, JAY LITE 1416IN 3535CM |         |
| 2B   | BJL1418P    | BASE, POS, JAY LITE 1418IN 3535CM |         |
| 2B   | BJL1515P    | BASE, POS, JAY LITE 1515IN        |         |
| 2B   | BJL1516P    | BASE, POS, JAY LITE 1516IN 3840CM |         |
| 2B   | BJL1517P    | BASE, POS, JAY LITE 1517IN 3844CM |         |
| 2B   | BJL1518P    | BASE, POS, JAY LITE 1518IN        |         |
| 2B   | BJL1615P    | BASE, POS, JAY LITE 1615IN 4038CM |         |
| 2B   | BJL1616P    | BASE, POS, JAY LITE 1616IN 4040CM |         |
| 2B   | BJL1617P    | BASE, POS, JAY LITE 1617IN 4043CM |         |
| 2B   | BJL1618P    | BASE, POS, JAY LITE 1618IN 4045CM |         |
| 2B   | BJL1620P    | BASE, POS, JAY LITE 1620IN 4050CM |         |
| 2B   | BJL1715P    | BASE, POS, JAY LITE 1715IN 4338CM |         |
| 2B   | BJL1716P    | BASE, POS, JAY LITE 1716IN 4340   |         |
| 2B   | BJL1717P    | BASE, POS, JAY LITE 1717IN 4343CM |         |
| 2B   | BJL1718P    | BASE, POS, JAY LITE 1718IN 4345CM |         |
| 2B   | BJL1816P    | BASE, POS, JAY LITE 1816IN 4540CM |         |
| 2B   | BJL1817P    | BASE, POS, JAY LITE 1817IN 4543CM |         |
| 2B   | BJL1818P    | BASE, POS, JAY LITE 1818IN 4545CM |         |
| 2B   | BJL1820P    | BASE, POS, JAY LITE 1820IN 4550CM |         |
| 2B   | BJL1916P    | BASE, POS, JAY LITE 1916IN 4840CM |         |
| 2B   | BJL1917P    | BASE, POS, JAY LITE 1917IN 4843CM |         |
| 2B   | BJL1918P    | BASE, POS, JAY LITE 1918IN 4845CM |         |
| 2B   | BJL1920P    | BASE, POS, JAY LITE 1920IN 4850CM |         |
| 2B   | BJL2016P    | BASE, POS, JAY LITE 2016IN 5040CM |         |
| 2B   | BJL2018P    | BASE, POS, JAY LITE 2018IN 5045CM |         |

| Pos. | Item Number | Description                       | Remarks       |
|------|-------------|-----------------------------------|---------------|
| 2B   | BJL2020P    | BASE, POS, JAY LITE 2020IN 5050CM |               |
| 2B   | BJL2218P    | BASE, POS, JAY LITE 2218IN 5545CM |               |
| 2B   | BJL2220P    | BASE, POS, JAY LITE 2220IN 5550CM |               |
| 2B   | BJL2222P    | BASE, POS, JAY LITE 2222IN 5555CM |               |
| 2B   | BJL2418P    | BASE, POS, JAY LITE 2418IN 6045CM |               |
| 2B   | BJL2420P    | BASE, POS, JAY LITE 2420IN 6050CM |               |
| 2B   | BJL2422P    | BASE, POS, JAY LITE 2422IN 6055CM |               |
| 2B   | BJL2424P    | BASE, POS, JAY LITE 2424IN 6060CM |               |
| 3    | PUIJLA      | PU FOAM INSERT JAY LITE PLA A     | 14" wide only |
| 3    | PUIJLB      | PU FOAM INSERT JAY LITE PLA B     | 15"-17" wide  |
| 3    | PUIJLC      | PU FOAM INSERT JAY LITE PLA C     | 18"-20" wide  |

**MEASURING GUIDE**

| PROFILE         | "A"   |
|-----------------|-------|
| STANDARD        | 4 1/2 |
| REDUCED PROFILE | 3 1/2 |

| Pos. | Item Number | Description                    | Pos. | Item Number | Description                    |
|------|-------------|--------------------------------|------|-------------|--------------------------------|
|      |             | STANDARD OUTER COVER           |      |             |                                |
| 1A   | CJU1414     | OUTER CVR JAY UNION CUSH 14X14 | 1A   | CJU1717     | OUTER CVR JAY UNION CUSH 17X17 |
| 1A   | CJU1415     | OUTER CVR JAY UNION CUSH 14X15 | 1A   | CJU1718     | OUTER CVR JAY UNION CUSH 17X18 |
| 1A   | CJU1416     | OUTER CVR JAY UNION CUSH 14X16 | 1A   | CJU1719     | OUTER CVR JAY UNION CUSH 17X19 |
| 1A   | CJU1417     | OUTER CVR JAY UNION CUSH 14X17 | 1A   | CJU1720     | OUTER CVR JAY UNION CUSH 17X20 |
| 1A   | CJU1418     | OUTER CVR JAY UNION CUSH 14X18 | 1A   | CJU1721     | OUTER CVR JAY UNION CUSH 17X21 |
| 1A   | CJU1514     | OUTER CVR JAY UNION CUSH 15X14 | 1A   | CJU1814     | OUTER CVR JAY UNION CUSH 18X14 |
| 1A   | CJU1515     | OUTER CVR JAY UNION CUSH 15X15 | 1A   | CJU1815     | OUTER CVR JAY UNION CUSH 18X15 |
| 1A   | CJU1516     | OUTER CVR JAY UNION CUSH 15X16 | 1A   | CJU1816     | OUTER CVR JAY UNION CUSH 18X16 |
| 1A   | CJU1517     | OUTER CVR JAY UNION CUSH 15X17 | 1A   | CJU1817     | OUTER CVR JAY UNION CUSH 18X17 |
| 1A   | CJU1518     | OUTER CVR JAY UNION CUSH 15X18 | 1A   | CJU1818     | OUTER CVR JAY UNION CUSH 18X18 |
| 1A   | CJU1519     | OUTER CVR JAY UNION CUSH 15X19 | 1A   | CJU1819     | OUTER CVR JAY UNION CUSH 18X19 |
| 1A   | CJU1520     | OUTER CVR JAY UNION CUSH 15X20 | 1A   | CJU1820     | OUTER CVR JAY UNION CUSH 18X20 |
| 1A   | CJU1614     | OUTER CVR JAY UNION CUSH 16X14 | 1A   | CJU1821     | OUTER CVR JAY UNION CUSH 18X21 |
| 1A   | CJU1615     | OUTER CVR JAY UNION CUSH 16X15 | 1A   | CJU1822     | OUTER CVR JAY UNION CUSH 18X22 |
| 1A   | CJU1616     | OUTER CVR JAY UNION CUSH 16X16 | 1A   | CJU1915     | OUTER CVR JAY UNION CUSH 19X15 |
| 1A   | CJU1617     | OUTER CVR JAY UNION CUSH 16X17 | 1A   | CJU1916     | OUTER CVR JAY UNION CUSH 19X16 |
| 1A   | CJU1618     | OUTER CVR JAY UNION CUSH 16X18 | 1A   | CJU1917     | OUTER CVR JAY UNION CUSH 19X17 |
| 1A   | CJU1619     | OUTER CVR JAY UNION CUSH 16X19 | 1A   | CJU1918     | OUTER CVR JAY UNION CUSH 19X18 |
| 1A   | CJU1620     | OUTER CVR JAY UNION CUSH 16X20 | 1A   | CJU1919     | OUTER CVR JAY UNION CUSH 19X19 |
| 1A   | CJU1714     | OUTER CVR JAY UNION CUSH 17X14 | 1A   | CJU1920     | OUTER CVR JAY UNION CUSH 19X20 |
| 1A   | CJU1715     | OUTER CVR JAY UNION CUSH 17X15 | 1A   | CJU1921     | OUTER CVR JAY UNION CUSH 19X21 |
| 1A   | CJU1716     | OUTER CVR JAY UNION CUSH 17X16 | 1A   | CJU1922     | OUTER CVR JAY UNION CUSH 19X22 |
|      |             |                                | 1A   | CJU1923     | OUTER CVR JAY UNION CUSH 19X23 |

| Pos. | Item Number | Description                    | Pos. | Item Number | Description                    |
|------|-------------|--------------------------------|------|-------------|--------------------------------|
| 1A   | CJU2015     | OUTER CVR JAY UNION CUSH 20X15 | 2A   | ICJU1619    | INNER CVR JAY UNION CUSH 16X19 |
| 1A   | CJU2016     | OUTER CVR JAY UNION CUSH 20X16 | 2A   | ICJU1620    | INNER CVR JAY UNION CUSH 16X20 |
| 1A   | CJU2017     | OUTER CVR JAY UNION CUSH 20X17 | 2A   | ICJU1714    | INNER CVR JAY UNION CUSH 17X14 |
| 1A   | CJU2018     | OUTER CVR JAY UNION CUSH 20X18 | 2A   | ICJU1715    | INNER CVR JAY UNION CUSH 17X15 |
| 1A   | CJU2019     | OUTER CVR JAY UNION CUSH 20X19 | 2A   | ICJU1716    | INNER CVR JAY UNION CUSH 17X16 |
| 1A   | CJU2020     | OUTER CVR JAY UNION CUSH 20X20 | 2A   | ICJU1717    | INNER CVR JAY UNION CUSH 17X17 |
| 1A   | CJU2021     | OUTER CVR JAY UNION CUSH 20X21 | 2A   | ICJU1718    | INNER CVR JAY UNION CUSH 17X18 |
| 1A   | CJU2022     | OUTER CVR JAY UNION CUSH 20X22 | 2A   | ICJU1719    | INNER CVR JAY UNION CUSH 17X19 |
| 1A   | CJU2023     | OUTER CVR JAY UNION CUSH 20X23 | 2A   | ICJU1720    | INNER CVR JAY UNION CUSH 17X20 |
| 1A   | CJU2024     | OUTER CVR JAY UNION CUSH 20X24 | 2A   | ICJU1721    | INNER CVR JAY UNION CUSH 17X21 |
| 1A   | CJU2117     | OUTER CVR JAY UNION CUSH 21X17 | 2A   | ICJU1814    | INNER CVR JAY UNION CUSH 18X14 |
| 1A   | CJU2118     | OUTER CVR JAY UNION CUSH 21X18 | 2A   | ICJU1815    | INNER CVR JAY UNION CUSH 18X15 |
| 1A   | CJU2119     | OUTER CVR JAY UNION CUSH 21X19 | 2A   | ICJU1816    | INNER CVR JAY UNION CUSH 18X16 |
| 1A   | CJU2120     | OUTER CVR JAY UNION CUSH 21X20 | 2A   | ICJU1817    | INNER CVR JAY UNION CUSH 18X17 |
| 1A   | CJU2121     | OUTER CVR JAY UNION CUSH 21X21 | 2A   | ICJU1818    | INNER CVR JAY UNION CUSH 18X18 |
| 1A   | CJU2122     | OUTER CVR JAY UNION CUSH 21X22 | 2A   | ICJU1819    | INNER CVR JAY UNION CUSH 18X19 |
| 1A   | CJU2123     | OUTER CVR JAY UNION CUSH 21X23 | 2A   | ICJU1820    | INNER CVR JAY UNION CUSH 18X20 |
| 1A   | CJU2124     | OUTER CVR JAY UNION CUSH 21X24 | 2A   | ICJU1821    | INNER CVR JAY UNION CUSH 18X21 |
| 1A   | CJU2217     | OUTER CVR JAY UNION CUSH 22X17 | 2A   | ICJU1822    | INNER CVR JAY UNION CUSH 18X22 |
| 1A   | CJU2218     | OUTER CVR JAY UNION CUSH 22X18 | 2A   | ICJU1915    | INNER CVR JAY UNION CUSH 19X15 |
| 1A   | CJU2219     | OUTER CVR JAY UNION CUSH 22X19 | 2A   | ICJU1916    | INNER CVR JAY UNION CUSH 19X16 |
| 1A   | CJU2220     | OUTER CVR JAY UNION CUSH 22X20 | 2A   | ICJU1917    | INNER CVR JAY UNION CUSH 19X17 |
| 1A   | CJU2221     | OUTER CVR JAY UNION CUSH 22X21 | 2A   | ICJU1918    | INNER CVR JAY UNION CUSH 19X18 |
| 1A   | CJU2222     | OUTER CVR JAY UNION CUSH 22X22 | 2A   | ICJU1919    | INNER CVR JAY UNION CUSH 19X19 |
| 1A   | CJU2223     | OUTER CVR JAY UNION CUSH 22X23 | 2A   | ICJU1920    | INNER CVR JAY UNION CUSH 19X20 |
| 1A   | CJU2224     | OUTER CVR JAY UNION CUSH 22X24 | 2A   | ICJU1921    | INNER CVR JAY UNION CUSH 19X21 |
| 1A   | CJU2319     | OUTER CVR JAY UNION CUSH 23X19 | 2A   | ICJU1922    | INNER CVR JAY UNION CUSH 19X22 |
| 1A   | CJU2320     | OUTER CVR JAY UNION CUSH 23X20 | 2A   | ICJU1923    | INNER CVR JAY UNION CUSH 19X23 |
| 1A   | CJU2321     | OUTER CVR JAY UNION CUSH 23X21 | 2A   | ICJU2015    | INNER CVR JAY UNION CUSH 20X15 |
| 1A   | CJU2322     | OUTER CVR JAY UNION CUSH 23X22 | 2A   | ICJU2016    | INNER CVR JAY UNION CUSH 20X16 |
| 1A   | CJU2323     | OUTER CVR JAY UNION CUSH 23X23 | 2A   | ICJU2017    | INNER CVR JAY UNION CUSH 20X17 |
| 1A   | CJU2324     | OUTER CVR JAY UNION CUSH 23X24 | 2A   | ICJU2018    | INNER CVR JAY UNION CUSH 20X18 |
| 1A   | CJU2418     | OUTER CVR JAY UNION CUSH 24X18 | 2A   | ICJU2019    | INNER CVR JAY UNION CUSH 20X19 |
| 1A   | CJU2419     | OUTER CVR JAY UNION CUSH 24X19 | 2A   | ICJU2020    | INNER CVR JAY UNION CUSH 20X20 |
| 1A   | CJU2420     | OUTER CVR JAY UNION CUSH 24X20 | 2A   | ICJU2021    | INNER CVR JAY UNION CUSH 20X21 |
| 1A   | CJU2421     | OUTER CVR JAY UNION CUSH 24X21 | 2A   | ICJU2022    | INNER CVR JAY UNION CUSH 20X22 |
| 1A   | CJU2422     | OUTER CVR JAY UNION CUSH 24X22 | 2A   | ICJU2023    | INNER CVR JAY UNION CUSH 20X23 |
| 1A   | CJU2423     | OUTER CVR JAY UNION CUSH 24X23 | 2A   | ICJU2024    | INNER CVR JAY UNION CUSH 20X24 |
| 1A   | CJU2424     | OUTER CVR JAY UNION CUSH 24X24 | 2A   | ICJU2117    | INNER CVR JAY UNION CUSH 21X17 |
|      |             | STANDARD INNER COVER           | 2A   | ICJU2118    | INNER CVR JAY UNION CUSH 21X18 |
| 2A   | ICJU1414    | INNER CVR JAY UNION CUSH 14X14 | 2A   | ICJU2119    | INNER CVR JAY UNION CUSH 21X19 |
| 2A   | ICJU1415    | INNER CVR JAY UNION CUSH 14X15 | 2A   | ICJU2120    | INNER CVR JAY UNION CUSH 21X20 |
| 2A   | ICJU1416    | INNER CVR JAY UNION CUSH 14X16 | 2A   | ICJU2121    | INNER CVR JAY UNION CUSH 21X21 |
| 2A   | ICJU1417    | INNER CVR JAY UNION CUSH 14X17 | 2A   | ICJU2122    | INNER CVR JAY UNION CUSH 21X22 |
| 2A   | ICJU1418    | INNER CVR JAY UNION CUSH 14X18 | 2A   | ICJU2123    | INNER CVR JAY UNION CUSH 21X23 |
| 2A   | ICJU1514    | INNER CVR JAY UNION CUSH 15X14 | 2A   | ICJU2124    | INNER CVR JAY UNION CUSH 21X24 |
| 2A   | ICJU1515    | INNER CVR JAY UNION CUSH 15X15 | 2A   | ICJU2217    | INNER CVR JAY UNION CUSH 22X17 |
| 2A   | ICJU1516    | INNER CVR JAY UNION CUSH 15X16 | 2A   | ICJU2218    | INNER CVR JAY UNION CUSH 22X18 |
| 2A   | ICJU1517    | INNER CVR JAY UNION CUSH 15X17 | 2A   | ICJU2219    | INNER CVR JAY UNION CUSH 22X19 |
| 2A   | ICJU1518    | INNER CVR JAY UNION CUSH 15X18 | 2A   | ICJU2220    | INNER CVR JAY UNION CUSH 22X20 |
| 2A   | ICJU1519    | INNER CVR JAY UNION CUSH 15X19 | 2A   | ICJU2221    | INNER CVR JAY UNION CUSH 22X21 |
| 2A   | ICJU1520    | INNER CVR JAY UNION CUSH 15X20 | 2A   | ICJU2222    | INNER CVR JAY UNION CUSH 22X22 |
| 2A   | ICJU1614    | INNER CVR JAY UNION CUSH 16X14 | 2A   | ICJU2223    | INNER CVR JAY UNION CUSH 22X23 |
| 2A   | ICJU1615    | INNER CVR JAY UNION CUSH 16X15 | 2A   | ICJU2224    | INNER CVR JAY UNION CUSH 22X24 |
| 2A   | ICJU1616    | INNER CVR JAY UNION CUSH 16X16 | 2A   | ICJU2319    | INNER CVR JAY UNION CUSH 23X19 |
| 2A   | ICJU1617    | INNER CVR JAY UNION CUSH 16X17 | 2A   | ICJU2320    | INNER CVR JAY UNION CUSH 23X20 |
| 2A   | ICJU1618    | INNER CVR JAY UNION CUSH 16X18 | 2A   | ICJU2321    | INNER CVR JAY UNION CUSH 23X21 |

| Pos. | Item Number | Description                    | Pos. | Item Number | Description                    |
|------|-------------|--------------------------------|------|-------------|--------------------------------|
| 2A   | ICJU2322    | INNER CVR JAY UNION CUSH 23X22 | 3A   | BJU2016     | BASE FOAM JAY UNION CUSH 20X16 |
| 2A   | ICJU2323    | INNER CVR JAY UNION CUSH 23X23 | 3A   | BJU2017     | BASE FOAM JAY UNION CUSH 20X17 |
| 2A   | ICJU2324    | INNER CVR JAY UNION CUSH 23X24 | 3A   | BJU2018     | BASE FOAM JAY UNION CUSH 20X18 |
| 2A   | ICJU2418    | INNER CVR JAY UNION CUSH 24X18 | 3A   | BJU2019     | BASE FOAM JAY UNION CUSH 20X19 |
| 2A   | ICJU2419    | INNER CVR JAY UNION CUSH 24X19 | 3A   | BJU2020     | BASE FOAM JAY UNION CUSH 20X20 |
| 2A   | ICJU2420    | INNER CVR JAY UNION CUSH 24X20 | 3A   | BJU2021     | BASE FOAM JAY UNION CUSH 20X21 |
| 2A   | ICJU2421    | INNER CVR JAY UNION CUSH 24X21 | 3A   | BJU2022     | BASE FOAM JAY UNION CUSH 20X22 |
| 2A   | ICJU2422    | INNER CVR JAY UNION CUSH 24X22 | 3A   | BJU2023     | BASE FOAM JAY UNION CUSH 20X23 |
| 2A   | ICJU2423    | INNER CVR JAY UNION CUSH 24X23 | 3A   | BJU2024     | BASE FOAM JAY UNION CUSH 20X24 |
| 2A   | ICJU2424    | INNER CVR JAY UNION CUSH 24X24 | 3A   | BJU2117     | BASE FOAM JAY UNION CUSH 21X17 |
|      |             | STANDARD FOAM BASE             | 3A   | BJU2118     | BASE FOAM JAY UNION CUSH 21X18 |
| 3A   | BJU1414     | BASE FOAM JAY UNION CUSH 14X14 | 3A   | BJU2119     | BASE FOAM JAY UNION CUSH 21X19 |
| 3A   | BJU1415     | BASE FOAM JAY UNION CUSH 14X15 | 3A   | BJU2120     | BASE FOAM JAY UNION CUSH 21X20 |
| 3A   | BJU1416     | BASE FOAM JAY UNION CUSH 14X16 | 3A   | BJU2121     | BASE FOAM JAY UNION CUSH 21X21 |
| 3A   | BJU1417     | BASE FOAM JAY UNION CUSH 14X17 | 3A   | BJU2122     | BASE FOAM JAY UNION CUSH 21X22 |
| 3A   | BJU1418     | BASE FOAM JAY UNION CUSH 14X18 | 3A   | BJU2123     | BASE FOAM JAY UNION CUSH 21X23 |
| 3A   | BJU1514     | BASE FOAM JAY UNION CUSH 15X14 | 3A   | BJU2124     | BASE FOAM JAY UNION CUSH 21X24 |
| 3A   | BJU1515     | BASE FOAM JAY UNION CUSH 15X15 | 3A   | BJU2217     | BASE FOAM JAY UNION CUSH 22X17 |
| 3A   | BJU1516     | BASE FOAM JAY UNION CUSH 15X16 | 3A   | BJU2218     | BASE FOAM JAY UNION CUSH 22X18 |
| 3A   | BJU1517     | BASE FOAM JAY UNION CUSH 15X17 | 3A   | BJU2219     | BASE FOAM JAY UNION CUSH 22X19 |
| 3A   | BJU1518     | BASE FOAM JAY UNION CUSH 15X18 | 3A   | BJU2220     | BASE FOAM JAY UNION CUSH 22X20 |
| 3A   | BJU1519     | BASE FOAM JAY UNION CUSH 15X19 | 3A   | BJU2221     | BASE FOAM JAY UNION CUSH 22X21 |
| 3A   | BJU1520     | BASE FOAM JAY UNION CUSH 15X20 | 3A   | BJU2222     | BASE FOAM JAY UNION CUSH 22X22 |
| 3A   | BJU1614     | BASE FOAM JAY UNION CUSH 16X14 | 3A   | BJU2223     | BASE FOAM JAY UNION CUSH 22X23 |
| 3A   | BJU1615     | BASE FOAM JAY UNION CUSH 16X15 | 3A   | BJU2224     | BASE FOAM JAY UNION CUSH 22X24 |
| 3A   | BJU1616     | BASE FOAM JAY UNION CUSH 16X16 | 3A   | BJU2319     | BASE FOAM JAY UNION CUSH 23X19 |
| 3A   | BJU1617     | BASE FOAM JAY UNION CUSH 16X17 | 3A   | BJU2320     | BASE FOAM JAY UNION CUSH 23X20 |
| 3A   | BJU1618     | BASE FOAM JAY UNION CUSH 16X18 | 3A   | BJU2321     | BASE FOAM JAY UNION CUSH 23X21 |
| 3A   | BJU1619     | BASE FOAM JAY UNION CUSH 16X19 | 3A   | BJU2322     | BASE FOAM JAY UNION CUSH 23X22 |
| 3A   | BJU1620     | BASE FOAM JAY UNION CUSH 16X20 | 3A   | BJU2323     | BASE FOAM JAY UNION CUSH 23X23 |
| 3A   | BJU1714     | BASE FOAM JAY UNION CUSH 17X14 | 3A   | BJU2324     | BASE FOAM JAY UNION CUSH 23X24 |
| 3A   | BJU1715     | BASE FOAM JAY UNION CUSH 17X15 | 3A   | BJU2418     | BASE FOAM JAY UNION CUSH 24X18 |
| 3A   | BJU1716     | BASE FOAM JAY UNION CUSH 17X16 | 3A   | BJU2419     | BASE FOAM JAY UNION CUSH 24X19 |
| 3A   | BJU1717     | BASE FOAM JAY UNION CUSH 17X17 | 3A   | BJU2420     | BASE FOAM JAY UNION CUSH 24X20 |
| 3A   | BJU1718     | BASE FOAM JAY UNION CUSH 17X18 | 3A   | BJU2421     | BASE FOAM JAY UNION CUSH 24X21 |
| 3A   | BJU1719     | BASE FOAM JAY UNION CUSH 17X19 | 3A   | BJU2422     | BASE FOAM JAY UNION CUSH 24X22 |
| 3A   | BJU1720     | BASE FOAM JAY UNION CUSH 17X20 | 3A   | BJU2423     | BASE FOAM JAY UNION CUSH 24X23 |
| 3A   | BJU1721     | BASE FOAM JAY UNION CUSH 17X21 | 3A   | BJU2424     | BASE FOAM JAY UNION CUSH 24X24 |
| 3A   | BJU1814     | BASE FOAM JAY UNION CUSH 18X14 |      |             | REDUCED PROFILE OUTER COVER    |
| 3A   | BJU1815     | BASE FOAM JAY UNION CUSH 18X15 | 1B   | CJU1414RP   | OUTER CVR UNION RP CUSH 14X14  |
| 3A   | BJU1816     | BASE FOAM JAY UNION CUSH 18X16 | 1B   | CJU1415RP   | OUTER CVR UNION RP CUSH 14X15  |
| 3A   | BJU1817     | BASE FOAM JAY UNION CUSH 18X17 | 1B   | CJU1416RP   | OUTER CVR UNION RP CUSH 14X16  |
| 3A   | BJU1818     | BASE FOAM JAY UNION CUSH 18X18 | 1B   | CJU1417RP   | OUTER CVR UNION RP CUSH 14X17  |
| 3A   | BJU1819     | BASE FOAM JAY UNION CUSH 18X19 | 1B   | CJU1418RP   | OUTER CVR UNION RP CUSH 14X18  |
| 3A   | BJU1820     | BASE FOAM JAY UNION CUSH 18X20 | 1B   | CJU1514RP   | OUTER CVR UNION RP CUSH 15X14  |
| 3A   | BJU1821     | BASE FOAM JAY UNION CUSH 18X21 | 1B   | CJU1515RP   | OUTER CVR UNION RP CUSH 15X15  |
| 3A   | BJU1822     | BASE FOAM JAY UNION CUSH 18X22 | 1B   | CJU1516RP   | OUTER CVR UNION RP CUSH 15X16  |
| 3A   | BJU1915     | BASE FOAM JAY UNION CUSH 19X15 | 1B   | CJU1517RP   | OUTER CVR UNION RP CUSH 15X17  |
| 3A   | BJU1916     | BASE FOAM JAY UNION CUSH 19X16 | 1B   | CJU1518RP   | OUTER CVR UNION RP CUSH 15X18  |
| 3A   | BJU1917     | BASE FOAM JAY UNION CUSH 19X17 | 1B   | CJU1519RP   | OUTER CVR UNION RP CUSH 15X19  |
| 3A   | BJU1918     | BASE FOAM JAY UNION CUSH 19X18 | 1B   | CJU1520RP   | OUTER CVR UNION RP CUSH 15X20  |
| 3A   | BJU1919     | BASE FOAM JAY UNION CUSH 19X19 | 1B   | CJU1614RP   | OUTER CVR UNION RP CUSH 16X14  |
| 3A   | BJU1920     | BASE FOAM JAY UNION CUSH 19X20 | 1B   | CJU1615RP   | OUTER CVR UNION RP CUSH 16X15  |
| 3A   | BJU1921     | BASE FOAM JAY UNION CUSH 19X21 | 1B   | CJU1616RP   | OUTER CVR UNION RP CUSH 16X16  |
| 3A   | BJU1922     | BASE FOAM JAY UNION CUSH 19X22 | 1B   | CJU1617RP   | OUTER CVR UNION RP CUSH 16X17  |
| 3A   | BJU1923     | BASE FOAM JAY UNION CUSH 19X23 | 1B   | CJU1618RP   | OUTER CVR UNION RP CUSH 16X18  |
| 3A   | BJU2015     | BASE FOAM JAY UNION CUSH 20X15 | 1B   | CJU1619RP   | OUTER CVR UNION RP CUSH 16X19  |


| Pos. | Item Number | Description                   | Pos. | Item Number | Description                     |
|------|-------------|-------------------------------|------|-------------|---------------------------------|
| 1B   | CJU1620RP   | OUTER CVR UNION RP CUSH 16X20 | 2B   | ICJU1817RP  | INNER CVR UNION RP CUSH 18X17   |
| 1B   | CJU1714RP   | OUTER CVR UNION RP CUSH 17X14 | 2B   | ICJU1818RP  | INNER CVR UNION RP CUSH 18X18   |
| 1B   | CJU1715RP   | OUTER CVR UNION RP CUSH 17X15 | 2B   | ICJU1819RP  | INNER CVR UNION RP CUSH 18X19   |
| 1B   | CJU1716RP   | OUTER CVR UNION RP CUSH 17X16 | 2B   | ICJU1820RP  | INNER CVR UNION RP CUSH 18X20   |
| 1B   | CJU1717RP   | OUTER CVR UNION RP CUSH 17X17 | 2B   | ICJU1915RP  | INNER CVR UNION RP CUSH 19X15   |
| 1B   | CJU1718RP   | OUTER CVR UNION RP CUSH 17X18 | 2B   | ICJU1916RP  | INNER CVR UNION RP CUSH 19X16   |
| 1B   | CJU1719RP   | OUTER CVR UNION RP CUSH 17X19 | 2B   | ICJU1917RP  | INNER CVR UNION RP CUSH 19X17   |
| 1B   | CJU1720RP   | OUTER CVR UNION RP CUSH 17X20 | 2B   | ICJU1918RP  | INNER CVR UNION RP CUSH 19X18   |
| 1B   | CJU1814RP   | OUTER CVR UNION RP CUSH 18X14 | 2B   | ICJU1919RP  | INNER CVR UNION RP CUSH 19X19   |
| 1B   | CJU1815RP   | OUTER CVR UNION RP CUSH 18X15 | 2B   | ICJU1920RP  | INNER CVR UNION RP CUSH 19X20   |
| 1B   | CJU1816RP   | OUTER CVR UNION RP CUSH 18X16 | 2B   | ICJU2015RP  | INNER CVR UNION RP CUSH 20X15   |
| 1B   | CJU1817RP   | OUTER CVR UNION RP CUSH 18X17 | 2B   | ICJU2016RP  | INNER CVR UNION RP CUSH 20X16   |
| 1B   | CJU1818RP   | OUTER CVR UNION RP CUSH 18X18 | 2B   | ICJU2017RP  | INNER CVR UNION RP CUSH 20X17   |
| 1B   | CJU1819RP   | OUTER CVR UNION RP CUSH 18X19 | 2B   | ICJU2018RP  | INNER CVR UNION RP CUSH 20X18   |
| 1B   | CJU1820RP   | OUTER CVR UNION RP CUSH 18X20 | 2B   | ICJU2019RP  | INNER CVR UNION RP CUSH 20X19   |
| 1B   | CJU1915RP   | OUTER CVR UNION RP CUSH 19X15 | 2B   | ICJU2020RP  | INNER CVR UNION RP CUSH 20X20   |
| 1B   | CJU1916RP   | OUTER CVR UNION RP CUSH 19X16 |      |             | REDUCED PROFILE FOAM BASE COVER |
| 1B   | CJU1917RP   | OUTER CVR UNION RP CUSH 19X17 | 3B   | BJU1414RP   | BASE FOAM UNION RP CUSH 14X14   |
| 1B   | CJU1918RP   | OUTER CVR UNION RP CUSH 19X18 | 3B   | BJU1415RP   | BASE FOAM UNION RP CUSH 14X15   |
| 1B   | CJU1919RP   | OUTER CVR UNION RP CUSH 19X19 | 3B   | BJU1416RP   | BASE FOAM UNION RP CUSH 14X16   |
| 1B   | CJU1920RP   | OUTER CVR UNION RP CUSH 19X20 | 3B   | BJU1417RP   | BASE FOAM UNION RP CUSH 14X17   |
| 1B   | CJU2015RP   | OUTER CVR UNION RP CUSH 20X15 | 3B   | BJU1418RP   | BASE FOAM UNION RP CUSH 14X18   |
| 1B   | CJU2016RP   | OUTER CVR UNION RP CUSH 20X16 | 3B   | BJU1514RP   | BASE FOAM UNION RP CUSH 15X14   |
| 1B   | CJU2017RP   | OUTER CVR UNION RP CUSH 20X17 | 3B   | BJU1515RP   | BASE FOAM UNION RP CUSH 15X15   |
| 1B   | CJU2018RP   | OUTER CVR UNION RP CUSH 20X18 | 3B   | BJU1516RP   | BASE FOAM UNION RP CUSH 15X16   |
| 1B   | CJU2019RP   | OUTER CVR UNION RP CUSH 20X19 | 3B   | BJU1517RP   | BASE FOAM UNION RP CUSH 15X17   |
| 1B   | CJU2020RP   | OUTER CVR UNION RP CUSH 20X20 | 3B   | BJU1518RP   | BASE FOAM UNION RP CUSH 15X18   |
|      |             | REDUCED PROFILE INNER COVER   | 3B   | BJU1519RP   | BASE FOAM UNION RP CUSH 15X19   |
| 2B   | ICJU1414RP  | INNER CVR UNION RP CUSH 14X14 | 3B   | BJU1520RP   | BASE FOAM UNION RP CUSH 15X20   |
| 2B   | ICJU1415RP  | INNER CVR UNION RP CUSH 14X15 | 3B   | BJU1614RP   | BASE FOAM UNION RP CUSH 16X14   |
| 2B   | ICJU1416RP  | INNER CVR UNION RP CUSH 14X16 | 3B   | BJU1615RP   | BASE FOAM UNION RP CUSH 16X15   |
| 2B   | ICJU1417RP  | INNER CVR UNION RP CUSH 14X17 | 3B   | BJU1616RP   | BASE FOAM UNION RP CUSH 16X16   |
| 2B   | ICJU1418RP  | INNER CVR UNION RP CUSH 14X18 | 3B   | BJU1617RP   | BASE FOAM UNION RP CUSH 16X17   |
| 2B   | ICJU1514RP  | INNER CVR UNION RP CUSH 15X14 | 3B   | BJU1618RP   | BASE FOAM UNION RP CUSH 16X18   |
| 2B   | ICJU1515RP  | INNER CVR UNION RP CUSH 15X15 | 3B   | BJU1619RP   | BASE FOAM UNION RP CUSH 16X19   |
| 2B   | ICJU1516RP  | INNER CVR UNION RP CUSH 15X16 | 3B   | BJU1620RP   | BASE FOAM UNION RP CUSH 16X20   |
| 2B   | ICJU1517RP  | INNER CVR UNION RP CUSH 15X17 | 3B   | BJU1714RP   | BASE FOAM UNION RP CUSH 17X14   |
| 2B   | ICJU1518RP  | INNER CVR UNION RP CUSH 15X18 | 3B   | BJU1715RP   | BASE FOAM UNION RP CUSH 17X15   |
| 2B   | ICJU1519RP  | INNER CVR UNION RP CUSH 15X19 | 3B   | BJU1716RP   | BASE FOAM UNION RP CUSH 17X16   |
| 2B   | ICJU1520RP  | INNER CVR UNION RP CUSH 15X20 | 3B   | BJU1717RP   | BASE FOAM UNION RP CUSH 17X17   |
| 2B   | ICJU1614RP  | INNER CVR UNION RP CUSH 16X14 | 3B   | BJU1718RP   | BASE FOAM UNION RP CUSH 17X18   |
| 2B   | ICJU1615RP  | INNER CVR UNION RP CUSH 16X15 | 3B   | BJU1719RP   | BASE FOAM UNION RP CUSH 17X19   |
| 2B   | ICJU1616RP  | INNER CVR UNION RP CUSH 16X16 | 3B   | BJU1720RP   | BASE FOAM UNION RP CUSH 17X20   |
| 2B   | ICJU1617RP  | INNER CVR UNION RP CUSH 16X17 | 3B   | BJU1814RP   | BASE FOAM UNION RP CUSH 18X14   |
| 2B   | ICJU1618RP  | INNER CVR UNION RP CUSH 16X18 | 3B   | BJU1815RP   | BASE FOAM UNION RP CUSH 18X15   |
| 2B   | ICJU1619RP  | INNER CVR UNION RP CUSH 16X19 | 3B   | BJU1816RP   | BASE FOAM UNION RP CUSH 18X16   |
| 2B   | ICJU1620RP  | INNER CVR UNION RP CUSH 16X20 | 3B   | BJU1817RP   | BASE FOAM UNION RP CUSH 18X17   |
| 2B   | ICJU1714RP  | INNER CVR UNION RP CUSH 17X14 | 3B   | BJU1818RP   | BASE FOAM UNION RP CUSH 18X18   |
| 2B   | ICJU1715RP  | INNER CVR UNION RP CUSH 17X15 | 3B   | BJU1819RP   | BASE FOAM UNION RP CUSH 18X19   |
| 2B   | ICJU1716RP  | INNER CVR UNION RP CUSH 17X16 | 3B   | BJU1820RP   | BASE FOAM UNION RP CUSH 18X20   |
| 2B   | ICJU1717RP  | INNER CVR UNION RP CUSH 17X17 | 3B   | BJU1915RP   | BASE FOAM UNION RP CUSH 19X15   |
| 2B   | ICJU1718RP  | INNER CVR UNION RP CUSH 17X18 | 3B   | BJU1916RP   | BASE FOAM UNION RP CUSH 19X16   |
| 2B   | ICJU1719RP  | INNER CVR UNION RP CUSH 17X19 | 3B   | BJU1917RP   | BASE FOAM UNION RP CUSH 19X17   |
| 2B   | ICJU1720RP  | INNER CVR UNION RP CUSH 17X20 | 3B   | BJU1918RP   | BASE FOAM UNION RP CUSH 19X18   |
| 2B   | ICJU1814RP  | INNER CVR UNION RP CUSH 18X14 | 3B   | BJU1919RP   | BASE FOAM UNION RP CUSH 19X19   |
| 2B   | ICJU1815RP  | INNER CVR UNION RP CUSH 18X15 | 3B   | BJU1920RP   | BASE FOAM UNION RP CUSH 19X20   |
| 2B   | ICJU1816RP  | INNER CVR UNION RP CUSH 18X16 | 3B   | BJU2015RP   | BASE FOAM UNION RP CUSH 20X15   |

| Pos. | Item Number | Description                   | Pos. | Item Number | Description                   |
|------|-------------|-------------------------------|------|-------------|-------------------------------|
| 3B   | BJU2016RP   | BASE FOAM UNION RP CUSH 20X16 | 3B   | BJU2019RP   | BASE FOAM UNION RP CUSH 20X19 |
| 3B   | BJU2017RP   | BASE FOAM UNION RP CUSH 20X17 | 3B   | BJU2020RP   | BASE FOAM UNION RP CUSH 20X20 |
| 3B   | BJU2018RP   | BASE FOAM UNION RP CUSH 20X18 |      |             |                               |



| Pos. | Item Number | Description                 | Remarks |
|------|-------------|-----------------------------|---------|
| 1    | CJX21414    | COVER, JAY X2, OUTER, 14X14 |         |
| 1    | CJX21416    | COVER, JAY X2, OUTER, 14X16 |         |
| 1    | CJX21515    | COVER, JAY X2, OUTER, 15X15 |         |
| 1    | CJX21516    | COVER, JAY X2, OUTER, 15X16 |         |
| 1    | CJX21616    | COVER, JAY X2, OUTER, 16X16 |         |
| 1    | CJX21617    | COVER, JAY X2, OUTER, 16X17 |         |
| 1    | CJX21618    | COVER, JAY X2, OUTER, 16X18 |         |
| 1    | CJX21620    | COVER, JAY X2, OUTER, 16X20 |         |
| 1    | CJX21716    | COVER, JAY X2, OUTER, 17X16 |         |
| 1    | CJX21717    | COVER, JAY X2, OUTER, 17X17 |         |
| 1    | CJX21816    | COVER, JAY X2, OUTER, 18X16 |         |
| 1    | CJX21818    | COVER, JAY X2, OUTER, 18X18 |         |
| 1    | CJX21820    | COVER, JAY X2, OUTER, 18X20 |         |
| 1    | CJX22020    | COVER, JAY X2, OUTER, 20X20 |         |
| 2A   | FJX2-1      | FL PD JAY X2 14X14 14X16    |         |

| Pos. | Item Number | Description                    | Remarks                                        |
|------|-------------|--------------------------------|------------------------------------------------|
| 2A   | FJX2-2      | FL PD JAY X2 15X15 15X16 16X16 |                                                |
| 2A   | FJX2-3      | FL PD JAY X2 16X17 16X18 16X20 |                                                |
| 2A   | FJX2-4      | FL PD JAY X2 17X16 18X16       |                                                |
| 2A   | FJX2-5      | FL PD JAY X2 17X17 18X18 18X20 |                                                |
| 2A   | FJX2-6      | FL PD JAY X2 20X20             |                                                |
| 2B   | F128        | FLUID PAD JA 14X16             | ACTIVE FLUID PAD FOR JAY YOUR WAY MODIFICATION |
| 2B   | F125        | FLUID PAD JA 16X16             | ACTIVE FLUID PAD FOR JAY YOUR WAY MODIFICATION |
| 2B   | F126        | FLUID PAD JA 16X18             | ACTIVE FLUID PAD FOR JAY YOUR WAY MODIFICATION |
| 2B   | F120        | FLUID PAD JA 18X16             | ACTIVE FLUID PAD FOR JAY YOUR WAY MODIFICATION |
| 2B   | F121        | FLUID PAD JA 18X18             | ACTIVE FLUID PAD FOR JAY YOUR WAY MODIFICATION |
| 2B   | F122        | FLUID PAD JA 20X18             | ACTIVE FLUID PAD FOR JAY YOUR WAY MODIFICATION |
| 3    | B119JA      | 1" ADDUCTOR WEDGES             | PAIR                                           |
| 3    | B113JA      | 2" ADDUCTOR WEDGES (PR)        | PAIR                                           |
| 4    | ICJX21414   | COVER, JAY X2, INNER, 14X14    |                                                |
| 4    | ICJX21416   | COVER, JAY X2, INNER, 14X16    |                                                |
| 4    | ICJX21515   | COVER, JAY X2, INNER, 15X15    |                                                |
| 4    | ICJX21516   | COVER, JAY X2, INNER, 15X16    |                                                |
| 4    | ICJX21616   | COVER, JAY X2, INNER, 16X16    |                                                |
| 4    | ICJX21617   | COVER, JAY X2, INNER, 16X17    |                                                |
| 4    | ICJX21618   | COVER, JAY X2, INNER, 16X18    |                                                |
| 4    | ICJX21620   | COVER, JAY X2, INNER, 16X20    |                                                |
| 4    | ICJX21716   | COVER, JAY X2, INNER, 17X16    |                                                |
| 4    | ICJX21717   | COVER, JAY X2, INNER, 17X17    |                                                |
| 4    | ICJX21816   | COVER, JAY X2, INNER, 18X16    |                                                |
| 4    | ICJX21818   | COVER, JAY X2, INNER, 18X18    |                                                |
| 4    | ICJX21820   | COVER, JAY X2, INNER, 18X20    |                                                |
| 4    | ICJX22020   | COVER, JAY X2, INNER, 20X20    |                                                |
| 5    | BJX21414    | BASE, JAY X2, 14X14            |                                                |
| 5    | BJX21416    | BASE, JAY X2, 14X16            |                                                |
| 5    | BJX21515    | BASE, JAY X2, 15X15            |                                                |
| 5    | BJX21516    | BASE, JAY X2, 15X16            |                                                |
| 5    | BJX21616    | BASE, JAY X2, 16X16            |                                                |
| 5    | BJX21617    | BASE, JAY X2, 16X17            |                                                |
| 5    | BJX21618    | BASE, JAY X2, 16X18            |                                                |
| 5    | BJX21620    | BASE, JAY X2, 16X20            |                                                |
| 5    | BJX21716    | BASE, JAY X2, 17X16            |                                                |
| 5    | BJX21717    | BASE, JAY X2, 17X17            |                                                |
| 5    | BJX21816    | BASE, JAY X2, 18X16            |                                                |
| 5    | BJX21818    | BASE, JAY X2, 18X18            |                                                |
| 5    | BJX21820    | BASE, JAY X2, 18X20            |                                                |
| 5    | BJX22020    | BASE, JAY X2, 20X20            |                                                |




| Pos. | Item Number | Description               | Pos. | Item Number | Description               |
|------|-------------|---------------------------|------|-------------|---------------------------|
|      |             | OUTER COVER               | 1    | CJZ1008     | JAY ZIP OUTER COVER 10X08 |
| 1    | CJZ0808     | JAY ZIP OUTER COVER 08X08 | 1    | CJZ1009     | JAY ZIP OUTER COVER 10X09 |
| 1    | CJZ0809     | JAY ZIP OUTER COVER 08X09 | 1    | CJZ1010     | JAY ZIP OUTER COVER 10X10 |
| 1    | CJZ0810     | JAY ZIP OUTER COVER 08X10 | 1    | CJZ1011     | JAY ZIP OUTER COVER 10X11 |
| 1    | CJZ0811     | JAY ZIP OUTER COVER 08X11 | 1    | CJZ1012     | JAY ZIP OUTER COVER 10X12 |
| 1    | CJZ0812     | JAY ZIP OUTER COVER 08X12 | 1    | CJZ1013     | JAY ZIP OUTER COVER 10X13 |
| 1    | CJZ0813     | JAY ZIP OUTER COVER 08X13 | 1    | CJZ1014     | JAY ZIP OUTER COVER 10X14 |
| 1    | CJZ0814     | JAY ZIP OUTER COVER 08X14 | 1    | CJZ1015     | JAY ZIP OUTER COVER 10X15 |
| 1    | CJZ0908     | JAY ZIP OUTER COVER 09X08 | 1    | CJZ1016     | JAY ZIP OUTER COVER 10X16 |
| 1    | CJZ0909     | JAY ZIP OUTER COVER 09X09 | 1    | CJZ1108     | JAY ZIP OUTER COVER 11X08 |
| 1    | CJZ0910     | JAY ZIP OUTER COVER 09X10 | 1    | CJZ1109     | JAY ZIP OUTER COVER 11X09 |
| 1    | CJZ0911     | JAY ZIP OUTER COVER 09X11 | 1    | CJZ1110     | JAY ZIP OUTER COVER 11X10 |
| 1    | CJZ0912     | JAY ZIP OUTER COVER 09X12 | 1    | CJZ1111     | JAY ZIP OUTER COVER 11X11 |
| 1    | CJZ0913     | JAY ZIP OUTER COVER 09X13 | 1    | CJZ1112     | JAY ZIP OUTER COVER 11X12 |
| 1    | CJZ0914     | JAY ZIP OUTER COVER 09X14 | 1    | CJZ1113     | JAY ZIP OUTER COVER 11X13 |
| 1    | CJZ0915     | JAY ZIP OUTER COVER 09X15 | 1    | CJZ1114     | JAY ZIP OUTER COVER 11X14 |

| Pos. | Item Number | Description               | Pos. | Item Number | Description               |
|------|-------------|---------------------------|------|-------------|---------------------------|
| 1    | CJZ1115     | JAY ZIP OUTER COVER 11X15 | 2    | ICJZ0914    | JAY ZIP INNER COVER 09X14 |
| 1    | CJZ1116     | JAY ZIP OUTER COVER 11X16 | 2    | ICJZ0915    | JAY ZIP INNER COVER 09X15 |
| 1    | CJZ1117     | JAY ZIP OUTER COVER 11X17 | 2    | ICJZ1008    | JAY ZIP INNER COVER 10X08 |
| 1    | CJZ1209     | JAY ZIP OUTER COVER 12X09 | 2    | ICJZ1009    | JAY ZIP INNER COVER 10X09 |
| 1    | CJZ1210     | JAY ZIP OUTER COVER 12X10 | 2    | ICJZ1010    | JAY ZIP INNER COVER 10X10 |
| 1    | CJZ1211     | JAY ZIP OUTER COVER 12X11 | 2    | ICJZ1011    | JAY ZIP INNER COVER 10X11 |
| 1    | CJZ1212     | JAY ZIP OUTER COVER 12X12 | 2    | ICJZ1012    | JAY ZIP INNER COVER 10X12 |
| 1    | CJZ1213     | JAY ZIP OUTER COVER 12X13 | 2    | ICJZ1013    | JAY ZIP INNER COVER 10X13 |
| 1    | CJZ1214     | JAY ZIP OUTER COVER 12X14 | 2    | ICJZ1014    | JAY ZIP INNER COVER 10X14 |
| 1    | CJZ1215     | JAY ZIP OUTER COVER 12X15 | 2    | ICJZ1015    | JAY ZIP INNER COVER 10X15 |
| 1    | CJZ1216     | JAY ZIP OUTER COVER 12X16 | 2    | ICJZ1016    | JAY ZIP INNER COVER 10X16 |
| 1    | CJZ1217     | JAY ZIP OUTER COVER 12X17 | 2    | ICJZ1108    | JAY ZIP INNER COVER 11X08 |
| 1    | CJZ1218     | JAY ZIP OUTER COVER 12X18 | 2    | ICJZ1109    | JAY ZIP INNER COVER 11X09 |
| 1    | CJZ1310     | JAY ZIP OUTER COVER 13X10 | 2    | ICJZ1110    | JAY ZIP INNER COVER 11X10 |
| 1    | CJZ1311     | JAY ZIP OUTER COVER 13X11 | 2    | ICJZ1111    | JAY ZIP INNER COVER 11X11 |
| 1    | CJZ1312     | JAY ZIP OUTER COVER 13X12 | 2    | ICJZ1112    | JAY ZIP INNER COVER 11X12 |
| 1    | CJZ1313     | JAY ZIP OUTER COVER 13X13 | 2    | ICJZ1113    | JAY ZIP INNER COVER 11X13 |
| 1    | CJZ1314     | JAY ZIP OUTER COVER 13X14 | 2    | ICJZ1114    | JAY ZIP INNER COVER 11X14 |
| 1    | CJZ1315     | JAY ZIP OUTER COVER 13X15 | 2    | ICJZ1115    | JAY ZIP INNER COVER 11X15 |
| 1    | CJZ1316     | JAY ZIP OUTER COVER 13X16 | 2    | ICJZ1116    | JAY ZIP INNER COVER 11X16 |
| 1    | CJZ1317     | JAY ZIP OUTER COVER 13X17 | 2    | ICJZ1117    | JAY ZIP INNER COVER 11X17 |
| 1    | CJZ1318     | JAY ZIP OUTER COVER 13X18 | 2    | ICJZ1209    | JAY ZIP INNER COVER 12X09 |
| 1    | CJZ1411     | JAY ZIP OUTER COVER 14X11 | 2    | ICJZ1210    | JAY ZIP INNER COVER 12X10 |
| 1    | CJZ1412     | JAY ZIP OUTER COVER 14X12 | 2    | ICJZ1211    | JAY ZIP INNER COVER 12X11 |
| 1    | CJZ1413     | JAY ZIP OUTER COVER 14X13 | 2    | ICJZ1212    | JAY ZIP INNER COVER 12X12 |
| 1    | CJZ1414     | JAY ZIP OUTER COVER 14X14 | 2    | ICJZ1213    | JAY ZIP INNER COVER 12X13 |
| 1    | CJZ1415     | JAY ZIP OUTER COVER 14X15 | 2    | ICJZ1214    | JAY ZIP INNER COVER 12X14 |
| 1    | CJZ1416     | JAY ZIP OUTER COVER 14X16 | 2    | ICJZ1215    | JAY ZIP INNER COVER 12X15 |
| 1    | CJZ1417     | JAY ZIP OUTER COVER 14X17 | 2    | ICJZ1216    | JAY ZIP INNER COVER 12X16 |
| 1    | CJZ1418     | JAY ZIP OUTER COVER 14X18 | 2    | ICJZ1217    | JAY ZIP INNER COVER 12X17 |
| 1    | CJZ1512     | JAY ZIP OUTER COVER 15X12 | 2    | ICJZ1218    | JAY ZIP INNER COVER 12X18 |
| 1    | CJZ1513     | JAY ZIP OUTER COVER 15X13 | 2    | ICJZ1310    | JAY ZIP INNER COVER 13X10 |
| 1    | CJZ1514     | JAY ZIP OUTER COVER 15X14 | 2    | ICJZ1311    | JAY ZIP INNER COVER 13X11 |
| 1    | CJZ1515     | JAY ZIP OUTER COVER 15X15 | 2    | ICJZ1312    | JAY ZIP INNER COVER 13X12 |
| 1    | CJZ1516     | JAY ZIP OUTER COVER 15X16 | 2    | ICJZ1313    | JAY ZIP INNER COVER 13X13 |
| 1    | CJZ1517     | JAY ZIP OUTER COVER 15X17 | 2    | ICJZ1314    | JAY ZIP INNER COVER 13X14 |
| 1    | CJZ1518     | JAY ZIP OUTER COVER 15X18 | 2    | ICJZ1315    | JAY ZIP INNER COVER 13X15 |
| 1    | CJZ1613     | JAY ZIP OUTER COVER 16X13 | 2    | ICJZ1316    | JAY ZIP INNER COVER 13X16 |
| 1    | CJZ1614     | JAY ZIP OUTER COVER 16X14 | 2    | ICJZ1317    | JAY ZIP INNER COVER 13X17 |
| 1    | CJZ1615     | JAY ZIP OUTER COVER 16X15 | 2    | ICJZ1318    | JAY ZIP INNER COVER 13X18 |
| 1    | CJZ1616     | JAY ZIP OUTER COVER 16X16 | 2    | ICJZ1411    | JAY ZIP INNER COVER 14X11 |
| 1    | CJZ1617     | JAY ZIP OUTER COVER 16X17 | 2    | ICJZ1412    | JAY ZIP INNER COVER 14X12 |
| 1    | CJZ1618     | JAY ZIP OUTER COVER 16X18 | 2    | ICJZ1413    | JAY ZIP INNER COVER 14X13 |
|      |             | INNER COVER               | 2    | ICJZ1414    | JAY ZIP INNER COVER 14X14 |
| 2    | ICJZ0808    | JAY ZIP INNER COVER 08X08 | 2    | ICJZ1415    | JAY ZIP INNER COVER 14X15 |
| 2    | ICJZ0809    | JAY ZIP INNER COVER 08X09 | 2    | ICJZ1416    | JAY ZIP INNER COVER 14X16 |
| 2    | ICJZ0810    | JAY ZIP INNER COVER 08X10 | 2    | ICJZ1417    | JAY ZIP INNER COVER 14X17 |
| 2    | ICJZ0811    | JAY ZIP INNER COVER 08X11 | 2    | ICJZ1418    | JAY ZIP INNER COVER 14X18 |
| 2    | ICJZ0812    | JAY ZIP INNER COVER 08X12 | 2    | ICJZ1512    | JAY ZIP INNER COVER 15X12 |
| 2    | ICJZ0813    | JAY ZIP INNER COVER 08X13 | 2    | ICJZ1513    | JAY ZIP INNER COVER 15X13 |
| 2    | ICJZ0814    | JAY ZIP INNER COVER 08X14 | 2    | ICJZ1514    | JAY ZIP INNER COVER 15X14 |
| 2    | ICJZ0908    | JAY ZIP INNER COVER 09X08 | 2    | ICJZ1515    | JAY ZIP INNER COVER 15X15 |
| 2    | ICJZ0909    | JAY ZIP INNER COVER 09X09 | 2    | ICJZ1516    | JAY ZIP INNER COVER 15X16 |
| 2    | ICJZ0910    | JAY ZIP INNER COVER 09X10 | 2    | ICJZ1517    | JAY ZIP INNER COVER 15X17 |
| 2    | ICJZ0911    | JAY ZIP INNER COVER 09X11 | 2    | ICJZ1518    | JAY ZIP INNER COVER 15X18 |
| 2    | ICJZ0912    | JAY ZIP INNER COVER 09X12 | 2    | ICJZ1613    | JAY ZIP INNER COVER 16X13 |
| 2    | ICJZ0913    | JAY ZIP INNER COVER 09X13 | 2    | ICJZ1614    | JAY ZIP INNER COVER 16X14 |

| Pos. | Item Number | Description               | Pos. | Item Number | Description        |
|------|-------------|---------------------------|------|-------------|--------------------|
| 2    | ICJZ1615    | JAY ZIP INNER COVER 16X15 | 3    | BJZ1210     | JAY ZIP BASE 12X10 |
| 2    | ICJZ1616    | JAY ZIP INNER COVER 16X16 | 3    | BJZ1211     | JAY ZIP BASE 12X11 |
| 2    | ICJZ1617    | JAY ZIP INNER COVER 16X17 | 3    | BJZ1212     | JAY ZIP BASE 12X12 |
| 2    | ICJZ1618    | JAY ZIP INNER COVER 16X18 | 3    | BJZ1213     | JAY ZIP BASE 12X13 |
|      |             | BASE                      | 3    | BJZ1214     | JAY ZIP BASE 12X14 |
| 3    | BJZ0808     | JAY ZIP BASE 08X08        | 3    | BJZ1215     | JAY ZIP BASE 12X15 |
| 3    | BJZ0809     | JAY ZIP BASE 08X09        | 3    | BJZ1216     | JAY ZIP BASE 12X16 |
| 3    | BJZ0810     | JAY ZIP BASE 08X10        | 3    | BJZ1217     | JAY ZIP BASE 12X17 |
| 3    | BJZ0811     | JAY ZIP BASE 08X11        | 3    | BJZ1218     | JAY ZIP BASE 12X18 |
| 3    | BJZ0812     | JAY ZIP BASE 08X12        | 3    | BJZ1310     | JAY ZIP BASE 13X10 |
| 3    | BJZ0813     | JAY ZIP BASE 08X13        | 3    | BJZ1311     | JAY ZIP BASE 13X11 |
| 3    | BJZ0814     | JAY ZIP BASE 08X14        | 3    | BJZ1312     | JAY ZIP BASE 13X12 |
| 3    | BJZ0908     | JAY ZIP BASE 09X08        | 3    | BJZ1313     | JAY ZIP BASE 13X13 |
| 3    | BJZ0909     | JAY ZIP BASE 09X09        | 3    | BJZ1314     | JAY ZIP BASE 13X14 |
| 3    | BJZ0910     | JAY ZIP BASE 09X10        | 3    | BJZ1315     | JAY ZIP BASE 13X15 |
| 3    | BJZ0911     | JAY ZIP BASE 09X11        | 3    | BJZ1316     | JAY ZIP BASE 13X16 |
| 3    | BJZ0912     | JAY ZIP BASE 09X12        | 3    | BJZ1317     | JAY ZIP BASE 13X17 |
| 3    | BJZ0913     | JAY ZIP BASE 09X13        | 3    | BJZ1318     | JAY ZIP BASE 13X18 |
| 3    | BJZ0914     | JAY ZIP BASE 09X14        | 3    | BJZ1411     | JAY ZIP BASE 14X11 |
| 3    | BJZ0915     | JAY ZIP BASE 09X15        | 3    | BJZ1412     | JAY ZIP BASE 14X12 |
| 3    | BJZ1008     | JAY ZIP BASE 10X08        | 3    | BJZ1413     | JAY ZIP BASE 14X13 |
| 3    | BJZ1009     | JAY ZIP BASE 10X09        | 3    | BJZ1414     | JAY ZIP BASE 14X14 |
| 3    | BJZ1010     | JAY ZIP BASE 10X10        | 3    | BJZ1415     | JAY ZIP BASE 14X15 |
| 3    | BJZ1011     | JAY ZIP BASE 10X11        | 3    | BJZ1416     | JAY ZIP BASE 14X16 |
| 3    | BJZ1012     | JAY ZIP BASE 10X12        | 3    | BJZ1417     | JAY ZIP BASE 14X17 |
| 3    | BJZ1013     | JAY ZIP BASE 10X13        | 3    | BJZ1418     | JAY ZIP BASE 14X18 |
| 3    | BJZ1014     | JAY ZIP BASE 10X14        | 3    | BJZ1512     | JAY ZIP BASE 15X12 |
| 3    | BJZ1015     | JAY ZIP BASE 10X15        | 3    | BJZ1513     | JAY ZIP BASE 15X13 |
| 3    | BJZ1016     | JAY ZIP BASE 10X16        | 3    | BJZ1514     | JAY ZIP BASE 15X14 |
| 3    | BJZ1108     | JAY ZIP BASE 11X08        | 3    | BJZ1515     | JAY ZIP BASE 15X15 |
| 3    | BJZ1109     | JAY ZIP BASE 11X09        | 3    | BJZ1516     | JAY ZIP BASE 15X16 |
| 3    | BJZ1110     | JAY ZIP BASE 11X10        | 3    | BJZ1517     | JAY ZIP BASE 15X17 |
| 3    | BJZ1111     | JAY ZIP BASE 11X11        | 3    | BJZ1518     | JAY ZIP BASE 15X18 |
| 3    | BJZ1112     | JAY ZIP BASE 11X12        | 3    | BJZ1613     | JAY ZIP BASE 16X13 |
| 3    | BJZ1113     | JAY ZIP BASE 11X13        | 3    | BJZ1614     | JAY ZIP BASE 16X14 |
| 3    | BJZ1114     | JAY ZIP BASE 11X14        | 3    | BJZ1615     | JAY ZIP BASE 16X15 |
| 3    | BJZ1115     | JAY ZIP BASE 11X15        | 3    | BJZ1616     | JAY ZIP BASE 16X16 |
| 3    | BJZ1116     | JAY ZIP BASE 11X16        | 3    | BJZ1617     | JAY ZIP BASE 16X17 |
| 3    | BJZ1117     | JAY ZIP BASE 11X17        | 3    | BJZ1618     | JAY ZIP BASE 16X18 |
| 3    | BJZ1209     | JAY ZIP BASE 12X09        |      |             |                    |

Note: Custom Drop Seat is only available for the Jay Soft Combi & Soft Combi P, Jay Gel, Jay Easy, Jay 2 & J2 DC, Jay GS and Jay Triad Cushions.



| Pos. | Item Number | Description                   | Remarks |
|------|-------------|-------------------------------|---------|
| 1    | SD2444      | SOLID DROP SEAT 14X14 W/HRDWR |         |
| 1    | SD2446      | SOLID DROP SEAT 14X16 W/HRDWR |         |
| 1    | SD2455      | SOLID DROP SEAT 15X15 W/HRDWR |         |
| 1    | SD2457      | SOLID DROP SEAT 15X17 W/HRDWR |         |
| 1    | SD2460      | SOLID DROP SEAT 16X20 W/HRDWR |         |
| 1    | SD2466      | SOLID DROP SEAT 16X16 W/HRDWR |         |
| 1    | SD2468      | SOLID DROP SEAT 16X18 W/HRDWR |         |
| 1    | SD2475      | SOLID DROP SEAT 17X15 W/HRDWR |         |
| 1    | SD2477      | SOLID DROP SEAT 17X17 W/HRDWR |         |
| 1    | SD2480      | SOLID DROP SEAT 18X20 W/HRDWR |         |
| 1    | SD2488      | SOLID DROP SEAT 18X18 W/HRDWR |         |
| 1    | SD2486      | SOLID DROP SEAT 18X16 W/HRDWR |         |
| 1    | SD2406      | SOLID DROP SEAT 20X16 W/HRDWR |         |
| 1    | SD2408      | SOLID DROP SEAT 20X18 W/HRDWR |         |

| Pos. | Item Number | Description                      | Remarks   |
|------|-------------|----------------------------------|-----------|
| 1    | SD2400      | SOLID DROP SEAT 20X20 W/HRDWR    |           |
| 1    | SD2426      | SOLID DROP SEAT 22X16 W/HRDWR    |           |
| 1    | SD2428      | SOLID DROP SEAT 22X18 W/HRDWR    |           |
| 1    | SD2448      | SOLID DROP SEAT 24X18 W/HRDWR    |           |
| 2-4  | M-24S       | SHORT 1/2" DROP HOOKS            | Pack of 4 |
| 2-4  | M-24M       | MEDIUM 1" DROP HOOKS             | Pack of 4 |
| 2-4  | M-24L       | LONG 1 1/2" DROP HOOKS           | Pack of 4 |
| 5    | 283         | HOOK STOPPER W/HARDWARE          |           |
| 6    | 282         | HOOK RETAINER/STOP WITH HARDWARE |           |



| Pos. | Item Number | Description                         | Remarks |
|------|-------------|-------------------------------------|---------|
| 1A   | C120A       | JAY ACTIVE COVER AIR EXCHANGE 18x16 |         |
| 1A   | C121A       | JAY ACTIVE COVER AIR EXCHANGE 18x18 |         |
| 1A   | C125A       | JAY ACTIVE COVER AIR EXCHANGE 16x16 |         |
| 1A   | C126A       | JAY ACTIVE COVER AIR EXCHANGE 16x18 |         |
| 1A   | C128A       | JAY ACTIVE COVER AIR EXCHANGE 14x16 |         |
| 1A   | C122A       | JAY ACTIVE COVER AIR EXCHANGE 20x18 |         |
| 1A   | C123A       | JAY ACTIVE COVER AIR EXCHANGE 18x20 |         |
| 1A   | C127A       | JAY ACTIVE COVER AIR EXCHANGE 16x20 |         |
| 1B   | C120N       | JAY ACTIVE COVER INCONTINENT 18x16  |         |
| 1B   | C121N       | JAY ACTIVE COVER INCONTINENT 18x18  |         |
| 1B   | C125N       | JAY ACTIVE COVER INCONTINENT 16x16  |         |
| 1B   | C126N       | JAY ACTIVE COVER INCONTINENT 16x18  |         |
| 1B   | C128N       | JAY ACTIVE COVER INCONTINENT 14x16  |         |
| 1B   | C122N       | JAY ACTIVE COVER INCONTINENT 20x18  |         |
| 1B   | C123N       | JAY ACTIVE COVER INCONTINENT 18x20  |         |
| 1B   | C127N       | JAY ACTIVE COVER INCONTINENT 16x20  |         |
| 2    | F120        | JAY ACTIVE FLUID PAD 18x16          |         |
| 2    | F121        | JAY ACTIVE FLUID PAD 18x18          |         |
| 2    | F125        | JAY ACTIVE FLUID PAD 16x16          |         |
| 2    | F126        | JAY ACTIVE FLUID PAD 16x18          |         |
| 2    | F128        | JAY ACTIVE FLUID PAD 14x16          |         |
| 2    | F122        | JAY ACTIVE FLUID PAD 20x18          |         |

| Pos. | Item Number | Description                 | Remarks                          |
|------|-------------|-----------------------------|----------------------------------|
| 3A   | B119JA      | 1" ADDUCTOR WEDGER (PR)     | Pair-ships standard with cushion |
| 3B   | B113JA      | 2" ADDUCTOR WEDGES (PR)     | Pair                             |
| 4    | CB120       | JAY ACTIVE BASE COVER 18x16 |                                  |
| 4    | CB121       | JAY ACTIVE BASE COVER 18x18 |                                  |
| 4    | CB125       | JAY ACTIVE BASE COVER 16x16 |                                  |
| 4    | CB126       | JAY ACTIVE BASE COVER 16x18 |                                  |
| 4    | CB128       | JAY ACTIVE BASE COVER 14x16 |                                  |
| 4    | CB122       | JAY ACTIVE BASE COVER 20x18 |                                  |
| 4    | CB123       | JAY ACTIVE BASE COVER 18x20 |                                  |
| 4    | CB127       | JAY ACTIVE BASE COVER 16x20 |                                  |
| 3A-5 | B120        | JAY ACTIVE BASE 18x16       |                                  |
| 3A-5 | B121        | JAY ACTIVE BASE 18x18       |                                  |
| 3A-5 | B125        | JAY ACTIVE BASE 16x16       |                                  |
| 3A-5 | B126        | JAY ACTIVE BASE 16x18       |                                  |
| 3A-5 | B128        | JAY ACTIVE BASE 14x16       |                                  |
| 3A-5 | B122        | JAY ACTIVE BASE 20x18       |                                  |
| 3A-5 | B123        | JAY ACTIVE BASE 18x20       |                                  |
| 3A-5 | B127        | JAY ACTIVE BASE 16x20       |                                  |



improving people's lives



Dear Valued Customer,

Effective April 6th, 2015, the Jay Cushion has been made obsolete and parts will no longer be offered by Sunrise Medical. This product has been discontinued and due to the diminished demand for its parts Sunrise Medical can no longer offer replacement parts.

Sincerely,

Sunrise Medical

2842 Business Park Avenue • Fresno, CA 93727-1328 • USA

800.333.4000 [www.SunriseMedical.com](http://www.SunriseMedical.com)

| Description      |
|------------------|
| OBSOLETE PRODUCT |



improving people's lives



Dear Valued Customer,

Effective April 6th, 2015, the Jay Triad Cushion has been made obsolete and parts will no longer be offered by Sunrise Medical. This product has been discontinued and due to the diminished demand for its parts Sunrise Medical can no longer offer replacement parts.

Sincerely,

Sunrise Medical

2842 Business Park Avenue • Fresno, CA 93727-1328 • USA

800.333.4000 [www.SunriseMedical.com](http://www.SunriseMedical.com)

| Description      |
|------------------|
| OBSOLETE PRODUCT |



improving people's lives



Dear Valued Customer,

Effective April 6th, 2015, the Jay Triad Cushion has been made obsolete and parts will no longer be offered by Sunrise Medical. This product has been discontinued and due to the diminished demand for its parts Sunrise Medical can no longer offer replacement parts.


Sincerely,

Sunrise Medical

2842 Business Park Avenue • Fresno, CA 93727-1328 • USA

800.333.4000 [www.SunriseMedical.com](http://www.SunriseMedical.com)

| Description      |
|------------------|
| OBSOLETE PRODUCT |



| Pos. | Item Number | Description                  | Remarks | Pos. | Item Number | Description                    | Remarks |
|------|-------------|------------------------------|---------|------|-------------|--------------------------------|---------|
| 1A   | C900T       | XTREME COVER TRITEX 20X20    |         | 1B   | C966LS      | XTREME COVER LOW SHEAR 16X16   |         |
| 1A   | C944T       | XTREME COVER TRITEX 14X14    |         | 1B   | C967LS      | XTREME COVER LOW SHEAR 16X17   |         |
| 1A   | C946T       | XTREME COVER TRITEX 14X16    |         | 1B   | C968LS      | XTREME COVER LOW SHEAR 16X18   |         |
| 1A   | C955T       | XTREME COVER TRITEX 15X15    |         | 1B   | C976LS      | XTREME COVER LOW SHEAR 17X16   |         |
| 1A   | C956T       | XTREME COVER TRITEX 15X16    |         | 1B   | C977LS      | XTREME COVER LOW SHEAR 17X17   |         |
| 1A   | C960T       | XTREME COVER TRITEX 16X20    |         | 1B   | C980LS      | XTREME COVER LOW SHEAR 18X20   |         |
| 1A   | C966T       | XTREME COVER TRITEX 16X16    |         | 1B   | C986LS      | XTREME COVER LOW SHEAR 18X16   |         |
| 1A   | C967T       | XTREME COVER TRITEX 16X17    |         | 1B   | C988LS      | XTREME COVER LOW SHEAR 18X18   |         |
| 1A   | C968T       | XTREME COVER TRITEX 16X18    |         | 1C   | C900T       | XTREME COVER INCONTINENT 20X20 |         |
| 1A   | C976T       | XTREME COVER TRITEX 17X16    |         | 1C   | C944T       | XTREME COVER INCONTINENT 14X14 |         |
| 1A   | C977T       | XTREME COVER TRITEX 17X17    |         | 1C   | C946T       | XTREME COVER INCONTINENT 14X16 |         |
| 1A   | C980T       | XTREME COVER TRITEX 18X20    |         | 1C   | C955T       | XTREME COVER INCONTINENT 15X15 |         |
| 1A   | C986T       | XTREME COVER TRITEX 18X16    |         | 1C   | C956T       | XTREME COVER INCONTINENT 15X16 |         |
| 1A   | C988T       | XTREME COVER TRITEX 18X18    |         | 1C   | C960T       | XTREME COVER INCONTINENT 16X20 |         |
| 1B   | C900LS      | XTREME COVER LOW SHEAR 20X20 |         | 1C   | C966T       | XTREME COVER INCONTINENT 16X16 |         |
| 1B   | C944LS      | XTREME COVER LOW SHEAR 14X14 |         | 1C   | C967T       | XTREME COVER INCONTINENT 16X17 |         |
| 1B   | C946LS      | XTREME COVER LOW SHEAR 14X16 |         | 1C   | C968T       | XTREME COVER INCONTINENT 16X18 |         |
| 1B   | C955LS      | XTREME COVER LOW SHEAR 15X15 |         | 1C   | C976T       | XTREME COVER INCONTINENT 17X16 |         |
| 1B   | C956LS      | XTREME COVER LOW SHEAR 15X16 |         | 1C   | C977T       | XTREME COVER INCONTINENT 17X17 |         |
| 1B   | C960LS      | XTREME COVER LOW SHEAR 16X20 |         | 1C   | C980T       | XTREME COVER INCONTINENT 18X20 |         |

| Pos. | Item Number | Description              | Remarks | Pos. | Item Number | Description             | Remarks |
|------|-------------|--------------------------|---------|------|-------------|-------------------------|---------|
| 1C   | C986T       | XTREME COVER INCONTINENT | 18X16   | 4    | CB988       | XTREME BASE COVER 18X18 |         |
| 1C   | C988T       | XTREME COVER INCONTINENT | 18X18   | 4,5  | B900        | XTREME BASE 20X20       |         |
| 2    | F900        | XTREME FLUID PAD 20X20   |         | 4,5  | B944        | XTREME BASE 14X14       |         |
| 2    | F944        | XTREME FLUID PAD 14X14   |         | 4,5  | B946        | XTREME BASE 14X16       |         |
| 2    | F946        | XTREME FLUID PAD 14X16   |         | 4,5  | B955        | XTREME BASE 15X15       |         |
| 2    | F955        | XTREME FLUID PAD 15X15   |         | 4,5  | B956        | XTREME BASE 15X16       |         |
| 2    | F956        | XTREME FLUID PAD 15X16   |         | 4,5  | B960        | XTREME BASE 16X20       |         |
| 2    | F960        | XTREME FLUID PAD 16X20   |         | 4,5  | B966        | XTREME BASE 16X16       |         |
| 2    | F966        | XTREME FLUID PAD 16X16   |         | 4,5  | B967        | XTREME BASE 16X17       |         |
| 2    | F967        | XTREME FLUID PAD 16X17   |         | 4,5  | B968        | XTREME BASE 16X18       |         |
| 2    | F968        | XTREME FLUID PAD 16X18   |         | 4,5  | B976        | XTREME BASE 17X16       |         |
| 2    | F976        | XTREME FLUID PAD 17X16   |         | 4,5  | B977        | XTREME BASE 17X17       |         |
| 2    | F977        | XTREME FLUID PAD 17X17   |         | 4,5  | B980        | XTREME BASE 18X20       |         |
| 2    | F980        | XTREME FLUID PAD 18X20   |         | 4,5  | B986        | XTREME BASE 18X16       |         |
| 2    | F986        | XTREME FLUID PAD 18X16   |         | 4,5  | B988        | XTREME BASE 18X18       |         |
| 2    | F988        | XTREME FLUID PAD 18X18   |         | 6    | S900        | SOLID SEAT XTM 20X20    |         |
| 3    | B119JA      | 1" ADDUCTOR WEDGES       | PAIR    | 6    | S944        | SOLID SEAT XTM 14X14    |         |
| 4    | CB900       | XTREME BASE COVER 20X20  |         | 6    | S946        | SOLID SEAT XTM 14X16    |         |
| 4    | CB944       | XTREME BASE COVER 14X14  |         | 6    | S955        | SOLID SEAT XTM 15X15    |         |
| 4    | CB946       | XTREME BASE COVER 14X16  |         | 6    | S956        | SOLID SEAT XTM 15X16    |         |
| 4    | CB955       | XTREME BASE COVER 15X15  |         | 6    | S960        | SOLID SEAT XTM 16X20    |         |
| 4    | CB956       | XTREME BASE COVER 15X16  |         | 6    | S966        | SOLID SEAT XTM 16X16    |         |
| 4    | CB960       | XTREME BASE COVER 16X20  |         | 6    | S967        | SOLID SEAT XTM 16X17    |         |
| 4    | CB966       | XTREME BASE COVER 16X16  |         | 6    | S968        | SOLID SEAT XTM 16X18    |         |
| 4    | CB967       | XTREME BASE COVER 16X17  |         | 6    | S976        | SOLID SEAT XTM 17X16    |         |
| 4    | CB968       | XTREME BASE COVER 16X18  |         | 6    | S977        | SOLID SEAT XTM 17X17    |         |
| 4    | CB976       | XTREME BASE COVER 17X16  |         | 6    | S980        | SOLID SEAT XTM 18X20    |         |
| 4    | CB977       | XTREME BASE COVER 17X17  |         | 6    | S986        | SOLID SEAT XTM 18X16    |         |
| 4    | CB980       | XTREME BASE COVER 18X20  |         | 6    | S988        | SOLID SEAT XTM 18X18    |         |
| 4    | CB986       | XTREME BASE COVER 18X16  |         |      |             |                         |         |



improving people's lives



Dear Valued Customer,

Effective April 6th, 2015, the Jay Xtreme Cushion has been made obsolete and parts will no longer be offered by Sunrise Medical. This product has been discontinued and due to the diminished demand for its parts Sunrise Medical can no longer offer replacement parts.

Sincerely,

Sunrise Medical

2842 Business Park Avenue • Fresno, CA 93727-1328 • USA

800.333.4000 [www.SunriseMedical.com](http://www.SunriseMedical.com)

| Description      |
|------------------|
| OBSOLETE PRODUCT |

improving people's lives



Dear Valued Customer,

As of **6/1/2015**, the sale of **Jay J-Gel Cushion** has been discontinued. In compliance with Sunrise Medical's policy, replacement parts will continue to be available until **6/1/2021**. Effective **6/1/2021**, replacement parts for the **Jay J-Gel Cushion** will be obsolete and therefore no longer available for sale.

Sincerely,


Sunrise Medical

6899 Winchester Circle, Suite 200 • Boulder, CO 80301 • USA

800.333.4000 [www.SunriseMedical.com](http://www.SunriseMedical.com)

| Description  |
|--------------|
| DISCONTINUED |

**Note:** For exact cushion dimensions and measuring guidelines please reference the Jay Backs and Cushions Brochure. (Literature part number 0110.)



| Pos. | Item Number | Description                     | Pos. | Item Number | Description                     |
|------|-------------|---------------------------------|------|-------------|---------------------------------|
| 1    | C4206LS     | COVER J-GEL LOW SHEAR 20x16     | 2    | G4248       | GEL PAD J-GEL REPLACEMENT 14x18 |
| 1    | C4208LS     | COVER J-GEL LOW SHEAR 20x18     | 2    | G4260       | GEL PAD J-GEL REPLACEMENT 16x20 |
| 1    | C4246LS     | COVER J-GEL LOW SHEAR 14x16     | 2    | G4266       | GEL PAD J-GEL REPLACEMENT 16x16 |
| 1    | C4248LS     | COVER J-GEL LOW SHEAR 14x18     | 2    | G4268       | GEL PAD J-GEL REPLACEMENT 16x18 |
| 1    | C4260LS     | COVER J-GEL LOW SHEAR 16x20     | 2    | G4270       | GEL PAD J-GEL REPLACEMENT 17x20 |
| 1    | C4266LS     | COVER J-GEL LOW SHEAR 16x16     | 2    | G4276       | GEL PAD J-GEL REPLACEMENT 17x16 |
| 1    | C4268LS     | COVER J-GEL LOW SHEAR 16x18     | 2    | G4278       | GEL PAD J-GEL REPLACEMENT 17x18 |
| 1    | C4270LS     | COVER J-GEL LOW SHEAR 17x20     | 3    | G4280       | GEL PAD J-GEL REPLACEMENT 18x20 |
| 1    | C4276LS     | COVER J-GEL LOW SHEAR 17x16     | 2    | G4286       | GEL PAD J-GEL REPLACEMENT 18x16 |
| 1    | C4278LS     | COVER J-GEL LOW SHEAR 17x18     | 2    | G4288       | GEL PAD J-GEL REPLACEMENT 18x18 |
| 1    | C4280LS     | COVER J-GEL LOW SHEAR 18x20     | 3    | B4206       | BASE J-GEL 20x16                |
| 1    | C4286LS     | COVER J-GEL LOW SHEAR 18x16     | 3    | B4208       | BASE J-GEL 20x18                |
| 1    | C4288LS     | COVER J-GEL LOW SHEAR 18x18     | 3    | B4246       | BASE J-GEL 14x16                |
| 2    | G4206       | GEL PAD J-GEL REPLACEMENT 20x16 | 3    | B4248       | BASE J-GEL 14x18                |
| 2    | G4208       | GEL PAD J-GEL REPLACEMENT 20x18 | 3    | B4260       | BASE J-GEL 16x20                |
| 2    | G4246       | GEL PAD J-GEL REPLACEMENT 14x16 | 3    | B4266       | BASE J-GEL 16x16                |

| Pos. | Item Number | Description            | Pos. | Item Number | Description            |
|------|-------------|------------------------|------|-------------|------------------------|
| 3    | B4268       | BASE J-GEL 16x18       | 4    | S4260       | SOLID SEAT J-GEL 16x20 |
| 3    | B4270       | BASE J-GEL 17x20       | 4    | S4266       | SOLID SEAT J-GEL 16x16 |
| 3    | B4276       | BASE J-GEL 17x16       | 4    | S4268       | SOLID SEAT J-GEL 16x18 |
| 3    | B4278       | BASE J-GEL 17x18       | 4    | S4270       | SOLID SEAT J-GEL 17x20 |
| 3    | B4280       | BASE J-GEL 18x20       | 4    | S4276       | SOLID SEAT J-GEL 17x16 |
| 3    | B4286       | BASE J-GEL 18x16       | 4    | S4278       | SOLID SEAT J-GEL 17x18 |
| 3    | B4288       | BASE J-GEL 18x18       | 4    | S4280       | SOLID SEAT J-GEL 18x20 |
| 4    | S4206       | SOLID SEAT J-GEL 20x16 | 4    | S4286       | SOLID SEAT J-GEL 18x16 |
| 4    | S4208       | SOLID SEAT J-GEL 20x18 | 4    | S4288       | SOLID SEAT J-GEL 18x18 |
| 4    | S4246       | SOLID SEAT J-GEL 14x16 | 5    | M-4204      | REPAIR KIT J-GEL       |
| 4    | S4248       | SOLID SEAT J-GEL 14x18 |      |             |                        |

improving people's lives



Dear Valued Customer,

As of **6/1/2015**, the sale of **Jay Soft Combi Cushion** has been discontinued. In compliance with Sunrise Medical's policy, replacement parts will continue to be available until **6/1/2021**. Effective **6/1/2021**, replacement parts for the **Jay Soft Combi Cushion** will be obsolete and therefore no longer available for sale.


Sincerely,

Sunrise Medical

6899 Winchester Circle, Suite 200 • Boulder, CO 80301 • USA

800.333.4000 [www.SunriseMedical.com](http://www.SunriseMedical.com)

| Description  |
|--------------|
| DISCONTINUED |



| Pos. | Item Number | Description                      | Remarks                 |
|------|-------------|----------------------------------|-------------------------|
| 1,2  | 2250        | LUMBAR SUPPORT SOFT COMBI        | Ships standard with 1-4 |
| 1,2  | 2251-MJ     | LUMBAR SUPPORT SOFT COMBI JUNIOR |                         |
| 3    | C2200-MJ    | SOFT COMBI COVER 18x16           |                         |
| 3    | C2201       | SOFT COMBI COVER 18x18           |                         |
| 3    | C2202       | SOFT COMBI COVER 20x18           |                         |
| 3    | C2203       | SOFT COMBI COVER 18x20           |                         |
| 3    | C2205       | SOFT COMBI COVER 15.5x16         |                         |
| 3    | C2206       | SOFT COMBI COVER 15.5x18         |                         |
| 3    | C2207       | SOFT COMBI COVER 16x20           |                         |
| 3    | C2210       | SOFT COMBI COVER 14x14           |                         |
| 3    | C2215       | SOFT COMBI COVER 12x12           |                         |
| 5    | S200        | COMBI SOLID SEAT 18x16           |                         |
| 5    | S201        | COMBI SOLID SEAT 18x18           |                         |
| 5    | S202        | COMBI SOLID SEAT 20x18           |                         |
| 5    | S203        | COMBI SOLID SEAT 18x20           |                         |
| 5    | S205        | COMBI SOLID SEAT 15.5x16         |                         |

| Pos. | Item Number | Description              | Remarks |
|------|-------------|--------------------------|---------|
| 5    | S206        | COMBI SOLID SEAT 15.5x18 |         |
| 5    | S207        | COMBI SOLID SEAT 16x20   |         |
| 5    | S210        | COMBI SOLID SEAT 14x14   |         |
| 5    | S215        | COMBI SOLID SEAT 12x12   |         |

improving people's lives



Dear Valued Customer,

As of **1/1/2011**, the sale of **Jay Combi Cushion** has been discontinued. In compliance with Sunrise Medical's policy, replacement parts will continue to be available until **1/1/2016**. Effective **1/1/2016**, replacement parts for the **Jay Combi Cushion** will be obsolete and therefore no longer available for sale.


Sincerely,

Sunrise Medical

6899 Winchester Circle, Suite 200 • Boulder, CO 80301 • USA

800.333.4000 [www.SunriseMedical.com](http://www.SunriseMedical.com)

| Description  |
|--------------|
| DISCONTINUED |




| Pos. | Item Number | Description                 | Remarks                                    |
|------|-------------|-----------------------------|--------------------------------------------|
| 1,2  | 250XB       | LUMBAR SUPPORT COMBI ADULT  | SIZE 15 X 9" Ships standard with 1,2A,3A,4 |
| 1,2  | 251XB       | LUMBAR SUPPORT COMBI JUNIOR | SIZE 12 X 7"                               |
| 3    | C200NB      | COMBI COVER 18x16           |                                            |
| 3    | C201NB      | COMBI COVER 18x18           |                                            |
| 3    | C202NB      | COMBI COVER 20x18           |                                            |
| 3    | C203NB      | COMBI COVER 18x20           |                                            |
| 3    | C205NB      | COMBI COVER 15.5x16         |                                            |
| 3    | C206NB      | COMBI COVER 15.5x18         |                                            |
| 3    | C207NB      | COMBI COVER 16x20           |                                            |
| 3    | C210NB      | COMBI COVER 14x14           |                                            |
| 3    | C215NB      | COMBI COVER 12x12           |                                            |
| 5    | S200        | COMBI SOLID SEAT 18x16      |                                            |
| 5    | S201        | COMBI SOLID SEAT 18x18      |                                            |
| 5    | S202        | COMBI SOLID SEAT 20x18      |                                            |
| 5    | S203        | COMBI SOLID SEAT 18x20      |                                            |
| 5    | S205        | COMBI SOLID SEAT 15.5x16    |                                            |

| Pos. | Item Number | Description              | Remarks |
|------|-------------|--------------------------|---------|
| 5    | S206        | COMBI SOLID SEAT 15.5x18 |         |
| 5    | S207        | COMBI SOLID SEAT 16x20   |         |
| 5    | S210        | COMBI SOLID SEAT 14x14   |         |
| 5    | S215        | COMBI SOLID SEAT 12x12   |         |

NOTE: MAY ALSO BE USED WITH JAY AND JAY J2 CUSHIONS. ITEMS 1 & 5 MAY BE USED WITH JAY J2 DEEP CONTOUR CUSHION.


NOTE: ITEMS 1" OR GREATER COMBINED WITH J2 DEEP CONTOUR CUSHION ARE RECOGNIZED BY THE SADMERC AS A K0108 ADJUSTABLE SKIN PROTECTION AND POSITIONING CUSHION.

NOTE: JAY GS ACCESSORY BOX REPLACED BY JAY BOX OF ACCESSORIES (P/N JBOXACCS).




| Pos. | Item Number | Description                | Remarks          |
|------|-------------|----------------------------|------------------|
| 1A   | F119        | FLUID SUPPLEMENT (8 OZ)    | Size 5" x 10.5"  |
| 1B   | F1119       | MINI FLUID SUPPLEMENT      | Size 3.5" x 7.5" |
| 2A   | B113J       | 2" ADDUCTOR WEDGES         | Pair, Height 2"  |
| 2B   | B1113       | WEDGES, 3" ADDUCTOR (PAIR) | Pair, Height 3"  |
| 3    | B118W       | 2" ABDUCTOR WEDGE          |                  |
| 4    | B118        | BLD UP, 2" ABDUCTOR, JAY   |                  |


| Pos.       | Item Number | Description                         | Remarks |
|------------|-------------|-------------------------------------|---------|
| 5          | B168        | BUILD-UP 2" JUNIOR ABDUCTOR         |         |
| 6          | B114        | BASE WEDGE 1 1/2"                   |         |
| 7          | 086         | FRONT PAINT 4 OZ + PAINT BRUSH      |         |
| 8          | S105        | SOLID SEAT JAY 15.5X16              |         |
| 8          | S101        | SOLID SEAT JAY 18 x 18              |         |
| 8          | S106        | SOLID SEAT JAY 15.5 x 18            |         |
| 8          | S108        | SOLID SEAT JAY 14 x 16              |         |
| 8          | S111        | SOLID SEAT JAY 10 x 10              |         |
| 8          | S115        | SOLID SEAT JAY 12 x 11              |         |
| 9          | B163        | 2" JUNIOR ADDUCTOR WEDGES           | Pair    |
| 10         | B1167       | HIP GUIDES,2",JAY & GS              | Pair    |
| 11         | B116R       | 1" PELVIC OBLIQUITY BLD UP          | Right** |
| 11         | B116L       | 1" PELVIC OBLIQUITY BLD UP          | Left**  |
| 12         | B119J       | 1" ADDUCTOR WEDGES -JAY             | Pair    |
| 13         | B2160L      | 1/2" PELVIC OBLIQUITY BUILD UP 2100 | Left**  |
| 13         | B2160R      | 1/2" PELVIC OBLIQUITY BUILD UP 2100 | Right** |
| 13         | B2161L      | 1/2" PELVIC OBLIQUITY BUILD UP 2101 | Left**  |
| 13         | B2161R      | 1/2" PELVIC OBLIQUITY BUILD UP 2101 | Right** |
| 13         | B2162L      | 1/2" PELVIC OBLIQUITY BUILD UP 2102 | Left**  |
| 13         | B2162R      | 1/2" PELVIC OBLIQUITY BUILD UP 2102 | Right** |
| 13         | B2163L      | 1/2" PELVIC OBLIQUITY BUILD UP 2103 | Left**  |
| 13         | B2163R      | 1/2" PELVIC OBLIQUITY BUILD UP 2103 | Right** |
| 13         | B2164L      | 1/2" PELVIC OBLIQUITY BUILD UP 2104 | Left**  |
| 13         | B2164R      | 1/2" PELVIC OBLIQUITY BUILD UP 2104 | Right** |
| 13         | B2165L      | 1/2" PELVIC OBLIQUITY BUILD UP 2105 | Left**  |
| 13         | B2165R      | 1/2" PELVIC OBLIQUITY BUILD UP 2105 | Right** |
| 13         | B2166L      | 1/2" PELVIC OBLIQUITY BUILD UP 2106 | Left**  |
| 13         | B2166R      | 1/2" PELVIC OBLIQUITY BUILD UP 2106 | Right** |
| 13         | B2167L      | 1/2" PELVIC OBLIQUITY BUILD UP 2107 | Left**  |
| 13         | B2167R      | 1/2" PELVIC OBLIQUITY BUILD UP 2107 | Right** |
| 13         | B2168L      | 1/2" PELVIC OBLIQUITY BUILD UP 2108 | Left**  |
| 13         | B2168R      | 1/2" PELVIC OBLIQUITY BUILD UP 2108 | Right** |
| 13         | B2169L      | 1/2" PELVIC OBLIQUITY BUILD UP 2116 | Left**  |
| 13         | B2169R      | 1/2" PELVIC OBLIQUITY BUILD UP 2116 | Right** |
| 13         | B2177L      | 1/2" PELVIC OBLIQUITY BUILD UP 2117 | Left**  |
| 13         | B2177R      | 1/2" PELVIC OBLIQUITY BUILD UP 2117 | Right** |
| 13         | B2170L      | 1/2" PELVIC OBLIQUITY BUILD UP 2118 | Left**  |
| 13         | B2170R      | 1/2" PELVIC OBLIQUITY BUILD UP 2118 | Right** |
| 14         | B1118       | 3" ABDUCTOR BUILD-UP                |         |
| 15         | K1100       | WELL BUILD-UP GS 18"                |         |
| 15         | K1105       | WELL BUILD-UP GS 16" 1              |         |
| 15         | K1110       | WELL BUILD-UP GS 14"                |         |
| 15         | K1115       | WELL BUILD-UP GS (12")              |         |
| 1-12,14,15 | JBOXACCS    | JAY BOX OF ACCESSORIES              |         |
| 6,10,11    | 156         | CHILDREN'S ACCESSORIES KIT          |         |



| Pos. | Item Number | Description                 |
|------|-------------|-----------------------------|
| 1    | DJ3SPNKA    | ASSESSMENT KIT, SPINE ALIGN |




| Pos.        | Item Number | Description              |
|-------------|-------------|--------------------------|
| 1,2A,3A,4-7 | J3CEKS      | J3 CUSH STD PLA EVAL KIT |
| 1,2B,3B,4-7 | J3CEKD      | J3 CUSH DP PLA EVAL KIT  |



| Pos. | Item Number | Description                    | Remarks        |
|------|-------------|--------------------------------|----------------|
| 2    | J3CALM1     | LATERAL THIGH 1" HEIGHT MEDIUM | 8.25" LONG     |
| 2    | J3CALL1     | LATERAL THIGH 1" HEIGHT LARGE  | 11" LONG       |
| 2    | J3CALS2     | LATERAL THIGH 2" HEIGHT SMALL  | 5.5" LONG      |
| 2    | J3CALM2     | LATERAL THIGH 2" HEIGHT MEDIUM | 8.25" LONG     |
| 2    | J3CALL2     | LATERAL THIGH 2" HEIGHT LARGE  | 11" LONG       |
| 3    | J3CMM1      | MED THIGH 1" HEIGHT MEDIUM     | 6" LONG        |
| 3    | J3CAML1     | MED THIGH 1" HEIGHT LARGE      | 8.5" LONG      |
| 3    | J3CAMS2     | MED THIGH 2" HEIGHT SMALL      | 3.5" LONG      |
| 3    | J3CMM2      | MED THIGH 2" HEIGHT MEDIUM     | 6" LONG        |
| 3    | J3CAML2     | MED THIGH 2" HEIGHT LARGE      | 8.5" LONG      |
| 4    | J3CAPS2     | LAT PEL 2" HEIGHT SMALL        | 4.75" LONG     |
| 4    | J3CAPM2     | LAT PEL 2" HEIGHT MEDIUM       | 6.25" LONG     |
| 5    | J3CSAOR     | OBLIQUITY PAD STD A RIGHT      | 1/2" in height |
| 5    | J3CSAOL     | OBLIQUITY PAD STD A LEFT       | 1/2" in height |
| 5    | J3CSBOR     | OBLIQUITY PAD STD B RIGHT      | 1/2" in height |
| 5    | J3CSBOL     | OBLIQUITY PAD STD B LEFT       | 1/2" in height |
| 5    | J3CSCOR     | OBLIQUITY PAD STD C RIGHT      | 1/2" in height |
| 5    | J3CSCOL     | OBLIQUITY PAD STD C LEFT       | 1/2" in height |
| 5    | J3CDAOR     | OBLIQUITY PAD DP A RIGHT       | 1/2" in height |
| 5    | J3CDAOL     | OBLIQUITY PAD DP A LEFT        | 1/2" in height |
| 5    | J3CDBOR     | OBLIQUITY PAD DP B RIGHT       | 1/2" in height |

| Pos. | Item Number | Description                     | Remarks        |
|------|-------------|---------------------------------|----------------|
| 5    | J3CDBOL     | OBLIQUITY PAD DP B LEFT         | 1/2" in height |
| 5    | J3CDCOR     | OBLIQUITY PAD DP C RIGHT        | 1/2" in height |
| 5    | J3CDCOL     | OBLIQUITY PAD DP C LEFT         | 1/2" in height |
| 6    | J3CPRT      | REPLACEMENT TAPE 30' ROLL       |                |
| 1-6  | J3CPEK      | J3 CUSHION POSITIONING EVAL KIT |                |


Note: Custom Drop Seat is only available for the Jay Soft Combi & Soft Combi P, Jay Gel, Jay Easy, Jay 2 & J2 DC, Jay GS and Jay Triad Cushions.



| Pos. | Item Number | Description                   | Remarks |
|------|-------------|-------------------------------|---------|
| 1    | SD2444      | SOLID DROP SEAT 14X14 W/HRDWR |         |
| 1    | SD2446      | SOLID DROP SEAT 14X16 W/HRDWR |         |
| 1    | SD2455      | SOLID DROP SEAT 15X15 W/HRDWR |         |
| 1    | SD2457      | SOLID DROP SEAT 15X17 W/HRDWR |         |
| 1    | SD2460      | SOLID DROP SEAT 16X20 W/HRDWR |         |
| 1    | SD2466      | SOLID DROP SEAT 16X16 W/HRDWR |         |
| 1    | SD2468      | SOLID DROP SEAT 16X18 W/HRDWR |         |
| 1    | SD2475      | SOLID DROP SEAT 17X15 W/HRDWR |         |
| 1    | SD2477      | SOLID DROP SEAT 17X17 W/HRDWR |         |
| 1    | SD2480      | SOLID DROP SEAT 18X20 W/HRDWR |         |
| 1    | SD2488      | SOLID DROP SEAT 18X18 W/HRDWR |         |
| 1    | SD2486      | SOLID DROP SEAT 18X16 W/HRDWR |         |
| 1    | SD2406      | SOLID DROP SEAT 20X16 W/HRDWR |         |
| 1    | SD2408      | SOLID DROP SEAT 20X18 W/HRDWR |         |

| Pos. | Item Number | Description                      | Remarks   |
|------|-------------|----------------------------------|-----------|
| 1    | SD2400      | SOLID DROP SEAT 20X20 W/HRDWR    |           |
| 1    | SD2426      | SOLID DROP SEAT 22X16 W/HRDWR    |           |
| 1    | SD2428      | SOLID DROP SEAT 22X18 W/HRDWR    |           |
| 1    | SD2448      | SOLID DROP SEAT 24X18 W/HRDWR    |           |
| 2-4  | M-24S       | SHORT 1/2" DROP HOOKS            | Pack of 4 |
| 2-4  | M-24M       | MEDIUM 1" DROP HOOKS             | Pack of 4 |
| 2-4  | M-24L       | LONG 1 1/2" DROP HOOKS           | Pack of 4 |
| 5    | 283         | HOOK STOPPER W/HARDWARE          |           |
| 6    | 282         | HOOK RETAINER/STOP WITH HARDWARE |           |

NOTE: ITEM 11 REQUIRES DRILLING OF (2) 7/32" HOLES FOR PROPER MOUNTING AND IS USED TO KEEP THE CUSHION FROM SLIDING OFF THE SEAT.



| Pos. | Item Number | Description                                | Remarks                   |
|------|-------------|--------------------------------------------|---------------------------|
| 2    | M-271       | 1/4"X20X1" BOLT                            |                           |
| 3    | M-272       | WHIZ NUT                                   |                           |
| 7    | 283         | HOOK STOPPER W/HARDWARE                    |                           |
| 8    | 284         | HOOK RETAINER W/HARDWARE                   |                           |
| 9    | 287         | HOOK&BRACKET MOD. 1 3/4"DROP               |                           |
| 9    | 298         | HOOK & BRACKET MODFID MAX DROP             |                           |
| 10   | 286-MJ      | HOOK & BRACKET MOD. FLUSH MOUNT            |                           |
| 10   | 296         | HOOK & "L" ASSM FOR 20" FLUSH              |                           |
| 11   | M-285       | LIP KIT - FOR ADJUSTABLE SOLID SEAT        |                           |
| 2-6  | 2HL         | HOOK & BRACKET - FASTENED                  |                           |
| 2-6  | 3HL         | LRG AD HOOK & BRACKET-FASTENED             |                           |
| 2-6  | 4HL         | HOOK & "L"-FASTENED FOR 10"                |                           |
| 2-8  | 285         | JAY SOLID SEAT HARDWARE PACKET             |                           |
| 1-8  | 270         | SOLID SEAT ADJUSTABLE 15-18 X 12-18        | INCLUDES HARDWARE P/N 2HL |
| 1-8  | 271         | SOLID SEAT ADJUSTABLE 10/11 X 10-12        | INCLUDES HARDWARE P/N 4HL |
| 1-8  | 272         | SOLID SEAT ADJUSTABLE 11.5-14.5 X 10-16    | INCLUDES HARDWARE P/N 2HL |
| 1-8  | 270.1       | SOLID SEAT ADJUSTABLE 15-18 W/2.5" DROP OP | INCLUDES HARDWARE P/N 3HL |

| Pos.   | Item Number | Description                                   | Remarks                   |
|--------|-------------|-----------------------------------------------|---------------------------|
| 1-8    | 271.1       | SOLID SEAT ADJUSTABLE 10-11 W/2.5" DROP OP    | INCLUDES HARDWARE P/N 6HL |
| 1-8    | 272.1       | SOLID SEAT ADJUSTABLE 11.5-14.5 W/2.5" DRP OP | INCLUDES HARDWARE P/N 3HL |
| 1-8    | 273.1       | SOLID SEAT ADJUSTABLE 18.5-22 X 12-20         | INCLUDES HARDWARE P/N 3HL |
| 1-8,11 | 270L        | SOLID SEAT ADJUSTABLE W/LIP KIT               |                           |
| 1-8,11 | 271L        | SOLID SEAT ADJUSTABLE 10" W/LIP KIT           |                           |
| 1-8,11 | 272L        | SOLID SEAT ADJUSTABLE W/LIP KIT               |                           |
| 1-8,11 | 270.1L      | SOLID SEAT ADJUSTABLE 15-18 W/2.5" DRP,LIP KT |                           |
| 1-8,11 | 271.1L      | SOLID SEAT ADJUSTABLE 10-11 W/2.5" DRP,LIP KT |                           |
| 1-8,11 | 272.1L      | SOLID SEAT ADJUSTABLE 11.5-14.5 /2.5 DRP,LIP  |                           |
| 1-8,11 | 273.1L      | SOLID SEAT ADJUSTABLE LARGE ADULT W/LIP KIT   |                           |



| Pos.  | Item Number | Description                      |
|-------|-------------|----------------------------------|
| 1     | C141        | COVER PROTECTOR CUSHION          |
| 1A    | C144        | JUNIOR COVER PROTECTOR CUSHION   |
| 2     | 141         | PROTECTOR CUSHION                |
| 2A    | 144         | JUNIOR PROTECTOR CUSHION         |
| 3     | 142         | PROTECTOR SLING                  |
| 3A    | 145         | JUNIOR PROTECTOR SLING           |
| 1-3   | 140         | PROTECTOR CUSHION AND SLING      |
| 1A-3A | 143-MJ      | JUNIOR PROTECTOR CUSHION & SLING |

Note: Spinal Fluid PAdS can be added to any Jay back with a hard shell.

Note: The Spinal Fluid Pad measures 12" tall but the pad can be cut down to accommodate shorter backs by cutting inbetween the cells to remove cells that are not needed



| PART # | HEIGHT |
|--------|--------|
| XF1505 | 20"    |
| XF506  | 23"    |

| Pos. | Item Number | Description             |
|------|-------------|-------------------------|
| 1    | MF506       | FLUID PAD TALL JAY BACK |



#### FOOT BOX WIDTH MATRIX

##### Flared Frame Formula:

(Seat Width) + 2" = Foot Box Width

##### Offset Hanger Receiver Formula:

(Seat Width) + 2" = Foot Box Width

##### Inset Formulas:

0" Inset (Seat Width) = Foot Box Width

1" Inset (Seat Width) - 2" = Foot Box Width

2" Inset (Seat Width) - 4" = Foot Box Width

| Pos. | Item Number | Description                    | Remarks              |
|------|-------------|--------------------------------|----------------------|
| 1    | CCFOOTS22   | FULL FOOTBOX 22W FOOT COVER    |                      |
| 1    | CCFOOTS20   | FULL FOOTBOX 20W FOOT COVER    |                      |
| 1    | CCFOOTS18   | FULL FOOTBOX 18W FOOT COVER    |                      |
| 1    | CCFOOTS16   | FULL FOOTBOX 16W FOOT COVER    |                      |
| 2    | CCFOOTS22   | FUL FOOTBX 22W SEP CALF CVR SD | Cover only           |
| 2    | CCFOOTS20   | FUL FOOTBX 20W SEP CALF CVR SD | Cover only           |
| 2    | CCFOOTS18   | FUL FOOTBX 18W SEP CALF CVR SD | Cover only           |
| 2    | CCFOOTS16   | FUL FOOTBX 16W SEP CALF CVR SD | Cover only           |
| 2    | CCFOOTS22   | FUL FOOTBX 22W SEP CALF CVR TL | Cover only           |
| 2    | CCFOOTS20   | FUL FOOTBX 20W SEP CALF CVR TL | Cover only           |
| 2    | CCFOOTS18   | FUL FOOTBX 18W SEP CALF CVR TL | Cover only           |
| 2    | CCFOOTS16   | FUL FOOTBX 16W SEP CALF CVR TL | Cover only           |
| 3,4  | HDWFOOT     | DUAL FOOTBOX HRDWR KIT         | 4 screws and washers |
| 5    | 132826      | JAY FOOTBOARD ASM 12W KIT      | Cover, foam & board  |
| 5    | 132827      | JAY FOOTBOARD ASM 14W KIT      | Cover, foam & board  |
| 5    | 132828      | JAY FOOTBOARD ASM 16W KIT      | Cover, foam & board  |
| 5    | 132829      | JAY FOOTBOARD ASM 18W KIT      | Cover, foam & board  |
| 5    | 132830      | JAY FOOTBOARD ASM 20W KIT      | Cover, foam & board  |

| Pos. | Item Number | Description               | Remarks             |
|------|-------------|---------------------------|---------------------|
| 5    | 132831      | JAY FOOTBOARD ASM 22W KIT | Cover, foam & board |




| Pos. | Item Number | Description                 | Remarks              |
|------|-------------|-----------------------------|----------------------|
| 1    | CCFOOT22    | FUL FOOTBX 22W CALF CVR STD |                      |
| 1    | CCFOOT20    | FUL FOOTBX 20W CALF CVR STD |                      |
| 1    | CCFOOT18    | FUL FOOTBX 18W CALF CVR STD |                      |
| 1    | CCFOOT16    | FUL FOOTBX 16W CALF CVR STD |                      |
| 1    | CCFOOT14    | FUL FOOTBX 14W CALF CVR STD |                      |
| 1    | CCFOOT12    | FUL FOOTBX 12W CALF CVR STD |                      |
| 1    | CCFOOTT22   | FUL FOOTBX 22W CALF CVR TAL |                      |
| 1    | CCFOOTT20   | FUL FOOTBX 20W CALF CVR TAL |                      |
| 1    | CCFOOTT18   | FUL FOOTBX 18W CALF CVR TAL |                      |
| 1    | CCFOOTT16   | FUL FOOTBX 16W CALF CVR TAL |                      |
| 2    | CCFOOTS22   | FULL FOOTBOX 22W FOOT COVER |                      |
| 2    | CCFOOTS20   | FULL FOOTBOX 20W FOOT COVER |                      |
| 2    | CCFOOTS18   | FULL FOOTBOX 18W FOOT COVER |                      |
| 2    | CCFOOTS16   | FULL FOOTBOX 16W FOOT COVER |                      |
| 2    | CCFOOTS14   | FULL FOOTBOX 14W FOOT COVER |                      |
| 2    | CCFOOTS12   | FULL FOOTBOX 12W FOOT COVER |                      |
| 3,4  | HDWFOOT     | DUAL FOOTBOX HRDWR KIT      | 4 screws and washers |



| FOOT BOX WIDTH MATRIX                  |                  |
|----------------------------------------|------------------|
| <b>Flared Frame Formula:</b>           |                  |
| (Seat Width) + 2"                      | = Foot Box Width |
| <b>Offset Hanger Receiver Formula:</b> |                  |
| (Seat Width) + 2"                      | = Foot Box Width |
| <b>Inset Formulas:</b>                 |                  |
| 0" Inset (Seat Width)                  | = Foot Box Width |
| 1" Inset (Seat Width) - 2"             | = Foot Box Width |
| 2" Inset (Seat Width) - 4"             | = Foot Box Width |


| Pos. | Item Number | Description                  | Remarks                |
|------|-------------|------------------------------|------------------------|
| 1    | CCFOOTSS    | SINGLE FOOTBOX CALF CVR STD  | Cover only             |
| 1    | CCFOOTST    | SINGLE FOOTBOX CALF CVR TAL  | Cover only             |
| 2    | FCFOOTS     | SINGLE FOOTBOX FOOT COVER    | Cover only             |
| 3,4  | HDWSFOOT    | SINGLE FOOTBOX HRDWR KIT     | 2 screws and 2 washers |
| 5    | 132824      | JAY FOOTBOARD ASM SINGLE STD | Cover, foam & board    |



| Pos. | Item Number   | Description                  |
|------|---------------|------------------------------|
| 1    | CJARM14       | ARM TROUGH 14" COVER         |
| 1    | CJARM18       | ARM TROUGH 18" COVER         |
| 1    | CJARM22       | ARM TROUGH 22" COVER         |
| 2    | JARMSTRAP14   | JAY ARMTROUGH 14IN CVR STRAP |
| 2    | JARMSTRAP1822 | JAY ARMTRGH 1418IN CVR STRAP |
| 5    | JARMLEVER     | ARM TROUGH ADJ CAM LEVER     |
| 3-9  | SAARMHRDW     | ARM TROUGH S/A HARDWARE      |
| 7-11 | FXARMHRDW     | ARM TROUGH FIXED HARDWARE    |



## MEASURING GUIDE



| Part Number | Dim A (Width) | Dim B (Height) | Dim C (Thickness) | Dim D (Right Side Strap*) | Dim E (Left Side Strap*) |
|-------------|---------------|----------------|-------------------|---------------------------|--------------------------|
| 2250        | 15"           | 8.75"          | 2"                | 7"                        | 7"                       |
| 2251-MJ     | 11"           | 7"             | 2"                | 6.5"                      | 15"                      |

\* Strap side is located as if you are resting against the back with the lumbar support on it.

| Pos. | Item Number | Description                   | Remarks |
|------|-------------|-------------------------------|---------|
| 1A   | 2250        | LUMBAR SUPPORT SOFT COMBI     | Adult   |
| 1B   | 2251-MJ     | LUMBAR SUPPORT SFT COMBI JUNR | Junior  |

# JAY<sup>®</sup>

